



ORDER FORM

Squiggles Seat

EILK03

Ultimate combination of postural seating and everyday comfort for young children while still maintaining a fun, colorful design.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Squiggles Seat

Weight Capacity: 48.4 lbs

<input checked="" type="checkbox"/> OPTION	DESCRIPTION	RETAIL
Seat Shell		\$
Selection Required.		
<input type="checkbox"/> 151FT72	Squiggles Seat Shell	\$2,082
	151BS28 Seat Base with Depth Adjustment	
	151BF38 Seat Back - Fixed	
	151BF39 Seat Back Reclining Mechanism	
	151GU5 Seat Base Hip Guides / Support	
Seat Cover		\$
Selection Required. Includes seat base and backrest cushion covers.		
<input type="checkbox"/> 151YJ14	Squiggles Seat Covers - Blue	\$540
<input type="checkbox"/> 151YJ15	Squiggles Seat Covers - Green	\$540
<input type="checkbox"/> 151YJ16	Squiggles Seat Covers - Orange	\$540
<input type="checkbox"/> 151YJ17	Squiggles Seat Covers - Pink	\$540
Base		\$
<input type="checkbox"/> 151NB116	Size 1 Hi Low Chassis (Foot Pedal Operated.)	\$2,082
	151PH35 Hi Low Chassis Push Handle	\$301
<input type="checkbox"/> 151NB114	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	\$2,265
	151JI04 Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis	\$210
<input type="checkbox"/> 151NB117	Mobility Base Interface Plate Allows mounting Squiggles Seat onto Sunrise Wheelchairs	\$644
Foot Positioning		\$
Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest		
<input type="checkbox"/> 151DI67	Squiggles Seat Adjustable Footrest	\$238
<input type="checkbox"/> 1518A107	Small Sandals w/ Straps (Pair)	\$269
<input type="checkbox"/> 1518A108	Medium Sandals w/ Straps (Pair)	\$287
<input type="checkbox"/> 1519A24	QTY Footplate Raiser 1" (Each)	\$96
Harness		\$
<input type="checkbox"/> 151NH109	Pelvic Harness Small - Blue	\$155
<input type="checkbox"/> 151NH110	Pelvic Harness Small - Green	\$155
<input type="checkbox"/> 151NH111	Pelvic Harness Small - Orange	\$155
<input type="checkbox"/> 151NH112	Pelvic Harness Small - Pink	\$155
<input type="checkbox"/> 151NH113	Pelvic Harness Medium - Blue	\$155
<input type="checkbox"/> 151NH114	Pelvic Harness Medium - Green	\$155
<input type="checkbox"/> 151NH115	Pelvic Harness Medium - Pink	\$155
<input type="checkbox"/> 151NH108	Pelvic Harness Spacer Pads - Black (Pair)	\$95
<input type="checkbox"/> 151NH105	Butterfly Harness Small	\$409
<input type="checkbox"/> 151NH106	Butterfly Harness Medium	\$409

Squiggles Seat

✓ OPTION	DESCRIPTION	RETAIL
Lateral Support		\$
Cover selection required with Lateral Support Hardware		
<input type="checkbox"/> 151AW26	Flip-Away Lateral Support Hardware (Pair)	\$338
<input type="checkbox"/> 151IU19	Lateral Support Covers - Blue (Pair)	\$155
<input type="checkbox"/> 151IU20	Lateral Support Covers - Green (Pair)	\$155
<input type="checkbox"/> 151IU21	Lateral Support Covers - Orange (Pair)	\$155
<input type="checkbox"/> 151IU22	Lateral Support Covers - Pink (Pair)	\$155
<input type="checkbox"/> 151AW27	Rigid Lateral Support Hardware (Pair)	\$207
<input type="checkbox"/> 151IU19	Lateral Support Covers - Blue (Pair)	\$155
<input type="checkbox"/> 151IU20	Lateral Support Covers - Green (Pair)	\$155
<input type="checkbox"/> 151IU21	Lateral Support Covers - Orange (Pair)	\$155
<input type="checkbox"/> 151IU22	Lateral Support Covers - Pink (Pair)	\$155
<input type="checkbox"/> 151KE54	Laterals w/ Protraction Pads - Black (Pair)	\$432
<input type="checkbox"/> 151ZY25	Chest Support - Black (Requires Laterals)	\$201
Tray		\$
<input type="checkbox"/> 151A230	Tray - Gray	\$387
<input type="checkbox"/> 151A231	Extra High Tray +2" Adjustment	\$558
<input type="checkbox"/> 151A216	Padded Tray Insert - Black	\$326
<input type="checkbox"/> 1518A106	Grab Rail / Toy Bar for Tray	\$321
Pommel		\$
<input type="checkbox"/> 1518A118	Pommel	\$119

Squiggles Seat

✓ OPTION	DESCRIPTION	RETAIL
Headrest		\$
Cushion selection required with hardware.		
Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.		
<input type="checkbox"/>	151AW24 Contoured Headrest Support Hardware	\$390
<input type="checkbox"/>	151Y040 Contoured Headrest Cushion - Blue	\$191
<input type="checkbox"/>	151Y042 Contoured Headrest Cushion - Green	\$191
<input type="checkbox"/>	151Y057 Contoured Headrest Cushion - Orange	\$191
<input type="checkbox"/>	151Y043 Contoured Headrest Cushion - Pink	\$191
<input type="checkbox"/>	151AW22 Flat Headrest Support Hardware	\$243
<input type="checkbox"/>	151Y048 Flat Headrest Cushion - Blue	\$191
<input type="checkbox"/>	151Y050 Flat Headrest Cushion - Green	\$191
<input type="checkbox"/>	151Y051 Flat Headrest Cushion - Pink	\$191
<input type="checkbox"/>	151AW23 Lateral Bracket Assembly Hardware (Pair) Requires Flat Headrest Support Hardware and Cushion	\$191
<input type="checkbox"/>	151Y048 Flat Headrest Lateral Cushions - Blue (Pair)	\$191
<input type="checkbox"/>	151Y050 Flat Headrest Lateral Cushions - Green (Pair)	\$191
<input type="checkbox"/>	151Y051 Flat Headrest Lateral Cushions - Pink (Pair)	\$191
<input type="checkbox"/>	151JG16 Removable Head Support Bracket for 1/2" Box Section	\$160
Accessories		\$
<input type="checkbox"/>	151KE55 1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	\$790
<input type="checkbox"/>	151KE56 1" Back Spacer Pad and Flip-Away Lateral Supports - Green	\$790
<input type="checkbox"/>	151KE57 1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	\$790
<input type="checkbox"/>	151AW28 Shoulder Support Hardware	\$304
<input type="checkbox"/>	151SS13 Shoulder Support Cover - Blue	\$174
<input type="checkbox"/>	151SS14 Shoulder Support Cover - Green	\$174
<input type="checkbox"/>	151SS15 Shoulder Support Cover - Pink	\$174

Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	1 - 5 years
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	10° / 25°
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	10° / 10°
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kgs
Tray Height	4.5" - 6.7"	115 mm x 170 mm
Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kgs
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, Ont., L4K 4P8



MK-129856-CAN © 7/2024

