



ORDER FORM

February 2024

Squiggles Seat

Ultimate combination of postural seating and everyday comfort for young children while still maintaining a fun, colorful design.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Squiggles Seat

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	RETAIL
Seat Shell		\$
Selection Required.		
<input type="checkbox"/> 120-610S	Squiggles Seat Shell	\$2,082
120-610-1	Seat Base with Depth Adjustment	
120-610-2	Seat Back - Fixed	
120-610-3	Seat Back Reclining Mechanism	
120-610-4	Seat Base Hip Guides / Support	
Seat Cover		\$
Selection Required. Includes seat base and backrest cushion covers.		
<input type="checkbox"/> 120-710-03	Squiggles Seat Covers - Blue	\$540
<input type="checkbox"/> 120-710-01	Squiggles Seat Covers - Green	\$540
<input type="checkbox"/> 120-710-02	Squiggles Seat Covers - Orange	\$540
<input type="checkbox"/> 120-710-04	Squiggles Seat Covers - Pink	\$540
Base		\$
<input type="checkbox"/> 120-800	Size 1 Hi Low Chassis (Foot Pedal Operated.)	\$2,082
122-710	Hi Low Chassis Push Handle	\$301
<input type="checkbox"/> 102-2753	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	\$2,265
120-818	Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis	\$210
<input type="checkbox"/> 194-600	Mobility Base Interface Plate Allows mounting Squiggles Seat onto Sunrise Wheelchairs	\$644
Foot Positioning		\$
Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest		
<input type="checkbox"/> 120-639	Squiggles Seat Adjustable Footrest	\$238
<input type="checkbox"/> 152-1600	Small Sandals w/ Straps (Pair)	\$269
<input type="checkbox"/> 152-2600	Medium Sandals w/ Straps (Pair)	\$287
<input type="checkbox"/> 120-816	QTY Footplate Raiser 1" (Each)	\$96
Harness		\$
<input type="checkbox"/> 120-844-03	Pelvic Harness Small - Blue	\$155
<input type="checkbox"/> 120-844-01	Pelvic Harness Small - Green	\$155
<input type="checkbox"/> 120-844-02	Pelvic Harness Small - Orange	\$155
<input type="checkbox"/> 120-844-04	Pelvic Harness Small - Pink	\$155
<input type="checkbox"/> 117-848-03	Pelvic Harness Medium - Blue	\$155
<input type="checkbox"/> 117-848-01	Pelvic Harness Medium - Green	\$155
<input type="checkbox"/> 117-848-04	Pelvic Harness Medium - Pink	\$155
<input type="checkbox"/> 120-789	Pelvic Harness Spacer Pads - Black (Pair)	\$95
<input type="checkbox"/> 155-N001-L600	Butterfly Harness Small	\$409
<input type="checkbox"/> 155-N001-L601	Butterfly Harness Medium	\$409

Squiggles Seat

✓ OPTION	DESCRIPTION	RETAIL
Lateral Support		\$
Cover selection required with Lateral Support Hardware		
<input type="checkbox"/> 120-855	Flip-Away Lateral Support Hardware (Pair)	\$338
<input type="checkbox"/> 120-838-03	Lateral Support Covers - Blue (Pair)	\$155
<input type="checkbox"/> 120-838-01	Lateral Support Covers - Green (Pair)	\$155
<input type="checkbox"/> 120-838-02	Lateral Support Covers - Orange (Pair)	\$155
<input type="checkbox"/> 120-838-04	Lateral Support Covers - Pink (Pair)	\$155
<input type="checkbox"/> 120-656	Rigid Lateral Support Hardware (Pair)	\$207
<input type="checkbox"/> 120-838-03	Lateral Support Covers - Blue (Pair)	\$155
<input type="checkbox"/> 120-838-01	Lateral Support Covers - Green (Pair)	\$155
<input type="checkbox"/> 120-838-02	Lateral Support Covers - Orange (Pair)	\$155
<input type="checkbox"/> 120-838-04	Lateral Support Covers - Pink (Pair)	\$155
<input type="checkbox"/> 120-954	Laterals w/ Protraction Pads - Black (Pair)	\$432
<input type="checkbox"/> 137-808-03	Chest Support - Blue (Requires Laterals)	\$201
<input type="checkbox"/> 137-808-01	Chest Support - Green (Requires Laterals)	\$201
<input type="checkbox"/> 137-808-02	Chest Support - Orange (Requires Laterals)	\$201
<input type="checkbox"/> 137-808-04	Chest Support - Pink (Requires Laterals)	\$201
Tray		\$
<input type="checkbox"/> 120-672	Tray - Gray	\$387
<input type="checkbox"/> 120-688	Extra High Tray +2" Adjustment	\$558
<input type="checkbox"/> 188-1854-05	Padded Tray Insert - Black	\$326
<input type="checkbox"/> 117-769	Grab Rail / Toy Bar for Tray	\$321
Pommel		\$
<input type="checkbox"/> 120-641	Pommel	\$119

Squiggles Stander

✓ OPTION	DESCRIPTION	RETAIL	
Headrest		\$	
Cushion selection required with hardware.			
Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.			
<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware	\$390
<input type="checkbox"/>	137-800-03	Contoured Headrest Cushion - Blue	\$191
<input type="checkbox"/>	137-800-01	Contoured Headrest Cushion - Green	\$191
<input type="checkbox"/>	137-800-02	Contoured Headrest Cushion - Orange	\$191
<input type="checkbox"/>	137-800-04	Contoured Headrest Cushion - Pink	\$191
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	\$243
<input type="checkbox"/>	117-850-03	Flat Headrest Cushion - Blue	\$122
<input type="checkbox"/>	117-850-01	Flat Headrest Cushion - Green	\$122
<input type="checkbox"/>	117-850-04	Flat Headrest Cushion - Pink	\$122
<input type="checkbox"/>	118-667	Lateral Bracket Assembly Hardware (Pair) Requires Flat Headrest Support Hardware and Cushion	\$191
<input type="checkbox"/>	137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	\$191
<input type="checkbox"/>	137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	\$191
<input type="checkbox"/>	137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	\$191
<input type="checkbox"/>	117-745	Removable Head Support Bracket for 1/2" Box Section	\$117
Accessories		\$	
<input type="checkbox"/>	120-779-03	1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	\$790
<input type="checkbox"/>	120-779-01	1" Back Spacer Pad and Flip-Away Lateral Supports - Green	\$790
<input type="checkbox"/>	120-779-04	1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	\$790
<input type="checkbox"/>	120-702	Shoulder Support Hardware	\$304
<input type="checkbox"/>	120-828-03	Shoulder Support Cover - Blue	\$174
<input type="checkbox"/>	120-828-01	Shoulder Support Cover - Green	\$174
<input type="checkbox"/>	120-828-04	Shoulder Support Cover - Pink	\$174

Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	1 - 5 years
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	10° / 25°
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	10° / 10°
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kgs
Tray Height	4.5" - 6.7"	115 mm x 170 mm
Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kgs
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm



Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada

Customer Service: (800) 263-3390 E-mail: cscanada@sunmed.com www.sunrisemedical.ca

Specifications subject to change without notice. © 2024 Sunrise Medical Canada Inc. MK-129856-CAN Rev 8

