



ORDER FORM

Squiggles Seat

EILK03

Ultimate combination of postural seating and everyday comfort for young children while still maintaining a fun, colorful design.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Squiggles Seat

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	RETAIL
Seat Shell		\$
Selection Required.		
<input type="checkbox"/>	151FT72 Squiggles Seat Shell	\$2,082
	151BS28 Seat Base with Depth Adjustment	
	151BF38 Seat Back - Fixed	
	151BF39 Seat Back Reclining Mechanism	
	151GU5 Seat Base Hip Guides / Support	
Seat Cover		\$
Selection Required. Includes seat base and backrest cushion covers.		
<input type="checkbox"/>	151YJ14 Squiggles Seat Covers - Blue	\$540
<input type="checkbox"/>	151YJ15 Squiggles Seat Covers - Green	\$540
<input type="checkbox"/>	151YJ17 Squiggles Seat Covers - Pink	\$540
Base		\$
<input type="checkbox"/>	151NB116 Size 1 Hi Low Chassis (Foot Pedal Operated.)	\$2,082
	151PH35 Hi Low Chassis Push Handle	\$301
<input type="checkbox"/>	151NB114 Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	\$2,265
	151JI04 Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis	\$210
<input type="checkbox"/>	151NB117 Mobility Base Interface Plate Allows mounting Squiggles Seat onto Sunrise Wheelchairs	\$644
Foot Positioning		\$
Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest		
<input type="checkbox"/>	151DI67 Squiggles Seat Adjustable Footrest	\$238
<input type="checkbox"/>	1518A107 Small Sandals w/ Straps (Pair)	\$269
<input type="checkbox"/>	1518A108 Medium Sandals w/ Straps (Pair)	\$287
<input type="checkbox"/>	1519A24 QTY Footplate Raiser 1" (Each)	\$96
Harness		\$
<input type="checkbox"/>	151NH109 Pelvic Harness Small - Blue	\$155
<input type="checkbox"/>	151NH110 Pelvic Harness Small - Green	\$155
<input type="checkbox"/>	151NH111 Pelvic Harness Small - Pink	\$155
<input type="checkbox"/>	151NH113 Pelvic Harness Medium - Blue	\$155
<input type="checkbox"/>	151NH114 Pelvic Harness Medium - Green	\$155
<input type="checkbox"/>	151NH115 Pelvic Harness Medium - Pink	\$155
<input type="checkbox"/>	151NH108 Pelvic Harness Spacer Pads - Black (Pair)	\$95
<input type="checkbox"/>	151NH105 Butterfly Harness Small	\$409
<input type="checkbox"/>	151NH106 Butterfly Harness Medium	\$409

Squiggles Seat

✓ OPTION	DESCRIPTION	RETAIL
Lateral Support		\$
Cover selection required with Lateral Support Hardware		
<input type="checkbox"/> 151AW26	Flip-Away Lateral Support Hardware (Pair)	\$338
<input type="checkbox"/>	151IU19 Lateral Support Covers - Blue (Pair)	\$155
<input type="checkbox"/>	151IU20 Lateral Support Covers - Green (Pair)	\$155
<input type="checkbox"/>	151IU22 Lateral Support Covers - Pink (Pair)	\$155
<input type="checkbox"/> 151AW27	Rigid Lateral Support Hardware (Pair)	\$207
<input type="checkbox"/>	151IU19 Lateral Support Covers - Blue (Pair)	\$155
<input type="checkbox"/>	151IU20 Lateral Support Covers - Green (Pair)	\$155
<input type="checkbox"/>	151IU22 Lateral Support Covers - Pink (Pair)	\$155
<input type="checkbox"/> 151KE54	Laterals w/ Protraction Pads - Black (Pair)	\$432
<input type="checkbox"/>	151ZY25 Chest Support - Black (Requires Laterals)	\$201
Tray		\$
<input type="checkbox"/>	151A230 Tray - Gray	\$387
<input type="checkbox"/>	151A231 Extra High Tray +2" Adjustment	\$558
<input type="checkbox"/>	151A216 Padded Tray Insert - Black	\$326
<input type="checkbox"/>	1518A106 Grab Rail / Toy Bar for Tray	\$321
Pommel		\$
<input type="checkbox"/>	1518A118 Pommel	\$119

Squiggles Seat

✓ OPTION	DESCRIPTION	RETAIL	
Headrest		\$	
Cushion selection required with hardware.			
Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.			
<input type="checkbox"/>	151AW24	Contoured Headrest Support Hardware	\$390
<input type="checkbox"/>	151Y040	Contoured Headrest Cushion - Blue	\$191
<input type="checkbox"/>	151Y042	Contoured Headrest Cushion - Green	\$191
<input type="checkbox"/>	151Y043	Contoured Headrest Cushion - Pink	\$191
<input type="checkbox"/>	151AW22	Flat Headrest Support Hardware	\$243
<input type="checkbox"/>	151Y048	Flat Headrest Cushion - Blue	\$191
<input type="checkbox"/>	151Y050	Flat Headrest Cushion - Green	\$191
<input type="checkbox"/>	151Y051	Flat Headrest Cushion - Pink	\$191
<input type="checkbox"/>	151AW23	Lateral Bracket Assembly Hardware (Pair) Requires Flat Headrest Support Hardware and Cushion	\$191
<input type="checkbox"/>	151Y048	Flat Headrest Lateral Cushions - Blue (Pair)	\$191
<input type="checkbox"/>	151Y050	Flat Headrest Lateral Cushions - Green (Pair)	\$191
<input type="checkbox"/>	151Y051	Flat Headrest Lateral Cushions - Pink (Pair)	\$191
<input type="checkbox"/>	151JG16	Removable Head Support Bracket for 1/2" Box Section	\$160
Accessories		\$	
<input type="checkbox"/>	151KE55	1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	\$790
<input type="checkbox"/>	151KE56	1" Back Spacer Pad and Flip-Away Lateral Supports - Green	\$790
<input type="checkbox"/>	151KE57	1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	\$790
<input type="checkbox"/>	151AW28	Shoulder Support Hardware	\$304
<input type="checkbox"/>	151SS13	Shoulder Support Cover - Blue	\$174
<input type="checkbox"/>	151SS14	Shoulder Support Cover - Green	\$174
<input type="checkbox"/>	151SS15	Shoulder Support Cover - Pink	\$174

Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	1 - 5 years
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	10° / 25°
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	10° / 10°
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kgs
Tray Height	4.5" - 6.7"	115 mm x 170 mm
Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kgs
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, Ont., L4K 4P8



MK-129856-CAN Rev 2 © 12/2024

