


VEYPR SUB



ORDER FORM

Customer Dealer Information		Shipping Address ☐ Identical To Billing Address	
Dealer		Client	
Client		Address	
Address			
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Information	
Purchase Order No.:		Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):		User Disability:	
		User Weight:	kg User Height
☐ Build4me			

RGK VEYPR Sub4

SBCXZ010

Veypr Sub4

\$22,647.00

Max User Weight 110kg (242lbs)

Please Note: All RGK chairs are hand built and have a manufacturing measurement tolerance of +/- 5mm

Product Style



XA070

Open Design

No Cost

CAD Approval



XW009

Quote only (No CAD For Approval)



XW010

Order CAD Drawing For Approval
(2 Customer Amendments ONLY)

\$300.00



XW011

Order (No CAD Required)



XW012

Quote & CAD For Approval
(2 Customer Amendments ONLY)

\$550.00

Key Measurements

XY010

Seat Width (Inside of sideguards)

mm (Min - 250mm / Max - 440mm)

XY020

Seat Length

mm (Min - 250mm / Max - 440mm)



XY030

Extra Seat Length (Provides additional leg support and hand grip)

mm (Max - 440mm Total Seat)

No Cost

XY050

Front Frame Height

mm (Min - 380mm / Max - 520mm)

XY060

Rear Frame Height

mm (Min - 380mm / Max - 500mm)

XY070

Frame Length

mm (Min = XY020 + XY030 +140mm / Max - 625mm)

XY072

Castor Position

mm



Optimise Measurements





XA020

Ergonomic Seat

mm Specify Length

\$125.00



Front Frame




XY080

Seat Taper Design

mm

(When combined with XY091, dimension is inbetween the centre of the main bends)



No Cost


XY091

Parallel Taper Design

mm

(When combined with XY080, dimension must be smaller by 20mm)


No Cost


Footrest Measurements



XY100

Round Design

(Castor wheel size 4", max height achievable 130mm)


Std

XY115

Front Footrest Height

mm

XY120


Rear Footrest Height

mm STD = 20mm lower than front height

XY125

Footrest Length

mm Maximum length = 120mm (certified weight)





Footrest Style



SBCXB400

Carbon Fibre Footrest with Carbon Plate

(Height adjustable +/- 10mm)

STD


Castor Wheels

☐	XC135	Aluminium Centre 4" x 1.25"	Std
☐		Silver Hub / Black Tyre (Std)	
☐		Black Hub / Black Tyre	
☐		Purple Hub / Black Tyre	
☐		Green Hub / Black Tyre	
☐		Red Hub / Black Tyre	
☐		Blue Hub / Black Tyre	
☐	SBCXC155	Forged Carbon Centre 4" x 1.4" ④	No Cost
☐	XC140	Aluminium Centre 5" x 1.4" (Silver Hub/Black Tyre) ④	\$239.00
☐	XC145	Aluminium Centre 6" x 1.4" (Silver Hub/Black Tyre) ④	\$319.00

Castor Forks

☐	XC190	Performance castor fork (Std)	Std
Choose colour below:			
☐	Black	☐	Red
☐	Silver	☐	Blue
☐	Purple		
☐	XC192	High Trail Castor Fork	No Cost
Note: Greater distance between contact patch and steering axis reduces flutter at higher speeds. Colour follows accent			

Castor Pin Width

XY330

Castor Pin Width




Optimised Measurements
(Minimum width for clear castor rotation and maximum wheelbase length)








Wheel Lock Options

☐	SBCXD036	Compact Scissor with Lightweight Fixing	No Cost
☐	XD010	Omit Wheel Locks	No Cost


Sideguards

☐	XE055	Screw on Lightweight Carbon Level with Tyre Individual Above Tyre (Follow Tyre - Max 25mm)	
☐	XE065	Screw on Lightweight Carbon with Fender (Fender centre of tyre) 10mm Above Tyre (Std) Individual Above Tyre (Min 10mm)	 

Backrest

☒	YG010	Fixed Backrest (Non Folding)		Std
	XY230	Backrest Height (Max 400mm)	mm Max 400mm	
	XY360	Backrest Angle	☐ -5 ☐ -3 ☐ 0 ☐ +3 ☐ +5 Individual (Max -10/+20 Deg) Deg	
☐	YG100	Ergonomic Backrest		\$150.00
		Vertical Height	mm	
		Rearward Position	mm	
☐	YG110	Inward Taper (Inside Of Backrest Posts - Specify Width Or Per Side)	mm	\$110.00
				
☐	YG120	Outward Flare (Inside Of Backrest Posts - Specify Width Or Per Side)	mm	\$150.00
				

Camber Bar

	XY390	Camber Angle	☐ 1 ☐ 2 ☐ 3	
☐	SBCXH100	Single Post Fixed Centre of Gravity (25-125mm)	mm	

Push Handles

☐	YG080	Flip Down Push Handles 	\$454.00
--------------------------	-------	--	----------

Anti-Tip

☐	XJ100	Single Sport Anti-Tip 	\$519.00
--------------------------	-------	---	----------



Frame Finish

- | | | | |
|--------------------------|----------|---------------------------|----------|
| ☐ | SBCXN300 | Carbon Fibre Matte Finish | No Cost |
| ☐ | SBCXN305 | Carbon Fibre Gloss Finish | \$922.00 |

Accent Colour

Accent colour includes decals only

- | | | | |
|--------------------------------|--------------------------------|---------------------------------|---------------------------------|
| ☐ Black | ☐ White | ☐ Purple | ☐ Yellow |
| ☐ Red | ☐ Blue | ☐ Orange | ☐ None |
| ☐ Pink | ☐ Green | ☐ Grey | |

Seating Options

- | | | | |
|--------------------------|-------|--|----------|
| ☐ | XP045 | Sub4 Seat Sling Strap Adjustable | Std |
| ☐ | XR050 | Carbon Fibre Seat Plate (With velcro to attach to seat sling - Not fixed to frame)  | \$494.00 |
| | | mm (STD XY020 - 20mm) | |

Sub4 Backrest Upholstery Options

- | | | | |
|--------------------------|-------|--|-----|
| ☐ | XP035 | Airtech Sub4 backrest strap adjustable | Std |
|--------------------------|-------|--|-----|

Non Sub4 Backrest Upholstery Options

- | | | | |
|--------------------------|-------|---|----------|
| ☐ | XP050 | Airtech backrest strap adjustable  | No Cost |
| ☐ | XP070 | Backrest upgrade to leather corners (Black)  | \$125.00 |
| ☐ | XP110 | Backrest with zip pocket  | \$131.00 |
| ☐ | XP170 | Individual backrest embroidery (Please supply JPeg)  | \$247.00 |

Backrest Upholstery Stitching Colour

- | | |
|-------------------------------------|---------------------------------|
| ☐ Silver | ☐ Red |
| ☐ Black | ☐ Orange |
| ☐ Blue | ☐ Purple |
| ☐ Individual | |

RGK Logo Colour

- | | |
|--------------------------------------|---------------------------------|
| ☐ White (Std) | ☐ Red |
| ☐ Black | ☐ Orange |
| ☐ Blue | ☐ Purple |
| ☐ Individual | |

Rear Wheel Size

- ☐ XY420 24" Wheels ☐ XY430 25" Wheels
- XY460 Rear Wheel Gap - 15 ☐ 20 ☐ 25 ☐

Rear Wheel Style

- ☐ XK075 Spinergy BLXL 18 Blade Spoke (N/A in Silver Hub) Std
- ☐ SBCXK150 Spinergy CLX 18 Carbon Blade Spoke (N/A in Silver Hub) \$1,354.00
- ☐ SBCXK040 Supply Without Wheels Wheel Type Tyre Type NC

Spinergy Spoke Colour

- ☐ Black ☐ Red ☐ Purple
- ☐ White ☐ Blue ☐ Orange
- ☐ Yellow ☐ Pink ☐ Green
- ☐ XK180 Individual Spoke Design Inside Outside \$219.00

Spinergy Hub / Rim Colour

- XK190 ☐ Black Hub (Std) Std
- XK200 ☐ Silver Hub \$131.00
- XK201 ☐ Red Hubs & Rims \$421.00
- XK202 ☐ Purple Hubs & Rims \$421.00
- XK203 ☐ Blue Hubs & Rims \$421.00
- XK204 ☐ Gold Hubs & Rims \$555.00

Handrims Style

- ☐ XL010 Aluminium hard anodised 6 tab Std
- ☐ XL160 Ellipse 6 Tab \$650.00
- ☐ XL020 Tetra Grip 6 Tab \$412.00
- ☐ XL100 Para Grip 6 Tab \$412.00
- ☐ XL040 Titanium 6 Tab \$459.00
- ☐ XL035 Natural Fit LT \$765.00
- ☐ XL170 Carbolife Curve L (BLK) \$1,435.00
- ☐ XL150 Carbolife Gekko - ☐ Black ☐ Silver \$1,448.00

Tyre Style

- ☐ XM010 Kenda - ☐ Grey ☐ Blue ☐ Red (Grey standard if no option selected) Std
- ☐ XM060 Schwalbe One (Lightest Option) \$137.00
- ☐ XM020 Schwalbe Marathon Plus \$101.00
- ☐ XM040 SHOX Solid (26" N/A with Fusion Wheels) \$165.00
- ☐ XM130 Kenda Black \$29.00
- ☐ XM050 Kenda K885 mountain bike tyre (24",25" only) \$41.00

Inner Tube

- ☐ XY500 Schraeder (Car Valve) - Standard if no option selected No Cost
- ☐ XY490 Presta (Sports Valve) No Cost

Wheel Build Up Options

- ☐ XY480 Narrow Set Handrim
- ☐ XY485 Wide Set Handrim
- ☐ XL090 Remove Excess Mounting Tabs \$150.00

Wheel Accessories

- ☐ XK145 Black Gloss Spokeguards ☐ Fitted ☐ Supply Non - Fitted \$655.00

Cushion

☐	XP080	Cushion Foam	\$125.00
		☐ Soft (CMS040)	
		☐ Medium (CMS050)	
		☐ 50mm	
		☐ 75mm	
		☐ 100mm	

Cushion Cover

☐	XP085	Nylon	Std
☐	XP090	Airlite	\$101.00
☐	XP100	Leather	\$262.00

Accessories

☐	XP120	Frame protectors Neoprene zip fasten (Black)	\$87.00
☐	XP420	Frame protectors Nylon velcro fasten (Black)	\$83.00
☐	XR010	Seat pouch velcro attach (Black)	\$61.00
☐	XR200	Rucksack with RGK logo (Black)	\$156.00
☒	XR090	Tool Kit	No Cost
☒	XS065	Sub4 Calf Strap	No Cost

Spec Notes



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8

Customer Service: (800) 263-3390 E-Mail: cscanada@sunmed.com www.sunrisemedical.ca
Specifications are subject to change without notice. © 2026 Sunrise Medical Canada Inc. CN-100067-CAN REV (6)