

# JAY<sup>®</sup> J3<sup>®</sup> PLUS BACKREST



Designed to specifically address the pelvic and trunk stability needs of the bariatric client.



**SUPERIOR CLINICAL SEATING**

# J3<sup>®</sup> PLUS

## BACKREST

Designed for the bariatric population, J3<sup>®</sup> Plus backs can accommodate up to 500 lbs. in widths of up to 26" and are available with either TH 2-point or FS 4-point hardware. The upper section of the backrest has a 2.2" deep contour to provide for non-restricted movement while still supporting the client.

Available on the Mid Thoracic shell is an optional lower vanity shell extension that can pivot to accommodate redundant tissues allowing proper positioning without interference.

### WIDTH

			22"	24"	26"
MT	<b>S</b>	<b>12"</b>	*	*	*
	<b>M</b>	<b>13.5"</b>			
	<b>T</b>	<b>15"</b>	*	*	*
UT	<b>S</b>	<b>16.5"</b>			
	<b>M</b>	<b>18"</b>			
	<b>T</b>	<b>19.5"</b>			
SH	<b>S</b>	<b>21"</b>			
	<b>M</b>	<b>22.5"</b>			
	<b>T</b>	<b>24"</b>			

\*Optional Vanity Shell Extension available (adds 8" in height)

#### JAY TECHNOLOGY IN THE J3 PLUS BACK



DENSITY



STABILITY



CONTOURING



4-WAY STRETCH



Sunrise Medical (US) LLC, Fresno, CA 93727  
MK-129562 Rev.A © 11.2019



### J3 PLUS BACK SPECIFICATIONS:

**HCPCS Code**  
E2614

**Width**  
22" to 26"

**Contour**  
2.2"

**Height**  
12", 15", 18"  
20", 21" 23"

**Cover Material**  
3DX™ Spacer Fabric  
with Lycra® Stretch

**TH**  
Contact: 2-Point  
Angle: +/- 20°  
Width: 2"  
Depth: 2.25"  
Weight Capacity: 500 lbs.

**FS**  
Contact: 4-Point  
Angle: +/- 10°  
Width: 2"  
Depth: 1.38"  
Weight Capacity: 500 lbs.

**JAY Your Way**  
Customizations Available



800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)