

**Zippie IRIS - AADL W145****October 2026****EIZ5A**

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ Prov: _____ PC: _____

Dealer Phone: () _____ Fax: () _____

Confirm Via: ☐ Fax ☐ Email

Mark For: _____

SUBMITTING FOR: ☐ Quote ☐ Order

PO#: _____

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ Prov: _____ PC: _____

Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

1. FRAME**Base Model**

Maximum user weight is 225 lbs.

Maximum weight capacity for Z-Finity footrests is 165 lbs.

☐	188WM48	Zippie IRIS Includes True Fit Growth	\$4,433
--------------------------	---------	---	---------

☐	188WM91	Zippie IRIS Folding Includes True Fit Growth. N/A with Mono Back or Rocker Cover	\$4,671
--------------------------	---------	--	---------

Frame Struts (Rigid Frames Only)

☐	188GY01	Fixed Struts Growth thru kits N/A w/ Folding Frame	STD
--------------------------	---------	--	-----

☐	188GY02	Width Adjusting System Built in growth of 2" N/A w/ Mono Back or Folding	\$70
--------------------------	---------	--	------

Frame Options

☐	188FT46	Contracture Seat Frame N/A w/10" Frame Depth or Folding Frame	NC
--------------------------	---------	--	----

☐	188FS29	Battery Tray N/A w/ 10" & 11" frame widths, foot wheel lock or Folding Frame	\$617
--------------------------	---------	--	-------

2. SIZE / SEATING**Frame Width (Personalized set-up on 17" & 18" folding frame thru B4Me)**

☐	188W1	10" - N/A w/ Mono Backrest	NC
☐	188W2	11" - N/A w/ Mono Backrest	NC
☐	188W3	12"	NC
☐	188W4	13"	NC
☐	188W5	14"	STD
☐	188W6	15"	NC
☐	188W7	16" - Folding Frame Grows to 18"	NC
☐	188W8	17" - Folding Frame Grows to 18"	NC
☐	188W9	18" - No Growth w/ Folding Frame	NC

Seats

☐	188ST35	Seat Pan-Growable Some limits may exist. See Page 9 for matrix of offering	\$247
☐	188ST36	Drop Seat Pan	\$247
☐	188ST54	Omit Seat Pan	STD

Transit

☐	188TO45	Transit Optional	NC
--------------------------	---------	------------------	----

Tilt Release Options

☐	188TT8	Foot Release N/A with Vent and Battery Tray. May limit folded width, folding frame defaults to right side ☐ Left ☐ Right (facing from rear)	NC
--------------------------	--------	---	----

☐	188TT12	Dual Trigger Cable Release	NC
--------------------------	---------	----------------------------	----

Pre-Set Tilt Angle (55° range) (Default in Bold, N/A w/ reclining back)

☐	188TA5	-5° to 50° N/A with 12.5" STFH	NC
☐	188TA6	0° to 55°	STD
☐	188TA7	+5° to 60°	NC

Pre-Set Tilt Angle (40° range) (Default in Bold)

☐	188TA8	-5° to 35° N/A with 12.5" STFH	NC
☐	188TA9	0° to 40°	STD
☐	188TA10	+5° to 45°	NC

Made to Order Seating

For best seating to wheelchair integration, utilize universal hardware.

☐	188QC01	JAY SureFit	
☐	188QC03	JAY Fit	

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth.

☐	188JC20	JAY Bundled Cushion		Price List
		Name / Part #		
		Width:	Depth:	

2. SIZE / SEATING CONT.

Frame Depth

(Choose Selections from A OR B)

Selection A: Frame Depth Set up By Zippie


Desired Seat Depth

- ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"
☐ 15" ☐ 16" ☐ 17" ☐ 18"

Back Type


- ☐ 3" ☐ 2" ☐ 1"
☐ Flush ☐ -1" ☐ -2"

This section determines frame length



Top of Cushion to Seat Rail

- ☐ 2 1/2" or Greater
☐ Less than 2 1/2" (STD)
☐ Below Seat Rail



Selection B: Frame Depth Range

- | | | |
|----------------------------------|---------------------------|----|
| ☐ 188SN01 | Short Seat Frame 10"-14" | NC |
| ☐ 188SN02 | Medium Seat Frame 12"-16" | NC |
| ☐ 188SN03 | Long Seat Frame 14"-18" | NC |

Selection B: Frame Depth

- | | | |
|---------------------------------|--|----|
| ☐ 188FD1 | 10" rec. w/ Short Seat Frame
N/A with Contracture Footrest System | NC |
| ☐ 188FD2 | 11" rec. w/ Short Seat Frame | NC |
| ☐ 188FD3 | 12" rec. w/ Short Seat Frame | NC |
| ☐ 188FD4 | 13" rec. w/ Med. Seat Frame | NC |
| ☐ 188FD5 | 14" rec. w/ Med. Seat Frame | NC |
| ☐ 188FD6 | 15" rec. w/ Med. Seat Frame | NC |
| ☐ 188FD7 | 16" rec. w/ Long Seat Frame | NC |
| ☐ 188FD8 | 17" rec. w/ Long Seat Frame | NC |
| ☐ 188FD9 | 18" rec. w/ Long Seat Frame | NC |

Select from A OR B

3. CASTERS / SEAT HEIGHT

Caster Wheel

☐ 188CW31	4" Micro Light Up - N/A with Folding	\$110
☐ 188CW73	4" x 1.25" Semi-Pneumatic -N/A w/ Folding	NC
☐ 188CW51	5" Micro Light Up	\$110
☐ 188CW5	5" Low-Profile Polyurethane	STD
☐ 188CW7	5" x 1.5" Semi-Pneumatic	\$110
☐ 188CW33	6" Pneumatic	\$110
☐ 188CW34	6" Polyurethane	NC
☐ 188CW35	6" x 1.5" Semi-Pneumatic	\$110
☐ 188CW21	8" Pneumatic	\$146
☐ 188CW32	8" Polyurethane	NC
☐ 188CW39	8" x 1.5" Semi-Pneumatic	\$146
☐ 188CW36	8" x 2" Pneumatic	\$247
☐ 188CW37	8" x 2" Pneumatic w/Airless Insert	\$393

Please call customer service for 13.25" on a 12" RW and 4" Caster (Rigid only)

Rear Wheel Size	CASTER SIZE					
	☐ 4"	☐ 5"	☐ 6"	☐ 8"x1"	☐ 8"x1.5"	☐ 8"x2"
☐ 12"	☐ 13.0"	☐ 16.0"	☐ 16.0"	N/A	N/A	N/A
	N/A	☐ 16.5"	☐ 16.5"		☐ 16.5"	
		☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"
		☐ 17.5"	N/A	☐ 17.5"	☐ 17.5"	
☐ 12" Low Poly	N/A	☐ 15.0"	☐ 15.0"	N/A	N/A	N/A
		☐ 15.5"	☐ 15.5"			
		☐ 16.0"	☐ 16.0"	N/A	N/A	N/A
		☐ 16.5"	N/A			
☐ 16"	N/A	☐ 13.5"	N/A	N/A	N/A	N/A
		☐ 14.0"	☐ 14.5"			
		☐ 14.5"	☐ 17.0"	☐ 17.5"	☐ 17.5"	N/A
		☐ 17.5"	☐ 17.5"	☐ 18.0"	☐ 18.0"	
		☐ 18.0"	☐ 18.0"	☐ 18.5"	☐ 18.5"	
		☐ 18.5"	☐ 18.5"	N/A	☐ 19.0"	
☐ 20"	N/A	☐ 15.5"	☐ 16.0"	N/A	N/A	N/A
		☐ 16.0"	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"
		☐ 16.5"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"
☐ 22"	N/A	☐ 17.0"	☐ 16.5"	N/A	☐ 16.5"	N/A
		☐ 17.5"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"
		☐ 18.0"	☐ 17.5"	☐ 17.5"	☐ 17.5"	N/A
		N/A	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"
☐ 24"	N/A	N/A	☐ 17.0"	N/A	☐ 17.5"	N/A
		☐ 17.5"	☐ 17.5"	☐ 18.0"	☐ 18.0"	☐ 18.0"
		☐ 18.0"	☐ 18.0"	☐ 18.5"	☐ 18.5"	N/A
		☐ 18.5"	☐ 18.5"	☐ 19.0"	☐ 19.0"	☐ 19.0"
N/A with Folding Frame						

4. REAR WHEELS / SEAT HEIGHT

One Arm Drive

Only available with 6 spoke Mag wheels. One Arm Drive N/A with folding frame

☐	188AL17	One Arm Drive Avail. w/ 20", 22", 24" Mag Rear Wheels	\$1,382
☐	188OA2	Standard Handrim	NC
☐	188HM1	Right	NC
☐	188HM2	Left	NC
☐	188OA1	Small Handrim	NC
☐	188HM1	Right	NC
☐	188HM2	Left	NC

Axle

☐	188AX26	One Arm Drive Axle Only avail. w/One Arm Drive Option	STD
☐	188AX1	Quick Release	NC
☐	188AX11	Quad Quick Release N/A with Lite Spokes	\$182

Rear Wheel Size

12"-16" N/A w/handrims. 12" w/attendant wheel lock N/A w/ +5 - 60° tilt angle.

☐	188RW13	12" Mag - N/A with handrims	STD
☐	188RW14	16" Mag - N/A with handrims	NC
☐	188RW3	20" Mag	NC
☐	188RW8	22" Mag	NC
☐	188RW1	24" Mag	NC
☐	188RW90	20" Lite Spoke	\$320
☐	188RW91	22" Lite Spoke	\$320
☐	188RW54	24" Lite Spoke	\$320
☐	188RW103	20" Mag Wheel 5 Spoke	NC
☐	188RW104	22" Mag Wheel 5 Spoke	NC
☐	188RW102	24" Mag Wheel 5 Spoke	NC

Rear Wheel Tire

☐	188RT1	Pneumatic	NC
☐	188RT2	Pneumatic Airless Insert	\$247
☐	188RT3	Full Polyurethane N/A w/ 12" & 16" rear wheels	NC
☐	188RT4	Low Polyurethane Only Available with 12" wheels	NC

Rear Tire & Wheel Matrix

Tire	Rear Wheel		
	5 Spoke Mag	Spoke	Mag Wheel
Pneumatic	20,22,24	24	12,16,20,22,24
Pneumatic w/ Airless Insert	20,22,24	24	12,16,20,22,24
Full Poly	20,22,24	24	20,22,24
Low Poly	N/A	N/A	12

12" = mm 16" = mm 20" = 451mm 22" = 501mm 24" = 540mm

Handrim

12" and 16" wheels not available with hand rims. See matrix for availability

☐	188HR1	Aluminum Anodized	NC
☐	188HR4	Plastic Coated	\$146
☐	188HR7	Omit N/A with wheel/axle credit	NC

Projections (Not Available on Plastic Coated)

Only available with 24" Mag wheels & not available with spoke guards.

☐	188HR8	8 Oblique	\$146
--------------------------	--------	-----------	-------

Rear Wheel & Handrim Matrix

Handrim	Rear Wheel		
	5 Spoke Mag	Spoke	Mag Wheel
Aluminum Anodized	20,22,24	24	20,22,24
Plastic Coat	20,22,24	24	20,22,24
Projection 8 Oblique	N/A	N/A	24

5. WHEEL LOCKS / HUB LOCK

Wheel Lock

Not available with Hub Lock

☐	188WL34	High Mount Push	NC
☐	188WL36	Foot Lock N/A w/ One Arm Drive, Folding frame defaults to left side, may limit folded width	\$89
☐	188WL10	Omit	NC

Extension Handle

☐	188EH1	E0961 6" Pair	\$73
--------------------------	--------	---------------	------

Hub Lock

Only available with 12" or 16" wheels. N/A with one arm drive.

☐	188LK01	12" Hub Lock	\$162
☐	188LK02	16" Hub Lock	\$162

6. FOOTRESTS (See Appendix on Page 9 for matrix of Leg Length ranges)

Hanger Mounted Footrest (Footplates Mount to Hanger Tubing)

Hanger Type (See Matrix on Page 9 for Lower Leg Length ranges)			
☐	188H140	70° Swing In/Out Front Mount	NC
☐	188H142	80° Swing In/Out Front Mount	NC
☐	188H158	90° Swing In/Out Front Mount	NC
☐	188H9	70° Heavy Duty Front Mount Lift-Off	NC
☐	188H10	80° Heavy Duty Front Mount Lift-Off	NC
☐	188H11	90° Heavy Duty Front Mount Lift-Off	NC
☐	188H161	Elev.Fnt Mount - N/A 10-11" wide	\$339

Extensions (may interfere with some seat heights & casters)			
☐	188E23	2" - not mounted	\$53
☐	188E24	4" - not mounted	\$53

Footplate			
☐	188F9	Kids Angle Adjustable N/A with 17"-18" widths	NC
☐	188F8	Adult Angle Adj. Available with 10"-18" widths	\$176
☐	188F62	Footboard Available with 10"-16" widths	\$230

Ext. Tube Mounted Footrests (Footplates Mount to Ext. Tube)

Hanger Type (See Matrix on Page 9 for Lower Leg Length ranges)			
☐	188H145	60° Swing In/Out	NC
☐	188H139	70° Swing In/Out	STD
☐	188H141	80° Swing In/Out - N/A with 8" Casters	STD
☐	188H159	Articulating - N/A 10"-13" widths	\$526
☐	188H160	Elevating N/A w/ platform footplates or < 14" width	\$339

Extensions (See Matrix on Page 9 for Lower Leg Length ranges)			
☐	188E1	Short	NC
☐	188E2	Medium	NC
☐	188E5	Long	NC
☐	188E16	High Mount	\$97

Footplate			
☐	188F1	Composite N/A In 10", 11" or 12" width	STD
☐	188F3	Adult Angle Adjustable Available with 10"-18" widths	\$176
☐	188F4	Kids Angle Adjustable N/A with 17"-18" widths	NC
☐	188F5	Aluminum Available with elevating 14"-18" wide	NC
☐	188F48	Platform Flip-Up - Ang Adj Available on 60° 70° and 80° ext. mount hangers only. N/A w/Folding Frame	\$230

Hanger Mounts

1" Offset N/A with -5-50° and -5-35° preset tilt angle and possible interference with 80° & 90° hangers

☐	188MZ02	Standard Hanger Mount	STD
☐	188MZ01	1" Offset Hanger Mount	NC

Articulating / Elevating Hanger Cover

☐	188DI47	Neoprene Cover	\$325
☐	188DI48	Gel Pad Cover	\$380

Footrest Hanger Options

☐	188IG2	Neoprene Covers Only available for ext. mount hangers	\$61
☐	188HL1	Heel Loops N/A with Z-Finity 10-14 widths	NC
☐	188HL5	Leg Strap	NC
☐	188A54	Extra Single Leg Strap	\$61
☐	188A110	Toe Loop Only avail. w/ adult ang adj. footplates, N/A with Z-Finity footrest	\$43
☐	188H32	Omit Hanger	NC
☐	188F7	Omit Footplate and Ext. Tube	NC

Z-Finity Hangers/Footrest System

Maximum weight capacity for Z-Finity footrests is 165 lbs.
Some interference may exist in certain configurations.

☐	188H162	Z-Finity Hanger	STD
--------------------------	---------	-----------------	-----

Z-Finity Footrest Mount (Must Chose see Illustration on Page 9)

☐	188DI22	Direct Hanger Mt/Fixed Frame - High	NC
☐	188DI23	Direct Hanger Mt/Fixed Frame - Low	NC
☐	188DI24	Extension Tube Mount	NC

Z-Finity Footrest Pivot Arm (Must Choose, use matrix on Page 9 for leg length)

☐	188DI20	Short 2"	NC
☐	188DI21	Long 3 1/4"	NC

Z-Finity Footplate

Some interference may exist in certain configurations.

☐	188F320	Multi Angle Adjustable - Pair N/A with 8"-9" widths	\$247
☐	188F321	Angle Adjustable Platform Footboard	\$231
☐	188F6	Omit Footplate	NC

Z-Finity Extension Tube Length (When Choosing Ext. Tube Mount)

☐	188E62	Hemi Short	NC
☐	188E1	Short	NC
☐	188E2	Medium	NC
☐	188E5	Long	NC

JAY FOOT BOXES

Only available with Swing Away, Lift-off Hangers, Adult Angle Adj Footplates and ALR/ELR. Single Foot Box pairs N/A with 10"-15" frame widths. Full Foot boxes N/A with Heel Loop or Leg Strap or in odd number sizes, Tall Height N/A w/ 14"-15" frame widths and ALR/ELR. N/A with Z-Finity.

Single Foot Box

☐	188JX01	Standard Height - Left	\$266
☐	188JX02	Tall Height - Left	\$302
☐	188JX03	Standard Height - Right	\$266
☐	188JX04	Tall Height - Right	\$302
☐	188JX05	Standard Height - Pair	\$533
☐	188JX06	Tall Height - Pair	\$563

OR

Full Foot Box

☐	188JX07	Full Foot box Std Height	\$533
☐	188JX08	Full Foot box Tall Height	\$604
☐	188JX09	Full Foot box w/ Divider Std Ht	\$533
☐	188JX10	Full Foot box w/ Divider Tall Ht	\$604

6. FOOTRESTS CONT. (See Appendix for matrix)

Contracture Footrest System (CFS)

N/A w/ Neoprene Impact Guards, Kids Angle Adj. Aluminum Footrests, Dual Post Armrest or Folding Frame.

☐	188H26	Contracture Footrest System	\$598
☐	188FP01	CFS - Center Mount	NC
☐	188FP02	CFS - Outside Mount	NC
☐	188E1	SHORT	NC
☐	188EP1	10.5" Extension Height	NC
☐	188EP2	11.0" Extension Height	NC
☐	188EP3	11.5" Extension Height	NC
☐	188EP4	12.0" Extension Height	NC
☐	188EP5	12.5" Extension Height	NC
☐	188EP6	13.0" Extension Height	NC

☐	188E2	REGULAR	NC
☐	188EP5	12.5" Extension Height	NC
☐	188EP6	13.0" Extension Height	NC
☐	188EP7	13.5" Extension Height	NC
☐	188EP8	14.0" Extension Height	NC
☐	188EP9	14.5" Extension Height	NC
☐	188EP10	15.0" Extension Height	NC

7. BACKREST / SEATING

DYNAMIC BACKREST

☐	188BT101	Dynamic Backrest N/A w/ Reclining or Folding Frame. N/A in some configurations.	\$1,558
--------------------------	----------	---	---------

JAY Backs (Optional)

Ordering JAY backs on-chair, See the JAY Price List. Requires Part #

☐	188JB01	JAY Bundled Backrest Backrest Name or Part # Width: Height:	Price List
--------------------------	---------	---	------------

STROLLER HANDLE BACK POSTS - BACK RESTS (Choose A, B, C or D)

A. Fixed Height

☐	188BP9	Low (19" Ht)	NC
☐	188BP10	Medium (22" Ht)	NC
☐	188BP11	Tall (25" Ht)	NC
☐	188BP13	Removable Extension	\$176
☐	188RG1	Folding Adjustable Rigidizer N/A w/ Mono Back	\$207

A

OR

B Height Adjustable

Rigid Back w/ Angle Adjust Push Handle

☐	188B170	Std. (16"-22" Height.)	NC
☐	188PH8	Std. Angle Adj. Push Handle Lgth	NC
☐	188PH9	Long Angle Adj. Push Handle Length	NC
☐	188B171	Tall (19"-25" Height.)	NC
☐	188PH8	Std. Angle Adj. Push Handle Lgth	NC
☐	188PH9	Long Angle Adj. Push Handle Length	NC

B

OR

C. Mono Back System

Must select JAY J3 when ordering Mono Back.
N/A w/ Folding Frame.

Select One

☐	188MO1	Mono Backrest System	\$609
☐	188MO2	Mono Backrest System with Dynamic Backrest	\$1,895

C

AND Select desired Push Handle

☐	188PH8	Std. Angle Adj. Push Handle Length	NC
☐	188PH9	Long Angle Adj. Push Handle Length	NC


OR

D. Reclining Backrest

Only available with 40° Tilt Range. N/A w/Folding Frame

☐	188BT21	Reclining Backrest	\$1,570
	☐	188B172 Medium (21" Height)	NC
	☐	188B173 Tall (24" Height)	NC
☐	188BP64	Angle Adj stroller handles Height adj range 19"- 23" N/A on 10" widths w/Dual Trigger tilt release	\$207
	☐	188PH8 Standard push handle	NC
	☐	188PH9 Long push handle	NC

D



7. BACKREST / SEATING CONT.

Mono Backrest Options ONLY

N/A on Folding Frame, Vent/Battery Tray Base(LTV only), Reclining Backrest, Desk Length Arm pads, Canopy and 10"-11" frame widths and depths.

JAY J3 Back for Mono Backrest System

☐	188JB04	JAY J3 Back Mono	See OF
☐	188JB05	JAY J3 via JAY Your Way	See OF

JAY J3 Backrest Contour

☐	188JP06	J3 PA Back - 2.2" of Contour Depth	\$923
☐	188JP08	J3 PL Back - 3" of Contour Depth	\$997
☐	188JP09	J3 PDL Back - 6" of Contour Depth N/A w/ 12" width	\$1,070
☐	188JP10	J3 PDC Back - 6" of Contour Depth N/A w/ 12" width	\$1,080

JAY J3 Backrest Height

☐	188JR07	UT Standard - 16.5" shell height	NC
☐	188JR08	UT Medium - 18" shell height N/A w/ 12" width	NC
☐	188JR09	UT Tall - 19.5" shell height N/A w/ 12" width	NC
☐	188JR10	SH Standard - 21" shell height N/A w/ 12" width	NC
☐	188JR11	SH Medium - 22.5" shell height N/A w/ 12" width	NC
☐	188JR12	SH Tall - 24" shell height N/A w/ 12" width	NC

JAY J3 Backrest Width

Some interference may exist w/ armrest when ordering a backrest width greater than 4" wider than the base.

☐	188JQ06	12" Backrest Width Available w/UT Standard 16.5" shell ht	NC
☐	188JQ01	14" Backrest Width	NC
☐	188JQ02	16" Backrest Width	NC
☐	188JQ03	18" Backrest Width	NC

JAY J3 Backrest Color

☐	188GZ01	Black	NC
☐	188GZ02	Match Frame- Standard Colors	NC
☐	188GZ03	Match Frame- Personalized Color	\$77

JAY J3 Back for Mono Backrest System Cont.

Select desired Headrest Mount - optional


☐	188MH10	Axys Mount & Bracket	\$362
☐	188MH15	Cobra Flip Back Mount and Bracket	\$664
☐	188MH16	Onyx Mount and Bracket	\$272
☐	188MH17	Omit	NC

Select desired Headrest Pad - optional

☐	188HE65	Plush Head Support Pad 8"	\$272
☐	188HE66	Plush Head Support Pad 10"	\$272
☐	188HE78	Plush Head Support Pad 14"	\$272
☐	188HE87	HeadsUp Lrg: 5.3" Lng Arms	\$362
☐	188HE88	HeadsUp Med: 4.1" Sht Arms	\$362
☐	188HE89	HeadsUp Sml: 3.5" Sht Arms	\$362
☐	188HE83	Adjust-A-Plush Standard	\$281
☐	188HE82	Adjust-A-Plush Narrow	\$281

Select desired Lateral - optional

☐	LAT 4X4 S/A MOUNTS includes 1" OFFSET - PR		\$795
	188ZV14	LATERAL 4X4 W/ 1" OFFSET-PAIR	\$470
	188ZS07	LATERAL SWING AWAY MOUNT-PAIR	\$325
☐	LAT 4X5 S/A MOUNTS includes 1" OFFSET - PR		\$795
	188ZV15	LATERAL 4X5 W/ 1" OFFSET-PAIR	\$470
	188ZS07	LATERAL SWING AWAY MOUNT-PAIR	\$325
☐	LAT 4X6 S/A MOUNTS includes 1" OFFSET - PR		\$795
	188ZV16	LATERAL 4X6 W/ 1" OFFSET-PAIR	\$470
	188ZS07	LATERAL SWING AWAY MOUNT-PAIR	\$325
☐	LAT 6X4 S/A MOUNTS includes 1" OFFSET - PR		\$795
	188ZV17	LATERAL 6X4 W/ 1" OFFSET-PAIR	\$470
	188ZS07	LATERAL SWING AWAY MOUNT-PAIR	\$325
☐	LAT 6X5 S/A MOUNTS includes 1" OFFSET - PR		\$795
	188ZV18	E0956 LATERAL 6X5 W/ 1" OFFSET-PAIR	\$470
	188ZS07	E1028 LATERAL SWING AWAY MOUNT-PAIR	\$325
☐	LAT 6X6 S/A MOUNTS includes 1" OFFSET - PR		\$795
	188ZV19	E0956 LATERAL 6X6 W/ 1" OFFSET-PAIR	\$470
	188ZS07	E1028 LATERAL SWING AWAY MOUNT-PAIR	\$325



8. ARMRESTS / SIDE GUARDS

Armrest - Height Adjustable

Dual Post not available with 0 - 13" deep or with 14 - 15" deep seat pan, Dual Post also not available with contracture Footrest. Some interference may exist in some configurations.

☐	188AA12	T-Post Regular w/Foam	NC
☐	188AA13	T-Post Tall w/Foam	NC
☐	188AA3	Flip-Up	NC
☐	188AA1	Adjustable Locking Flip-Up	\$339
☐	188AA7	Length Adjust. Locking Flip-Up N/A with Contour Armpad	\$339
☐	188AA62	Dual Post	NC
☐	188AA63	Dual Post Height Adjustable	\$362
☐	188AA43	Dual Post Flip Back	\$168
☐	188AA48	Dual Post Flip Back Height Adj.	\$339
☐	188AA11	Single-Post Height Adj. Std.	\$339
☐	188AA29	Single-Post Height Adj. Low	\$339
☐	188AA154	Height Adj. Cantilever Armrest Avail. w/ Mono Back only. N/A w/ 20", 22" and 24" rear wheels.	\$308

Armpad

Contour & Ergo only w/ Single Post Height Adj. - Foam only Length Adj.

☐	188AM23	Desk Classic	NC
☐	188AM21	Desk Contour	NC
☐	188AM24	Full Classic	NC
☐	188AM22	Full Contour	NC
☐	188AM9	Foam Short	NC
☐	188AM10	Foam Long	NC
☐	188AM27	Ergo armpad	NC

Armrest Receiver for Single Post Height Adjustable Armrest

☐	188AR17	Low	NC
☐	188AR18	Std	NC

JAY Arm Troughs

N/A with Length Adj. Locking Flip Up and Adj Locking Flip Up armrests
and Single Post Height Adj Armrests. Full Length Classic Armpads will automatically be selected when choosing only one side Arm Trough.

☐	188JL01	Arm Trough - 14" long - Left	\$184
☐	188JL04	Arm Trough - 14" long - Right	\$184
☐	188JL05	Arm Trough - 14" long - Pair	\$367
☐	188JL02	Arm Trough - 18" long - Left	\$184
☐	188JL06	Arm Trough - 18" long - Right	\$184
☐	188JL07	Arm Trough - 18" long - Pair	\$367
☐	188JL03	Arm Trough - 22" long - Left	\$184
☐	188JL08	Arm Trough - 22" long - Right	\$184
☐	188JL09	Arm Trough - 22" long - Pair	\$367

JAY Arm Trough - Hardware


☐	188JN01	Fixed - Each	NC
☐	188JN03	Fixed - Pair	NC
☐	188JN02	Swing Away - Rotational - Each	\$290
☐	188JN04	Swing Away - Rotational - Pair	\$574

Side Guard

Not available with height adjustable armrests, T Posts or
Cantilever armrests. N/A with Mono Back.

☐	188SG33	Plastic Std	\$188
☐	188SG29	Plastic Low	\$188

Armrest Matrix for Single Post


Armrest Matrix for Single Post					
Armrest Size	Transfer Tube Length D	Receiver	Classic Pad		
			H Range	Pad length L	
				Desk	Full
Std	8"	Std	10-3/4" to 15-3/4"	10-1/4"	14"
		Low	8-3/4" to 13-3/4"		
Low	5-1/4"	Std	9" to 13"		
		Low	7-1/4" to 11-1/4"		
		Receiver	Contour Pad		
			H Range	Pad length L	
				Desk	Full
		Std	10" to 15"	10-1/2"	14-1/4"
		Low	8-1/4" to 13-1/4"		
		Std	8-1/2" to 12-1/2"		
		Low	6-3/4" to 10-3/4"		
		Receiver	Ergo Pad		
			H Range	Pad length L	
				11-3/4"	
		Std	10-1/2" to 15-1/2"		
		Low	8-3/4" to 13-3/4"		
		Std	9" to 13"		
		Low	7-1/4" to 11-1/4"		

9. ACCESSORIES			
Back Pack (Kids Backpack included)			
All Zippie chairs ship with one Kids Backpack at NC			
☐	188A1	Utility Backpack - Black	\$112
☐	188A2	Additional Zippie Backpack	\$105
Seat Pouch			
☐	188A8	Black	\$83
Miscellaneous			
☐	188A83	Tool Kit	\$43
Positioning Belt (w/clamps)			
☐	188A112	1" Dual Pull Side Squeeze Padded	\$132
☐	188A117	Belt 2" Aircraft Buckle	\$139
☐	188A118	Belt 2" Aircraft Buckle Padded	\$204
☐	188A814	JAY 2 Pt 1" Push Btn - 4.5" Pad	\$146
☐	188ZZ25	JAY 2 Pt 1.5" Push Btn - No Pad	\$64
☐	188ZZ28	JAY 2 Pt 1.5" Push Btn - 7.5" Pad	\$126
☐	188ZZ27	JAY 2 Pt 2" Push Btn - No Pad	\$64
☐	188ZZ26	JAY 2 Pt 2" Push Btn - 9.5" Pad	\$126
Vertical Mount Vent Tray-Back cane Mount			
☐	188LF12	Dovetail Vertical Frame Mounted Hardware for LTV Vent (limit 15lbs) includes 1 Utility Hook. - <i>Must remove for folding.</i>	\$522
☐	188LF13	Dovetail Vertical Frame Mt Hardware for Trilogy Vent (limit 15lbs) includes 1 Utility Hook - <i>Must remove for folding.</i>	\$522
☐	188LF7	Utility Hook	\$31
Upper Extremity Support Tray			
When ordering a tray with full armrest or desk armrest pads, add 1" to seat width to determine best tray size. N/A with JAY Arm Trough. Must remove for folding.			
☐	188A29	X-Sm (Suggested 10"-12" Width)	\$381
☐	188A30	Sm (Suggested 13"-14" Width)	\$381
☐	188A31	Med (Suggested 15"-17" Width)	\$381
☐	188A32	Lrg (Suggested 18"-20" Width)	\$381
10. FRAME COLOR			
Standard Colors			
☐	188C211	Celestial Sapphire	NC
☐	188C212	Desert Sage	NC
☐	188C213	Desert Tan	NC
☐	188C214	Gloss Black	NC
☐	188C215	Matte Black	NC
☐	188C216	Metallic Shadow	NC
☐	188C217	Midnight Blue	NC
☐	188C218	Mystic Purple	NC
☐	188C219	Petal Pink	NC
☐	188C220	Pink Power	NC
☐	188C221	Red Pepper	NC
☐	188C222	Royal Blue	NC
☐	188C223	Salute Green	NC
☐	188C224	Scarlet Noir	NC
☐	188C225	Silver Ore	NC
☐	188C226	Solar Orange	NC
☐	188C227	Starlight White	NC
☐	188C228	Storm Grey	NC
☐	188C229	Turbo Green	NC
☐	188C230	Velocity Yellow	NC
Canopy (SunShade)			
Need 2" back post clearance to mount canopy hardware, may require next height based on seating; please note some interference may occur with some headrest configurations. N/A with Mono Back. Must remove for folding.			
☐	1888A06	Canopy (Sunshade) 10"-14" Width	\$428
☐	1888A07	Canopy (Sunshade) 15"-20" Width	\$428
JAY Anterior Thoracic Support			
Style			
☐	188JY01	Classic Dynamic	\$245
☐	188JY02	Classic Structure	\$245
☐	188JY03	Contour Dynamic	\$245
☐	188JY04	Contour Structure	\$245
☐	188JY05	Center Opening with Buckle	\$245
Size			
☐	188JZ01	Small	NC
☐	188JZ02	Medium	NC
☐	188JZ03	Large	NC
☐	188JZ04	X-Large	NC
JAY Padded Trunk Strap			
☐	188JA01	Padded Trunk Strap - Small	\$245
☐	188JA02	Padded Trunk Strap - Medium	\$245
☐	188JA03	Padded Trunk Strap - Large	\$245
☐	188JA04	Padded Trunk Strap - XL	\$245
Spoke Guard			
Not available with 12" & 16" rear wheels & Projections Handrims			
☐	188A23	Black	\$157
☐	188A303	Clear Zippie	\$157
Anti-Tip			
☐	188AT1	Rear Anti-Tip	\$121


APPENDIX

Z-Finity Matrix

Direct Hanger Mount



Extension Tube Mount



Direct Hanger Mounted


Footrest Adaptor	Hanger Mounting Position	Min	Max
Short	High Mount	-0.9"	3.4"
	Low Mount	2.6"	6.7"
Long	High Mount	-2.1"	4.6"
	Low Mount	1.3"	7.9"

Extension Tube Mounted

Footrest Adaptor	Extension Tube Length	Min	Max
Short	Hemi Short	3.9"	12.4"
	Short	4.9"	13.4"
	Medium	6.6"	15.4"
	Long	8.6"	17.4"
Long	Hemi Short	2.6"	13.6"
	Short	4.1"	14.6"
	Medium	5.6"	16.6"
	Long	7.6"	18.6"

Lower Leg Length Matrix


Hanger	Extension	Seat Pan to Top of Footrest		Min Allowable FSH
		Min	Max	
70° Front Mount	No Ext.	3 1/2"	9"	11 1/2"
	2"	3 1/2"	11 1/2"	13 1/2"
	4"	3 1/2"	13 1/2"	15 1/2"
80° Front Mount	No Ext.	3 1/2"	9 1/2"	11 1/2"
	2"	3 1/2"	12"	14"
	4"	3 1/2"	14"	16"
90° Front Mount	No Ext.	3"	9"	11 1/2"
ELR Front Mount	No Ext.	3"	7"	10 1/2"
60° Swing In/Out	Short	9"	12 1/2"	12 1/2"
	Medium	9"	14 1/2"	12 1/2"
	Long	11 1/2"	16"	14 1/2"
70° Swing In/Out	Short	10 1/2"	15 1/2"	12 1/2"
	Medium	10 1/2"	17 1/2"	13"
	Long	13 1/2"	19"	15"
80° Swing In/Out	Short	11"	14 1/2"	13"
	Medium	11"	16 1/2"	13"
	Long	13"	18"	15"
ALR	Short	14"	17 1/2"	12 1/2"
	Medium	15"	19 1/2"	13 1/2"
	Long	16 1/2"	19"	14 1/2"
ELR Ext Mount	Short	13"	17"	12 1/2"
	Medium	13"	18 1/2"	13 1/2"
	Long	14 1/2"	20"	14 1/2"



Growing Seat Pan Matrix - Rigid and Folding Frame Options


Rigid Frame - Frame Depth Growth Range Matrix

Seat Pan Size	Dim A Range
10" - 16"	11-1/4" - 15-1/4"
14" - 18"	15" - 19"



Folding Frame - Frame Depth Growth Range Matrix

Seat Depth Pan Size	Seat Depth											L Range	
	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	Min	Max
Short	x	x	x	x	x							10"	14"
Medium					x	x						11.75"	16.25"
Long							x	x				13.75"	18.25"
X-Long									x	x	x	15.75"	20.25"



Depth

Additional Options Through Built-4-Me

To inquire about any unique personalization for Zippie products, contact:
Phone: (800) 263-3390, Email: cscanada@sunmed.com



Commonly Requested Options Include:

- ✓ Reverse Configuration
- ✓ Angle Adjustable Handle for Reclining Backrest
- ✓ Dual Mount Wheel Locks - Attendant Foot Lock
- ✓ Dynamic Back on Frame Depth up to 20"
- ✓ O2 Holder with Transit
- ✓ IV Pole with Transit
- ✓ Solid Seat w/ Quick Release Hardware for folding frame



Sunrise Medical Canada Inc.

1000 Creditstone Road, Unit 2 • Concord, ON L4K 4P8 • Canada