

Mark For:

Submitting for:

Quote Order

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ Prov: _____ PC: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ Prov: _____ PC: _____
 Ship To Phone: () _____ Fax: () _____

JAY ZIP CUSHION

JAY Zip Cushion			
<input type="checkbox"/>	JZXXYY	JAY Zip Cushion 8-16 wide	\$700.00
Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)			
Includes: Ionic+ [™] fabric outer cover, inner cover, and foam base			

2 Outer Covers included with each Cushion!



		JAY Zip Size Matrix										
		WIDTH										
		8	9	10	11	12	13	14	15	16		
DEPTH	8					N/A	N/A	N/A	N/A	N/A		
	9						N/A	N/A	N/A	N/A		
	10							N/A	N/A	N/A		
	11								N/A	N/A		
	12									N/A		
	13										N/A	
	14											N/A
	15	N/A										
	16	N/A	N/A									
17	N/A	N/A	N/A									
18	N/A	N/A	N/A	N/A								


	=	Standard Offering
JYW	=	Offered Via Jay Your Way
N/A	=	Not available

JAY Zip Measuring Guide

Measurements made from completed cushion with covers.

Outer Cover
 The JAY Zip outer cover is designed to promote comfort, convenience and breathability by the use of "Ionic+[™]" material and 3DX Spacer Fabric.

Inner Cover
 The water resistant inner cover, combined with the Aquaguard[™] Zipper and anti wicking thread provides the ultimate protection against odor.



Width
 X = Width
 Note: Measurement made from bottom outside edge to outside edge of cushion.
 Ex: JZ1214 X = 12"

Depth
 Y = Depth
 Note: Measurement made from bottom rear edge to front lower edge of cushion.
 Ex: JZ1214 Y = 14"

Bottom View

Note: To order via JAY Your Way add a "M" to the front of the part number and follow the prompts. For additional modifications see Page 2 or visit www.sunrisemedical.ca to view even more!

Please visit www.sunparts.us for spare parts information.

JAY Zip Cushion

✓	Option #	Definition	Price	✓	Option #	Definition	Price
<input type="checkbox"/>	JZ0808	JAY ZIP 8" X 8"	\$700.00	<input type="checkbox"/>	JZ1212	JAY ZIP 12" X 12"	\$700.00
<input type="checkbox"/>	JZ0809	JAY ZIP 8" X 9"	\$700.00	<input type="checkbox"/>	JZ1213	JAY ZIP 12" X 13"	\$700.00
<input type="checkbox"/>	JZ0810	JAY ZIP 8" X 10"	\$700.00	<input type="checkbox"/>	JZ1214	JAY ZIP 12" X 14"	\$700.00
<input type="checkbox"/>	JZ0811	JAY ZIP 8" X 11"	\$700.00	<input type="checkbox"/>	JZ1215	JAY ZIP 12" X 15"	\$700.00
<input type="checkbox"/>	JZ0812	JAY ZIP 8" X 12"	\$700.00	<input type="checkbox"/>	JZ1216	JAY ZIP 12" X 16"	\$700.00
<input type="checkbox"/>	JZ0813	JAY ZIP 8" X 13"	\$700.00	<input type="checkbox"/>	JZ1217	JAY ZIP 12" X 17"	\$700.00
<input type="checkbox"/>	JZ0814	JAY ZIP 8" X 14"	\$700.00	<input type="checkbox"/>	JZ1218	JAY ZIP 12" X 18"	\$700.00
<input type="checkbox"/>	JZ0908	JAY ZIP 9" X 8"	\$700.00	<input type="checkbox"/>	JZ1310	JAY ZIP 13" X 10"	\$700.00
<input type="checkbox"/>	JZ0909	JAY ZIP 9" X 9"	\$700.00	<input type="checkbox"/>	JZ1311	JAY ZIP 13" X 11"	\$700.00
<input type="checkbox"/>	JZ0910	JAY ZIP 9" X 10"	\$700.00	<input type="checkbox"/>	JZ1312	JAY ZIP 13" X 12"	\$700.00
<input type="checkbox"/>	JZ0911	JAY ZIP 9" X 11"	\$700.00	<input type="checkbox"/>	JZ1313	JAY ZIP 13" X 13"	\$700.00
<input type="checkbox"/>	JZ0912	JAY ZIP 9" X 12"	\$700.00	<input type="checkbox"/>	JZ1314	JAY ZIP 13" X 14"	\$700.00
<input type="checkbox"/>	JZ0913	JAY ZIP 9" X 13"	\$700.00	<input type="checkbox"/>	JZ1315	JAY ZIP 13" X 15"	\$700.00
<input type="checkbox"/>	JZ0914	JAY ZIP 9" X 14"	\$700.00	<input type="checkbox"/>	JZ1316	JAY ZIP 13" X 16"	\$700.00
<input type="checkbox"/>	JZ0915	JAY ZIP 9" X 15"	\$700.00	<input type="checkbox"/>	JZ1317	JAY ZIP 13" X 17"	\$700.00
<input type="checkbox"/>	JZ1008	JAY ZIP 10" X 8"	\$700.00	<input type="checkbox"/>	JZ1318	JAY ZIP 13" X 18"	\$700.00
<input type="checkbox"/>	JZ1009	JAY ZIP 10" X 9"	\$700.00	<input type="checkbox"/>	JZ1411	JAY ZIP 14" X 11"	\$700.00
<input type="checkbox"/>	JZ1010	JAY ZIP 10" X 10"	\$700.00	<input type="checkbox"/>	JZ1412	JAY ZIP 14" X 12"	\$700.00
<input type="checkbox"/>	JZ1011	JAY ZIP 10" X 11"	\$700.00	<input type="checkbox"/>	JZ1413	JAY ZIP 14" X 13"	\$700.00
<input type="checkbox"/>	JZ1012	JAY ZIP 10" X 12"	\$700.00	<input type="checkbox"/>	JZ1414	JAY ZIP 14" X 14"	\$700.00
<input type="checkbox"/>	JZ1013	JAY ZIP 10" X 13"	\$700.00	<input type="checkbox"/>	JZ1415	JAY ZIP 14" X 15"	\$700.00
<input type="checkbox"/>	JZ1014	JAY ZIP 10" X 14"	\$700.00	<input type="checkbox"/>	JZ1416	JAY ZIP 14" X 16"	\$700.00
<input type="checkbox"/>	JZ1015	JAY ZIP 10" X 15"	\$700.00	<input type="checkbox"/>	JZ1417	JAY ZIP 14" X 17"	\$700.00
<input type="checkbox"/>	JZ1016	JAY ZIP 10" X 16"	\$700.00	<input type="checkbox"/>	JZ1418	JAY ZIP 14" X 18"	\$700.00
<input type="checkbox"/>	JZ1108	JAY ZIP 11" X 8"	\$700.00	<input type="checkbox"/>	JZ1512	JAY ZIP 15" X 12"	\$700.00
<input type="checkbox"/>	JZ1109	JAY ZIP 11" X 9"	\$700.00	<input type="checkbox"/>	JZ1513	JAY ZIP 15" X 13"	\$700.00
<input type="checkbox"/>	JZ1110	JAY ZIP 11" X 10"	\$700.00	<input type="checkbox"/>	JZ1514	JAY ZIP 15" X 14"	\$700.00
<input type="checkbox"/>	JZ1111	JAY ZIP 11" X 11"	\$700.00	<input type="checkbox"/>	JZ1515	JAY ZIP 15" X 15"	\$700.00
<input type="checkbox"/>	JZ1112	JAY ZIP 11" X 12"	\$700.00	<input type="checkbox"/>	JZ1516	JAY ZIP 15" X 16"	\$700.00
<input type="checkbox"/>	JZ1113	JAY ZIP 11" X 13"	\$700.00	<input type="checkbox"/>	JZ1517	JAY ZIP 15" X 17"	\$700.00
<input type="checkbox"/>	JZ1114	JAY ZIP 11" X 14"	\$700.00	<input type="checkbox"/>	JZ1518	JAY ZIP 15" X 18"	\$700.00
<input type="checkbox"/>	JZ1115	JAY ZIP 11" X 15"	\$700.00	<input type="checkbox"/>	JZ1613	JAY ZIP 16" X 13"	\$700.00
<input type="checkbox"/>	JZ1116	JAY ZIP 11" X 16"	\$700.00	<input type="checkbox"/>	JZ1614	JAY ZIP 16" X 14"	\$700.00
<input type="checkbox"/>	JZ1117	JAY ZIP 11" X 17"	\$700.00	<input type="checkbox"/>	JZ1615	JAY ZIP 16" X 15"	\$700.00
<input type="checkbox"/>	JZ1209	JAY ZIP 12" X 9"	\$700.00	<input type="checkbox"/>	JZ1616	JAY ZIP 16" X 16"	\$700.00
<input type="checkbox"/>	JZ1210	JAY ZIP 12" X 10"	\$700.00	<input type="checkbox"/>	JZ1617	JAY ZIP 16" X 17"	\$700.00
<input type="checkbox"/>	JZ1211	JAY ZIP 12" X 11"	\$700.00	<input type="checkbox"/>	JZ1618	JAY ZIP 16" X 18"	\$700.00

JAY Your Way - Modifications

REDUCED HEIGHT			
<input type="checkbox"/>	RH	Remove 1/2" of foam from bottom	NC
NOTCH FOR BACK POST			
<input type="checkbox"/>	NP	Notch for Back Post	NC
Specify distance from back of cushion:			
<input type="checkbox"/>	1"	<input type="checkbox"/>	2"
<input type="checkbox"/>	3"		
NOTCH FOR SEAT RAILS			
<input type="checkbox"/>	NR	Notch for Seat Rails	NC
COVER OPTIONS			
<input type="checkbox"/>	CVR	Cover Options	NC
<input type="checkbox"/>	Dartex	<input type="checkbox"/>	Incontinent
<input type="checkbox"/>	Microclimatic		

ANTERIOR OR POSTERIOR SLOPE			
<input type="checkbox"/>	WGP	Posterior Slope - Wedge	NC
Wedge glued underneath cushion (Cover sewn to match)			
Select desired height increase at front of cushion			
<input type="checkbox"/>	.5"	<input type="checkbox"/>	1"
<input type="checkbox"/>	1.5"		
<input type="checkbox"/>	WGA	Anterior Slope - Wedge	NC
Wedge glued underneath cushion (Cover sewn to match)			
Select desired height increase at rear of cushion			
<input type="checkbox"/>	.5"	<input type="checkbox"/>	1"
<input type="checkbox"/>	1.5"		

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:



Sunrise Medical Canada Inc

237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada

Customer Service (800) 263-3390 | Fax (800) 561-5834 | Email cscanada@sunmed.com | www.sunrisemedical.ca
 Specifications subject to change without notice. © 02/2024 Sunrise Medical Canada Inc. MK-100077-CAN REV (4).