



ORDER FORM

Advance Bath Chair

EILK04

Versatile bath chair for all ages
ensuring safety and comfort for your child

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Advance Bath Chair

Weight Capacity: 159 lbs

<input checked="" type="checkbox"/> OPTION	DESCRIPTION	RETAIL
Packages		\$
Includes Chair Frame, Chair Cover, Hip Belt, and Chest Belt. Substitutions are not available for packages.		
<input type="checkbox"/> 152QE23	Size 1 Advance Bath Chair Package - Pink	\$1064
<input type="checkbox"/> 152QE24	Size 1 Advance Bath Chair Package - Royal Blue	\$1064
<input type="checkbox"/> 152QE25	Size 2 Advance Bath Chair Package - Pink	\$1064
<input type="checkbox"/> 152QE26	Size 2 Advance Bath Chair Package - Royal Blue	\$1064
<input type="checkbox"/> 152QE27	Size 3 Advance Bath Chair Package - Pink	\$1139
<input type="checkbox"/> 152QE28	Size 3 Advance Bath Chair Package - Royal Blue	\$1139
<input type="checkbox"/> 152QE29	Size 4 Advance Bath Chair Package - Pink	\$1169
<input type="checkbox"/> 152QE30	Size 4 Advance Bath Chair Package - Royal Blue	\$1169
Chair Frame		\$
Selection Required if an Advance Bath Chair Package was not selected.		
<input type="checkbox"/> 152BS29	Size 1 Chair Frame	\$619
<input type="checkbox"/> 152BS30	Size 2 Chair Frame	\$619
<input type="checkbox"/> 152BS31	Size 3 Chair Frame	\$678
<input type="checkbox"/> 152BS32	Size 4 Chair Frame	\$708
Chair Cover		\$
Selection Required. Must match chair frame size.		
<input type="checkbox"/> 152IU23	Size 1 Cover - Pink	\$208
<input type="checkbox"/> 152IU24	Size 1 Cover - Royal Blue	\$208
<input type="checkbox"/> 152IU25	Size 2 Cover - Pink	\$208
<input type="checkbox"/> 152IU26	Size 2 Cover - Royal Blue	\$208
<input type="checkbox"/> 152IU27	Size 3 Cover - Pink	\$224
<input type="checkbox"/> 152IU28	Size 3 Cover - Royal Blue	\$224
<input type="checkbox"/> 152IU29	Size 4 Cover - Pink	\$224
<input type="checkbox"/> 152IU30	Size 4 Cover - Royal Blue	\$224
Accessories		\$
Abduction Belts must match chair frame size.		
<input type="checkbox"/> 152A876	Hip Belt - Pink	\$116
<input type="checkbox"/> 152A877	Hip Belt - Royal Blue	\$116
<input type="checkbox"/> 152EJ07	Chest Belt - Pink	\$121
<input type="checkbox"/> 152EJ08	Chest Belt - Royal Blue	\$121
<input type="checkbox"/> 152A878	Size 1 and 2 Abduction Belt - Pink	\$121
<input type="checkbox"/> 152A879	Size 1 and 2 Abduction Belt - Royal Blue	\$121
<input type="checkbox"/> 152A880	Size 3 and 4 Abduction Belt - Pink	\$121
<input type="checkbox"/> 152A881	Size 3 and 4 Abduction Belt - Royal Blue	\$121
<input type="checkbox"/> 152HQ08	Head Laterals - Pink	\$121
<input type="checkbox"/> 152HQ09	Head Laterals - Royal Blue	\$121
<input type="checkbox"/> 1528A119	Shower Trolley	\$1307

Technical Specifications - Imperial

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 5 years	4 - 9 years	8 - 14 years	12 - 18 years
User Height (Max / Min)	30" / 41"	37" / 49"	45" / 57"	53" / 65"
Head Support	6"	6"	8.5"	8.5"
Back Support Height	14"	16"	18"	22"
Seat Depth	9"	11"	13"	16"
Calf Support	10"	10"	13"	13"
Seat Height	2" - 17"	2" - 17"	2" - 17"	2" - 17"
Maximum Weight	159 lbs	159 lbs	159 lbs	159 lbs
Total Width	17"	17"	17"	17"

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 5 years	4 - 9 years	8 - 14 years	12 - 18 years
User Height (Max / Min)	750 mm / 1050 mm	950 mm / 1250 mm	1150 mm / 1450 mm	1350 mm / 1650 mm
Head Support	150 mm	150 mm	215 mm	215 mm
Back Support Height	350 mm	400 mm	460 mm	550 mm
Seat Depth	225 mm	275 mm	335 mm	400 mm
Calf Support	260 mm	260 mm	340 mm	340 mm
Seat Height	50 mm - 425 mm	50 mm - 425 mm	50 mm - 425 mm	50 mm - 425 mm
Maximum Weight	72 kgs	72 kgs	72 kgs	72 kgs
Total Width	430 mm	430 mm	430 mm	430 mm

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, Ontario L4K 4P8



MK-129844-CAN © 07/2024

