

GEMINO 30
 GEMINO 30 M
 GEMINO 30 S
 September 2026



Product sold with backrest in Canada.

ORDER FORM

Invoice address or customer number:

Customer number:

Name:

Street:

Town:

Postal code:

Country:

E-mail:

Tel:

Delivery address:

☐ Deliver to Invoice address

Customer number:

Name:

Street:

Town:

Postal code:

Country:

E-mail:

Tel:

Order information

☐ Order

☐ Quote

Order date:


Your reference:

Reimbursement code: ADP#: AA0001082

Customer Service: cscanada@sunmed.com

Base price \$ 743.00

DIMENSIONS



	Gemino 30	Gemino 30 M	Gemino 30 S
A Total height	78-100 cm	69-88 cm	65-77 cm
B Total length	65 cm	65 cm	68 cm
C Total width	60 cm	60 cm	55 cm
D Width between push handles	47 cm	47 cm	42 cm
E Seat height	62 cm	55 cm	48 cm
Height folded	80 cm	76 cm	69 cm
Length folded	65 cm	65 cm	68 cm
Width folded	23 cm	23 cm	23 cm
Turning circle	84 cm	84 cm	79 cm
F Front / rear wheels	20 x 3.6 cm	20 x 3.6 cm	20 x 3.6 cm
Total weight without basket	7.1 kg	7 kg	6.6 kg
G Weight basket, standard	0.3 kg	0.3 kg	0.25 kg
Max. user weight	150 kg	130 kg	125 kg
Max. weight in basket	5 kg	5 kg	5 kg
Recommended user height	150-200 cm	135-170 cm	125-165 cm

ORDER FORM GEMINO 30 / GEMINO 30 M / GEMINO 30 S

SIZES AND COLOURS		Price \$	Quantity	Total
Gemino 30 + HC-9012222 Adj. Comfort Back				
☐ HC-7100000	Gemino 30, silver grey	743.00		
☐ HC-7100001	Gemino 30, champagne	743.00		
☐ HC-7100002	Gemino 30, sapphire blue	743.00		
☐ HC-7100003	Gemino 30, pink	743.00		
Gemino 30 M + HC-9012222 Adj. Comfort Back				
☐ HC-7100200	Gemino 30 M, silver grey	743.00		
☐ HC-7100201	Gemino 30 M, champagne	743.00		
☐ HC-7100202	Gemino 30 M, sapphire blue	743.00		
Gemino 30 S + HC-9012222 Adj. Comfort Back				
☐ HC-7100100	Gemino 30 S, silver grey	743.00		
☐ HC-K7100101	Gemino 30 S, sapphire blue	743.00		
☐ HC-7100103	Gemino 30 S, pink	743.00		

ACCESSORIES		Price \$	Quantity	Total
Bag for Gemino 30 S. <i>NA for Gemino 30 and 30 M</i>				
☐ HC-7104007	Bag, advanced with lid and handles, black/black	129.00		
Other accessories				
☐ HC-7109011	Crutch holder/ Cane holder	106.00		
☐ HC-7109140	Table tray	106.00		
☐ HC-7109060	Umbrella, black	264.00		
☐ HC-7109050	Basket for oxygen cylinder - Gemino 30, Gemino 30 M	264.00		
☐ HC-7109051	Basket for oxygen cylinder - Gemino 30 S	264.00		
☐ HC-7109080	Slowdown brake, set	252.00		
☐ HC-7109095	Light	129.00		
☐ HC-7109100	Infusion holder	264.00		

GEMINO 30 / GEMINO 30 M / GEMINO 30 S		Price \$
Base price		
Option price		
Total price		

GEMINO ACCESSORIES

GEMINO

BAG - GEMINO 20 S & GEMINO 30 S



HC-7104007 Bag,
advanced, black/black **\$129.00**

OTHER ACCESSORIES



HC-7109011 **\$106.00**
Crutch holder



HC-7109140 **\$106.00**
Table tray



HC-7109095 **\$129.00**
Light



HC-7109060 **\$264.00**
Umbrella, black



HC-7109050 & HC-7109051 **\$264.00**
Basket for oxygen cylinder*



HC-7109100 **\$264.00**
Infusion holder



Sunrise Medical
1000 Creditstone Road, Unit 2
Concord, ON
L4K 4P8, Canada
Tel. +1.800.263.3390
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca