



ORDER FORM

# Squiggles + Stander

# EILK02

Prone, upright, or supine standing with up to 60° abduction.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

## Squiggles+ Stander

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

### Multi-Position Packages \$

Package includes Squiggles+ Standard Model Multi-Position, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, and Small Sandals with Straps.

<input type="checkbox"/>	141QE15	Squiggles Plus Stander Package - Blue	\$6,642
<input type="checkbox"/>	141QE16	Squiggles Plus Stander Package - Purple	\$6,642
<input type="checkbox"/>	141QE17	Squiggles Plus Stander Package - Red	\$6,642
<input type="checkbox"/>	141QE18	Squiggles Plus Stander Package - Red Vinyl	\$6,642

### Standard Model \$

Selection Required.

Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplates and tray attachment.

<input type="checkbox"/>	141KJ14	Standard Support Shell - Multi-Position (Prone, Supine, and Upright)	\$3,137
<input type="checkbox"/>	141FC6	Squiggles+ Pivot Base	\$1,688

### Positioning Pads \$

Selection Required. Includes chest and hip pad covers and knee support covers.

<input type="checkbox"/>	141ZI11	Squiggles Plus Stander Positioning Pads - Blue	\$699
<input type="checkbox"/>	141ZI12	Squiggles Plus Stander Positioning Pads - Purple	\$699
<input type="checkbox"/>	141ZI13	Squiggles Plus Stander Positioning Pads - Red	\$699
<input type="checkbox"/>	141ZI14	Squiggles Plus Stander Positioning Pads - Red Vinyl	\$699

### Accessories \$

<input type="checkbox"/>	141A208	Tray - Gray	\$396
<input type="checkbox"/>	141A216	Padded Tray Insert - Black	\$326
<input type="checkbox"/>	1418A106	Grab Rail / Toy Bar for Tray	\$321
<input type="checkbox"/>	1418A107	Small Sandals w/ Straps (Pair)	\$269
<input type="checkbox"/>	1418A108	Medium Sandals w/ Straps (Pair)	\$287
<input type="checkbox"/>	1418A110	QTY Sandal Raiser 1" (Each) Requires Sandal Selection.	\$105

### Headrest \$

Cushion selection required with hardware.

Headrest selection NOT available with Prone Configuration.

<input type="checkbox"/>	141AW22	<b>Flat Headrest Support Hardware</b>	\$243
		<b>Flat Headrest</b>	
<input type="checkbox"/>	141Y032	Flat Headrest Cushions - Blue	\$124
<input type="checkbox"/>	141Y033	Flat Headrest Cushions - Purple	\$124
<input type="checkbox"/>	141Y034	Flat Headrest Cushions - Red	\$124
<input type="checkbox"/>	141Y035	Flat Headrest Cushions - Red Vinyl (Pair)	\$124
		<b>Flat Headrest with Laterals</b>	
<input type="checkbox"/>	141Y036	Flat Headrest and Lateral Cushions - Blue	\$210
<input type="checkbox"/>	141Y037	Flat Headrest and Lateral Cushions - Purple	\$210
<input type="checkbox"/>	141Y038	Flat Headrest and Lateral Cushions - Red	\$210
<input type="checkbox"/>	141Y039	Flat Headrest and Lateral Cushions - Red Vinyl (Pair)	\$210

## Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Distance Between Midline of Knees	5.5" - 8.4"	140 mm - 215 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.5" - 31.5"	480 mm - 800 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 mm L x 390 mm W x 390 mm H
Pivot Base Weight	14.3 lbs	6.5 kgs
Pivot Base Dimensions (Unfolded)	37.4" L x 23.2" W x 20.3" H	950 mm L x 590 mm W x 515 mm H
Pivot Base Dimensions (Folded)	26" L x 23.2" W x 20.3" H	660 mm L x 590 mm W x 515 mm H
Pivot Base Angle Range	Prone 90-70° / Supine 90-160°	Prone 90-70° / Supine 90-160°

Sunrise Medical Canada Inc.  
1000 Creditstone Road, Unit 2  
Concord, Ont., L4K 4P8



MK-129855-CAN © 7/2024

