



QUICKIE

HIGHLY ADJUSTABLE  
PERSONALIZING  
YOUR MOBILITY

5<sup>R</sup>

QUICKIE®

5<sup>R</sup>

The QUICKIE® 5<sup>R</sup> is designed to deliver superior rigid wheelchair performance day after day. As you and your environment evolve, or as you refine your riding capabilities, the 5<sup>R</sup>'s frame and components adjust to elevate your mobility.

## OPEN FRAME DESIGN

The QUICKIE 5<sup>R</sup>'s open frame takes your ride quality and geometry to a new level. The swept-back hydro-formed caster design not only provides a smooth ride but is stylish.

The 5<sup>R</sup> offers two frame options. The **Classic** features a traditional frame length for maximum stability. The shortened **Performance** frame length offers maximum responsiveness.

## INNOVATIVE AXLE PLATE DESIGN

The 5<sup>R</sup>'s swept-back axle plate design isn't just for looks. It is also incredibly quick to adjust and is available in two lengths to ensure you are getting the adjustments you need without adding extra weight.






### QUICKIE Xtender Power Assist

Developed by QUICKIE and Yamaha®, the Xtender features motorized rear wheels that increase the force you apply to the handrims by *four times*.



### XTR Suspension

The XTR Suspension Frame helps you fearlessly conquer curbs, bumps, and tough terrain. Add Mountain wheels, Knobby tires, wide casters, and Frog Legs suspension caster forks for the ultimate off-road wheelchair experience.



\*Black Zebra Kolorfusion™ pattern is available on any base color. Shown here on Candy Blue base color.

### Available Colors

The color options shown are for reference only. They are not intended as accurate representations. Refer to the QUICKIE Color Brochure or Quickie Color Ring for representations of Quickie colors.



### Matching-Colored JAY® J3™ Backs

Add a pop of matching or contrasting color to your 5<sup>R</sup> with a painted JAY J3 back. Choose from any of the standard or matte frame colors available on QUICKIE wheelchairs, and we will paint your J3 back at no additional charge.

THE 5<sup>R</sup>  
GIVES YOU  
ACCESS TO THE  
PREMIUM  
OPTIONS  
YOU WANT TO  
PERSONALIZE  
YOUR MOBILITY

5<sup>R</sup>

# COMPLETE OPTION MIX FOR FUNCTION & STYLE

The 5R offers a great mix of rigid wheelchair options, including HCPCS-coded options, handrims, tension-adjustable upholstery, single-sided forks, and colored components.



Plastic Coated



Aluminum Anodized



Natural Fit LT\*\*



The Surge



Tension-Adjustable Upholstery



Single-Sided Fork

\*\*Available with Standard or Super Grip



Our Built-4-Me Team is dedicated to creating product personalized for you, beyond what is listed on our product order forms. Experts in the engineering and manufacturing of personalized rehab products will listen to your specific requests and build a chair as individual as your fingerprint!

#### Personalized design includes:

- Hot pink graduated to white painted frame
- Chrome rims
- Hot pink accents with white sideguards  
*(Shown on Quickie 5R to the left)*

#### 5<sup>R</sup> SPECS

HCPCS Code: K0005

Overall Weight: 21.5 lbs

Transport Weight: 12.9 lbs

Weight Capacity: 265 lbs

Seat Width: 12" to 20"

Seat Depth: 12" to 20"

Center of Gravity Adjustment:

-1.5" to 4"

Back Angle Adjustment:

Fore (-) 17° to Aft (+) 8°

Front Seat Height:

16" to 21"

Rear Seat Height:

13" to 20"

Transit Approved: Yes\*

Warranty:

Limited/Lifetime

\*Conforms with ANSI/RESNA WC-19 standards. Please refer to QUICKIE® & ZIPPIE® WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at [www.SunriseMedical.com](http://www.SunriseMedical.com).



Sunrise Medical (US) LLC, Fresno, CA 93727  
MK-100136 Rev.D © 09.2021

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

