

ELITE^X

SERIAL NO. _____

PROD. DATE: _____

April 2025



Customer Contact Details

Customer Name: _____

Advisor: _____

Address: _____

Phone: _____

E-mail: _____

Shipping Name: Shipping _____

Address: _____

Phone: _____

E-mail: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Guide:

☒ Base product

☐ Required option

☐ Optional Extra



Critical information



Information needed

PRODUCT STYLE

☒ EXAXZ010 ELITE X
\$ 7,637.00

☒ XA600 Heat treated frame


☐ EXAXA062 Standard

☐ EXAXA063 Standard with 25mm seatrail

☐ EXAXA280 Heavy Duty
\$ 1,179.00
(Spinergy R10 wheels also required)

Max User Weight 125kg (275 lbs.)

☐ Yes, I am within the max user weight



CAD APPROVAL

☐ XW012 QUOTE CAD drawing for approval
\$550 FLAT FEE (includes two free revisions)

☐ XW009 QUOTE only no CAD \$0

☐ XW010 ORDER CAD drawing for approval
\$300 FLAT FEE (includes two free revisions)

☐ XW011 ORDER only no CAD \$0

KEY MEASUREMENTS

☒ XY010 Seat width (Inside to inside of sideguards)



mm

(Max - 500mm)

☒ XY020 Seat length



mm

(Max - 500mm)

☐ XY030 Extra seat length (Provides additional leg support & hand grip)



mm

☒ XY050 Front frame height



mm

(Max - 630mm)
(Min - 360mm)

☒ XY060 Rear frame height



mm

(Max - 630mm)
(Min - 360mm)

☒ XY070 Frame length



mm


(Min = XY020 + XY030 + 140mm)

☐ XA020 Ergonomic seat
\$ 121.00



Specify length
mm


(Standard = 120mm)



Front Frame Options


☐ XY080 Seat taper design

i Specify taper width
 mm



☐ XA130 V-seat design
 \$ 307.00


i Specify taper width
 mm



Advanced Front Frame Options

☐ XY082 Seat taper position
 \$ 90.00


i Specify measurement
 mm



☐ XA131 V-seat and taper design
 \$ 331.00


i Specify inside width at front of sideguard (A)
 mm

i Specify taper width (B)
 mm



☐ XA136 Parallel taper design
 \$ 90.00 (Not with Seat Taper XY080)


i Specify taper width
 mm



☐ XA135 Asymmetrical taper design
 \$ 265.00

i Specify taper width (A)
 mm

i Specify offset (B)
 (Inside of sideguard to inside of taper width)
 mm




Seat Options

☐ XR050 Carbon fibre seat plate
 \$ 480.00 (riveted to frame)

i Specify length of plate
 mm (Std XY020 - 20mm)

☐ XA260 Drop seat design
 \$ 219.00 (only with either XR040 or 050)

i Specify drop
 mm (Std = 20mm)



Footrest Type	
☒ EXAXB080 Fixed position footrest	☐ EXAXB145 Multi position footrest tubes \$ 336.00
Footrest Measurements	
☐ Front footrest height mm (min = 35mm)	☐ Footrest width mm (Std = 200mm)
☐ Rear footrest height mm (min = 35mm)	☐ Footrest length mm (Std = 250mm)
☐ Footrest position (front of footplate to castor pot) mm	
Footrest Plate	
☐ XB190 Carbon fibre footplate \$560.00	☐ XB100 ABS moulded foot tray \$ 183.00
☐ EXAXB140 Spacer kit \$45.00 (Allows for 25mm increase in height)	Specify length mm Max = 380mm (Full) Std = 330mm
Footrest Plate Options	
☐ XB240 Aluminium raised edges \$188.00 (Cost per Edge)	
☐ Front ☐ Back ☐ Left ☐ Right	
Specify height (Max = 60mm)	
Caster Pin Width	
☒ XY330 Caster pin width mm	
☐ Optimize measurements	
CAMBER BAR, WING AND ANTI-TIP	
Centre of Gravity and Camber Angle	
☒ EXAXY510 Fixed centre of gravity Specify measurement mm (Note: dimension is taken from front of backrest to centre of axle)	☒ XY390 Camber angle
Wing Options	
☐ EXAXA110 Round design basketball wing 	
☐ EXAXA120 Straight sided design basketball wing \$ 376.00 	
Anti-tip options	
☐ EXAXJ040 Single fixed anti tip	☐ EXAXJ060 Double fixed anti tip \$ 564.00

PROD DATE: _____ / _____
CASTER AND REAR WHEELS

Caster Options

- | | |
|--|---|
| ☐ XC160 Sport 61 mm caster
(Black) | ☐ XC181 Bones Reds bearings
(All caster bearings) |
| ☐ XC200 Sport micro fork | ☐ XC182 Bones Swiss Ceramic bearings
(All caster bearings)
\$ 601.00 |

Rear Wheel Size

- | | | | |
|---|---|---|--|
| ☐ XY420 24" Wheels | ☐ XY430 25" Wheels | ☐ XY440 26" Wheels | ☐ XY450 700C Wheels |
|---|---|---|--|

Rear Wheel Gap

☐ XY460 Rear wheel gap  15202530 Individual _____ mm	OR	☐ XY470 Inside to inside of wheel  _____ mm
---	----	---

Rear Wheel Options

- | | |
|---|--|
| ☐ XK080 Spinergy SLX 24 spoke
\$ 1006.00 | ☐ XK330 Spinergy SLX R10 heavy duty 24 spoke
(18 deg camber angle only)
\$ 1,859.00 |
| ☐ XK101 Spinergy XSLX 24 spoke + Caps
\$ 1,367.00 | ☐ XK521 Spinergy Ceramic bearing upgrade
\$ 551.00 |
- Spinergy Colour Options:** 
- | | | | |
|--|--|--|---|
| ☐ Black spokes | ☐ White spokes | ☐ Yellow spokes | ☐ XK180 Individual spoke design
(Specify colour)
\$ 213.00 |
| ☐ Red spokes | ☐ Blue spokes | ☐ Pink spokes | Inside: _____ |
| ☐ Purple spokes | ☐ Orange spokes | ☐ Green spokes | Outside: _____ |
- Spinergy Anodize Options**
- | | | | |
|--|--|--|--|
| ☐ XK190 Black hub (Std) | ☐ XK200 Silver hub
\$ 127.00 | ☐ XK201 Red hub/rim
\$ 409.00 | ☐ XK205 Individual hub/rim
\$ 568.00
Select colour: _____ |
| ☐ XK202 Purple hub/rim
\$ 409.00 | ☐ XK203 Blue hub/rim
\$ 409.00 | ☐ XK204 Custom hub/rim
\$ 539.00 | |

Supply Without Wheels

- | | | |
|--|--|---|
| ☐ XK040 Supply without Wheels |  Wheel Type: <div style="border: 1px solid black; width: 150px; height: 20px;"></div> | Tire Type: <div style="border: 1px solid black; width: 150px; height: 20px;"></div> |
|--|--|---|

Handrim Position / Options

- | | |
|---|--|
| ☐ XY480 Narrow set handrim | ☐ XY485 Wide set handrim |
| ☐ XL090 Remove excess mounting tabs
(recommended for small wheel gaps)
\$ 144.00 | ☐ XL080 Super narrow set handrim
(only with XL040 and includes XL090)
\$ 290.00 |

Handrim Style

- | | |
|---|---|
| ☐ XL010 Aluminium hard anodised 6 tab | ☐ XL040 Titanium 6 tab
\$ 446.00
(Not with 700C) |
| ☐ XL070 Aluminium hard anodised 12 tab
\$ 146.00
(24" not available) | ☐ XL050 Stainless steel 6 tab
\$ 248.00 |

Tire Options

- | | | | |
|--|------------------------------|-------------------------------|---|
| ☐ XM070 Kenda Kaliente | ☐ Red | ☐ Blue | ☐ XM080 Tufo MS3 tubeless
\$81.00
(Blue only) |
| ☐ XM150 IRC Exeracer Pro
\$ 188.00
(Red Only) | | | ☐ XM090 Tufo WS3 tubeless
\$81.00
(24" only and Blue only) |

Spare Wheel

- | | |
|---|--|
| ☐ EXAXK500 Spare wheel Spinergy
\$ 2,258.00
(All inclusive + axle pin) | ☐ XK520 Spare wheel Titanium / Super narrow upcharge
\$ 634.00 |
|---|--|

Inner Tube Valve Options

- | | |
|---|--|
| ☐ XY500 Schraeder inner tube valve | ☐ XY490 Presta inner tube valve
(Standard if no option selected) |
|---|--|

FINISH OPTIONS

Frame finish

☐ XN010 Powder coated frame

Frame Gloss

- ☐ Black gloss 103
- ☐ White gloss 100
- ☐ Grey gloss 132
- ☐ Cream gloss 135
- ☐ Red gloss 108
- ☐ Yellow gloss 120
- ☐ Forest green gloss 143
- ☐ Green gloss 107
- ☐ Dark blue gloss 151
- ☐ Royal blue gloss 154
- ☐ Blue gloss 174
- ☐ Duckegg blue gloss 141
- ☐ Pink gloss 188
- ☐ Orange gloss 180
- ☐ Purple gloss 157

Frame Matt

- ☐ Black matt 105
- ☐ White matt 101
- ☐ Grey matt 131
- ☐ Light grey matt 133
- ☐ Navy blue matt 152

Frame Special

- ☐ Chrome Mirror 300
- ☐ Sparkle silver 170
- ☐ Metallic black matt 104
- ☐ Metallic grey matt 130
- ☐ Black sandpaper 201
- ☐ Grey sandpaper 202
- ☐ White sandpaper 200

☐ RGK-XN200 Individual paint colour (State Colour/RAL code)
\$ 565.00 **INDIVIDUAL**

☐ EXAXN030 Brushed frame
\$ 741.00

☐ EXAXN050 Polished frame
\$ 1,938.00

Accent Colour

Note: Accent colour includes, forks, caster caps, decals, ELITE logo and stitching. Cushion, toe strap, frame protectors and RGK logo are white.

- ☐ Black accent

☐ Red accent

☐ Purple accent

☐ Silver accent

☐ Blue accent

☐ Carbon accent (caps only, black for others)

☐ RGK-XN250 Individual accent
\$ 315.00

Select Colour:

UPHOLSTERY (N/A WITH CX SEAT)

Seat Upholstery

☐ EXAXP040 Seat sling strap adjustable
(Not required with seat plate)

Backrest Upholstery

☒ EXAXP055 Airtech sport backrest strap adjustable
(comes with Ti buckles)

☐ XP070 Backrest upgrade to leather corners
\$ 121.00 (Black)

Backrest Stitching Options

- | | | | |
|---------------------------------|-------------------------------------|---------------------------------|--------------------------------------|
| ☐ Silver | ☐ Black | ☐ Blue | ☐ Individual: |
| ☐ Red | ☐ Light blue | ☐ Purple | _____ |

Backrest RGK Logo Colour

- | | | | |
|--------------------------------------|-------------------------------------|---------------------------------|--------------------------------------|
| ☐ White (Std) | ☐ Green | ☐ Blue | ☐ Individual: |
| ☐ Red | ☐ Light blue | ☐ Purple | _____ |

Backrest Accessories

☐ XP170 Individual backrest embroidery
\$ 240.00 (Please provide .jpeg file and positioning)

CUSHION

Cushion foam options

- ☐ XP080 Cushion foam
\$ 121.00
(select density and size below)
- ☐ Soft (CMS040) ☐ Medium (CMS050) ☐ Hard (9lb recon)
- ☐ 50mm ☐ 75mm ☐ 100mm ☐ Other
- ☐ XP340 Plastazote cushion upgrade
\$ 115.00
(Indicate max height)
- ☐ 580mm ☐ 630mm

Cushion Cover Upgrades

- ☐ XP081 Nylon cushion cover
- ☐ XP090 Air Lite cushion cover (Black)
\$ 98.00
- ☐ XP380 Neoprene cushion cover (Black)
\$ 127.00


STRAPS

Velcro straps


- ☐ XS050 100mm Pro adjustable velcro strap
\$ 173.00
- ☐ XS055 50mm Pro adjustable velcro strap
\$ 135.00
- ☐ XS057 50mm Pro adjustable velcro **Knee** strap
\$ 250.00
(measurement from tube to tube over knees)

Velcro strap position

☐ Strap 1




Specify strap length




Specify distance from back post

☐ Strap 1



Specify strap length




Specify distance from back post

Ratchet Straps

- ☐ XS020 Sport ratchet strap (Click strap)
\$ 140.00
- ☐ XS025 Pelvi Loc 2DA ratchet strap
\$ 406.00
(For thighs and knees)
- ☐ XS030 Pelvi Loc 3DA ratchet strap
\$ 419.00
(For hips)
- ☐ XS035 M2 ratchet strap
\$ 283.00
(For thighs and knees)
- Select size: ☐ Medium (8") ☐ Large (10")
- Select size: ☐ Small (8") ☐ Medium (9") ☐ Large (10") ☐ Extra Large (11")
- ☐ 8" ☐ 10" ☐ 12"


Ratchet strap position

☐ Strap 1



Specify distance from back post

☐ Strap 2



Specify distance from back post

☐ XS040 Fixed mounting points
\$ 158.00
(one set per click strap)

☐ x1 ☐ x2 ☐ x3

FRAME OPTIONS

- ☐ XR060 Calf strap mounting loops
\$ 160.00
- ☐ XR065 Strap location loops
\$ 118.00

ACCESSORIES

- ☐ XR030 Triple wheel bag (Black)
\$ 201.00
- ☐ XP130 Frame protectors full length
\$ 127.00
- ☐ XA150 Plastic tube fitted to offensive wing
\$ 93.00
- ☐ XS070 Foot strap (select if required)
- ☐ XS060 Calf strap (select if required)
- ☐ XP400 Back bar protector
- ☐ XR090 Tool kit