

**ELITE X**

SERIAL NO. \_\_\_\_\_

PROD. DATE: \_\_\_\_\_

February 2024



**Customer Contact Details**

<p>Customer Name: _____</p> <p>Advisor: _____</p> <p>Address: _____</p> <p>_____</p> <p>_____</p> <p>Phone: _____</p> <p>E-mail: _____</p>	<p>Shipping Name: Shipping _____</p> <p>Address: _____</p> <p>_____</p> <p>_____</p> <p>Phone: _____</p> <p>E-mail: _____</p>
--	---


**Additional Information**

<p>Notes: _____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Current Wheelchair: _____</p> <p>Customer Weight: _____</p> <p>Customer Height: _____</p> <p>Customer Disability: _____</p>
--	--

Guide:  Base product     Required option     Optional Extra    Critical information    Information needed

### PRODUCT STYLE

- EXAXZ010 ELITE X**  
\$ **7,524.00**
- XA600 Heat treated frame
- EXAXA062 Standard
- EXAXA063 Standard with 25mm seatrail
- EXAXA280 Heavy Duty  
\$ 1,179.00 (Spinergy R10 wheels also required)



**Max User Weight 125kg (275 lbs.)**


Yes, I am within the max user weight

### CAD APPROVAL

- XW012 QUOTE CAD drawing for approval \$550 FLAT FEE (includes two free revisions)
- XW009 QUOTE only no CAD \$0
- XW010 ORDER CAD drawing for approval \$300 FLAT FEE (includes two free revisions)
- XW011 ORDER only no CAD \$0

### KEY MEASUREMENTS

- XY010 Seat width (Inside to inside of sideguards)  
  (Max - 500mm)
- XY020 Seat length  
  (Max - 500mm)
- XY030 Extra seat length (Provides additional leg support & hand grip)
- XY050 Front frame height  
  (Max - 630mm) (Min - 360mm)
- XY060 Rear frame height  
  (Max - 630mm) (Min - 360mm)
- XY070 Frame length  
  (Min = XY020 + XY030 + 140mm)
- XA020 Ergonomic seat \$ 119.00  
 Specify length  (Standard = 120mm)




PROD DATE: \_\_\_\_ / \_\_\_\_ / \_\_\_\_  
**FRONT FRAME AND SEAT**

Front Frame Options


XY080 Seat taper design

**i** Specify taper width  
 mm



XA130 V-seat design  
 \$ 302.00


**i** Specify taper width  
 mm



Advanced Front Frame Options

XY082 Seat taper position  
 \$ 89.00


**i** Specify measurement  
 mm



XA131 V-seat and taper design  
 \$ 326.00


**i** Specify inside width at front of sideguard (A)  
 mm

**i** Specify taper width (B)  
 mm



XA136 Parallel taper design  
 \$ 89.00 (Not with Seat Taper XY080)


**i** Specify taper width  
 mm



XA135 Asymmetrical taper design  
 \$ 261.00

**i** Specify taper width (A)  
 mm

**i** Specify offset (B)  
 (Inside of sideguard to inside of taper width)  
 mm




Seat Options

XR050 Carbon fibre seat plate  
 \$ 473.00 (riveted to frame)

**i** Specify length of plate  
 mm (Std XY020 - 20mm)

XA260 Drop seat design  
 \$ 216.00 (only with either XR040 or 050)

**i** Specify drop  
 mm (Std = 20mm)



Footrest Type	
<input checked="" type="checkbox"/> <b>EXAXB080</b> Fixed position footrest	<input type="checkbox"/> <b>EXAXB145</b> Multi position footrest tubes \$ 331.00
Footrest Measurements	
<input type="checkbox"/> Front footrest height <input style="width: 100px;" type="text"/> mm (min = 35mm)	<input type="checkbox"/> Footrest width <input style="width: 100px;" type="text"/> mm (Std = 200mm)
<input type="checkbox"/> Rear footrest height <input style="width: 100px;" type="text"/> mm (min = 35mm)	<input type="checkbox"/> Footrest length <input style="width: 100px;" type="text"/> mm (Std = 250mm)
<input type="checkbox"/> Footrest position (front of footplate to castor pot) <input style="width: 100px;" type="text"/> mm	
Footrest Plate	
<input type="checkbox"/> <b>XB190</b> Carbon fibre footplate \$552.00	<input type="checkbox"/> <b>XB100</b> ABS moulded foot tray \$ 180.00 Specify length <input style="width: 100px;" type="text"/> mm (Max = 380mm (Full) Std = 330mm)
<input type="checkbox"/> <b>EXAXB140</b> Spacer kit \$44.00 (Allows for 25mm increase in height)	
Footrest Plate Options	
<input type="checkbox"/> <b>XB240</b> Aluminium raised edges \$185.00 (Cost per Edge) <input type="checkbox"/> Front <input type="checkbox"/> Back <input type="checkbox"/> Left <input type="checkbox"/> Right Specify height <input style="width: 100px;" type="text"/> (Max = 60mm)	
Caster Pin Width	
<input checked="" type="checkbox"/> <b>XY330</b> Caster pin width <input style="width: 100px;" type="text"/> mm <input type="checkbox"/> Optimize measurements	
<b>CAMBER BAR, WING AND ANTI-TIP</b>	
Centre of Gravity and Camber Angle	
<input checked="" type="checkbox"/> <b>EXAXY510</b> Fixed centre of gravity Specify measurement <input style="width: 100px;" type="text"/> mm (Note: dimension is taken from front of backrest to centre of axle)	<input checked="" type="checkbox"/> <b>XY390</b> Camber angle <input type="checkbox"/> 16 <input type="checkbox"/> 18 <input type="checkbox"/> 20 <input type="checkbox"/> 22
Wing Options	
<input type="checkbox"/> <b>EXAXA110</b> Round design basketball wing 	<input type="checkbox"/> <b>EXAXA120</b> Straight sided design basketball wing \$ 370.00 
Anti-tip options	
<input type="checkbox"/> <b>EXAXJ040</b> Single fixed anti tip	<input type="checkbox"/> <b>EXAXJ060</b> Double fixed anti tip \$ 556.00

PROD DATE: \_\_\_\_\_ / \_\_\_\_\_ / \_\_\_\_\_  
**SIDEGUARDS**  
 Sideguard Style

EXAXE230 Welded aluminium sport  
 \$ 914.00

XE310 High design upgrade  
 \$ 326.00 *(Only with option ELAXE230)*




EXAXE250 Welded aluminium sport with Carbon Fibre  
 \$ 1,607.00

Sideguard Measurements

Sideguard length

mm



Sideguard height (back)

mm



Sideguard height (front)  
*(optional)*

mm



Follow Tire Design

25mm above tire


Individual above tire

mm

Level with tire



Individual Shape



**BACKREST**


Backrest Height

XY230 Fixed backrest height



Specify height

mm




Ergonomic Options

XY360 Backrest angle



Individual  Deg  
 (Max -10/+10 Deg)





XG100 Ergonomic backrest  
 \$ 144.00



Vertical height  mm




Rearward position  mm



XG110 Inward taper  
 \$ 144.00  
*(inside of backrest posts)*




Specify width / per side  mm



XG120 Outward flare  
 \$ 144.00  
*(inside of backrest posts)*



Specify width / per side  mm



PROD DATE: \_\_\_\_\_ / \_\_\_\_\_ / \_\_\_\_\_  
**CASTER AND REAR WHEELS**

Caster Options			
<input type="checkbox"/>	XC160 Sport 61mm caster <i>(Black)</i>	<input type="checkbox"/>	XC181 Bones Reds bearings <i>(All caster bearings)</i>
<input type="checkbox"/>	XC200 Sport micro fork	<input type="checkbox"/>	XC182 Bones Swiss Ceramic bearings <i>(All caster bearings)</i> \$ 592.00
Rear Wheel Size			
<input type="checkbox"/>	XY420 24" Wheels	<input type="checkbox"/>	XY430 25" Wheels
		<input type="checkbox"/>	XY440 26" Wheels
		<input type="checkbox"/>	XY450 700C Wheels
Rear Wheel Gap			
<input type="checkbox"/>	XY460 Rear wheel gap		<input type="checkbox"/>
		OR	<input type="checkbox"/>
	Individual <input style="width: 100px;" type="text"/> mm		XY470 Inside to inside of wheel
Rear Wheel Options			
<input type="checkbox"/>	KX080 Spinerger SLX 24 spoke \$ 991.00	<input type="checkbox"/>	KX330 Spinerger SLX R10 heavy duty 24 spoke <i>(18 deg camber angle only)</i> \$ 1,832.00
<input type="checkbox"/>	KX101 Spinerger XSLX 24 spoke + Caps \$ 1,347.00	<input type="checkbox"/>	KX521 Spinerger Ceramic bearing upgrade \$ 543.00
<b>Spinerger Colour Options:</b>			
<input type="checkbox"/>	Black spokes	<input type="checkbox"/>	White spokes
<input type="checkbox"/>	Red spokes	<input type="checkbox"/>	Blue spokes
<input type="checkbox"/>	Purple spokes	<input type="checkbox"/>	Orange spokes
		<input type="checkbox"/>	Yellow spokes
		<input type="checkbox"/>	Pink spokes
		<input type="checkbox"/>	Green spokes
<b>Spinerger Anodize Options</b>			
<input type="checkbox"/>	KX190 Black hub (Std)	<input type="checkbox"/>	KX200 Silver hub \$ 125.00
<input type="checkbox"/>	KX202 Purple hub/ribs \$ 403.00	<input type="checkbox"/>	KX203 Blue hub/ribs \$ 403.00
		<input type="checkbox"/>	KX201 Red hub/ribs \$ 403.00
		<input type="checkbox"/>	KX204 Custom hub/ribs \$ 531.00
		<input type="checkbox"/>	KX205 Individual hub/ribs Select colour: \$ 560.00
<b>Spinerger Inside/Outside:</b> Inside: _____ Outside: _____			
Supply Without Wheels			
<input type="checkbox"/>	KX040 Supply without Wheels		Wheel Type: <input style="width: 150px;" type="text"/>
			Tire Type: <input style="width: 150px;" type="text"/>
Handrim Position / Options			
<input type="checkbox"/>	XY480 Narrow set handrim	<input type="checkbox"/>	XY485 Wide set handrim
<input type="checkbox"/>	XL090 Remove excess mounting tabs <i>(recommended for small wheel gaps)</i> \$ 144.00	<input type="checkbox"/>	XL080 Super narrow set handrim <i>(only with XL040 and includes XL090)</i> \$ 290.00
Handrim Style			
<input type="checkbox"/>	XL010 Aluminium hard anodised 6 tab	<input type="checkbox"/>	XL040 Titanium 6 tab <i>(Not with 700C)</i> \$ 439.00
<input type="checkbox"/>	XL070 Aluminium hard anodised 12 tab \$ 144.00	<input type="checkbox"/>	XL050 Stainless steel 6 tab \$ 244.00
			<i>(24" not available)</i>
Tire Options			
<input type="checkbox"/>	XM070 Kenda Kaliente <input type="checkbox"/> Red <input type="checkbox"/> Blue	<input type="checkbox"/>	XM080 Tufo MS3 tubeless \$80.00 <i>(Blue only)</i>
<input type="checkbox"/>	XM150 IRC Exeracer Pro <i>(Red Only)</i> \$ 185.00	<input type="checkbox"/>	XM090 Tufo WS3 tubeless \$80.00 <i>(24" only and Blue only)</i>
Spare Wheel			
<input type="checkbox"/>	EXAXK500 Spare wheel Spinerger \$ 2,225.00	<i>(All inclusive + axle pin)</i>	<input type="checkbox"/>
			KX520 Spare wheel Titanium / Super narrow upcharge \$ 625.00
Inner Tube Valve Options			
<input type="checkbox"/>	XY500 Schraeder inner tube valve		<input type="checkbox"/>
			XY490 Presta inner tube valve <i>(Standard if no option selected)</i> 

## FINISH OPTIONS

Frame finish

 XN010 Powder coated frame

 Frame Gloss

- Black gloss 103
- White gloss 100
- Grey gloss 132
- Cream gloss 135
- Red gloss 108
- Yellow gloss 120
- Forest green gloss 143
- Green gloss 107
- Dark blue gloss 151
- Royal blue gloss 154
- Blue gloss 174
- Duckegg blue gloss 141
- Pink gloss 188
- Orange gloss 180
- Purple gloss 157

 Frame Matt

- Black matt 105
- White matt 101
- Grey matt 131
- Light grey matt 133
- Navy blue matt 152

 Frame Special

- Chrome Mirror 300
- Sparkle silver 170
- Metallic black matt 104
- Metallic grey matt 130
- Black sandpaper 201
- Grey sandpaper 202
- White sandpaper 200

 RGK-XN200 Individual paint colour (State Colour/RAL code)  
 \$ 557.00 **INDIVIDUAL**

\_\_\_\_\_

 EXAXN030 Brushed frame  
 \$ 730.00

 EXAXN050 Polished frame  
 \$ 1,909.00

### Accent Colour

**Note:** Accent colour includes, forks, caster caps, decals, ELITE logo and stitching. Cushion, toe strap, frame protectors and RGK logo are white.

- |  |  |   |
|--|--|---|
| <input type="checkbox"/> Black accent  | <input type="checkbox"/> Silver accent                               | <input type="checkbox"/> RGK-XN250 Individual accent<br>\$ 310.00 |
| <input type="checkbox"/> Red accent    | <input type="checkbox"/> Blue accent                                 | Select Colour:  |
| <input type="checkbox"/> Purple accent | <input type="checkbox"/> Carbon accent (caps only, black for others) |   |

## UPHOLSTERY (N/A WITH CX SEAT)

### Seat Upholstery

 EXAXP040 Seat sling strap adjustable  
 (Not required with seat plate)

### Backrest Upholstery

 EXAXP055 Airtech sport backrest strap adjustable  
 (comes with Ti buckles)

 XP070 Backrest upgrade to leather corners  
 \$ 119.00 (Black)

### Backrest Stitching Options

- |                                 |                                     |                                 |  |
|---------------------------------|-------------------------------------|---------------------------------|--|
| <input type="checkbox"/> Silver | <input type="checkbox"/> Black      | <input type="checkbox"/> Blue   | <input type="checkbox"/> Individual: _____ |
| <input type="checkbox"/> Red    | <input type="checkbox"/> Light blue | <input type="checkbox"/> Purple |  |

### Backrest RGK Logo Colour

- |                                      |                                     |                                 |  |
|--------------------------------------|-------------------------------------|---------------------------------|--|
| <input type="checkbox"/> White (Std) | <input type="checkbox"/> Green      | <input type="checkbox"/> Blue   | <input type="checkbox"/> Individual: _____ |
| <input type="checkbox"/> Red         | <input type="checkbox"/> Light blue | <input type="checkbox"/> Purple |  |

### Backrest Accessories

 XP170 Individual backrest embroidery  
 \$ 236.00 (Please provide .jpeg file and positioning)

## CUSHION

Cushion foam options

- XP080** Cushion foam  
 \$ 119.00  
*(select density and size below)*
- Soft (CMS040)     Medium (CMS050)     Hard (9lb recon)  
 50mm     75mm     100mm     Other  mm
- XP340** Plastazote cushion upgrade  
 \$ 113.00  
*(Indicate max height)*
- 580mm     630mm

## Cushion Cover Upgrades

- XP081** Nylon cushion cover
- XP090** Air Lite cushion cover (Black)  
 \$ 97.00
- XP380** Neoprene cushion cover (Black)  
 \$ 125.00

## STRAPS

### Velcro straps

- XS050** 100mm Pro adjustable velcro strap  
 \$ 170.00
- XS055** 50mm Pro adjustable velcro strap  
 \$ 133.00
- XS057** 50mm Pro adjustable velcro **Knee** strap  
 \$ 246.00  
*(measurement from tube to tube over knees)*

### Velcro strap position

**Strap 1**

Specify strap length

 mm

Specify distance from back post

 mm

**Strap 1**

Specify strap length

 mm

Specify distance from back post

 mm

### Ratchet Straps

- XS020** Sport ratchet strap (Click strap)  
 \$ 138.00
- XS025** Pelvi Loc 2DA ratchet strap  
 \$ 400.00  
*(For thighs and knees)*
- Select size:*  
 Medium (8")     Large (10")
- XS030** Pelvi Loc 3DA ratchet strap  
 \$ 413.00  
*(For hips)*
- Select size:*  
 Small (8")     Medium (9")  
 Large (10")     Extra Large (11")
- XS035** M2 ratchet strap  
 \$ 279.00  
*(For thighs and knees)*
- 8"     10"     12"

### Ratchet strap position

**Strap 1**

Specify distance from back post

 mm

**Strap 2**

Specify distance from back post

 mm

- XS040** Fixed mounting points  
 \$ 158.00  
*(one set per click strap)*
- x1     x2     x3

## FRAME OPTIONS

- XR060** Calf strap mounting loops  
 \$ 152.00
- XR065** Strap location loops  
 \$ 116.00

## ACCESSORIES

- XR030** Triple wheel bag (Black)  
 \$ 198.00
- XP130** Frame protectors full length  
 \$ 125.00
- XA150** Plastic tube fitted to offensive wing  
 \$ 92.00
- XS070** Foot strap (select if required)
- XS060** Calf strap (select if required)
- XP400** Back bar protector
- XR090** Tool kit