



Serial No:



ORDER FORM			
Customer Dealer Information		Shipping Address	Identical To Billing Address
Dealer		Client	
Client		Address	
Address			
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Information	
Purchase Order No.:		Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):		User Disability:	
		User Weight:	kg User height
build4me			

Product Style

EXAXZ010	Elite X	\$7,866.00
EXAXA062	Standard Frame	Std
EXAXA063	Standard Frame with 25mm seat rail	\$134.00
XA600	Heat treated frame	No Cost

Max User Weight 125kg (275lbs)

PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM

CAD Approval

XW009	Quote only (No CAD For Approval)	
XW010	Order CAD Drawing For Approval (Includes 2 CAD Amendments)	\$300.00
XW011	Order (No CAD Required)	
XW012	Quote & CAD For Approval (Includes 2 CAD Amendments)	\$550.00

Key Measurements

XY010	Seat Width (Inside to inside of sideguards)	
	mm (Max - 500mm)	
XY020	Seat Length	
	mm (Max - 500mm)	
XY030	Extra Seat Length (Provides additional leg support & hand grip)	
	mm	
XY050	Front frame height	(Min - 360mm)
	mm	(Max - 630mm)
XY060	Rear frame height	(Min - 360mm)
	mm	(Max - 630mm)
XY070	Frame Length	
	mm	(Min = XY020 + XY030 + 140mm)
XA020	Ergonomic seat	\$125.00
	Specify Length	
	mm	(Standard = 120mm)

Front Frame and Seat

Front Frame Options

XY080	Seat taper design	XA130	V-Seat design	\$316.00
Specify taper width		Specify taper width		
mm		mm		
		NOTE: Inside to inside wheel (XY470) required		

Advanced Front Frame Options

XY082	Seat taper position			\$93.00
Specify measurement				
mm				
XA131	V-seat and taper design			\$341.00
Specify inside width at front of sideguard (A)	Specify taper width (B)			
mm	mm			
	NOTE: Inside to inside wheel (XY470) required			
XA136	Parallel taper design (Not with seat taper XY080)			\$93.00
Specify taper width				
mm				
XA135	Asymmetrical taper design			\$265.00
Specify taper width (A)	Specify offset (B) (inside of sideguard to inside of taper width)			
mm	mm			
Specify offset taper side				
Left Right				

Seat Options

XR050	Carbon fibre seat plate (riveted to frame)	\$494.00
Specify length of plate	(Std XY020 - 20mm)	
mm		
XA260	Drop seat design (Only with either XR040 or 050)	\$226.00
Specify drop		
mm	(Std = 20mm)	

Footrest

EXAXB080	Fixed position footrest	Std
EXAXB145	Multi position footrest tubes	\$346.00

Front footrest height

mm (min = 35mm)

Rear footrest height

mm (min = 35mm)

Footrest position

(front of footplate to castor pot)

Footrest width

mm (Std= 200mm)

Footrest length

mm (Std= 250mm)

Fixed Footrest Plate Options

XB190	Carbon fibre footplate	\$577.00
XB100	ABS Moulded foot tray	\$188.00
	Specify length	
	Max = 380mm (Full)	
	mm Std = 330mm	
EXAXB140	Spacer kit	\$46.00
	(Allows for 25mm increase in height)	

Additional Footrest Plate Options

XB240	Aluminium raised edges	\$194.00
	(Cost per edge)	
	Specify height	Specify Depth
	Front	(Max = 60mm)
	Back	
	Left	
	Right	
XB250	Carbon Fibre raised edge	\$1,391.00
	(1 side only)	
	Specify height	
	Front	(Max = 45mm)
	Back	

Castor pin width

XY330	Castor pin width
	mm MAX 320mm
	Optimise measurements

Camber Bar, Wing and Anti-Tip


EXAXY510	Fixed centre of gravity	No cost
	Specify measurement	(Note: Dimension is taken from front of backrest to centre of axle)
	mm	

Specify Camber Angle

XY390	Camber angle
	16 18 20 22

Wing Options

EXAXA110	Round design basketball wing	No cost
EXAXA120	Straight sided design basketball wing	\$387.00
EXAXA125	Arrow design basketball wing	\$320.00
EXAXA129	Individual design basketball wing	\$415.00



Anti-Tip options

EXAXJ040	Single fixed anti-tip	No cost
EXAXJ060	Double fixed anti-tip	\$581.00

Sideguards

EXAXE230	Welded aluminium sport	\$956.00
	Please Note: *if the sideguard is 50mm or above the tyre, rear wheel gap will be taken from the plate (unless specified differently). *if the sideguard is below 50mm, rear wheel gap will be taken from the tube of welded sideguard	
EXAXE250	Welded aluminium sport with carbon fibre	\$1,680.00
XE310	High design upgrade	\$341.00

Sideguard Design (Choose one)

Straight design

Round design

Follow tyre

Sideguard Measurement

Sideguard length

mm

Follow tyre design

25mm above tyre

Sideguard height (back)

mm

Individual above tyre

mm

Sideguard height (front)

(Optional for max height players)

mm

Level with tyre

Backrest Height

XY230 Fixed backrest height

Specify height

mm

Ergonomic Options

XY360 Backrest angle

-5 -3 0 +3 +5

Individual

Deg

(Max -10/+10 Deg)

XG100

Ergonomic backrest

\$150.00

Vertical height

mm

Rearward position

mm

XG110 Inward taper

\$150.00

Specify width or per side

mm

XG120 Outward flare

\$150.00

Specify width or per side

mm

Castor Options

XC160	Sport 61mm Castor (Black)	XC181	Bones Reds bearings (All castor bearings)	
XC200	Sport micro fork	XC182	Bones Swiss Ceramic bearings (All castor bearings)	\$619.00

Rear Wheels

Rear Wheel Size

XY420	24" Wheels	XY430	25" Wheels	XY440	26" Wheels	XY450	700C Wheels
-------	------------	-------	------------	-------	------------	-------	-------------

Rear Wheel Gap

XY460	Rear Wheel Gap	XY470	Inside to inside of wheel
15	20	25	30
Individual	mm		mm

Rear Wheel Options

XK080	Spinergy SLX 24 Spoke						\$1,036.00
XK101	Spinergy XSLX 24 Spoke + Caps						\$1,408.00
XK330	Spinergy SLX R10 heavy duty 24 Spoke (18 deg camber angle only)						\$1,915.00
XK521	Spinergy Ceramic bearing upgrade						\$568.00
XK040	Supply without wheels & axle pins	Wheel type:		Tyre type:			NC

Spinergy Colour Options

Black Spokes	White Spokes	Yellow Spokes	Red Spokes	Green Spokes
Blue Spokes	Pink Spokes	Purple Spokes	Orange Spokes	
XK180	Individual spoke design (Specify colour)	Inside: Outside:		
				\$219.00

Spinergy Anodise Options

XK190	Black hub (Std option)	XK200	Silver hub	\$421.00
XK201	Red hub/rims			\$421.00
XK202	Purple hub/rims			\$421.00
XK203	Blue hub/rims			\$421.00
XK204	Gold hub/rims			\$555.00

Spare Wheel

EXAXK500	Spare wheel spinergy (Built up + axle pin)	\$2,225.00
XK520	Spare wheel upcharge - Titanium rims	\$634.00
XK522	Spare wheel upcharge - Remove excess tabs	\$123.00
XK526	Spare wheel upcharge - R10 Hub	\$635.00
XK524	Spare wheel upcharge - Individual spoke colour	\$209.00

Handrim Position/Options

XY480	Narrow set handrim	
XY485	Wide set handrim	
XL090	Remove excess mounting tabs (recommended for small wheel gaps)	\$150.00
XL080	Super narrow set handrim (only with XL010 & XL070 and includes XL090)	\$294.00

Handrim Style

XL010	Aluminium hard anodised 6 tab	Std
XL040	Titanium 6 tab (Not with 700C)	\$459.00
XL070	Aluminium hard anodised 12 tab (24" not available)	\$150.00
XL050	Stainless steel 6 tab	\$255.00

Tyre Options

XM070	Kenda Kaliente	Red	Blue	Std
XM150	IRC Exeracer Pro	(Red only)		\$194.00
XM080	Tufo MS3 tubeless	(Blue only)		\$83.00
XM095	Kenda Wide	(Grey only)		\$80.00
XM096	Schwable Dunk King	(Not available in 700c)		\$250.00

Inner Tube Valve Options

XY500	Schraeder inner tube valve	No Cost
XY490	Presta inner tube valve (Standard if no option selected)	Std

Finish Options

XN010 Powder coated frame

No Cost

Frame

Gloss

Black gloss 103
White gloss 100
Grey gloss 132
Cream gloss 135
Red gloss 108
Yellow gloss 120
Forest green gloss 143
Green gloss 107
Dark blue gloss 151
Royal blue gloss 154
Blue gloss 174
Duckegg Blue gloss 141
Pink gloss 188
Orange gloss 180
Purple gloss 157

Frame

Matt

Black matt 105
White matt 101
Grey matt 131
Light grey matt 133
Navy blue matt 152

Frame

Special

Chrome mirror 300
Sparkle silver 170
Metallic black matt 104
Metallic grey matt 130
Black sandpaper 201
Grey sandpaper 202
White sandpaper 200

XN150 Powder coated components

(Frame colour if not specified)

No Cost

Specify component colour code from above list

EXAXN030

Brushed Frame

(All parts black unless XN160 chosen)

\$763.00

XN160

Brushed Components

\$125.00

EXAXN050

Polished Frame

(All parts black unless XN170 chosen)

\$1,99.00

XN170

Polished Components

\$219.00

Accent Colour

Note: Accent colour includes forks, castor caps, decals, ELITE logo & stitching. Cushion, toe strap, frame protectors and RGK Logo

Black accent

Silver accent

Red accent

Blue accent

Purple accent

Carbon accent

(caps only, black for others)

Upholstery

Seat upholstery

EXAXP040	Seat sling strap adjustable (Std)	No seat sling (Seat Plate)	No Cost
----------	-----------------------------------	----------------------------	---------

Backrest upholstery

EXAXP055	Airtech sport backrest strap adjustable (comes with Ti buckles)	XP070 Backrest upgrade Leather Corners	\$125.00
----------	---	--	----------

Backrest Stitching options

Silver	Black	Blue	Individual:
Red	Orange	Purple	

RGK Logo colour

White (std)	Green	Blue	Individual:
Red	Orange	Purple	

Model Logo colour

White (std)	Green	Blue	Individual:
Red	Orange	Purple	

XP170	Individual backrest embroidery (PLEASE PROVIDE LOGO AT TIME OF ORDER)	\$247.00
-------	--	----------

Cushion

XP080	Cushion Foam (Select density and size below)	\$280.00
-------	---	----------

Soft (CMS040)	Medium (CMS050)	
50mm	75mm	100mm
	Other	mm

XP340	Hard Plastazote cushion upgrade (Indicate max height)	\$118.00
	580mm	630mm

Cushion Cover Options

XP085	Nylon cushion cover (Black)	Std
XP380	Neoprene cushion cover (Black)	\$131.00
XP090	Air Lite cushion cover (Black)	\$101.00

Straps

XS050	100mm Pro adjustable velcro strap	\$178.00
XS055	50mm Pro adjustable velcro strap	\$139.00
XS057	50mm Pro knee strap (measurement from tube to tube over knees)	\$258.00

Strap 1 Type Measurement

Strap 2 Type Measurement

Strap 3 (Additional) Type Measurement Specify strap length Distance from back post
mm mm

XS020 Sport ratchet strap (Click strap)
(11" x 2" pad) Quantity: \$144.00
1 2

XS021 Sport ratchet strap (Single) Quantity: \$185.00
1 2

XS025 Pelvi Loc 2DA ratchet strap
(For thighs & knees) Quantity: \$418.00
1 2

Select size:

Medium (9") Large (11")

XS030 Pelvi Loc 3DA ratchet strap
(For hips) Quantity: \$432.00
1 2

Select size:

XS (8") S (12") M (14") L (17")

XS035 M2 Ratchet strap
(For thighs and knees) Quantity: \$291.00
6" 8" 10" 12" 1 2

XS036 M2 Ratchet strap (Single)
(For Feet) Quantity: \$308.00
6" 8" 10" 12" 1 2

Strap 1 Type Strap 2 Type

Strap 3 Type

Specify distance from back post
mm

Specify distance from back post
mm

Specify distance from back post
mm

XS040 Fixed mounting points x1 x2 x3 \$165.00
(One set per click strap)

Frame Options

XR060	Calf strap mounting loops	\$159.00
XR065	Strap location loops	\$122.00

Accessories

XR030	Triple wheel bag (Black)	\$207.00
XP130	Frame protectors full length	\$131.00
XA150	Plastic tube fitted to offensive wing	\$96.00
XS070	Foot strap (select if required)	No Cost
XS060	Calf strap (select if required)	No Cost
XP400	Back bar protector	No Cost
XR090	Tool kit	No Cost



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8

Customer Service: (800) 263-3390 E-Mail: cscanada@sunmed.com www.sunrisemedical.ca
Specifications are subject to change without notice. © 9/2026 Sunrise Medical Canada Inc. CN-100062-CAN
REV (8)