



Hoyer® Classics Product Overview

redefining patient handling



For more than 50 years, Hoyer lifters have defined the lifter category. From Ted Hoyer's early designs in the 1950's to more efficient models built with modern designs and materials, the Hoyer name has always stood for innovation. Hoyer's full line of truly functional and versatile lifts offers outstanding performance and style beyond compare!



HPL700/HPL700WSC

- Power operated base
- Base clearance of 4.5"
- 6-point cradle design minimizes bariatric patient discomfort
- Oversized handles offer various grip choices while encouraging proper body mechanics
- Floor pick-up capability
- Emergency stop, power manual lowering
- Utilizes HoyerPro Bariatric Sling
- 700 lb. patient capacity



HPL402

- Adjustable U-base
- Reverse-mounted actuator for additional knee clearance
- 6-point non-sway cradle offers maximum sling options
- Built-in emergency stop button
- Manual release for emergency lowering
- Lift leg guards and caster bumpers
- 400 lb. patient capacity



P-C-HLA-2

- Powered lifting provides more focus on the patient being lifted
- Adjustable U-base with 5.5" clearance
- 6-point cradle works with an array of Hoyer slings
- Easy cleaning, durable chrome finish
- 400 lb. patient capacity



C-HLA-2, C-HLA-2T, C-CB-L2, C-HLA

- Innovative 6-point cradle works with an array of Hoyer slings, including 2, 4, and 6-point slings
- Extended reach for floor pick-up capabilities
- Adjustable legs (does not apply to C-CB-L2)
- 400 lb. resident capacity
- Long-wearing chrome finish maintains appearance
- Mast and base separate for easy travel or storage
- Standard 4" casters are offset for greater stability and smooth operation



HML400

- Contemporary painted design
- Ergonomically placed hydraulic pump handle is easy to use
- Pump handle rotates, allowing operation from either side of lifter
- Adjustable U-legs for easier seating or toileting
- Disassembles for easy storage
- 400 lb. resident capacity
- Priced for reimbursement

Hoyer® Classics Product Overview

Hoyer Classic Lifts	
HPL700/HPL700WSC 700 lb. patient capacity	Power operated base, Base clearance of 4.5", 6-point bariatric cradle design minimizes patient sway
HPL402 400 lb. patient capacity	Adjustable U-base, Reverse-mounted actuator for additional knee clearance, 6-point cradle works with an array of Hoyer slings
P-C-HLA-2 400 lb. patient capacity	Adjustable U-base with 5.5" clearance, 6-point cradle works with an array of Hoyer slings
C-HLA-2, C-HLA-2T, C-HLA, C-CB-L2 400 lb. patient capacity	Adjustable legs, Cradle design minimizes sway, 6-point non-sway cradle offers extra stability
HML400 400 lb. patient capacity	Adjustable U-base

Hoyer® Classics Sling Overview



One Piece Commode Sling with Positioning Strap

- Durable, polyester with closed-cell foam padding for comfort
- Aperture for toileting without removing sling
- Weight Capacity: 600 lbs.



One Piece Sling with Positioning Strap

- Durable, polyester with closed-cell foam padding for comfort
- Provides easy position adjustment of patient
- Weight Capacity: 600 lbs.



U-Sling, Padded with Head Support

- Durable, polyester with closed-cell foam padding for comfort
- Nylon head support
- Weight Capacity: 600 lbs.



U-Sling, Padded

- Durable, polyester with closed-cell foam padding for comfort
- Weight Capacity: 600 lbs.

Hoyer Classic Slings	X-Large	Large	Medium	Small
One Piece Commode Sling with Positioning Strap	—	70056	70057	70058
One Piece Sling with Positioning Strap	—	70051	70052	70053
U-Sling, Padded with Head Support	70010	70011	70012	70013
U-Sling, Padded	70000	70001	70002	70003

HPL700 Slings	XX-Large	X-Large	Large	Medium
Full Back	NA1075	NA1074	NA1073	NA1072
Quick Fit	NA1012P	NA1011P	NA1010P	NA1009P



Full Back

- Weight Capacity: 850 lbs.



Quick Fit

- Weight Capacity: 850 lbs.



www.joerns.com • email: info@joerns.com

Post Acute, Acute, HomeCare
2430 WhiteHall Park Drive
Charlotte, NC 28273
(P) 800.826.0270
(F) 800.457.8827

VA/Government
19748 Dearborn Street
Chatsworth, CA 91311
(P) 800.966.6662
(F) 800.232.9796

Canadian Office
1000 Clarke Road, Suite 6
London, ON Canada N5V 3A9
(P) 866.546.1151
(F) 519.451.8662

© 2012 Joerns Healthcare • 6190034 RevC • 12-0350