

OPEN FRAME GAIT TRAINER supports children in an upright position, improving stepping and maximizing opportunities for exploration, interaction, and participation.

CORSET-LIKE HARNESS for maximized trunk support and upright, hands-free walking.

- 3. RE-SIZED FOR EARLY INTERVENTION Size 1 Frame fits users as young as 12 months old.
- 4. 7" INDOOR/OUTDOOR CASTERS with patent-pending four-way directional lock allow easy maneuverability over different surfaces. The only gait trainer that offers SIDE-STEPPING feature.
- NARROW SADDLE & DROP ATTACHMENT reduces abduction while addressing scissoring gait.
- 6. 3-WAY ADJUSTABLE HEADREST with adjustable laterals for additional head support.
- 3-WAY ADJUSTABLE USER HANDLES for additional forearm stability.
- MyWay Pedal provides elliptical motion in standing that promotes a natural gait pattern and strengthens anti-gravity muscles.











TECHNICAL DATA	
AGE RANGE	1 to 16 years
USER SHOULDER HEIGHT	23" to 59"
USER SADDLE HEIGHT	9.5" to 32"
FRAME MAX USER WEIGHT	Size 1 55 lbs. Size 2 110 lbs. Size 3 176 lbs.
UPPER SUPPORT UNIT WEIGHT	Size 1 N/A Size 2 23 lbs. Size 3 24 lbs.
FRAME WEIGHT	Size 1 26.5 lbs. Size 2 14.5 lbs. Size 3 21 lbs.
HARNESS MAX USER WEIGHT	Size 1 to 3 110 lbs. Size 4 to 6 176 lbs.
HARNESS WAIST MEASUREMENT	18" to 35.5"

Refer to order form as specifications may vary based on configuration.







Video & Resources













