



CLUB SPORT

»» A must in your team

Our aim is to inspire participation by building chairs that are suitable for each individual and every team to win at their discipline, regardless of age or aspiration.

We are proud of our history and strong links to universities and sporting federations alike, helping them to achieve the best results for their teams' performances with our chairs.



Club Sport - The more versatile sport chair

An entry level sports chair that is economical and ready to go when you are.

Feeling inspired to take up a wheelchair sport but don't want to make a major commitment just yet? Club Sport is the chair for you! An entry level sports chair that is economical and ready to go when you are. Customisable with your own embroidery or logo and compatible with any court sport. Its range of sizes and adjustments also makes it ideal for sport clubs and schools.

Ready to play with offensive wing, sideguards, cushion and complete with an anti tip feature to keep you safe at speed.

Key Features

- » Chrome moly steel frame
- » 4 configurable frame options
- » Aluminium sideguards
- » Integrated offensive wing
- » Height and angle adjustable footrest
- » Quick release height adjustable anti tip
- » Anodized aluminium handrims
- » Tension adjustable seat and back upholsteries
- » Height adjustable backrest (250-300 mm)
- » 18° camber
- » Adjustable COG
- » 3" castors with integrated fork
- » Leg strap and foot strap included standard
- » Wide range of colours and upholstery customisations
- » Max. user weight 125 kg



Aluminium sideguards



Integrated offensive wing



Height and angle adjustable footrest



Adjustable COG and quick release height adjustable anti tip