



Customer Contact Details

Customer Name: _____	Ship to Customer Address: <input type="checkbox"/>
Advisor: _____	Shipping Name: _____
Address: _____	Shipping Address: _____
_____	_____
_____	_____
Phone: _____	Phone: _____
Email: _____	Email: _____

Additional Information

Notes: _____	Current Wheelchair: _____
_____	Customer Weight: _____
_____	Customer Height: _____
_____	Customer Disability: _____

Guide:		<input checked="" type="checkbox"/> Base product	<input type="checkbox"/> Required option	<input type="checkbox"/> Optional Extra	Critical Information	Information Needed
PRODUCT STYLE						
<input checked="" type="checkbox"/>	FXAXZ010	TIGA FX				
	\$7,930.00					
<input checked="" type="checkbox"/>	XA070	Open design				
Max User Weight 110kg (241lbs)						
<input type="checkbox"/>	Yes, I am within the max user weight					
CAD APPROVAL						
<input type="checkbox"/>	XW012	QUOTE CAD drawing for approval \$550 FLAT FEE (Includes two free revisions)	<input type="checkbox"/>	XW010	ORDER CAD drawing for approval \$300 FLAT FEE (Includes two free revisions)	
<input type="checkbox"/>	XW009	QUOTE only no CAD \$0	<input type="checkbox"/>	XW011	ORDER only no CAD \$0	
KEY MEASUREMENTS						
<input checked="" type="checkbox"/>	XY010	Seat width (Inside of sideguards)				
	<input style="width: 100px;" type="text"/>	mm	(Max - 460mm) (Min - 250mm)			
<input checked="" type="checkbox"/>	XY020	Seat length				
	<input style="width: 100px;" type="text"/>	mm	(Max - 460mm) (Min - 250mm)			
<input type="checkbox"/>	XY030	Extra seat length (Provides additional leg support & hand grip)				
	<input style="width: 100px;" type="text"/>	mm				
<input checked="" type="checkbox"/>	XY050	Front frame height				
	<input style="width: 100px;" type="text"/>	mm	(Max - 550mm) (Min - 480mm)			
<input checked="" type="checkbox"/>	XY060	Rear frame height				
	<input style="width: 100px;" type="text"/>	mm	(Max - 550mm) (Min - 360mm)			
<input checked="" type="checkbox"/>	XY070	Frame length				
	<input style="width: 100px;" type="text"/>	mm	(Max 625mm) (Min XY020 + XY030 + 140mm)			
<input type="checkbox"/>	XY072	Caster position				
	<input style="width: 100px;" type="text"/>	mm				
	<input type="checkbox"/>	Optimize measurements				
<input type="checkbox"/>	XA020	Ergonomic seat				
	\$121.00					
	<input style="width: 100px;" type="text"/>	mm	Specify Length (Standard = 120mm)			


FRONT FRAME AND FOOTREST

Front Frame Options

 XY080 Seat taper design

Specify width



!
 mm

 Optimize measurements (XY010 - 100mm)




Footrest Design (Choose One)

 XY100 Round design

Caster Size	Max Height	With XB200
3"	115mm	165mm
4"	125mm	175mm
5"	135mm	185mm


 XB200 High level footrest (Only with XB010, XB080 & XB150)
\$201.00

 XY105 Square design

Caster Size	Max Height	With XB200
3"	135mm	185mm
4"	145mm	195mm
5"	155mm	205mm


 XB200 High level footrest (Only with XB010, XB080 & XB150)
\$201.00


Footrest Individual Measurements

 XY115 Front footrest height

!
 mm

 XY120 Rear footrest height


i
 mm

(Std = 20 lower than front height)

 XY125 Footrest length

i
 mm

(Std = 100mm)
(160mm and over requires heavy duty upgrade XB150)



Footrest Style and Cover

 FXAXB080 Tubular footrest welded w/ Aluminium plate
(Supplied with spacer for +20mm adjustment)

 FXAXB120 Clamp-on aluminium plate footrest
(Angle/height adjustable)

 FXAXB110 Tubular footrest welded w/ Carbon Fibre plate
(Supplied with spacer for +20mm adjustment)

Footrest Options

 XB165 Rubber grip tape (3M safety walk)
\$86.00

CASTER WHEELS, FORKS & WHEEL LOCKS													
Caster Wheels													
75mm (3")													
<input type="checkbox"/> XC150	Aluminium centre 3"x1.4" <i>(Silver hub and black urethane tire)</i>												
100mm (4")													
<input type="checkbox"/> XC060	Soft roll 4"x1.3" <i>(Black)</i>												
<input type="checkbox"/> XC110	Aluminium centre 4"x1" <i>(Silver hub and black urethane tire)</i>												
<input type="checkbox"/> XC130	Aluminium centre 4"x1.4" <i>(choose colour below)</i>												
i	<table style="width: 100%; border: none;"> <tr> <td><input type="checkbox"/> Silver hub & black tire</td> <td><input type="checkbox"/> Black hub & blue tire</td> <td><input type="checkbox"/> Purple hub & black tire</td> <td><input type="checkbox"/> Silver hub & blue tire</td> </tr> <tr> <td><input type="checkbox"/> Black hub & black tire</td> <td><input type="checkbox"/> Red hub & black tire</td> <td><input type="checkbox"/> Gold hub & black tire</td> <td><input type="checkbox"/> Silver hub & red tire</td> </tr> <tr> <td><input type="checkbox"/> Black hub & red tire</td> <td><input type="checkbox"/> Blue hub & black tire</td> <td><input type="checkbox"/> Green hub & black tire</td> <td></td> </tr> </table>	<input type="checkbox"/> Silver hub & black tire	<input type="checkbox"/> Black hub & blue tire	<input type="checkbox"/> Purple hub & black tire	<input type="checkbox"/> Silver hub & blue tire	<input type="checkbox"/> Black hub & black tire	<input type="checkbox"/> Red hub & black tire	<input type="checkbox"/> Gold hub & black tire	<input type="checkbox"/> Silver hub & red tire	<input type="checkbox"/> Black hub & red tire	<input type="checkbox"/> Blue hub & black tire	<input type="checkbox"/> Green hub & black tire	
<input type="checkbox"/> Silver hub & black tire	<input type="checkbox"/> Black hub & blue tire	<input type="checkbox"/> Purple hub & black tire	<input type="checkbox"/> Silver hub & blue tire										
<input type="checkbox"/> Black hub & black tire	<input type="checkbox"/> Red hub & black tire	<input type="checkbox"/> Gold hub & black tire	<input type="checkbox"/> Silver hub & red tire										
<input type="checkbox"/> Black hub & red tire	<input type="checkbox"/> Blue hub & black tire	<input type="checkbox"/> Green hub & black tire											
125mm (5")													
<input type="checkbox"/> XC070	Soft roll 5"x1.3" <i>(Black)</i>												
<input type="checkbox"/> XC100	Pneumatic Aluminium centre 5"x1" <i>(Grey)</i>												
<input type="checkbox"/> XC120	Aluminium centre 5"x1" \$232.00 <i>(Silver hub and black urethane tire)</i>												
<input type="checkbox"/> XC140	Aluminium centre 5"x1.4" \$232.00 <i>(Silver hub and black urethane tire)</i>												
150mm (6")													
<input type="checkbox"/> XC080	Soft roll 6"x1.3"												
<input type="checkbox"/> XC145	Aluminium centre 6"x1.4" \$310.00 <i>(Silver hub and black urethane tire)</i>												
Caster Forks													
<input type="checkbox"/> XC190	Low ride caster fork <i>(Std when no option selected)</i>												
<input type="checkbox"/> XC210	Froglegs suspension caster fork \$633.00												
<i>Note: Colour follows accent</i>	i <input type="checkbox"/> Black <input type="checkbox"/> Silver												
Caster Pin Width													
<input checked="" type="checkbox"/> XY330	Caster pin width (XY010 plus 70mm)												
Wheel Lock Options (Choose one)													
<input type="checkbox"/> XD020	Underseat Scissor												
<input type="checkbox"/> XD030	Compact scissor \$266.00												
<input type="checkbox"/> XD035	Compact scissor with lightweight fixing \$352.00 <i>(Only with fixed centre of gravity)</i>												
<input type="checkbox"/> XD045	Out Front Compact scissor \$409.00												
<input type="checkbox"/> XD050	Push to lock												
<input type="checkbox"/> XD046	Out Front push to lock \$409.00												
<input type="checkbox"/> XD010	Omit wheel locks												

SIDEGUARDS AND ARMRESTS

Sideguard		Sideguard with Fender	
Aluminium Options		Aluminium Options	
<input type="checkbox"/> XE080	Lift out aluminium \$532.00	<input type="checkbox"/> XE090	Lift out aluminium with fender \$570.00
<input type="checkbox"/> FXAXE130	Fold flat aluminium \$491.00	<input type="checkbox"/> FXAXE140	Fold flat aluminium with fender \$547.00
Carbon Fibre Options		Carbon Fibre Options	
<input type="checkbox"/> XE100	Lift out carbon fibre \$766.00	<input type="checkbox"/> XE110	Lift out carbon fibre with fender \$1,132.00
<input type="checkbox"/> FXAXE150	Fold flat carbon fibre \$748.00	<input type="checkbox"/> FXAXE160	Fold flat carbon fibre with fender \$1000.00
ABS Plastic Options		Sideguard Measurements	
<input type="checkbox"/> XE120	Lift out ABS plastic \$487.00	Fender height	
<input type="checkbox"/> FXAXE170	Fold flat ABS plastic \$491.00	<input type="checkbox"/> 10mm above tire <i>(Standard)</i>	
		<input type="checkbox"/> Individual above tire <i>(Min = 10mm)</i>	
Armrest Options		<input type="checkbox"/> <input type="text" value=""/> mm	
<input type="checkbox"/> XF010	Lift out armrest with plastic sideguard <i>(Indicate fixed height)</i>	Fender position	
<input type="checkbox"/> XF020	Lift out armrest with aluminium sideguard <i>(Indicate fixed height)</i>	<input type="checkbox"/> Centre of tire <i>(Standard if no option selected)</i>	
		<input type="checkbox"/> Outside of tire	
Sideguard Measurements			
Optimized measurements			
<input type="checkbox"/> 25mm above tire <i>(Follow tire)</i>			
<input type="checkbox"/> Individual above tire <i>(Follow tire)</i>			
<input type="checkbox"/> Level with tire <i>(Follow tire)</i>			
<input type="checkbox"/> Sideguard length	<input type="text" value=""/> mm		
<input type="checkbox"/> Sideguard height	<input type="text" value=""/> mm		

BACKREST	
Backrest Height	
<input type="checkbox"/> XY230 Fixed height backrest (Maximum 400mm) <div style="margin-top: 10px;"> Specify height <input style="width: 100px;" type="text"/> mm </div>	<input type="checkbox"/> XY245 Adjustable height backrest \$248.00 (Maximum height 400mm) (+/-25mm) <div style="margin-top: 10px;"> Specify height <input style="width: 100px;" type="text"/> mm </div>
Backrest Design	
<input checked="" type="checkbox"/> FXAXG060 Folding backrest with Q Lock (Black)	
Ergonomic Options	
<input checked="" type="checkbox"/> XY360 Backrest angle <div style="margin-top: 10px;"> <input type="checkbox"/> -5 <input type="checkbox"/> -3 <input type="checkbox"/> 0 <input type="checkbox"/> +3 <input type="checkbox"/> +5 </div> <div style="margin-top: 10px;"> Individual <input style="width: 100px;" type="text"/> Deg (Max -7/+7 Deg) </div>	<input type="checkbox"/> XG100 Ergonomic backrest \$146.00 <div style="margin-top: 10px;"> Vertical height <input style="width: 100px;" type="text"/> mm </div> <div style="margin-top: 10px;"> Rearward position <input style="width: 100px;" type="text"/> mm </div>
<input type="checkbox"/> XG110 Inward taper \$146.00 (inside of backrest posts) <div style="margin-top: 10px;"> Specify width or per side <input style="width: 100px;" type="text"/> mm </div>	<input type="checkbox"/> XG120 Outward flare \$146.00 (inside of backrest posts) <div style="margin-top: 10px;"> Specify width or per side <input style="width: 100px;" type="text"/> mm </div>
Push handle Options	
<input type="checkbox"/> XG070 Screw in push handles \$225.00 (N/A with XY245)	
<input type="checkbox"/> XG080 Flip down push handles \$441.00 (N/A with XY245)	
<input type="checkbox"/> XG085 Flip down push handles height adjustable \$720.00 (N/A with XY245)	
<input type="checkbox"/> XG097 Integrated push handle \$676.00 (Quickie design)	

CAMBER BAR	
<p style="text-align: center; background-color: #eee; margin-bottom: 5px;">Adjustable Centre of Gravity</p> <p><input type="checkbox"/> XH110 Square back adjustable centre of gravity (25-150mm) (Max seat width 500mm) Specify measurement <input style="width: 100px;" type="text"/> mm</p> <p><input type="checkbox"/> XH120 Round back adjustable centre of gravity (60-150mm) Specify measurement <input style="width: 100px;" type="text"/> mm</p>	<p style="text-align: center; background-color: #eee; margin-bottom: 5px;">Fixed Centre of Gravity</p> <p><input type="checkbox"/> FXAXH090 V-style fixed center of gravity \$613.00 (25-125mm) Specify measurement <input style="width: 100px;" type="text"/> mm</p>
Specify Camber Angle	Specify Camber Angle
<p><input checked="" type="checkbox"/> XY390 Camber angle</p> <p><input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 5</p>	<p><input checked="" type="checkbox"/> XY390 Camber angle</p> <p><input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 5 <input style="width: 50px;" type="text"/> Custom</p>
Camber Bar Options	
<p><input type="checkbox"/> FXAXK227 Standard Axle pin (120 or 130mm)</p> <p><input type="checkbox"/> XK240 Quad release axle nuts (pair) \$241.00 (N/A with FXAXK226)</p>	<p><input type="checkbox"/> XK220 Titanium 130mm quick release axle pins (pair) \$160.00 (N/A with XK010, XK070 or FXAXK226)</p>
Anti-Tip Options	
<p><input type="checkbox"/> XJ010 Swingaway anti tip LEFT side (Only with square back frame. Not with XH040) \$167.00</p> <p><input type="checkbox"/> XJ020 Swingaway anti tip RIGHT side (Only with square back frame. Not with XH040) \$167.00</p>	<p><input type="checkbox"/> XJ100 Single sport anti tip (centre mounted) \$504.00</p>

REAR WHEELS	
Rear Wheel Size	
<input type="checkbox"/> XY420 24" Wheels	<input type="checkbox"/> XY430 25" Wheels
Rear Wheel Gap	
<div style="display: flex; justify-content: space-between;"> <div style="width: 30%;"> <input checked="" type="checkbox"/> XY460 Rear wheel gap 15 <input type="checkbox"/> 20 <input type="checkbox"/> 25 Individual <input style="width: 100px;" type="text"/> mm </div> <div style="width: 30%; text-align: center;"> </div> </div>	
Rear Wheel Options	
<input type="checkbox"/> XK010 Fusion 16 spoke <i>(Silver Hub, black spokes)</i> <input type="checkbox"/> XK050 Spinerigy LX 12 spoke \$1000.00 <input type="checkbox"/> XK150 Spinerigy CLX 18 carbon blade spoke \$3,139.00 <i>(N/A in silver hub)</i>	<input type="checkbox"/> XK320 Spinerigy XLX 12 spoke, cross/radial spoked \$1,067.00 <input type="checkbox"/> XK060 Spinerigy SPOX 18 spoke \$665.00 <input type="checkbox"/> XK070 Spinerigy BLXL 18 blade spoke \$1,051.00 <i>(N/A in silver hub)</i>
Spinerigy Colour Options:	
<input type="checkbox"/> Black spokes <input type="checkbox"/> White spokes <input type="checkbox"/> Yellow spokes <input type="checkbox"/> Red spokes <input type="checkbox"/> Blue spokes <input type="checkbox"/> Pink spokes <input type="checkbox"/> Purple spokes <input type="checkbox"/> Orange spokes <input type="checkbox"/> Green spokes	<input type="checkbox"/> XK180 Individual spoke design \$213.00 Inside: _____ Outside: _____
Spinerigy Anodise Options	
<input type="checkbox"/> XK190 Black hub <i>(std option)</i> <input type="checkbox"/> XK200 Silver hub \$127.00 <i>(n/a with BLXL/CLX)</i>	<input type="checkbox"/> XK201 Red hub/rims \$409.00 <i>(LX, XLX & BLXL only)</i> <input type="checkbox"/> XK205 Individual hub/rims \$568.00 <i>Select colour:</i> _____
<input type="checkbox"/> XK202 Purple hub/rims \$409.00 <i>(LX, XLX & BLXL only)</i> <input type="checkbox"/> XK203 Blue hub/rims \$409.00 <i>(LX, XLX & BLXL only)</i>	<input type="checkbox"/> XK204 Gold hub/rims \$539.00 <i>(LX, XLX & BLXL only)</i>
Supply Without Wheels	
<input type="checkbox"/> XK040 Supply without Wheels	Wheel Type: Tire Type:
Handrim Position / Options	
<input type="checkbox"/> XY480 Narrow set handrim	<input type="checkbox"/> XY485 Wide set handrim
<input type="checkbox"/> XL090 Remove excess mounting tabs \$144.00 <i>(recommended for small wheel gaps)</i>	
Handrim Styles	
<input type="checkbox"/> XL010 Aluminium hard anodized 6 tab <i>(Standard if no option selected)</i>	<input type="checkbox"/> XL040 Titanium 6 tab \$446.00
<input type="checkbox"/> XL020 Tetra grip 6 tab \$400.00	<input type="checkbox"/> XL050 Stainless steel 6 tab \$248.00
<input type="checkbox"/> XL100 Para grip 6 tab \$400.00	
Tire Options	
<input type="checkbox"/> XM010 Kenda	<input type="checkbox"/> Grey <input type="checkbox"/> Blue <input type="checkbox"/> Red
<input type="checkbox"/> XM130 Kenda black \$28.00	
<input type="checkbox"/> XM020 Schwalbe Marathon Plus \$98.00	<input type="checkbox"/> XM060 Schwalbe One (24",25" only) \$133.00
<input type="checkbox"/> XM040 Shox solid airless (Black) \$160.00	<input type="checkbox"/> XM050 Kenda K885 mountain bike tire (24",25" only) \$40.00
Inner Tube Valve Options	
<input type="checkbox"/> XY500 Schraeder inner tube valve <i>(Standard if no option selected)</i>	<input type="checkbox"/> XY490 Presta inner tube valve
Wheel Accessories	
<input type="checkbox"/> XK310 Offroad wheel set 24" rim with handrim, MTB tire and axle pins (Replaces 25" wheel with XM010) \$933.00	<input type="checkbox"/> XK300 Offroad wheel set 22" rim with handrim, MTB tire and axle pins (Replaces 24" wheel with XM010) \$933.00

FINISH OPTIONS			
Frame & Component Finish			
<input type="checkbox"/> XN010	Powder coated frame	<input type="checkbox"/> XN150	Powder coated components <i>(Frame colour if not specified)</i>
Frame Comp. Gloss		Frame Comp. Matt	
<input type="checkbox"/>	<input type="checkbox"/> Black gloss 103	<input type="checkbox"/>	<input type="checkbox"/> Black matt 105
<input type="checkbox"/>	<input type="checkbox"/> White gloss 100	<input type="checkbox"/>	<input type="checkbox"/> White matt 101
<input type="checkbox"/>	<input type="checkbox"/> Grey gloss 132	<input type="checkbox"/>	<input type="checkbox"/> Grey matt 131
<input type="checkbox"/>	<input type="checkbox"/> Cream gloss 135	<input type="checkbox"/>	<input type="checkbox"/> Light grey matt 133
<input type="checkbox"/>	<input type="checkbox"/> Red gloss 108	<input type="checkbox"/>	<input type="checkbox"/> Navy blue matt 152
<input type="checkbox"/>	<input type="checkbox"/> Yellow gloss 120		
<input type="checkbox"/>	<input type="checkbox"/> Forest green gloss 143		
<input type="checkbox"/>	<input type="checkbox"/> Green gloss 107		
<input type="checkbox"/>	<input type="checkbox"/> Dark blue gloss 151		
<input type="checkbox"/>	<input type="checkbox"/> Royal blue gloss 154		
<input type="checkbox"/>	<input type="checkbox"/> Blue gloss 174		
<input type="checkbox"/>	<input type="checkbox"/> Duckegg blue gloss 141		
<input type="checkbox"/>	<input type="checkbox"/> Pink gloss 188		
<input type="checkbox"/>	<input type="checkbox"/> Orange gloss 180		
<input type="checkbox"/>	<input type="checkbox"/> Purple gloss 157		
<input type="checkbox"/>	FXAXN030 Brushed frame \$565.00 <i>(All parts black unless XN160 chosen)</i>	<input type="checkbox"/>	RGK-XN200 Individual paint colour (State Colour/RAL code) \$565.00
<input type="checkbox"/>	FXAXN050 Polished frame \$1,311.00 <i>(All parts black unless XN170 chosen)</i>	<input type="checkbox"/>	XN160 Brushed components \$121.00
		<input type="checkbox"/>	XN170 Polished components \$213.00
		<input type="checkbox"/>	Frame Comp
		<input type="checkbox"/>	<input type="checkbox"/>
Accent Colour			
Note: Accent colour includes, forks & castor caps.			
<input type="checkbox"/>	Black accent	<input type="checkbox"/>	Silver accent
<input type="checkbox"/>	Red accent	<input type="checkbox"/>	Blue accent
<input type="checkbox"/>	Purple accent	<input type="checkbox"/>	Carbon accent <i>(Caster caps only, black for rest)</i>
		<input type="checkbox"/>	RGK-XN250 Individual Accent Colour \$315.00
UPHOLSTERY			
Seat Upholstery			
<input type="checkbox"/>	XP040 Seat sling adjustable tension with straps	<input type="checkbox"/>	XR050 Carbon fibre seat plate \$480.00
<input type="checkbox"/>	FXAXP045 Sub4 Seat sling strap adjustable		
Backrest Upholstery			
<input type="checkbox"/>	XP050 Airtech backrest strap adjustable	<input type="checkbox"/>	XP060 Backrest upgrade to full leather \$500.00 <i>(Black)</i>
<input type="checkbox"/>	FXAXP035 Airtech Sub4 backrest strap adjustable	<input type="checkbox"/>	XP070 Backrest upgrade to leather corners \$121.00 <i>(Black)</i>
Backrest Stitching Colour			
<input type="checkbox"/>	Silver	<input type="checkbox"/>	Black
<input type="checkbox"/>	Red	<input type="checkbox"/>	Light blue
<input type="checkbox"/>		<input type="checkbox"/>	Blue
<input type="checkbox"/>		<input type="checkbox"/>	Purple
<input type="checkbox"/>		<input type="checkbox"/>	Individual: _____
Backrest RGK Logo Colour			
<input type="checkbox"/>	White (Std)	<input type="checkbox"/>	Black
<input type="checkbox"/>	Red	<input type="checkbox"/>	Light blue
<input type="checkbox"/>		<input type="checkbox"/>	Blue
<input type="checkbox"/>		<input type="checkbox"/>	Purple
<input type="checkbox"/>		<input type="checkbox"/>	Individual: _____
Backrest Accessories			
<input type="checkbox"/>	XP110 Backrest with zip pocket \$127.00	<input type="checkbox"/>	XP170 Individual backrest embroidery \$240.00 <i>(Please provide .jpeg file and positioning)</i>

CUSHION		
Cushion foam options		
<input type="checkbox"/> XP080 \$121.00	Cushion foam <i>(select density and size below)</i>	
<input type="checkbox"/> Soft (CMS040)	<input type="checkbox"/> Medium (CMS050)	
<input type="checkbox"/> 50mm	<input type="checkbox"/> 75mm	
<input type="checkbox"/> Hard (9lb recon)		
<input type="checkbox"/> 100mm		
Cushion cover Upgrades		
<input type="checkbox"/> XP081 	Nylon cushion cover <i>(Black)</i>	
<input type="checkbox"/> XP090 \$98.00	Air Lite cushion cover <i>(Black)</i>	
<input type="checkbox"/> XP100 \$254.00	Leather cushion cover <i>(Black)</i>	
Jay Products		
<input type="checkbox"/> JAY Cushion <i>(Price dependant on size / option)</i>	<input type="checkbox"/> Jay Order form attached	Notes: <div style="border: 1px solid black; height: 100px; width: 100%;"></div>
<input type="checkbox"/> JAY Backrest <i>(Price dependant on size / option)</i>	<input type="checkbox"/> Jay Order form attached	
<input type="checkbox"/> JAY Whitmyer Headrest <i>(Price dependant on size / option)</i>	<input type="checkbox"/> Jay Order form attached	
ACCESSORIES		
<input type="checkbox"/> XR080 \$1,053.00	Frontwheel assembly <i>(Only with tubular style footrest)</i>	
<input type="checkbox"/> XP120 \$84.00	Frame protectors Neoprene zip fasten <i>(Black)</i>	
<input type="checkbox"/> XP420 \$81.00	Frame protectors Nylon velcro fasten <i>(Black)</i>	
<input type="checkbox"/> XR010 \$59.00	Seat pouch velcro attach <i>(Black)</i>	
<input type="checkbox"/> XR200 \$151.00	Rucksack with RGK logo (Black)	
<input type="checkbox"/> XR100 \$334.00	FX carry bag	
<input checked="" type="checkbox"/> XR090	Tool kit	
<input checked="" type="checkbox"/> XS060	Calf strap	

NEW!

NEW!

NOTES

