



Serial No:



ORDER FORM			
Customer Dealer Information		Shipping Address	Identical To Billing Address
Dealer		Client	
Client		Address	
Address			
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Information	
Purchase Order No.:		Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):		User Disability:	
		User Weight:	kg User Height
Build4me			

Serial No:

Product Style

FXAXZ010 TIGA FX

\$8,047.00

XA070 Open design

Max User Weight 110kg (242lbs)

Please Note: All RGK chairs are hand built and have a manufacturing measurement tolerance of +/- 5mm

Std

CAD Approval

XW009 Quote only (No CAD For Approval)

XW010 Order CAD Drawing For Approval *(Includes 2 CAD Amendments)*

\$300.00

XW011 Order (No CAD Required)

XW012 Quote & CAD For Approval *(Includes 2 CAD Amendments)*

\$550.00

XY010 Seat Width (Inside of sideguards)

mm (Min - 250mm / Max - 460mm)

XY020 Seat Length

mm (Min - 250mm / Max - 460mm)

XY030 Extra Seat Length (Provides additional leg support and hand grip)

mm

XY050 Front Frame Height

mm (Min - 480mm / Max - 550mm)

XY060 Rear Frame Height

mm (Min - 360mm / Max - 550mm)

XY070 Frame Length

mm (Min = XY020 + XY030 +140mm / Max - 625mm)



XY072 Castor Position

mm Optimise Measurements

XA020 Ergonomic Seat

\$125.00

mm Specify Length - (Standard = 120mm)



Front Frame

XY080 Seat Taper Design

mm

Optimise Measurements
(XY010 Seat Width - 100mm)



Footrest Design

XY100 Round Design

(Castor wheel size 4", max height achievable 110mm)

No Cost

XY105 Square Design

(Castor wheel size 4", max height achievable 110mm)

No Cost

XB200 High Level Footrest

(Only with XB010, XB080 & XB150)

\$207.00

Footrest Individual Measurements

XY115 Front Footrest Height (Minimum 50mm)

mm

XY120 Rear Footrest Height (Minimum 50mm)

mm STD = 20mm lower than front height

XY125 Footrest Length

(Std = 100mm)

mm (160mm and over requires heavy duty upgrade)

Footrest Style and Cover

FXAXB080 Tubular Footrest welded w/Aluminium Plate

No Cost

FXAXB110 Tubular Footrest welded w/Carbon Fibre Plate

(Supplied with spacer for +20mm adjustment)

£108.00

FXAXB120 Clamp-on Aluminium Plate Footrest

(Angle/Height adjustable)

No Cost

Footrest Options

XB165 Rubber Grip Tape (3m safety walk)

\$89.00

XB150 Heavy Duty Tubular Footrest

\$185.00

Castor Wheels

XC060	Soft roll 4" x 1.3" (Black only)	No Cost
XC110	Aluminium centre 4" x 1" (Silver hub and black urethane tyre only)	\$239.00
XC130	Aluminium centre 4" x 1.4" (Indicate Colour)	\$239.00
	Silver hub & black tyre	
	Black hub & black tyre	
	Red hub & black tyre	
	Gold Hub & black tyre	
	Green hub & black tyre	
	Blue hub & black tyre	
	Purple hub & black tyre	
XC070	Soft roll 5" x 1.3" (Black only)	No Cost
XC100	Pneumatic Aluminium Centre 5" x 1" (Grey only)	\$255.00
XC120	Aluminium Centre 5" x 1" (Silver hub and black urethane tyre only)	\$239.00
XC140	Aluminium Centre 5" x 1.4" (Silver hub and black urethane tyre only)	\$239.00
XC080	Soft roll 6" x 1.3" (Black hub / grey tyre only)	No Cost
XC145	Aluminium centre 6" x 1.4" (Silver hub and black urethane tyre only)	\$319.00

Castor Forks

XC190	Performance Castor Fork (Std) Note: Colour follows accent	Std
XC192	High Trail Castor Fork Note: Greater distance between contact patch and steering axis reduces flutter at higher speeds. Colour follows accent	No Cost
XC210	Froglegs Suspension Castor Fork - Black - Silver	\$652.00

Castor Pin Width

XY330	Castor Pin Width (XY010 plus 70mm)
-------	---------------------------------------

Wheel Lock Options

XD020	Underseat Scissor	No Cost
XD050	Push To Lock	No Cost
XD030	Compact Scissor	\$274.00
XD046	Out front push to lock	\$421.00
XD035	Compact scissor with lightweight fixing (Only with fixed centre of gravity)	\$363.00
XD085	Eagle wheel locks (Indicate colour)	\$482.00
XD086	Eagle WL with L/W fixing (Indicate colour) (Only with fixed centre of gravity)	\$523.00
	Black Purple Red Blue	
XD010	Omit wheel locks	

Sideguards

Aluminium options

XE080	Lift out aluminium	\$548.00
FXAXE130	Fold flat aluminium	\$506.00
XE090	Lift out aluminium with fender	\$596.00
FXAXE140	Fold flat aluminium with fender	\$563.00

Carbon Fibre options

XE100	Lift out carbon fibre	\$789.00
FXTXE150	Fold flat carbon fibre	\$1,526.00
XE110	Lift out carbon fibre with fender	\$523.00
FXTXE160	Fold flat carbon fibre with fender	\$1,680.00

ABS Plastic options

XE120	Lift out ABS plastic	\$502.00
FXTXE170	Fold flat ABS plastic	\$741.00

Sideguard measurements

Optimised Measurements

25mm above tyre (Follow tyre)	Individual above tyre (Follow tyre)	mm
Level with tyre (Follow tyre)		
Sideguard length	Sideguard height	mm

Sideguard Fender Measurements

Fender Height

10mm above tyre (Standard)	Individual above tyre (Min = 10mm)	mm
Centre of tyre (Standard if no option selected)		
Outside of tyre		

Armrest Options

XF010	Lift out armrest with plastic sideguard (Indicate fixed height)	\$610.00
XF020	Lift out armrest with aluminium sideguard (Indicate fixed height)	\$481.00

Height (Top of seat rail to top of armrest)

mm

Backrest Height

XY230	Fixed height backrest (Std) (Maximum 400mm)	XY245	Adjustable height backrest (Maximum 400mm) (+/- 25MM)	\$255.00
	Specify height mm		Specify height mm	

Backrest Design

FXAXG060 Folding backrest with Q Lock
(Black)

Ergonomic Options

XY360	Backrest angle -5 -3 0 +3 +5	XG100	Ergonomic backrest	\$150.00
	Individual (Max -7/+7 Deg)		Vertical height mm	
	Deg		Rearward position mm	

XG110 Inward taper
(Inside of backrest posts)
Specify width or per side
mm

XG120 Outward flare
(Inside of backrest posts)
Specify width or per side
mm

Push Handle Options

XG070	Screw in push handles (N/A with XY245)	\$232.00
XG080	Flip down push handles (N/A with XY245)	\$454.00
XG085	Flip down push handles height adjustable (N/A with XY245)	\$742.00
XG097	Integrated push handles (Quickie design)	\$696.00

Camber Bar

Adjustable Centre Of Gravity

XH110	Square back adjustable centre of gravity (25-150mm)	No Cost
-------	--	---------

(Max seat width 500mm)

Specify measurement

mm

XH120	Round back adjustable centre of gravity (60-150mm)	No Cost
-------	---	---------

Specify measurement

mm

Fixed Centre Of Gravity

FXAXH090	V Style Fixed centre of gravity (25-150mm) Specify measurement	\$631.00
----------	--	----------

mm

Camber Angle

XY390	Camber angle 0 1 2 3 5
-------	---------------------------

Axle Pin Options

FXAXK227	Standard Axle Pin (120 or 130mm)	Std
FXAXK226	100mm Titanium Axle pin (Fixed Axle only)	\$416.00
XK220	Titanium 130mm quick release axle pins (pair) (N/A with XK010, XK070 or FXAXK226)	\$165.00
XK240	Quad release axle nuts (pair) (N/A with FXAXK226)	\$248.00

Anti-Tip Options

XJ010	Swingaway anti tip LEFT side (Only with square back frame. Not with XH040)	\$172.00
XJ020	Swingaway anti tip RIGHT side (Only with square back frame. Not with XH040)	\$172.00
XJ100	Single sport anti tip (Centre mounted as standard)	\$519.00

Rear Wheel Size

XY420 24" Wheels XY430 25" Wheels

Rear Wheel Gap

XY460 Rear Wheel Gap
15 20 25 Individual mm



Rear Wheel Options

XY460	Rear Wheel Gap	
15	20	25
Individual	mm	

XK010	Fusion 16 Spoke (Silver hub, black spokes)	Std
XK320	Spinergy XLX 12 Spoke, cross/radial spoked	\$1,099.00
XK050	Spinergy LX 12 Spoke	\$1,030.00
XK060	Spinergy Spox 18 Spoke	\$685.00
XK150	Spinergy CLX 18 Carbon blade spoke (N/A in silver hub)	\$3,233.00
XK070	Spinergy BLXL 18 blade spoke (N/A in silver hub)	\$1,099.00

Spinergy Spoke Colour options

Black spokes White spokes Yellow spokes
Red spokes Blue spokes Pink spokes
Purple spokes Orange spokes Green spokes

XK180	Individual spoke design (Specify colour)	\$219.00
	Inside Outside	

Spinergy Hub Colour options

XK190	Black Hub (Std option)	Std
XK200	Silver Hub (n/a with BLXL/CLX)	\$131.00
XK201	Red Hub/Rims (LX, XLX & BLXL only)	\$421.00
XK202	Purple Hub/Rims (LX, XLX & BLXL only)	\$421.00
XK203	Blue Hub/Rims (LX, XLX & BLXL only)	\$421.00
XK204	Gold Hub/Rims (LX, XLX & BLXL only)	\$555.00

Supply Without Wheels

XK040	Supply without Wheels	Wheel Type	NC
		Tyre Type	

Handrim Position / Options

XY480	Narrow set handrim	No Cost
XY485	Wide set handrim	No Cost
XL090	Remove excess mounting tabs (Recommended for small wheel gaps)	\$150.00

Handrim Styles

XL010	Aluminium hard anodised 6 tab (Standard if no option selected)	Std
XL040	Titanium 6 tab	\$459.00
XL160	Elipse 6 tab (replacement for surge)	\$650.00
XL035	Natural Fit LT	\$765.00
XL020	Tetra grip 6 tab	\$412.00
XL150	Carbolife Gekko BLACK SILVER	\$1,448.00
XL100	Para grip 6 tab	\$412.00
XL170	Carbolife Curve L (Black)	\$1,435.00

Tyre Options

XM010	Kenda Grey Blue Red (Grey standard if no option selected)	Std
XM060	Schwalbe one (24", 25" only)	\$137.00
XM130	Kenda black	\$29.00
XM050	Kenda K885 mountain bike tyre (24", 25" only)	\$41.00
XM020	Schwalbe marathon plus	\$101.00
XM040	Shox solid airless (Black) (26" N/A with Fusion Wheels)	\$165.00

Inner Tube Valve Options

XY500	Schraeder inner tube valve (Standard if no option selected)	No Cost
XY490	Presta inner tube valve	No Cost

Wheel Accessories

XK310	Offroad wheel set 24" rim with handrim, MTB tyre & axle pins (Replaces 25" wheel with XM010)	\$961.00
XK300	Offroad wheel set 22" rim with handrim, MTB tyre & axle pins (Replaces 24" wheel with XM010)	\$961.00

Finish Options

XN010 Powder coated frame

No Cost

Frame

Gloss

Black gloss 103
White gloss 100
Grey gloss 132
Cream gloss 135
Red gloss 108
Yellow gloss 120
Forest green gloss 143
Green gloss 107
Dark blue gloss 151
Royal blue gloss 154
Blue gloss 174
Duckegg Blue gloss 141
Pink gloss 188
Orange gloss 180
Purple gloss 157

Frame

Matt

Black matt 105
White matt 101
Grey matt 131
Light grey matt 133
Navy blue matt 152

Frame

Special

Chrome mirror 300
Sparkle silver 170
Metallic black matt 104
Metallic grey matt 130
Black sandpaper 201
Grey sandpaper 202
White sandpaper 200

XN150

Powder coated components

No Cost

(Frame colour if not specified)

Specify component colour code from above list

FXAXN030

Brushed Frame

\$582.00

(All parts black unless XN165 chosen)

XN160

Brushed Components

\$125.00

FXAXN050

Polished Frame

\$1,350.00

(All parts black unless XN170 chosen)

XN170

Polished Components

\$219.00

Accent Colour

Note: Accent colour includes forks only

Black accent

Silver accent

Red accent

Blue accent

Purple accent

Upholstery

Seat upholstery

XP040 Seat sling strap adjustable (Std)

XR050 Carbon fibre seat plate

\$494.00

Backrest upholstery

mm (Std XY020 - 20mm)

XP050 Airtech backrest strap adjustable

XP060 Backrest upgrade to full leather
(Black)

\$515.00

XP070 Backrest upgrade to leather corners
(Black)

\$125.00

RGK Logo colour

White (std)

Green

Blue

Individual:

Red

Orange

Purple

Model Logo colour

White (std)

Green

Blue

Individual:

Red

Orange

Purple

Backrest Accessories

XP110 Backrest with zip pocket (N/A with sub4 backrest)

\$131.00

XP170 Individual backrest embroidery (Please provide JPEG)

\$247.00

Cushion

XP080	Cushion foam (Select density and size below)	\$125.00
	Soft (CMS040)	Medium (CMS050)
	50mm	75mm
		100mm

Cushion cover upgrades

XP085	Nylon cushion cover (Black)	Std
XP090	Air lite cushion cover (Black)	\$101.00
XP100	Leather cushion cover (Black)	\$262.00

Jay Products

Jay Cushion
(Price dependant on size/option)

Jay order form attached

Whitmyer Headrest
(Price dependant on size/option)

Whitmyer order form attached

Accessories

XR080	Frontwheel assembly (only with tubular and angle adjustable style footrest)	\$1,085.00
XP120	Frame protectors neoprene zip fasten (Black)	\$87.00
XP420	Frame protectors nylon velcro fasten (Black)	\$83.00
XR010	Seat pouch velcro attach (Black)	\$61.00
XR200	Rucksack with RGK logo (Black)	\$156.00
XR100	FX Carry bag	\$344.00
XR090	Tool kit	Std
XS060	Calf strap	Std



1000 Creditstone Road, Unit 2, Concord ON L4K4P8

Customer Service: (800) 263-3390 E-mail: cscanada@sunmed.com www.sunrisemedical.ca
Specifications subject to change without notice. © 2026 Sunrise Medical Canada Inc. CN-100055-CAN REV (12).