



Quickie QS5 X - AADL

October 2026

EIF7



Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ Prov: _____ PC: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ Prov: _____ PC: _____

Ship To Phone: () _____ Fax: () _____

a Option # Definition Price a Option # Definition Price

1. FRAME

Base model

☐ 170WM112 Quickie QS5 X \$3,344.00

Frame Style

See figure 1. Measure from front frame to floor

☐	170FS31	Swing Away Frame	NC
☐	170FS43	Fixed Front 75°	NC
☐	170FS44	Fixed Front 85°	NC

Frame Inset (Fixed Front Only)

See FIGURE 2 Frame Inset is the distance from the outside of the seat tube to the outside of the footrest. Footrest width = frame width minus 2. N/A for 14" and 15" Wide using composite angle adj footrest. 75° only available with 0" inset.

☐	170FF1	0"	NC
☐	170FF2	1"	NC

Transit Option

ISO 7176-19 Transit Approved

☐ 170TO49 Transit System NC

Frame Options

Ergo Touch Points is a rubber cover on the knee bend area of the frame tubes

☐	1709A11	Ergo Touch Points Swing Away hangers or Fixed front	\$198.00
☐	1709A12	Frame Fold-N-Lock mechanism	NC
☐	1709A06	Twin LED lights	\$310.00



QS5 X - Swing Away



QS5 X - Fixed Front

2. SIZE / SEATING

Frame Width

See FIGURE 2 Measured from outside of seat tube to the outside of opposite seat tube.

☐	170W4	13"	NC
☐	170W5	14"	NC
☐	170W6	15"	NC
☐	170W7	16"	NC
☐	170W8	17"	NC
☐	170W9	18"	NC
☐	170W10	19"	NC
☐	170W11	20"	\$308.00

Seat Depth

See Seat Squeeze Matrix in Appendix section. Measured from back canes to front edge of upholstery. Frame Fold-N-Lock mechanism NA with 13"-14".

☐	170SD6	13"	NC
☐	170SD7	14"	NC
☐	170SD8	15"	NC
☐	170SD9	16"	NC
☐	170SD10	17"	NC
☐	170SD11	18"	NC
☐	170SD12	19"	NC
☐	170SD13	20"	\$308.00

a	Option #	Definition	Price	a	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

2. SIZE / SEATING (CONT)

Seat Sling

Standard seat sling includes utility bag sewn at the bottom of upholstery.

o	170SL1	Standard Seat Sling	NC
o	170SL6	Omit Seat Sling	NC




Figure 1



Figure 2

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

o	170JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width:	Depth:

3. CASTERS / SEAT HEIGHT

Front Seat Height

o	170SH41	14.0"	NC	o	170SH18	18.0"	NC
o	170SH42	14.5"	NC	o	170SH20	18.5"	NC
o	170SH12	15.0"	NC	o	170SH21	19.0"	NC
o	170SH13	15.5"	NC	o	170SH23	19.5"	NC
o	170SH14	16.0"	NC	o	170SH24	20.0"	NC
o	170SH15	16.5"	NC	o	170SH26	20.5"	NC
o	170SH16	17.0"	NC	o	170SH27	21.0"	NC
o	170SH17	17.5"	NC	o	170SH29	21.5"	NC
				o	170SH36	22.0"	NC

Caster Wheels (option A)

o	170CW3	4" Low Profile Polyurethane	NC
o	170CW44	4" x 1.25" Semi Pneumatic	\$110.00
o	170CW5	5" Low Profile Polyurethane	NC
o	170CW7	5" x 1.5" Semi-Pneumatic	\$110.00
o	170CW34	6" Polyurethane	NC
o	170CW33	6" Pneumatic	\$110.00
o	170CW35	6" x 1.5" Semi Pneumatic	\$110.00
o	170CW16	7" Polyurethane	NC
o	170CW32	8" Polyurethane	NC

Caster Fork

o	170CF5	Multiple Position	NC
o	170CF69	Multiple Position Single Sided Fork	\$176.00
o	170CF53	Frog Legs 0 weight to 100 lbs.	\$574.00
o	170CF54	Frog Legs 101 lbs. weight to 190 lbs.	\$574.00
o	170CF55	Frog Legs 191 lbs. weight to 265 lbs.	\$574.00

Caster Wheels colors (option B)


Casters are silver color unless you choose color option. Refer to Anodize package section on Page 7 for Anodize color selection. Match Ano color pkg pick is N/A with Black option.

o	170CW64	4" x 1.5" Aluminum Soft Roll	\$110.00
o	170CW66	5" x 1.5" Aluminum Soft Roll	\$110.00
o	170CW68	6" x 1.5" Aluminum Soft Roll	\$110.00
	Colors:	o Match anodize color package	Package price
		o Match frame color	NC
		o Black anodize caster only	\$68.00

Multi Position

Single side

Frog Leg



Seat Squeeze Matrix

Seat Depth	Max Seat Squeeze
13"-15"	2.0"
16"-19"	2.5"
20"	3.0"

<i>a</i>	Option #	Definition	Price	<i>a</i>	Option #	Definition	Price
4. REAR WHEELS / SEAT HEIGHT							

Rear Seat Height							
o	170RH13	14.0"	NC	o	170RH19	17.0"	NC
o	170RH14	14.5"	NC	o	170RH26	17.5"	NC
o	170RH15	15.0"	NC	o	170RH20	18.0"	NC
o	170RH16	15.5"	NC	o	170RH27	18.5"	NC
o	170RH17	16.0"	NC	o	170RH21	19.0"	NC
o	170RH18	16.5"	NC	o	170RH37	19.5"	NC
				o	170RH22	20.0"	NC

Axle			
o	170AX1	Quick Release Steel	NC
o	170AX11	Quad Quick Release	\$182.00

Rear Wheel Type			
Max weight capacity on Proton wheels is 220 lbs.			
o	170RW103	20" Mag Wheel 5 Spoke	NC
o	170RW104	22" Mag Wheel 5 Spoke	NC
o	170RW102	24" Mag Wheel 5 Spoke	NC
o	170RW90	20" Lite Spoke Wheel	\$320.00
o	170RW91	22" Lite Spoke Wheel	\$320.00
o	170RW92	24" Lite Spoke Wheel	NC
o	170RW93	25" Lite Spoke Wheel	NC
o	170RW94	26" Lite Spoke Wheel	\$126.00
o	170RW95	24" Spinergy SPOX Black (<i>silver hub</i>)	\$864.00
o	170RW96	25" Spinergy SPOX Black (<i>silver hub</i>)	\$864.00
o	170RW98	24" Spinergy LX Black (<i>silver hub</i>)	\$1,111.00
o	170RW99	25" Spinergy LX Black (<i>silver hub</i>)	\$1,111.00
o	170RW148	24" Proton Wheel	\$218.00
o	170RW149	25" Proton Wheel	\$218.00
o	170RW25	Omit Rear Wheel and Axle *	NC

Spinergy Spoke Colors (Optional)				
<i>Please select one of the following Spinergy wheel options first. Includes black hub.</i>				
Spinergy Spox spokes color:				
☐ Black	☐ Blue	☐ Red	☐ White	☐ Yellow
Spinergy LX spokes color:				
☐ Black	☐ Blue	☐ Green	☐ Orange	☐ Pink
☐ Purple	☐ Red	☐ White	☐ Yellow	

Rear Wheel Tire			
o	170RT1	Pneumatic	NC
o	170RT24	Pneumatic Airless Insert - Infinergy® Lightweight - <i>E-TPU material</i>	\$247.00
o	170RT3	Full Poly	NC
o	170RT38	SHOX Performance Poly	\$146.00
o	170RT25	Primo Orion Poly	\$271.00
o	170RT39	High-Pressure Clincher Black	\$105.00
o	170RT21	Schwalbe Marathon Plus Evolution	\$157.00

Handrim	
---------	--

Natural Fits available in 22", 24", and 25". Natural Fit w/ 5 spoke MAG only available in 24". Not all handrims available on all wheels. See page 6 for matrices.

o	170HR1	Aluminum Anodized	NC
o	170HR29	Aluminum Powder coated colored	\$95.00
o	170HR4	Plastic Coated	\$146.00
o	170HR22	Natural Fit - Standard Grip	\$374.00
o	170HR23	Natural Fit - Super Grip	\$434.00
o	170HR39	Ellipse 3R	\$550.00
o	170HR40	Carbolife - Gekko	\$575.00
o	170HR41	Carbolife - Curve L Tetra Grip	\$625.00
o	170HR28	Q-Grip Neoprene Coated	\$471.00
o	170HR7	Omit Handrim <i>Option may not be available on all wheels</i>	NC

Camber			
o	170CI1	0°	NC
o	170CI13	2° (<i>default</i>)	NC
o	170CI3	4°	NC

Center of Gravity Preset			
o	170CC20	1.0" (<i>default</i>)	NC
o	170CC15	1.5"	NC
o	170CC17	2.0"	NC
o	170CC18	2.5"	NC
o	170CC19	3.0"	NC

5. WHEEL LOCKS

Wheel Lock			
No Gos may exist depending on wheel size and position. Ergo & compact N/A w/ 13"-14" frame depth			
o	170WL1	High Mount Push	NC
o	170WL2	High Mount Pull	NC
o	170WL45	Ergo Scissor	NC
o	170WL88	Fold Down Push (Knee lever)	\$110.00
o	170WL30	Compact	\$162.00
o	170WL10	Omit	NC

Wheel Lock Option			
o	170WL12	Do Not Mount	NC

Extension Handle			
N/A w/ ergo scissor or compact wheel locks			
o	170EH7	Extension Handle (Pair)	\$73.00

a	Option #	Definition	Price	a	Option #	Definition	Price			
6. HANGER & FOOTREST										
Swing Away Hangers				Footrest Extension Tubes						
Swing Awat frame only				See Extension matrix for leg length range						
o	170H145	60° Swing In/Out Hanger	NC	o	170E1	Short	NC			
o	170H139	70° Swing In/Out Hanger	NC	o	170E2	Medium	NC			
o	170H141	80° Swing In/Out Hanger	NC	o	170E5	Long	NC			
o	170H160	Elevating Hanger - Ext Mount.	\$339.00							
o	170H159	Articulating Hanger - Ext Mount	\$526.00							
o	170H32	Omit Hanger and Footplate	NC							
Footplates										
Heel loop or leg strap not included, See Footplate Options Section.										
Platform footplates flip up (sideways).										
o	170F1	Composite	NC							
o	170F271	Composite Angle Adj-(16"-22"w)	\$176.00							
o	170F5	Aluminum (ELR & ALR)	NC							
o	170F3	Alum Adult Angle Adjustable - Ext Mount	\$176.00							
o	170F259	Multi Angle Adustable o Black NC o Std Frame color\$51	\$176.00							
o	170F50	Platform Footrest o Black NC o Std Frame color\$51	\$152.00							
o	107F335	Platform Ultra Lite Footrest	\$176.00							
o	107F336	Platform Carbon Fiber footrest	\$151.00							
o	170F7	Omit Footplate and Ext. Tube	NC							
Lower Leg Length (LLL)										
Lower Leg Length is measured from seat sling to the top rear of the footrest. See Appendix on Page 7. A minimum 2"ground clearance is required below the footrest. If the selected configuration results in less than 2"of ground clearance the footrest will be raised to 2"ground clearance, which may alter the lower leg length measurement. Available lower leg length is dependent on frame selection including adapter spacers. Reference LLL matrix in page 7										
o	10"	o 11"	o 12"	o 13"	o 14"	o 15"	o 16"	o 17"	o 18"	o 19"
7. BACKREST / SEATING										
Back Type					Backrest Upholstery					
o	107BT123	Standard Fix (90°)	NC	EXO standard and Tension adjustable upholstery binding is always black.						
o	170BT47	Angle Adjustable Back (80° - 105°)	NC	o	170BU95	Standard Upholstery EXO	NC			
Select desired setting			o 80° o 85° o 90° o 100°	Color Binding:	o Silver	o Orange	o Red	o Blue		
o	107BT124	Fixed Half Fold down (to the front)	\$246.00	o	170BU70	Backrest Tension Adjust Ballistic	\$190.00			
o	170BT43	Lumbar curve bend	NC	o	170BU89	Backrest Tension Adj EXO	\$449.00			
Select desired setting			o 5° forward o 5° backward	Color Binding:	o Silver	o Orange	o Red	o Blue		
Back Height Setting					o	170BU90	Backrest Tension Adj EXO Pro 3DX	\$449.00		
Half Fold Down N/A with the 10" - 13" back heights.					Color Binding:	o Silver	o Orange	o Red	o Blue	
Lumbar curve backrest with integral push handles available with 12" thru 20" back heights. Lumbar curve backrest with Fold Down or Omit push handles available with 14" thru 20" back heights.					o	170BU14	Omit Back Upholstery	NC		
o	170B161	Setting	10 "- 12" Back Height o 10" o 11" o 12"	NC						
o	170B162	Setting	12 "- 14" Back Height o 12" o 13" o 14"	NC						
o	170B163	Setting	14 "- 16" Back Height o 14" o 15" o 16"	NC						
o	170B164	Setting	16 "- 18" Back Height o 16" o 17" o 18"	NC						
o	170B165	Setting	18 "- 20" Back Height o 18" o 19" o 20"	NC						
Push Handles					JAY Back (Optional)					
Std EXO Uph N/A with Fold down. Height Adjustable push handles available with EXO / EXO Pro Tension Adjustab					To order a J3 Back please select Omit Back Upholstery. See J3 Back Order Form for Part #'s					
o	170PH6	Integral Push Handles	NC	o	170JB01	J3 Back (See J3 Order Form)	Price List			
o	170PH5	Fold Down Push Handles Available w/12"-20" back height	\$168.00							
o	170PH21	Height adjustable push Handles Available w/ 12"-20" back height. N/A w/Half fold down back	\$201.00							
o	170PH25	Height adjustable push Handles for JAY Backs N/A w/Half fold down back	\$207.00							
o	170PH22	Omit Push Handles - Straight Post	NC							

a	Option #	Definition	Price	a	Option #	Definition	Price
8. ARMREST / SIDE GUARDS							

Padded Swing Away

Some limitations when used w/ angle adj back in the rear most setting

Please select a setting, will default to 11" if no selection is made.

o	170AA162	Padded Swing-Away Fixed height (14")	NC
o	170AA18	Padded Swing-Away STD (9" to 14")	\$278.00
Setting o 9" o 10" o 11" o 12" o 13" o 14"			

Armrest - Height Adjustable

o	170AA11	Single Post Hgt. Adj. Std (8.5"-13.5")	\$339.00
o	170AA29	Single Post Hgt. Adj. Low (7"-11")	\$339.00
o	170AA65	Flip Back Height Adj LITE (8.5"-13.5")	\$381.00

Armpad - For Height Adjustable

Ergo armpad only available with Single post Height adj & Flip back armrest

o	170AM23	Desk Length - Classic (10")	NC
o	170AM24	Full Length - Classic (14")	NC
o	170AM27	Ergo Pad (9.5"-13.5")	NC

Omit Armrest

o	170AA14	Omit Armrest	NC
---	---------	--------------	----

9. KEY OPTIONS & ACCESSORIES

Anti-Tip

o	170AT1	Anti Tips (Pair)	\$121.00
o	170AT1	Anti Tip swing away	
		o 170AT88 Right Mount	\$152.00
		o 170AT89 Left Mount	\$152.00
		o 170AT90 Pair	\$305.00

Additional Accessories

o	170A1	Backpack - Adult	\$112.00
o	170A2	Backpack - Kids	\$105.00
o	1701A01	Seat Pouch Black	\$83.00

Additional Cover Accessories

o	1708A02	Gel Cover Hanger/Frame	\$162.00
o	1708A01	Gel Cover Swing Away Arm	\$162.00

Positioning Belts

o	170ZZ15	JAY 2PT Belt 1.5" Push Btn 7.5" Pad	\$80.00
o	170ZZ26	JAY 2PT Belt 2" Push Btn 9.5" Pad	\$126.00
o	170ZZ14	JAY 2PT Belt 1.5" Push Btn No Pad	\$51.00
o	170A118	Belt 2" Aircraft Buckle Pad w/Clamp	\$204.00
o	170A120	Positioning Belt Slide Clamps	\$61.00

Side Guard

Side guards N/A w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist

o	170SG33	Plastic Standard	\$188.00
o	170SG29	Plastic Low	\$188.00
o	170SG35	Aluminum	\$315.00
o	170SG36	Aluminum w/ Fender	\$315.00
o	170SG12	Carbon Fiber	\$399.00
o	170SG41	Carbon Fiber Long w/LG Fender	\$602.00

Side Guard - Color Option

Available with Aluminum or Aluminum with Fender. Not applicable to carbon fiber

o	170GE01	Black	NC
o	170GE02	Matching Frame- Standard Colors	NC
o	170GE03	Matching Frame- Personalized Color	\$68.00

Spoke Guards

Only available w/ Lite Spoke wheels. Not available with 5 spoke MAG wheels

o	1704A01	Quickie Black Spoke Guards	\$157.00
o	1704A02	Quickie Clear Spoke Guards	\$157.00

Miscellaneous

o	1708A03	Caddy	\$89.00
o	170A2100	Crutch Holder	\$115.00
o	1708A21	Tool Kit	\$43.00
o	170A999	Quickie Water Bottle & Holder	\$81.00
o	1708A47	Mobile phone pocket	\$34.00
o	170FO4	Transit wheels (Pair)	NC
o	170AT	Tip Assist	
		o 170AT83 Right Mount	\$20.00
		o 170AT84 Left Mount	\$20.00
		o 170AT91 Pair	\$41.00

a	Option #	Definition	Price	a	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

10. COLOR

Personalized X-Brace Color

Option has an extended leadtime. Frame colors only, NA with Personalized frame colors.
If no color pick crossbrace will be Gloss Black.

o	170GS15	Crossbrace Color package	NC
---	---------	--------------------------	----

Standard Frame Colors

Handrim color option only available if Aluminum Coated colored option was selected.

		Frame	X-brace	Handrim
170C211	Celestial Sapphire	o	o	o
170C212	Desert Sage	o	o	o
170C213	Desert Tan	o	o	o
170C214	Gloss Black	o	STD	o
170C215	Matte Black	o	o	o
170C216	Metallic Shadow	o	o	o
170C217	Midnight Blue	o	o	o
170C218	Mystic Purple	o	o	o
170C219	Petal Pink	o	o	o
170C220	Pink Power	o	o	o
170C221	Red Pepper	o	o	o
170C222	Royal Blue	o	o	o
170C223	Salute Green	o	o	o
170C224	Scarlet Noir	o	o	o
170C225	Silver Ore	o	o	o
170C226	Solar Orange	o	o	o
170C227	Starlight White	o	o	o
170C228	Storm Gray	o	o	o
170C229	Turbo Green	o	o	o
170C230	Velocity Yellow	o	o	o

NEW!

Legacy - Frame Colors

o	170C141	Black Cherry	NC
o	170C10	Black Opal	NC
o	170C14	Candy Purple	NC
o	170C142	Electric Blue	NC

Personalized Frame Color

Only frame will have the personalized color; Crossbrace, handrim, footplate and Aluminum side guard not included (unless selected).

o	170C112	Desert Camouflage	\$368.00
o	170C43	Mossy Oak Camouflage	\$368.00
o	170C42	Stars & Stripes	\$368.00
o	170C108	Zebra (Choose base color). N/A w/ Gloss Black and Matte Black.	\$368.00

Anodized Color Package

Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog suspension Legs.

o	170EB08	Blue	\$212.00
o	170EB13	Gold	\$212.00
o	170EB09	Lime Green	\$212.00
o	170EB10	Purple	\$212.00
o	170EB11	Red	\$212.00
o	170EB12	Ti Grey	\$212.00

Frame Decal Colors

Applies to front caster link decals and aluminum & carbon fiber sideguards.

o	170D16	Black	NC
o	170D17	Blue	NC
o	170D18	Orange	NC
o	170D19	Red	NC
o	170D20	Silver	NC

Built For Me

To inquire about any unique personalized for Quickie products, contact:
Phone: (800) 263-3390, Email: cscanada@sunmed.com



11. APPENDIX

Wheel - Handrims Matrix

Handrim	Rear Wheel			
	5 Spoke Mag	Lite Spoke	Spinergy	Proton
Aluminum Ano / Painted	20,22,24	20,22,24,25,26	24,25	24,25
Plastic Coat	20,22,24	22,24	24,25	-
Natural Fit - Standard Grip	24	22,24,25,26	24,25	-
Natural Fit - Super Grip	24	22,24,25,26	24,25	-
Natural Fit - No Thumb Grip	24	22,24,25,26	24,25	-
Surge - Oval w/Gription Strip	-	22,24,25,26	24,25	-
Surge LT - Oval w/Gription Strip	-	22,24,25,26	24,25	-
Q Grip Neoprene Coated	-	24,25	24,25	-
Omit Handrim	20,22,24	-	24,25	-

Wheel - Tire Matrix

Tire	Rear Wheel			
	5 Spoke Mag	Lite Spoke	Spinergy	Proton
Pneumatic	20,22,24	20,22,24,26	24	24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	24	-
Full Poly	20,22,24	20,22,24	24	-
SHOX Performance Poly	20,22,24	20,22,24,25	24,25	-
Primo Orion Poly	20,22,24	20,22,24	24	-
High-Pressure Clincher Bk	24	22,24,26	24	24
Schwalbe Marathon	22, 24	24,25,26	24,25	24,25
20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm				

11. APPENDIX

Note 1: Fork Hole 1+ indicates that a 10mm spacer is used to achieve additional Front Seat Height (FSH).

Note 2: FSH values are rounded estimates and assume flat seat heights (0 degree seat angle). Increasing seat angle can affect final FSH.

FSH for 4" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	14.0	14.5	15.0
	Multi 115	-	14.5	15.0	15.5	16.0
	Multi 131	-	15.0	15.5	16.0	16.5
	Multi 158	-	16.0	16.5	17.0	17.5
	SS Low	-	-	14.0	14.5	15.0
	SS Std	-	14.5	15.0	15.5	16.0
	Frog Legs	-	15.5	16.0	16.5	-
Short	Multi 90	-	-	15.5	16.0	16.0
	Multi 115	-	16.0	16.5	17.0	17.0
	Multi 131	-	16.5	17.0	17.5	18.0
	Multi 158	-	17.5	18.0	18.5	19.0
	SS Low	-	-	15.5	16.0	16.0
	SS Std	-	16.0	16.5	17.0	17.0
	Frog Legs	-	17.0	17.5	18.0	-
Tail	Multi 90	-	-	16.5	17.0	17.5
	Multi 115	-	17.0	17.5	18.0	18.5
	Multi 131	-	17.5	18.0	18.5	19.0
	Multi 158	-	18.5	19.0	19.5	20.0
	SS Low	-	-	16.5	17.0	17.5
	SS Std	-	17.0	17.5	18.0	18.5
	Frog Legs	-	18.0	18.5	19.0	-

FSH for 5" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	15.5	16.0	16.0
	Multi 131	15.0	15.5	16.0	16.5	17.0
	Multi 158	16.0	16.5	17.0	17.5	18.0
	SS Low	-	-	-	-	-
	SS Std	-	15.0	15.5	16.0	16.0
	Frog Legs	-	-	16.5	17.0	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	16.5	17.0	17.5
	Multi 131	16.5	17.0	17.5	18.0	18.0
	Multi 158	17.5	18.0	18.5	19.0	19.5
	SS Low	-	-	-	-	-
	SS Std	-	16.0	16.5	17.0	17.5
	Frog Legs	-	-	17.5	18.0	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	18.0	18.5	19.0
	Multi 131	17.5	18.0	18.5	19.0	19.5
	Multi 158	18.5	19.0	19.5	20.0	20.5
	SS Low	-	-	-	-	-
	SS Std	-	17.5	18.0	18.5	19.0
	Frog Legs	-	-	19.0	19.5	-

FSH for 6" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	16.0	16.5
	Multi 131	-	16.0	16.5	17.0	17.0
	Multi 158	16.5	17.0	17.5	18.0	18.5
	SS Low	-	-	-	-	-
	SS Std	-	-	15.5	16.0	16.5
	Frog Legs	-	-	-	17.0	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	17.5	18.0
	Multi 131	-	17.5	18.0	18.5	18.5
	Multi 158	18.0	18.5	19.0	19.5	19.5
	SS Low	-	-	-	-	-
	SS Std	-	-	17.0	17.5	18.0
	Frog Legs	-	-	-	18.5	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	19.0	19.5
	Multi 131	-	18.5	19.0	19.5	20.0
	Multi 158	19.0	19.5	20.0	20.5	21.0
	SS Low	-	-	-	-	-
	SS Std	-	-	18.5	19.0	19.5
	Frog Legs	-	-	-	20.0	-

FSH for 7" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	17.0	17.5	18.0
	Multi 158	17.0	17.5	18.0	18.5	19.0
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	18.5	19.0	19.0
	Multi 158	18.5	19.0	19.5	20.0	20.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	19.5	20.0	20.5
	Multi 158	19.5	20.0	20.5	21.0	21.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-

FSH for 8" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	17.5	18.0	18.5
	Multi 158	-	18.0	18.5	19.0	19.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	19.0	19.5	19.5
	Multi 158	-	19.5	20.0	20.5	20.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	20.0	20.5	21.0
	Multi 158	-	20.5	21.0	21.5	22.0
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-

Lower Leg Length Matrix				
Note for Swing Away Hangers only: subtract 1" from the max achievable LLL for long extension tubes if the Multi Angle Adjustable footplate is selected.				
Hanger	Extension	LLL (inches) (measured from top of seat sling to top of footrest)		Min Allowable FSH (due to 2" ground clearance requirement)
		Min	Max	
60° Swing In/Out	Short	10	14	No lower FSH limit
	Med	11	16	
	Long	12	18	
70° Swing In/Out	Short	10	16	No lower FSH limit
	Med	12	17	
	Long	13	19	
80° Swing In/Out	Short	11	15	No lower FSH limit
	Med	12	17	
	Long	14	19	
ELR	Short	11	16	No lower FSH limit
	Med	11	18	
	Long	13	19	
ALR	Short	11	16	No lower FSH limit
	Med	13	18	
	Long	15	19	
75° Fixed Front	Short	13	18	15.5
	Med	13	20	
	Long	14	20	
85° Fixed Front 0" inset	Short	13	18	15.5
	Med	13	20	
	Long	15	20	
85° Fixed Front 1" inset	Short	13	18	15.5
	Med	15	20	17.0
	Long	17	20	19.0



1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8