

LIVE
YOUR
LIFE



Because the world's not flat™





INNOVATION AND OUR THINKING

Magic Mobility is a leading Australian power wheelchair designer, developer and manufacturer. Our focus is on creating wheelchairs tailored to individual needs. We're determined to push the boundaries of what power wheelchairs can achieve and above all, provide riders with choice.

Our all terrain wheelchairs help bring more aspirational goals within reach because they take you where you need to go and help you stay comfortable for longer.

The modular design enables easy adaption with seating and accessories as your goals and needs change. This makes a Magic Mobility power wheelchair cost less to own and operate over its lifetime.





The Magic Mobility Extreme X8 provides the freedom you've always wanted. The X8 is your partner to almost any destination. It can handle sand, snow and every terrain you can throw at it. With unrivaled off-road performance and reliability, the X8 is perfect for outdoor adventurers, explorers and rural workers. If the great outdoors is where you want to be, the Extreme X8 is waiting to take you there.

The right wheelchair for: *

- Extreme off-road adventures
- Everyday off-road country living, working and rural productivity
- Active outer urban living



"WE HAVE A SIX ACRE PROPERTY AND THE EXTREME X8 ALLOWS ME TO GET TO ALL THE UNEVEN TERRAINS AND SEE PARTS OF THE FARM THAT I HAVE NOT BEEN ABLE TO ACCESS IN FOUR YEARS! THANK YOU FOR GIVING ME MY FREEDOM AND INDEPENDENCE."

Michael Trost
Maryland, USA



Our most popular powered wheelchair worldwide, the Magic Mobility Frontier V6 All Terrain is ideal for the outdoors, yet its size makes it a capable indoor performer. If performance and versatility is what you are after, the Frontier V6 All Terrain will impress those who prefer trails over sidewalks. Choose Off-Road or Urban wheels.

The right wheelchair for: *

- Adventurous off-road recreation
- Everyday off-road country living, working and rural productivity
- Urban outdoor living and working
- Roomy indoor spaces



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE TIME I'VE HAD IT. IT HAS NEVER LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays
Texas, USA



The Magic Mobility Magic 360 is one of the few compact powerchairs that is robust enough for off-road, maneuverable enough for indoors, and versatile enough for everything in between. The Magic 360 is compact enough to get around your home and is full of all-terrain technology to take you where you want to go. The chassis, suspension, motors, drive wheels, and batteries were all developed with off-road performance in mind. Pair with Off-Road wheels to enjoy the great outdoors or Urban tires to take on city living.

The right wheelchair for: *

- Adventurous off-road recreation
- Everyday off-road country living, working, and rural productivity
- Urban outdoor living and working
- Indoor spaces



"MY MAGIC 360 HAS GIVEN ME MORE CONFIDENCE BECAUSE IT'S VERY STABLE AND IT'S MUCH SAFER OFF-ROAD, WHICH MEANS IT'S GOOD FOR OUTDOOR MUSIC CONCERTS AND FESTIVALS. I CAN GET INTO SOME SHOPS AND GET AROUND WITHOUT KNOCKING THINGS OFF THE SHELVES."

Aidan Magee
Australia



*Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks or the sea.

"WE'VE BEEN CREATING MADE-TO-ORDER WHEELCHAIRS FOR OVER 25 YEARS, SO THAT YOU CAN LIVE THE LIFE YOU CHOOSE"

Ashley Daff
R&D Director



ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

BUILT TO LAST

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

PROUDLY AUSTRALIAN DESIGNED AND ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

SPECIFICATIONS	EXTREME X8	FRONTIER V6	MAGIC 360
Max User Weight	400 lbs., 340 lbs. with Elevate or Tilt		350 lbs.
Seat-to-Floor Height	18.5", 19.5", 20.5"		17", 18", 19"
Length x Width	40.5" x 28"	Off-Road 41.5" x 28" Urban 41.5" x 27.3"	Off-Road 39" x 26" Urban 39" x 24"
Drive Wheel Position	4 Wheel Drive	Mid Wheel Drive	Mid Wheel Drive
Max Speed	6 mph		
Battery Size	Group 24	Group 24 or 34	Group 24
Overall Product Weight	Starting at 320 lbs.		
Drive Wheels	14" All Terrain	14" Off-Road or 14" Urban	
Caster Wheels	N/A	Pneumatic or Solid Inserts	Solid Inserts
Power Seat Options	50° of CG Tilt/Power Recline/12" Seat Elevator/Power Elevating Legrest		
Colors	Midnight Black, Ocean Blue, Fire Engine Red, Clover Green, Tangerine and Camouflage	Black, Red, Mossy Green, and Royal Blue	Lucent Lime, Midnight Black, Horizon Blue, Dusk Purple, Vivid Pink, Desert Orange, Woods Camouflage, Desert Camouflage, Carbon Fiber, and American Flag

PUNCTURE RESISTANT TIRE LINERS

Puncture resistant tire liners provide peace of mind for those concerned about leaks and punctures when going off-road. An Australian rubber manufacturer developed an aramid fiber reinforced rubber compound for puncture resistance. The liners are included on all tires.



FULL LIGHTING PACKAGES

Mounted on flexible rubber brackets that resist bumps, our lights use an extremely bright and highly efficient 24V LED lighting system that runs directly from the wheelchair batteries.



RAM® MOUNT-PHONE HOLDER

The X-Grip® provides simple mounting and a secure hold for phones and tablets of all sizes. In extreme environments, it ensures added peace of mind.



ACCESSORY CHARGERS


Our accessory chargers can be fitted to all Magic Mobility power wheelchairs which use the R-Net system. Two chargers are available:

- 1 x 12V & 2 x USB ports
- 2 x USB ports



SEATING OPTIONS

From MPS seating to a traditional rehab seat, a whole range of powered functions are available to support your needs.



See more online.



DESIGNED IN AUSTRALIA



ASSEMBLED IN THE USA



Sunrise Medical (US) LLC, Fresno, CA 93727
MK-129945 Rev.C © 03.2024

800.333.4000 www.SunriseMedical.com

