



# JAY<sup>®</sup> Zip Pediatric Backrest

April 2025

## Mark For:

Date: \_\_\_\_\_  
Dealer Acct #: \_\_\_\_\_  
Dealer: \_\_\_\_\_  
Dealer Contact: \_\_\_\_\_  
Dealer Address: \_\_\_\_\_  
Dealer City: \_\_\_\_\_ Prov: \_\_\_\_\_ PC: \_\_\_\_\_  
Dealer Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_  
Confirmation Email: \_\_\_\_\_  
Confirm Via: ☐ Fax ☐ Email

## Submitting for:

☐ Quote ☐ Order

PO#: \_\_\_\_\_

## ADDITIONAL SHIPPING INFORMATION

Ship To: \_\_\_\_\_  
Attention: \_\_\_\_\_  
Address: \_\_\_\_\_  
Address: \_\_\_\_\_  
Ship To City: \_\_\_\_\_ Prov: \_\_\_\_\_ PC: \_\_\_\_\_  
Ship To Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

## JAY ZIP BACK

☐ JZIPBXXY JAY Zip Backrest All Sizes \$1,015.00





NOTE: Select appropriate assembly number below.  
165 lbs Weight Capacity



BACKREST WIDTH	8"	10"	12"	14"
FITS BACK CANE WIDTHS	8"-9"	10"-11"	12"-13"	14"-15"
MAX. USER TRUNK WIDTH	7"	9"	11"	13"
ACTUAL SHELL WIDTH	5.6"	7.6"	9.6"	11.6"
DEPTH OF CONTOUR	1"	1.5"	2"	2.5"

## EACH BACK SHIPS WITH THE FOLLOWING:

Ultra lightweight mounting hardware  
2" of anterior depth adjustment, 2.75" with reversing mounting pins  
10mm Combo Wrench  
4mm Hex Key  
Neoprene Backrest Hardware Covers (Allows for quick release)  
Water resistant inner cover & Ionic+™ outer cover  
Adjustable lumbar support

		Backrest Width			
		8"	10"	12"	14"
Backrest Height	 6" (S)	<input type="checkbox"/> JZIPB08S	<input type="checkbox"/> JZIPB10S	<input type="checkbox"/> JZIPB12S	<input type="checkbox"/> JZIPB14S
	 8" (M)	<input type="checkbox"/> JZIPB08M	<input type="checkbox"/> JZIPB10M	<input type="checkbox"/> JZIPB12M	<input type="checkbox"/> JZIPB14M
	 10" (T)	<input type="checkbox"/> JZIPB08T	<input type="checkbox"/> JZIPB10T	<input type="checkbox"/> JZIPB12T	<input type="checkbox"/> JZIPB14T
	 12" (X)	<input type="checkbox"/> JZIPB08X	<input type="checkbox"/> JZIPB10X	<input type="checkbox"/> JZIPB12X	<input type="checkbox"/> JZIPB14X

The JAY ZIP Backrest has been dynamically tested and meets all criteria for use in motor vehicles according to Section 5.1 of ISO 16840-4 and Section 20, ANSI/RESNA WC/Volume 4.  
Please visit [www.sunparts.us](http://www.sunparts.us) for spare parts information.

## JAY ZIP BACK ACCESSORIES

JAY ZIP Back accessories are optional. All accessories will be factory mounted as much as possible to ensure proper shipping.

### LATERAL SUPPORTS

Use this section for lateral supports.

#### FIXED LATERAL THORACIC SUPPORT

Includes fixed mounting hardware.


<input type="checkbox"/>	JZFX34F	LAT ZIP FIXED 3X4 FLAT	\$335.00	<input type="checkbox"/>	JZFX34C	LAT ZIP FIXED 3X4 CURVE	\$335.00
<input type="checkbox"/>	JZFX45F	LAT ZIP FIXED 4X5 FLAT	\$335.00	<input type="checkbox"/>	JZFX45C	LAT ZIP FIXED 4X5 CURVE	\$335.00

#### SWING-AWAY LATERAL THORACIC SUPPORT

Includes swing-away mounting hardware and all laterals are sold in pairs.

<input type="checkbox"/>	<b>JZSA34F</b>	<b>LAT 3X4 FLAT W/ S/A MOUNT</b>	<b>\$667.00</b>	<input type="checkbox"/>	<b>JZSA34CO</b>	<b>LAT 3X4 CRV 1" OFF W/ S/A MNT</b>	<b>\$667.00</b>
	JZ34F	LAT ZIP 3X4 FLAT	\$298.00		JZ34CO	LAT 3X4 CURVE W/ 1" OFFSET	\$298.00
	JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00		JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00
<input type="checkbox"/>	<b>JZSA45F</b>	<b>LAT 4X5 FLAT W/ S/A MOUNT</b>	<b>\$667.00</b>	<input type="checkbox"/>	<b>JZSA45CO</b>	<b>LAT 4X5 CRV 1" OFF W/ S/A MNT</b>	<b>\$667.00</b>
	JZ45F	LAT ZIP 4X5 FLAT	\$298.00		JZ45CO	LAT 4X5 CURVE W/ 1IN OFFSET	\$298.00
	JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00		JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00
<input type="checkbox"/>	<b>JZSA34C</b>	<b>LAT 3X4 CURVE W/ S/A MOUNT</b>	<b>\$667.00</b>	<input type="checkbox"/>	<b>JZSA34FOA</b>	<b>LAT 3X4 FLAT ADJ 1" OFF W/ S/A</b>	<b>\$667.00</b>
	JZ34C	LAT ZIP 3X4 CURVE	\$298.00		JZ34FOA	LAT 3X4 FLAT W/ ADJ 1" OFFSET	\$298.00
	JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00		JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00
<input type="checkbox"/>	<b>JZSA45C</b>	<b>LAT 4X5 CURVE W/ S/A MOUNT</b>	<b>\$667.00</b>	<input type="checkbox"/>	<b>JZSA45FOA</b>	<b>LAT 4X5 FLAT ADJ 1" OFF W/ S/A</b>	<b>\$667.00</b>
	JZ45C	LAT ZIP 4X5 CURVE	\$298.00		JZ45FOA	LAT 4X5 FLAT W/ ADJ 1" OFFSET	\$298.00
	JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00		JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00
<input type="checkbox"/>	<b>JZSA34FO</b>	<b>LAT 3X4 FLT 1" OFF W/ S/A MNT</b>	<b>\$667.00</b>	<input type="checkbox"/>	<b>JZSA34COA</b>	<b>LAT 3X4 CRV ADJ 1" OFF W/ S/A</b>	<b>\$667.00</b>
	JZ34FO	LAT ZIP 3X4 FLAT W/ 1" OFFSET	\$298.00		JZ34COA	LAT 3X4 CRV W/ ADJ 1" OFFSET	\$298.00
	JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00		JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00
<input type="checkbox"/>	<b>JZSA45FO</b>	<b>LAT 4X5 FLT 1" OFF W/ S/A MNT</b>	<b>\$667.00</b>	<input type="checkbox"/>	<b>JZSA45COA</b>	<b>LAT 4X5 CRV ADJ 1" OFF W/ S/A</b>	<b>\$667.00</b>
	JZ45FO	LAT ZIP 4X5 FLAT W/ 1" OFFSET	\$298.00		JZ45COA	LAT 4X5 CRV W/ ADJ 1" OFFSET	\$298.00
	JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00		JZSAMNT	LATERAL SWING AWAY MOUNT	\$369.00


#### LATERAL MEASURING GUIDE (Fixed and Swing Away)



		Back Width			
DIM	Description	8" Wide	10" Wide	12" Wide	14" Wide
A	<b>Lateral Adjustment Range</b> (Min Spacing - Max Spacing)	3.8" - 4.8"	5.1" - 6.1"	6.8" - 7.8"	8.4" - 9.4"
B	<b>Lateral Height Range</b> (Measured from bottom edge of shell to top of lateral)	6" Tall Back	Not Available		
		8" Tall Back	7"		
		10" Tall Back	9"		
		12" Tall Back	10"		

## HEADRESTS

<input type="checkbox"/>	<b>JZHED6</b>	<b>Headrest w/ Pad - 6" Plush</b>	<b>\$497.00</b>
	JPLSH06	Pad - 6" Plush	\$330.00
	JZIPBHH	Onyx post w/ Axys Hardware	\$167.00
<input type="checkbox"/>	<b>JZHED8</b>	<b>Headrest w/ Pad - 8" Plush</b>	<b>\$497.00</b>
	JPLSH08	Pad - 8" Plush	\$330.00
	JZIPBHH	Onyx post w/ Axys Hardware	\$167.00
<input type="checkbox"/>	<b>JZHEDCINF</b>	<b>Headrest w/ Pad - 6" Contoured</b>	<b>\$458.00</b>
	JCNTR06	Pad - 6" Contoured Cradle	\$291.00
	JZIPBHH	Onyx post w/ Axys Hardware	\$167.00
<input type="checkbox"/>	<b>JZHEDCS</b>	<b>Headrest w/ Pad - 8" Contoured</b>	<b>\$473.00</b>
	JCNTR08	Pad - 8" Contoured Cradle	\$306.00
	JZIPBHH	Onyx post w/ Axys Hardware	\$167.00



ONYX Post w/ AXYS Hdw & Adptr



Plush Pad



Contoured Pad


#### Vanity Flap

Vanity flaps are NOT INCLUDED as a standard option with the JAY ZIP Back.

<input type="checkbox"/>	VF08	8" to 9" Wide	\$88.00	<input type="checkbox"/>	VF12	12" to 13" Wide	\$88.00
<input type="checkbox"/>	VF10	10" to 11" Wide	\$88.00	<input type="checkbox"/>	VF14	14" to 15" Wide	\$88.00

### Spine Align Positioning Components


<input type="checkbox"/>	J3SPNBL	Large Bean Bag	\$42.00
<input type="checkbox"/>	J3SPNLS	Lateral Shim	\$67.00
<input type="checkbox"/>	J3SPNLP	Posterior/Lateral Pelvic	\$67.00



## CHEST STRAPS & HARNESSSES

Chest Straps & Harnesses include guide brackets and hardware.

□	JZ9107		Contoured - Dynamic Pull Down (EI)	\$239.00
□	JZ9100		Contoured - Dynamic Pull Down (XS)	\$239.00
□	JZ9101		Contoured - Dynamic Pull Down (S)	\$239.00
□	JZ9102		Contoured - Dynamic Pull Down (M)	\$239.00
□	JZ9300		Contoured - Strctrd Pull Down (XS)	\$239.00
□	JZ9301		Contoured - Strctrd Pull Down (S)	\$239.00
□	JZ9302		Contoured - Strctrd Pull Down (M)	\$239.00
□	JZ4050		Padded Anterior Trunk Straps (XS)	\$239.00
□	JZ4000		Padded Anterior Trunk Straps (S)	\$239.00
□	JZ4100		Padded Anterior Trunk Straps (M)	\$239.00



## JAY YOUR WAY MODIFICATIONS

*Note: To order options offered via JAY Your Way add an "M" to the front of the part number and follow the prompts.*

For additional modifications visit [www.sunrisemedical.com](http://www.sunrisemedical.com) to view even more!

### COVER COLOR OPTIONS - Rear Center Flap

<input type="checkbox"/>	COCPI		Pink Semi Stretch Cover	NC
<input type="checkbox"/>	COCB		Blue Semi Stretch Cover	NC
<input type="checkbox"/>	COCPU		Purple Cover	NC
<input type="checkbox"/>	COCR		Red Cover	NC
<input type="checkbox"/>	COCB		Toxic Green Cover	NC

## COLOR BACK SHELLS

*Paint your Zip back shell to match any of the Zippie or Quickie Frame Colors!*

<input type="checkbox"/>	XXCBS	Zippie/Quickie Back Shell Paint	NC
--------------------------	-------	---------------------------------	----

Write Desired Back Shell Color Here:

## EMBROIDERY

<input type="checkbox"/> EMBROIDERY	Embroidery (32 Character Max)	NC
-------------------------------------	-------------------------------	----

Location: ☐ Front

Font: ☐ Block ☐ Script

Embroidery Text:

## SPINAL FLUID PAD

<input type="checkbox"/>	SFPS	Small Spinal Fluid Pad	\$61.00
<input type="checkbox"/>	SFPM	Medium Spinal Fluid Pad	\$61.00

*Don't see the option you want? NO PROBLEM!*

*Write your request below and we can evaluate it and provide a quote.*

**Notes:**



Sunrise Medical Canada Inc.

1000 Creditstone Road, Unit 2 · Concord, ON L4K 4P8