

GEMINO 30 PARKINSON GEMINO 30 M PARKINSON

September 2026



Product sold with backrest in Canada.

ORDER FORM

Invoice address or customer number:

Customer number:

Name:

Street:

Town:

Postal code:

Country:

E-mail:

Tel:

Delivery address:

☐ Deliver to Invoice address

Customer number:

Name:

Street:

Town:

Postal code:

Country:

E-mail:

Tel:

Order information

☐ Order

☐ Quote


Order date:

Your reference:

Reimbursement code:

 Customer Service: cscanada@sunmed.com
Base price \$ 1,809.00

DIMENSIONS



	Gemino 30	Gemino 30 M
A Total height	78-100 cm	69-88 cm
B Total length	65 cm	65 cm
C Total width	60 cm	60 cm
D Width between push handles	47 cm	47 cm
E Seat height	62 cm	55 cm
Height folded	80 cm	76 cm
Length folded	65 cm	65 cm
Width folded	23 cm	23 cm
Turning circle	84 cm	84 cm
F Front / rear wheels	20 x 3.6 cm	20 x 3.6 cm
Total weight without basket	7.5 kg	7.4 kg
G Weight basket, standard	0.3 kg	0.3 kg
Max. user weight	130 kg	130 kg
Max. weight in basket	5 kg	5 kg

ORDER FORM GEMINO 30 PARKINSON / GEMINO 30 M PARKINSON

SIZES AND COLOURS		Price \$	Quantity	Total
Gemino 30 Parkinson + HC-9012222 Adj. Comfort Back				
☐ HC-7100070	Gemino 30 Parkinson, silver grey ***NOT AVAILABLE WITH BACKSTRAP***	1,809.00		
Gemino 30 M Parkinson + HC-9012222 Adj. Comfort Back				
☐ HC-7100270	Gemino 30 M Parkinson, silver grey ***NOT AVAILABLE WITH BACKSTRAP***	1,809.00		

ACCESSORIES		Price \$	Quantity	Total
☐ HC-7109011	Crutch holder	106.00		
☐ HC-7109140	Table tray	106.00		
☐ HC-7109060	Umbrella, black	264.00		
☐ HC-7109050	Basket for oxygen cylinder	264.00		
☐ HC-7109095	Light	129.00		
☐ HC-7109100	Infusion holder	264.00		

GEMINO 30 PARKINSON / GEMINO 30 M PARKINSON		Price \$
Base price		
Option price		
Total price		

GEMINO ACCESSORIES

GEMINO

BAG - GEMINO 20 S & GEMINO 30 S



HC-7104007 Bag,
advanced, black/black **\$129.00**

OTHER ACCESSORIES



HC-7109011 **\$106.00**
Crutch holder



HC-7109140 **\$106.00**
Table tray



HC-7109095 **\$129.00**
Light



HC-7109060 **\$264.00**
Umbrella, black



HC-7109050 & HC-7109051 **\$264.00**
Basket for oxygen cylinder*



HC-7109100 **\$264.00**
Infusion holder



Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2 Concord, ON
L4K 4P8, Canada
Tel. +1.800.263.3390
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca