

# GEMINO 30



**GEMINO 30 · GEMINO 30 M · GEMINO 30 S**  
Rollator with unique design, functionality and flexibility

# GEMINO 30 GEMINO 30 M GEMINO 30 S

The Gemino range features a variety of models and accessories, so you can always be sure of finding a Gemino model that suits you and your needs.



The rollator folds up in one easy step and locks automatically into the folded position



IV pole



Cane holder



Light



Tray



Umbrella



Slow-down brake



Bag for oxygen bottle



Bag with lid and shoulder strap

## Safety, support and stability

The unique and very stable design ensures that you can walk securely in all settings. When you need a rest, Gemino has a comfortable and convenient seat with backstrap.

## Easy to take with you

The rollator folds up in one easy step and locks automatically into the folded position.

## Comfortable, correct position for walking

It is easy to adjust Gemino to the correct height. The brake and pushing handles have been carefully designed to be comfortable and easy to use. A memory feature ensures handle heights can always return to the proper setting.

## Safe and easy to use

High-quality integrated reflectors and a reflective band on the bag make Gemino clearly visible to road users. The unique brake design ensures that the brakes stop the wheels very effectively. The design of the curb climber makes it very easy to negotiate obstacles such as curbs and door sills. A unique side deflector makes it easier to navigate doorways, table legs and other furniture, without the risk of the back wheels getting caught.

## Gemino – the award winner

The items with which we surround ourselves help define how we are perceived. The elegant and sleek design of Gemino appeals to both young and old alike. Sunrise is delighted and proud that Gemino has won two prestigious international design awards, Reddot and iF. It also won the German "Öko-test" consumer test.

### Specifications

	Gemino 30	Gemino 30 M	Gemino 30 S
Total height	30.7 - 39.4"	27.2 - 34.6"	25.6 - 30.3"
Total length	25.6"	25.6"	25.6"
Total width	23.6"	23.6"	21.7"
Distance between pushing handles	18"	18"	16"
Seat height	24.4"	21.7"	18.9"
Total weight	16.1 lbs	15.9 lbs	15 lbs
Max. user weight	330 lbs	287 lbs	275 lbs
Max. weight* bag	11 lbs	11 lbs	11 lbs
Recommended user height	4'11" - 6'7"	4'5" - 5'7"	4'1" - 5'5"

Gemino is available\* a number of different models; Gemino 30 walker and

Gemino 30 M walker – with support for the forearms.

Gemino 30 S – model for children. Gemino 20 – a simpler version of the rollator

