

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ Prov: _____ PC: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for:

Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ Prov: _____ PC: _____
 Ship To Phone: () _____ Fax: () _____

JAY® GO Cushion

JAY GO Cushion

✓	Option #	Definition	Price
<input type="checkbox"/>	JGOXXYY	JAY GO Cushion 14" - 21" wide	\$298
<i>Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)</i>			
<input type="checkbox"/>	JGOXXYY	JAY GO Cushion 22" - 24" wide	\$335
<i>Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)</i>			

JAY GO Size Matrix

		WIDTH													
		12	13	14	15	16	17	18	19	20	21	22	23	24	
DEPTH	12	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	13	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	14	N/A	N/A						N/A	N/A	N/A	N/A	N/A	N/A	
	15	N/A	N/A								N/A	N/A	N/A	N/A	
	16	N/A	N/A									N/A	N/A	N/A	
	17	N/A	N/A										N/A	N/A	
	18	N/A	N/A											N/A	N/A
	19	N/A	N/A												
	20	N/A	N/A												
	21	N/A	N/A	N/A											
	22	N/A	N/A	N/A	N/A										
	23	N/A	N/A	N/A	N/A	N/A									
	24	N/A	N/A	N/A	N/A	N/A	N/A								

	=	Standard Offering
JYW	=	Offered Via Jay Your Way
N/A	=	Not available

JAY GO Measuring Guide

Measurements made from completed cushion with covers.

Width

X = Width

Note: Measurement made from outside edge to outside edge of cushion.

Ex JGO1618

X = 16"


Depth

Y = Depth

Note: Measurement made from front edge of rear cutout and front lower edge of cushion.

Ex JGO1618

Y = 18"



JAY GO Cushion

✓	Option #	Definition	Price
<input type="checkbox"/>	JGO1414	JAY GO Cushion 14" x 14"	\$298
<input type="checkbox"/>	JGO1415	JAY GO Cushion 14" x 15"	\$298
<input type="checkbox"/>	JGO1416	JAY GO Cushion 14" x 16"	\$298
<input type="checkbox"/>	JGO1417	JAY GO Cushion 14" x 17"	\$298
<input type="checkbox"/>	JGO1418	JAY GO Cushion 14" x 18"	\$298
<input type="checkbox"/>	JGO1419	JAY GO Cushion 14" x 19"	\$298
<input type="checkbox"/>	JGO1420	JAY GO Cushion 14" x 20"	\$298
<input type="checkbox"/>	JGO1514	JAY GO Cushion 15" x 14"	\$298
<input type="checkbox"/>	JGO1515	JAY GO Cushion 15" x 15"	\$298
<input type="checkbox"/>	JGO1516	JAY GO Cushion 15" x 16"	\$298
<input type="checkbox"/>	JGO1517	JAY GO Cushion 15" x 17"	\$298
<input type="checkbox"/>	JGO1518	JAY GO Cushion 15" x 18"	\$298
<input type="checkbox"/>	JGO1519	JAY GO Cushion 15" x 19"	\$298
<input type="checkbox"/>	JGO1520	JAY GO Cushion 15" x 20"	\$298
<input type="checkbox"/>	JGO1521	JAY GO Cushion 15" x 21"	\$298
<input type="checkbox"/>	JGO1614	JAY GO Cushion 16" x 14"	\$298
<input type="checkbox"/>	JGO1615	JAY GO Cushion 16" x 15"	\$298
<input type="checkbox"/>	JGO1616	JAY GO Cushion 16" x 16"	\$298
<input type="checkbox"/>	JGO1617	JAY GO Cushion 16" x 17"	\$298
<input type="checkbox"/>	JGO1618	JAY GO Cushion 16" x 18"	\$298
<input type="checkbox"/>	JGO1619	JAY GO Cushion 16" x 19"	\$298
<input type="checkbox"/>	JGO1620	JAY GO Cushion 16" x 20"	\$298
<input type="checkbox"/>	JGO1621	JAY GO Cushion 16" x 21"	\$298
<input type="checkbox"/>	JGO1622	JAY GO Cushion 16" x 22"	\$298
<input type="checkbox"/>	JGO1714	JAY GO Cushion 17" x 14"	\$298
<input type="checkbox"/>	JGO1715	JAY GO Cushion 17" x 15"	\$298
<input type="checkbox"/>	JGO1716	JAY GO Cushion 17" x 16"	\$298
<input type="checkbox"/>	JGO1717	JAY GO Cushion 17" x 17"	\$298
<input type="checkbox"/>	JGO1718	JAY GO Cushion 17" x 18"	\$298
<input type="checkbox"/>	JGO1719	JAY GO Cushion 17" x 19"	\$298
<input type="checkbox"/>	JGO1720	JAY GO Cushion 17" x 20"	\$298
<input type="checkbox"/>	JGO1721	JAY GO Cushion 17" x 21"	\$298
<input type="checkbox"/>	JGO1722	JAY GO Cushion 17" x 22"	\$298
<input type="checkbox"/>	JGO1723	JAY GO Cushion 17" x 23"	\$298
<input type="checkbox"/>	JGO1814	JAY GO Cushion 18" x 14"	\$298
<input type="checkbox"/>	JGO1815	JAY GO Cushion 18" x 15"	\$298
<input type="checkbox"/>	JGO1816	JAY GO Cushion 18" x 16"	\$298
<input type="checkbox"/>	JGO1817	JAY GO Cushion 18" x 17"	\$298
<input type="checkbox"/>	JGO1818	JAY GO Cushion 18" x 18"	\$298
<input type="checkbox"/>	JGO1819	JAY GO Cushion 18" x 19"	\$298
<input type="checkbox"/>	JGO1820	JAY GO Cushion 18" x 20"	\$298
<input type="checkbox"/>	JGO1821	JAY GO Cushion 18" x 21"	\$298
<input type="checkbox"/>	JGO1822	JAY GO Cushion 18" x 22"	\$298
<input type="checkbox"/>	JGO1823	JAY GO Cushion 18" x 23"	\$298
<input type="checkbox"/>	JGO1824	JAY GO Cushion 18" x 24"	\$298
<input type="checkbox"/>	JGO1915	JAY GO Cushion 19" x 15"	\$298
<input type="checkbox"/>	JGO1916	JAY GO Cushion 19" x 16"	\$298
<input type="checkbox"/>	JGO1917	JAY GO Cushion 19" x 17"	\$298

✓	Option #	Definition	Price
<input type="checkbox"/>	JGO1918	JAY GO Cushion 19" x 18"	\$298
<input type="checkbox"/>	JGO1919	JAY GO Cushion 19" x 19"	\$298
<input type="checkbox"/>	JGO1920	JAY GO Cushion 19" x 20"	\$298
<input type="checkbox"/>	JGO1921	JAY GO Cushion 19" x 21"	\$298
<input type="checkbox"/>	JGO1922	JAY GO Cushion 19" x 22"	\$298
<input type="checkbox"/>	JGO1923	JAY GO Cushion 19" x 23"	\$298
<input type="checkbox"/>	JGO1924	JAY GO Cushion 19" x 24"	\$298
<input type="checkbox"/>	JGO2015	JAY GO Cushion 20" x 15"	\$298
<input type="checkbox"/>	JGO2016	JAY GO Cushion 20" x 16"	\$298
<input type="checkbox"/>	JGO2017	JAY GO Cushion 20" x 17"	\$298
<input type="checkbox"/>	JGO2018	JAY GO Cushion 20" x 18"	\$298
<input type="checkbox"/>	JGO2019	JAY GO Cushion 20" x 19"	\$298
<input type="checkbox"/>	JGO2020	JAY GO Cushion 20" x 20"	\$298
<input type="checkbox"/>	JGO2021	JAY GO Cushion 20" x 21"	\$298
<input type="checkbox"/>	JGO2022	JAY GO Cushion 20" x 22"	\$298
<input type="checkbox"/>	JGO2023	JAY GO Cushion 20" x 23"	\$298
<input type="checkbox"/>	JGO2024	JAY GO Cushion 20" x 24"	\$298
<input type="checkbox"/>	JGO2117	JAY GO Cushion 21" x 17"	\$298
<input type="checkbox"/>	JGO2118	JAY GO Cushion 21" x 18"	\$298
<input type="checkbox"/>	JGO2119	JAY GO Cushion 21" x 19"	\$298
<input type="checkbox"/>	JGO2120	JAY GO Cushion 21" x 20"	\$298
<input type="checkbox"/>	JGO2121	JAY GO Cushion 21" x 21"	\$298
<input type="checkbox"/>	JGO2122	JAY GO Cushion 21" x 22"	\$298
<input type="checkbox"/>	JGO2123	JAY GO Cushion 21" x 23"	\$298
<input type="checkbox"/>	JGO2124	JAY GO Cushion 21" x 24"	\$298

22" & Greater in Width

<input type="checkbox"/>	JGO2217	JAY GO Cushion 22" x 17"	\$335
<input type="checkbox"/>	JGO2218	JAY GO Cushion 22" x 18"	\$335
<input type="checkbox"/>	JGO2219	JAY GO Cushion 22" x 19"	\$335
<input type="checkbox"/>	JGO2220	JAY GO Cushion 22" x 20"	\$335
<input type="checkbox"/>	JGO2221	JAY GO Cushion 22" x 21"	\$335
<input type="checkbox"/>	JGO2222	JAY GO Cushion 22" x 22"	\$335
<input type="checkbox"/>	JGO2223	JAY GO Cushion 22" x 23"	\$335
<input type="checkbox"/>	JGO2224	JAY GO Cushion 22" x 24"	\$335
<input type="checkbox"/>	JGO2319	JAY GO Cushion 23" x 19"	\$335
<input type="checkbox"/>	JGO2320	JAY GO Cushion 23" x 20"	\$335
<input type="checkbox"/>	JGO2321	JAY GO Cushion 23" x 21"	\$335
<input type="checkbox"/>	JGO2322	JAY GO Cushion 23" x 22"	\$335
<input type="checkbox"/>	JGO2323	JAY GO Cushion 23" x 23"	\$335
<input type="checkbox"/>	JGO2324	JAY GO Cushion 23" x 24"	\$335
<input type="checkbox"/>	JGO2419	JAY GO Cushion 24" x 19"	\$335
<input type="checkbox"/>	JGO2420	JAY GO Cushion 24" x 20"	\$335
<input type="checkbox"/>	JGO2421	JAY GO Cushion 24" x 21"	\$335
<input type="checkbox"/>	JGO2422	JAY GO Cushion 24" x 22"	\$335
<input type="checkbox"/>	JGO2423	JAY GO Cushion 24" x 23"	\$335
<input type="checkbox"/>	JGO2424	JAY GO Cushion 24" x 24"	\$335


JAY Your Way - Modifications

*Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.
Please refer to the "JAY YOUR WAY" product catalog for complete listing of modifications*

Notch Adjustments

<input type="checkbox"/>	NR	Notch for Seat Rails	NC
--------------------------	----	----------------------	----


1.25" from bottom and 1.25" from side of cushion removed.



<input type="checkbox"/>	NP	Notch for Posts (Back Canes)	NC
--------------------------	----	------------------------------	----

From Back of Cushion: 1" 2" 3"

1" from side of cushion removed. No Modification to inner or outer cover.



Anterior or Posterior Slope

<input type="checkbox"/>	WGP	Anterior Slope	NC
<input type="checkbox"/>	WGA	Posterior Slope	NC

Height Increase:

Foam wedge glued underneath cushion, cover sewn to match.

Adjustments

<input type="checkbox"/>	LLCR	Leg Length Cut - Right	NC
<input type="checkbox"/>	LLCL	Leg Length Cut - Left	NC
<input type="checkbox"/>	CCT	Contracture Cut (1" x 45°)	NC
<input type="checkbox"/>	LTC	Leg Trough Cut (1" Standard)	NC
<input type="checkbox"/>	MSC	Medial Support Cut	NC

Cover Options

<input type="checkbox"/>	CC	Cover Color	\$42
--------------------------	----	-------------	------

Color:

Embroidery

<input type="checkbox"/>	EMBROIDERY	Embroidery (32 Character Max)	NC
--------------------------	------------	-------------------------------	----

Font: Block Script

Location:

Embroidery Text:

Embroidery Color:

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:

Please visit www.sunparts.us for spare parts information.