



Serial No:



ORDER FORM

Customer Dealer Information	Shipping Address	Identical To Billing Address
Dealer	Client	
Client	Address	
Address		
E-Mail:	E-Mail:	
Phone:	Phone:	
Order Information	User Information	
Purchase Order No.:	Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):	User Disability:	
	User Weight:	kg User Height
Build4me		

Serial No:

Product Style

TMLXZ010 MAXLITE

\$7,662.00

XA080 Closed design

Closed Design

Max User Weight 100kg (220lbs)

Please Note: All RGK chairs are hand built and have a manufacturing measurement tolerance of +/- 5mm

CAD Approval

XW009 Quote only (No CAD For Approval)

XW010 Order CAD Drawing For Approval *(Includes 2 CAD Amendments)*

\$300.00

XW011 Order (No CAD Required)

XW012 Quote & CAD For Approval *(Includes 2 CAD Amendments)*

\$550.00

XY010 Seat Width (Inside of sideguards)

mm (Min - 250mm / Max - 440mm)

XY020 Seat Length

mm (Min - 250mm / Max - 440mm)

XY030 Extra Seat Length (Provides additional leg support and hand grip)

mm

XY050 Front Frame Height

mm (Min - 380mm / Max - 550mm)

XY060 Rear Frame Height

mm (Min - 360mm / Max - 550mm)

XY070 Frame Length

mm (Min = XY020 + XY030 +110mm) (Max - 650mm)



XY072 Castor Position

mm Optimise Measurements

XA020 Ergonomic Seat

\$125.00

mm Specify Length - (Standard = 120mm)



Front Frame

XY080 Seat Taper Design

Specify width

mm

(When combined with XY091, dimension is inbetween the centre of the main bends)

XY091 Parallel Taper Design

Specify width

mm

(When combined with XY080, dimensions must be smaller by 20mm)

Footrest Design

XY100 Round Design

No Cost

XY105 Square Design

No Cost

XB200 High Level Footrest

(Only with XB010, XB080 & XB150)

\$207.00

Footrest Individual Measurements

XY115 Front Footrest Height

mm

XY120 Rear Footrest Height

mm STD = 20mm lower than front height

XY125 Footrest Length

(Std = 100mm)

mm (160mm and over requires heavy duty upgrade XB150)

Footrest Style and Cover

XB020 Tubular Titanium footrest (Height adjustable) No Cost

XB030 Flip up plastic footrest (Angle/height adjustable, with closed frame only) \$363.00

XB070 Angle adjustable aluminium footrest (Angle/height adjustable) \$213.00

XB150 Heavy Duty Tubular Footrest (Height adjustable, for Freewheel use) \$185.00

TMLXB080 Tubular footrest welded into frame (No height adjustment) \$173.00

Footrest Cover

XB050 Carbon fibre cover (Not with XB070 & XB030) \$225.00

XB170 Aluminium cover (Not with XB070 & XB030) \$101.00

XB040 Plastic cover (Not with XB070 & XB030) \$101.00

XB165 Rubber grip tape (3M safety walk) - (not with XB030) \$89.00

Castor Wheels

XC060	Soft roll 4" x 1.3" (Black only)	No Cost
XC110	Aluminium centre 4" x 1" (Silver hub and black urethane tyre only)	\$239.00
XC130	Aluminium centre 4" x 1.4" (Indicate Colour) Silver hub & black tyre Black hub & black tyre Red hub & black tyre Gold Hub & black tyre Green hub & black tyre Blue hub & black tyre Purple hub & black tyre	\$239.00
XC070	Soft roll 5" x 1.3" (Black only)	No Cost
XC100	Pneumatic Aluminium Centre 5" x 1" (Grey only)	\$255.00
XC120	Aluminium Centre 5" x 1" (Silver hub and black urethane tyre only)	\$239.00
XC140	Aluminium Centre 5" x 1.4" (Silver hub and black urethane tyre only)	\$239.00
XC080	Soft roll 6" x 1.3" (Black hub / grey tyre only)	No Cost
XC145	Aluminium centre 6" x 1.4" (Silver hub and black urethane tyre only)	\$319.00

Castor Forks

XC190	Performance Castor Fork (Std) Note: Colour follows accent	Std
XC192	High Trail Castor Fork Note: Greater distance between contact patch and steering axis reduces flutter at higher speeds. Colour follows accent	No Cost
XC210	Froglegs Suspension Castor Fork - Black - Silver	\$652.00

Castor Pin Width

XY330	Castor Pin Width (XY010 plus 70mm)
-------	---------------------------------------

Wheel Lock Options

XD020	Underseat Scissor	No Cost
XD050	Push To Lock	No Cost
XD030	Compact Scissor	\$274.00
XD046	Out front push to lock	\$421.00
XD035	Compact scissor with lightweight fixing (Only with fixed centre of gravity)	\$363.00
XD085	Eagle wheel locks (Indicate colour)	\$482.00
XD086	Eagle WL with L/W fixing (Indicate colour) (Only with fixed centre of gravity)	\$523.00
	Black Purple Red Blue	
XD010	Omit wheel locks	

Sideguards

Aluminium options

XE030	Screw on aluminium	\$337.00
XE080	Lift out aluminium	\$548.00
XE180	Flip back aluminium	\$405.00
XE040	Screw on aluminium with fender	\$287.00
XE090	Lift out aluminium with fender	\$596.00
XE190	Flip back aluminium with fender	\$344.00

Carbon Fibre options

XE050	Screw on carbon fibre	\$748.00
XE100	Lift out carbon fibre	\$789.00
XE200	Flip back carbon fibre	\$618.00
XE060	Screw on carbon fibre with fender	\$981.00
XE110	Lift out carbon fibre with fender	\$1,166.00
XE210	Flip back carbon fibre with fender	\$768.00

ABS Plastic options

XE070	Screw on ABS plastic	\$172.00
XE120	Lift out ABS plastic	\$502.00

Welded options

TMLXE220	Welded with 6mm bead	\$514.00
----------	----------------------	----------

Sideguard measurements

Optimised Measurements

25mm above tyre
(Follow tyre)

Individual above tyre
(Follow tyre) mm

Level with tyre
(Follow tyre)

Sideguard height

mm

Sideguard length
mm

Sideguard measurements

Fender Measurements

10mm above tyre
(Standard)

Individual above tyre
(Min = 10mm) mm

Centre of tyre
(Standard if no option selected)

Outside of tyre

Armrest options

XF010	Lift out armrest with plastic sideguard (Indicate fixed height)	\$610.00
XF020	Lift out armrest with aluminium sideguard (Indicate fixed height)	\$481.00

Height (Top of seat rail to top of armrest)

mm

Backrest Height

XY230 Fixed height backrest
(Maximum 400mm)

Specify height

mm

XY245 Adjustable height backrest **\$255.00**
(Maximum 400mm)
(+/- 25MM)

Specify height

mm

Backrest Design

XG010 Fixed backrest
(Non-folding)

TMLXG020 Folding backrest with 2 position snap lock and angle adjust (+/- 5 Deg) **\$506.00**

Ergonomic Options

XY360 Backrest angle

-5 -3 0 +3 +5

Individual Deg
(Max -10/+20 Deg)

XG100 Ergonomic backrest **\$150.00**

Vertical height

mm

Rearward position

mm

XG110 Inward taper **\$150.00**
(Inside of backrest posts)

Specify width or per side

mm

XG120 Outward flare **\$150.00**
(Inside of backrest posts)

Specify width or per side

mm

Push Handle Options

XG070 Screw in push handles **\$232.00**

XG080 Flip down push handles (N/A with XY245) **\$454.00**

XG085 Flip down push handles height adjustable (N/A with XY245) **\$742.00**

XG097 Integrated push handle (Quickie design) (N/A with XY245) **\$696.00**

XG090 Stroller push handles height adjustable **\$454.00**

Height from floor to top of handles

mm

Optimise (1000mm from floor)

Backrest Options

XG040 Heavy duty backrest upgrade **\$213.00**
(Recommended for hand cycling) (Not with XY245)

XG130 Extra deep square brace bar **\$118.00**

Depth

mm

(Recommended for Jay backs) (STD-90mm)

Camber Bar

Adjustable Centre Of Gravity

XH110	Square back adjustable centre of gravity (25-150mm)	No Cost
-------	--	---------

Specify measurement
mm

XH120	Round back adjustable centre of gravity (60-150mm)	No Cost
-------	---	---------

Specify measurement
mm

Fixed Centre Of Gravity

TMLXH090	V Style Fixed centre of gravity (25-125mm) Specify measurement mm	\$454.00
----------	--	----------

Camber Angle

XY390	Camber angle 0 1 2 3 5
-------	---------------------------

Camber Bar Options

XH040	Handbike camber bar fixed Camber angle (Handbike camber bar) 0 5 Castor height from floor for handbike fitting (Standard 30mm if not included) 30 50	\$412.00
-------	---	----------

Anti-Tip Options

XJ010	Swingaway anti tip LEFT side (Only with square back frame. Not with XH040)	\$172.00
XJ020	Swingaway anti tip RIGHT side (Only with square back frame. Not with XH040)	\$172.00
XJ100	Single sport anti tip (Centre mounted as standard)	\$519.00

Rear Wheel Size										
XY420		24" Wheels		XY430		25" Wheels				
Rear Wheel Gap										
XY460		Rear Wheel Gap			Individual		mm			
		15	20	25						
Rear Wheel Options										
XK010		Fusion 16 Spoke (Silver hub, black spokes)						Std		
XK320		Spinergy XLX 12 Spoke, cross/radial spoked						\$1,099.00		
XK050		Spinergy LX 12 Spoke						\$1,030.00		
XK060		Spinergy Spox 18 Spoke						\$685.00		
XK150		Spinergy CLX 18 Carbon blade spoke (N/A in silver hub)						\$3,233.00		
XK070		Spinergy BLXL 18 blade spoke (N/A in silver hub)						\$1,099.00		
XK040		Supply without Wheels		Wheel Type		Tyre Type				
Spinergy Spoke Colour options										
Black spokes		White spokes		Yellow spokes						
Red spokes		Blue spokes		Pink spokes						
Purple spokes		Orange spokes		Green spokes						
XK180		Individual spoke design (Specify colour)				Inside		Outside		\$219.00
Spinergy Hub Colour options										
XK190		Black Hub (Std option)						Std		
XK200		Silver Hub (n/a with BLXL/CLX)						\$131.00		
XK201		Red Hub/Rims (LX, XLX & BLXL only)						\$421.00		
XK202		Purple Hub/Rims (LX, XLX & BLXL only)						\$421.00		
XK203		Blue Hub/Rims (LX, XLX & BLXL only)						\$421.00		
XK204		Gold Hub/Rims (LX, XLX & BLXL only)						\$555.00		
Offroad Wheels										
XK300		Offroad wheel set 22" with handrim, MTB tyre & axle pins (Replaces 24" & 25" wheel with XM010)						\$961.00		
XK310		Offroad wheel set 24" with handrim, MTB tyre & axle pins (Replaces 26" wheel with XM010)						\$961.00		

Handrim Position / Options

XY480	Narrow set handrim	No Cost
XY485	Wide set handrim	No Cost
XL090	Remove excess mounting tabs (Recommended for small wheel gaps)	\$150.00

Handrim Styles

XL010	Aluminium hard anodised 6 tab (Standard if no option selected)	Std
XL040	Titanium 6 tab	\$459.00
XL160	Elipse 6 tab (replacement for surge)	\$650.00
XL035	Natural Fit LT	\$765.00
XL020	Tetra grip 6 tab	\$412.00
XL150	Carbolife Gekko BLACK SILVER	\$1,448.00
XL100	Para grip 6 tab	\$412.00
XL170	Carbolife Curve L (Black)	\$1,435.00

Tyre Options

XM010	Kenda Grey Blue Red (Grey standard if no option selected)	Std
XM060	Schwalbe one (24", 25" only)	\$137.00
XM130	Kenda black	\$29.00
XM050	Kenda K885 mountain bike tyre (24", 25" only)	\$41.00
XM020	Schwalbe marathon plus	\$101.00
XM040	Shox solid airless (Black) (26" N/A with Fusion Wheels)	\$165.00

Inner Tube Valve Options

XY500	Schraeder inner tube valve (Standard if no option selected)	No Cost
XY490	Presta inner tube valve	No Cost

Wheel Accessories

XK220	Titanium 130mm quick release axle pins (Pair) (N/A with XK010 & XK070)	\$165.00
XK240	Quad release axle nuts (Pair)	\$248.00
XK145	Black gloss spokeguards	Fitted Supply non - fitted \$655.00

Finish Options

XN010 Powder coated frame

No Cost

Frame

Gloss

Black gloss 103
White gloss 100
Grey gloss 132
Cream gloss 135
Red gloss 108
Yellow gloss 120
Forest green gloss 143
Green gloss 107
Dark blue gloss 151
Royal blue gloss 154
Blue gloss 174
Duckegg Blue gloss 141
Pink gloss 188
Orange gloss 180
Purple gloss 157

Frame

Matt

Black matt 105
White matt 101
Grey matt 131
Light grey matt 133
Navy blue matt 152

Frame

Special

Chrome mirror 300
Sparkle silver 170
Metallic black matt 104
Metallic grey matt 130
Black sandpaper 201
Grey sandpaper 202
White sandpaper 200

XN150

Powder coated components

No Cost

(Frame colour if not specified)

Specify component colour code from above list

TMLXN030

Bead Blast

(All parts black 105 unless specified)

\$582.00

XN165

Bead Blast components

\$258.00

TMLXN050

Polished Frame

(All parts black 103 unless specified)

\$1,198.00

XN170

Polished Components

\$219.00

Accent Colour

Note: Accent colour includes forks, clamps & castor caps.

Black accent

Silver accent

Red accent

Blue accent

Purple accent

Upholstery

Seat upholstery

XP040 Seat sling strap adjustable (Std)

XR050

Carbon fibre seat plate

\$494.00

mm (Std XY020 - 20mm)

Backrest upholstery

XP050 Airtech backrest strap adjustable

XP070

Backrest upgrade to leather corners

\$125.00

(Black)

Backrest stitching colour

Silver

Black

Blue

Individual:

Red

Orange

Purple

RGK Logo colour

White (std)

Green

Blue

Individual:

Red

Orange

Purple

Model Logo colour

White (std)

Green

Blue

Individual:

Red

Orange

Purple

Backrest Accessories

XP110

Backrest with zip pocket (N/A with sub4 backrest)

\$131.00

XP170

Individual backrest embroidery (Please provide JPEG)

\$247.00

Cushion

XP080	Cushion foam (Select density and size below)				\$125.00
	Soft (CMS040)	Medium (CMS050)			
	50mm	75mm	100mm		

Cushion cover upgrades

XP085	Nylon cushion cover (Black)	Std
XP090	Air lite cushion cover (Black)	\$101.00
XP100	Leather cushion cover (Black)	\$262.00

Frame Options

XR060	Welded calf strap mounting loops (To locate calf strap on front frame) Specify height from floor	\$159.00
	Optimise measurements	

Accessories

XR080	Frontwheel assembly - Specify colour: (Only with tubular and angle adjustable style footrest)	Black	Blue		\$1,085.00
		Silver	Purple	Red	
XP120	Frame protectors neoprene zip fasten (Black)				\$87.00
XP420	Frame protectors nylon velcro fasten (Black)				\$83.00
XR010	Seat pouch velcro attach (Black)				\$61.00
XR200	Rucksack with RGK logo (Black)				\$156.00
XR090	Tool kit				Std
XS060	Calf strap				Std

Spec Notes



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8

Customer Service: (800) 263-3390 E-mail: cscanada@sunmed.com www.sunrisemedical.ca
Specifications subject to change without notice. © 2026 Sunrise Medical Canada Inc. CN-100069 REV (3)