




Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		Province	
Buyer				PC	Phone		
Mark For				Email			

Please submit orders to: cscanada@sunmed.com

Seating Specifications			
Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options		
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	19", 20", 21"
Drive Base Specifications			
Drive Wheel Position	4 Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	400 lbs. 350 lbs. with Elevate or Tilt	Drive Motors	4 x 700W 4-pole motors
XT4 Base Size	Width 27.5" Length 39" Turning Radius 55.5"	Speed	6.2mph 3.1mph at full seat elevation
Min chair weight (No power options)	408lbs.	Tires	All Terrain Tires
Std chair weight (Lift, Tilt, Recline)	485lbs.	Charger	Dual Mode
Battery Size	Group 31	Ground Clearance	5"

Client Measurements

a)	Top of head to seat	_____	in
b)	Top of shoulder to seat	_____	in
c)	Underarm to seat	_____	in
d)	Elbow to seat	_____	in
e)	Top of leg to seat	_____	in
f)	Shoulder width	_____	in
g)	Chest width	_____	in
i)	Hip width	_____	in
j)	Actual sitting depth	_____	in
m)	Foot length	_____	in
n)	Knee to foot	_____	in
o)	Seat width	_____	in



Notes

✓ OPTION				DESCRIPTION	RETAIL	
Seating Options						\$
☐ 600SF4				REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC	
☐ 600SF5				MPS Seat and Back - Captain style seating	NC	
☐ 600SF6				MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC	
Base						
XT4 Base		End Item Number: EIPW60			\$	
All XT4 bases include four 700 watt 4 pole motors, LED lights and indicators, and ISO-7176 approved transit brackets as standard.						
☐ 600MW96				Magic Mobility - XT4 No Power Seating Options	\$ 28,835.00	
☐ 600MW97				Magic Mobility - XT4 Single Power Seating Options	\$ 28,835.00	
☐ 600MW98				Magic Mobility - XT4 Multi Power Seating Options	\$ 28,835.00	
☐ 600MW99				Magic Mobility - XT4 No Power Seating Options - HD	\$ 28,835.00	
☐ 600MW100				Magic Mobility - XT4 Single Power Seating Options - HD	\$ 28,835.00	
☐ 600MW101				Magic Mobility - XT4 Multi Power Seating Options - HD	\$ 28,835.00	
☐ 600MW102				Magic Mobility - XT4 No Power Seating Options - MPS (Captain)	\$ 28,835.00	
☐ 600MW103				Magic Mobility - XT4 Single Power Seating Options - MPS (Captain)	\$ 28,835.00	
☐ 600MW104				Magic Mobility - XT4 Multi Power Seating Options - MPS (Captain)	\$ 28,835.00	
☐ 600MW105				Magic Mobility - XT4 No Power Seating Options - HD MPS (Captain)	\$ 28,835.00	
☐ 600MW106				Magic Mobility - XT4 Single Power Seating Options - HD MPS (Captain)	\$ 28,835.00	
☐ 600MW107				Magic Mobility - XT4 Multi Power Seating Options - HD MPS (Captain)	\$ 28,835.00	
Base Accessories					\$	
☐ 6008A91				Jack for Drive Tire	\$ 280.00	
Optional Drive Tire Options					\$	
☐ 6008A101				Tire Sealant for 4 Pneumatic Drive Tires	\$ 106.00	
Shroud Color					\$	
☐ 600C185	Midnight Black	NC	☐ 600C207	Camouflage - Digital Desert	\$ 670.00	
☐ 600C188	Horizon Blue	NC	☐ 600C199	Camouflage - Woods	\$ 670.00	
☐ 600C187	Sunset Red	NC	☐ 600C209	Carbon Fiber	\$ 670.00	
			☐ 600C208	American Flag	\$ 670.00	
Seat To Floor Height					\$	
XT4/XT2 Seat to Floor Heights start at 18.75" +/- 0.25"						
☐ 600HT44				Lowest Seat To Floor Height	NC	
☐ 600HT45				Seat To Floor Height +1"	NC	
☐ 600HT46				Seat To Floor Height +2"	NC	
Batteries					\$	
☐ 600AG46				Group 31 Gel Battery (Pair) - Installed (Available with XT4)	\$ 1,494.00	
Battery Charger					\$	
☐ 600BJ28				10A Battery Charger	NC	
☐ 600BJ38				12A Battery Charger	\$ 716.00	

Seating: Seat Functions

✓ OPTION	DESCRIPTION	RETAIL
No Power Seat Options		\$
☐ 600ST87	No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, or future power seat functions	\$ 9,608.00
Power Tilt		\$
Utilized without Power Recline.		
☐ 600PA64	50° Power Tilt in Space - 350 lb. weight capacity	\$ 10,543.00
Power Seat Elevate		\$
Only available when Power Tilt or Power Tilt and Recline is selected.		
☐ 600PI16	12" Adjustable Seat Elevator - 350 lb. weight capacity	\$ 9,185.00
Power Recline		\$
Utilized without Power Tilt. Requires No Power Seating with Power Functionality.		
☐ 600PA65	150° Power Recline Only - 350 lb. weight capacity	\$ 7,505.00
Power Tilt and Recline		\$
☐ 600PA66	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 16,650.00

Seating: Seat Sizes

MPS (Captain) Seat Cushion (width x depth)		\$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.		
☐ 600SI24	16" x 18" Vinyl Seat	NC
☐ 600SI25	18" x 18" Vinyl Seat	NC
☐ 600SI26	18" x 20" Vinyl Seat	NC
☐ 600SI27	20" x 18" Vinyl Seat	NC
☐ 600SI28	20" x 20" Vinyl Seat	NC
☐ 600SI29	22" x 20" Vinyl Seat	\$ 636.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes		\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

☐ 600PN13	Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
☐ 600PN14	Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 640.00
☐ 600PN15	Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 640.00

REHAB Seat Width

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options		\$
-------------------------------	--	----

Cushion size will be determined by Rehab Seat size selected above

☐ 600CU59	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 677.00
☐ 600CU60	JAY® Ion® (Size 22" Wide)	\$ 719.00
☐ 600CU61	JAY® Union® (Sizes 16" - 21" Wide)	\$ 787.00
☐ 600CU62	JAY® Union® (Size 22" Wide)	\$ 830.00
☐ 600CU64	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 1,076.00
☐ 600CU65	JAY® Fusion® (Size 22" Wide)	\$ 1,224.00

Seating: Seat Sizes			
✓ OPTION		DESCRIPTION	RETAIL
Support Pads and Hardware (must pick pad and mount)			\$
Lateral Thigh Support Pads			
☐ 600KF22		Lateral Thigh Support Pad Short - Right (5x5)	\$ 77.00
☐ 600KF23		Lateral Thigh Support Pad Short - Left (5x5)	\$ 77.00
Lateral Thigh Support Mounts			
☐ 600KF26		Lateral Thigh Removable Mount - Right	\$ 335.00
☐ 600KF27		Lateral Thigh Removable Mount - Left	\$ 335.00
Seating: Seat Back Options			
MPS Backrest Options - Comes standard with Manual Recline			\$
Backrest Width			
☐ 600BN46		Standard Backrest Width - 16-20" Wide	NC
☐ 600BN47		Wide Backrest Width - 22-23" Wide	\$ 466.00
Backrest Height			
☐ 600B80		Standard Backrest Height - 22" Tall	NC
☐ 600B81		Tall Backrest Height - 26" Tall	\$ 766.00
Backrest Material			
☐ 600BU106		Cloth Backrest	NC
☐ 600BU107		Vinyl Backrest	\$ 422.00
MPS Back Accessories			\$
☐ 600PH33		Push Handle - Needed for O2 holder	\$ 217.00
☐ 6009A13		Accessory Mount - Needed for aftermarket headrest	\$ 164.00
REHAB Seating Backrest Options			\$
Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
JAY J3 Back Contour			
☐ 600JP01		J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 939.00
☐ 600JP02		J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 1,999.00
☐ 600JP11		J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 1,016.00
☐ 600JP13		J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 1,100.00
JAY® J3® Back Shell Options and Height Ranges			
☐ 600JR41		15" tall shell (Upper Thoracic Med) for 22" wide	NC
☐ 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
☐ 600JR15		18" tall shell (Upper Thoracic Med)	NC
☐ 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
☐ 600JR18		18" tall shell (Shoulder Height) for 22" wide	NC
☐ 600JR17		21" tall shell (Shoulder Height)	NC
☐ 600JR27		22.5" tall shell (Shoulder Height)	NC
☐ 600JR28		24" tall shell (Shoulder Height)	NC
REHAB Flip Forward Backrest Option			\$
☐ 600OD23		REHAB Flip Forward Backrest Option	\$ 220.00

Seating: Seat Back Options			
✓ OPTION		DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
JAY® Lateral Thoracic Support Pads			
☐ 600JW35		Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$ 77.00
☐ 600JW38		Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$ 77.00
☐ 600JW37		Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$ 77.00
☐ 600JW40		Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$ 77.00
JAY® Lateral Thoracic Swing-Away Mount			
☐ 600JW41		Lateral Thoracic Swing Away Bracket - Right	\$ 344.00
☐ 600JW42		Lateral Thoracic Swing Away Bracket - Left	\$ 344.00
JAY® Anterior Thoracic Support			\$
☐ 600ZR10		Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 234.00
☐ 600ZR11		Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 234.00
☐ 600K10		JAY® Anterior Thoracic Support - Medium ATS8302	\$ 332.00
☐ 600K11		JAY® Anterior Thoracic Support - Large ATS8303	\$ 332.00
Pelvic Positioning Options - Select One			\$
Pelvic Positioning			
☐ 600A848		Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
☐ 600A849		Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$ 106.00
☐ 600A850		Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 164.00
☐ 600A851		Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 280.00
Headrest Options - MPS Style			\$
☐ 600NU12		Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 569.00
☐ 600NU13		Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 569.00
Headrest Options - REHAB Style			\$
Headrest Pad			
☐ 600HE39		Whitmyer® Plush Headrest Pad, 10"	\$ 388.00
☐ 600HE46		Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 458.00
Headrest Mount			
☐ 600MH33		Whitmyer® Linx Headrest Mount with Quick Release	\$ 254.00
☐ 600MH34		Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 307.00
Armrest Style Options			\$
☐ 600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 1,842.00
☐ 600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 2,615.00
☐ 600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 3,262.00
Armrest Pad Options			\$
☐ 600AB81		10" Waterfall Armrest Pad - Right	NC
☐ 600AE81		10" Waterfall Armrest Pad - Left	NC
☐ 600AB82		14" Waterfall Armrest Pad - Right	NC
☐ 600AE82		14" Waterfall Armrest Pad - Left	NC

Legrest Options			
✓ OPTION	DESCRIPTION	RETAIL	
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
☐ 600H210	Center Mount Flip-Up Footplate	\$	628.00
Center Mount Fixed Legrest - Accessories			
☐ 600DI59	Fixed Center Mount Calf Pads	\$	718.00
☐ 600LS16	Foot Strap - 1.5" wide strap	\$	124.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest (Includes Calf Pads)			
☐ 600H215	Center Mount Power Elevating - 1 Piece footplate - Includes calf pads (11" x 11.5")	\$	5,750.00
☐ 600H216	Center Mount Power Elevating - 2 Piece footplate - Includes calf pads (5" x 10.5")	\$	5,580.00
Swing-Away Legrest			\$
Swing-Away Legrest			
☐ 600H211	Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$	903.00
Swing-Away Power Elevating Legrest			\$
Swing-Away Power Elevating Legrest			
☐ 600H214	Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate	\$	6,277.00
S/A Power Legrest Positioning Options			
☐ 600DI58	2" Metal Heel Guards for 2 Piece Footplates - May not fit with all seat sizes	\$	1,186.00
☐ 600HL15	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$	176.00
Accessories			\$
☐ 6008A93	Waterproof Wheelchair Cover	\$	229.00
☐ 6008A94	Waterproof Joystick Cover	\$	71.00
☐ 6008A95	Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$	53.00
☐ 600A221	Deluxe Backpack	\$	148.00
☐ 600A615	Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$	296.00
☐ 600A719	Articulating Oxygen Holder	\$	1,537.00
☐ 600A612	Cane Holder - Includes push handles	\$	986.00
☐ 6009A30	Luggage Rack for XT - Attaches to base behind the chair (10.5" x 13")	\$	514.00
☐ 6008A97	Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$	827.00
☐ 600A613	Cell Phone Holder - Adjustable to fit most cell phones	\$	458.00
☐ 600A614	Fishing Rod Holder	\$	458.00
☐ 6008A73	Armrest Bag	\$	113.00
☐ 6009A10	Utility Hooks	\$	88.00
☐ 6008A59	Transfer Bars	\$	541.00
Be Adaptive Products			
☐ 6008A98	Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$	2,309.00
☐ 6008A99	Be Adaptive Trigger Mechanism	\$	1,279.00
☐ 6008A100	Be Adaptive Limited Mobility Fishing Mount	\$	2,309.00
Electronics			
Joystick			\$
☐ 600J121	Standard LED Joystick (w/o Screen Display)	NC	
☐ 600J116	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC	
☐ 600J118	High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display	\$	780.00

Electronics									
✓ OPTION		DESCRIPTION						RETAIL	
Joystick Mount									\$
☐ 600JM34		Fixed on Armrest - Depth adjustable 8.25"						NC	
☐ 600JM45		Swing-Away Joystick Mount - Includes latch to lock in the forward position						\$ 874.00	
☐ 600JM46		Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)						\$ 2,112.00	
Joystick Mount Location									\$
☐ 600HM16		Right Hand Armrest						NC	
☐ 600HM17		Left Hand Armrest						NC	
Display Language Selection R-net									\$
☐ 600L1	English	NC	☐ 600L3	French	NC	☐ 600L7	Spanish	NC	
Joystick Handle Type									\$
Comes standard with PG Cone-Shaped Joystick Handle									
☐ 600JH13		Ball Handle						\$ 113.00	
☐ 600JH26		Bodypoint - Grooved Mushroom Joystick Handle						\$ 247.00	
☐ 600JH24		Bodypoint – Rubber Dome Joystick Handle						\$ 264.00	
☐ 600JH25		Bodypoint – U-Shaped Joystick Handle (4 inch)						\$ 264.00	
Expandable Power Module									\$
Must select both items									
☐ 600EX05		Expandable 120A Controller - R-net						\$ 1,132.00	
☐ 600EX06		Harness for 120A Expandable Controller						\$ 791.00	
Power Seat Function									\$
☐ 600EO37		Thru-Drive Control for 1 Actuator Must be selected if ordering 1 actuator function						\$ 1,486.00	
☐ 600EO38		Thru-Drive Control for 2+ Actuators Must be selected if ordering 2 or more actuator functions						\$ 2,990.00	
Power Seat Control Option									\$
☐ 600QG11		Ctrl+5 Switch Box with Free Assignable Functions						\$ 1,048.00	
☐ 600QG12		Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)						NC	
Ctrl+5 Switch Type									\$
☐ 600EY22		Push Button						NC	
☐ 600EY23		Toggle Switches						NC	
Power Seat Control Option - Ctrl+5 Mounting location									\$
☐ 600EY19	Right Armrest	NC	☐ 600EY20	Left Armrest					NC
Attendant Controls									\$
☐ 6000A05		Dual Control Module for Attendant - Includes accessory mount and push handle						\$ 2,534.00	
Attendant Control Mounting Location									\$
☐ 600DN5	Right Back Cane	NC	☐ 600DN6	Left Back Cane					NC
R-net Omni-2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth									\$
☐ 600ON14		R-net Omni-2 Only Input Control Module						\$ 3,719.00	
☐ 600ON15		R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices						NC	
R-net Omni-2 -Mounting Location									\$
☐ 600DS23	Right Hand Armrest	NC	☐ 600DS24	Left Hand Armrest					NC
Head Control Packages									\$
☐ 600EU25		Switch-It DualPro® Proportional Head Array - Facial Pads						\$ 7,269.00	
☐ 600EU26		Switch-It DualPro® Proportional Head Array - Spot Pads						\$ 7,269.00	