

LECKEY



Squiggles Seat

Giving little ones the support
they need to do big things

Supporting Children

Opening up a world of possibilities for little ones

It's easy to see why **Squiggles** is a firm favorite with little ones and their families when it comes to a starter seating system. **Squiggles** provides comfort and support in equal measure, meaning kids can get on and do the big things in life—eating dinner independently, watching a movie with siblings, using a tablet for homework, enjoying tabletop play, taking part in reading group at school.... the possibilities are endless with **Squiggles** supporting you!



PERFECTLY DESIGNED FOR THE UNIQUE NEEDS OF YOUNG CHILDREN

- **Squiggles** is not a scaled down version of another seat—importantly, kids are not over-supported in **Squiggles**, helping them to reach their full developmental potential. In **Squiggles**, children can achieve big things and still strive for more.
- The one-piece seat base and adjustable sacral support can be used to accommodate flexion or encourage extension, pertinent at this early stage.
- The optional contoured shoulder support is angle adjustable and helps bring hands to midline.
- **Squiggles** offers a fun, colourful, age-appropriate design, popular with both users and their families.

LOW TO THE GROUND POSITIONING POSSIBLE, IDEAL FOR CHILD-TO-CHILD INTERACTION

- **Squiggles** maximizes opportunities for inclusive play, learning and socializing at eye level, supporting kids in being kids.
- Perfect for use in nursery/school and home environments.



PROVIDES A CLOSE PROXIMAL FIT FOR COMFORT THAT LITTLE ONES LOVE

- A snug, comfortable fit is provided via the shape of the seat and back, the laterals with their fabric tip which wrap securely around the child and the way in which the 4-point pelvic harness integrates with the hip pads.
- The fit of **Squiggles** to the child gives a streamlined look, letting the child and their personality shine through.



DID YOU KNOW?

- Squiggles can be interfaced with a range of Mygo bases, offering more wheel type and chassis height adjustment options

Technical Information

Seat Dimensions

Age (approx.)	1-5
User Height	Min 750mm/29.5" Max 1100mm/43.3"
Max User Weight	22kg/48.5 lbs.
Seat Depth (backrest to back of knee)	Min 190mm/7.5" Max 290mm/11.4"
Seat Width	Min 160mm/6.3" Max 260mm/10.2"
Hip Width (distance between hip guides)	Min 160mm/6.3" Max 220mm/8.6"
Backrest Height (sitting shoulder height)	Min 320mm/12.6" Max 380mm/15"
Backrest Height (with shoulder support)	Min 420mm/16.5" Max 480mm/18.9"
Chest Width (distance between laterals)	Min 160mm/6.3" Max 220mm/8.6"
Knee Width (distance between pommel & side of seat)	90mm/3.5"
Seat to Footplate (sole of foot to back of knee)	Min 122mm/4.8" Max 270mm/10.6"
Top of Seat to Floor (Hi-Low Chassis)	Min 270mm/10.6" Max 650mm/25.6"
Backrest Angle	Prone 10° Recline 25°
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°
Tray Size	460mm x 560mm (18.1" x 22")