

elite

» Their loss is your game!

The logo for RGK, featuring the letters 'RGK' in a bold, white, sans-serif font with a stylized underline for the 'K'.

We are passionate about sporting life.

Just as you constantly push yourself to be on top of your game, we are driven in developing the highest quality sports chairs to support you and your goals.

Every athlete is unique and we understand the importance of having a chair made to fit you perfectly and to support your individual needs on the court.

With the confidence that your chair is individually created to maximise your potential, you just need to concentrate on the game, on your goals, on your victory.

Elite - maximum rigidity, absolute performance

Designed and created to fit you as a second skin..

The Elite has been designed based on in-depth talks with top athletes in the game.

We listened carefully to their requirements and focused on the most important aspects that a basketball chair needs. This has ensured that we have developed the lightest and most responsive basketball chair ever. It leads the way in performance and efficiency.

Integrated camber bar and footrest, offensive wing and anti tip make this professional, top of the range chair the best in the competition.

Key Features

- » Bespoke rigid frame Aluminum 7020
- » Modern look
- » Camber bar, offensive wing and anti tip integrated into the frame
- » Integrated footrest
- » Optional adjustable COG and footrests
- » Ergonomic seat and back
- » Camber options: 16°, 18°, 20°, 22°
- » Multiple offensive wing designs
- » Tension adjustable seat and back upholsteries
- » Aluminium and carbon rigid frame options
- » Max. user weight 125 kg



Integrated height adjustable anti tip (single or double)



Integrated camber bar



Integrated aluminium fork. Light and resistant



Integrated footrest with foot strap