

September 2026



Customer Contact Details

Customer Name: _____

Advisor: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address: ☐

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Guide:		☐ Base product		☒ Required option		☐ Optional Extra		Critical Information		Information Needed	
☒ SQCXZ010		CLUB SPORT		Max User Weight 125kg (275lbs)		☐ Yes, I am within the max user weight					
\$4,129.00											
Standard Specification Includes: - Basketball or Tennis (option 3 only) - Aluminium screw on sideguards - Adjustable position footrest - Adjustable centre of gravity (18 degree) - Removable single sport anti-tip - Adjustable height backrest 250-300mm - Cross laced wire spoke wheels - Aluminium hard anodized push rim/kenda tire 72mm caster wheels with extruded forks - Adjustable seat and backrest upholstery 5cm medium density foam cushion (*option 4) - Calf strap & foot strap											
Key Product Measurements											
Frame Option (Choose one)											
☐ XY520 Option 1 (Junior size)			☐ XY530 Option 2 (Low Point)			☐ XY540 Option 3 (Low/Mid Point)			☐ XY550 Option 4 (Max height)		
Sport			Sport			Sport			Sport		
☒ Basketball (with wing)			☒ Basketball (with wing)			☐ Basketball (with wing) ☐ Tennis (no wing) (5mm A/T)			☒ Basketball (with wing)		
Frame height			Frame height			Frame height			Frame height		
☒ Front frame height: 480mm			☒ Front frame height: 500mm			☒ Front frame height: 500mm			☒ Front frame height: 530mm		
☒ Rear frame height: 420mm			☒ Rear frame height: 400mm			☒ Rear frame height: 440mm			☒ Rear frame height: 530mm		
Wheel size			Wheel size			Wheel size			Wheel size		
☐ 22" ☐ 24" ☐ 25"			☐ 24" ☐ 25"			☐ 24" ☐ 25"			☐ 26" ☐ 700c		
XY415 XY420 XY430			XY420 XY430			XY420 XY430			XY440 XY450		
Seat width (inside of sides)			Seat width (inside of sides)			Seat width (inside of sides)			Seat width (inside of sides)		
☐ 340mm ☐ 370mm			☐ 370mm ☐ 400mm ☐ 440mm			☐ 370mm ☐ 400mm ☐ 440mm			☐ 370mm ☐ 400mm ☐ 440mm		
Seat length			Seat length			Seat length			Seat length		
☐ 300mm ☐ 330mm ☐ 360mm			☐ 360mm ☐ 390mm ☐ 430mm			☐ 360mm ☐ 390mm ☐ 430mm			☐ 360mm ☐ 390mm ☐ 430mm		
Frame length			Frame length			Frame length			Frame length		
☒ 520mm			☒ 580mm			☒ 580mm			☒ 580mm (SW = 370/400) ☒ 630mm SW = 440		
Options											
Finish Options						Accessories					
☒ XN010 Powder coated frame						☐ XS020 Sport ratchet strap (click strap) (not fitted) \$144.00					
Specify Code						☐ XS010 100mm velcro adjustable strap \$89.00					
☒ XN150 Powder coated components (Black)						☐ XK080 Spinergy SLX 24 spoke (Wheel upgrade) \$1036.00					
						☐ XP170 Individual backrest embroidery (Please provide .jpeg file and positioning) \$247.00					