



SUNRISE MEDICAL CANADA INC. SCOOTER MATRIX

Weight Capacity	300 lbs. (TA), 350 lbs. (DT)
Wheel Options	10" wheels all around
Total Width	24.5"
Total Length	47.5"
Turning Radius	42" (3-wheel), 57.5" (4-wheel)
Max Speed	8.8 km/h (3-wheel), 9.7 km/h (4-wheel)
Seat-to-Floor Height	17", 18", or 19" from floorpan (17" done by removal of seat slider)
Seating	Width: 18" or 19.75" (19.75" not available on 3wh models) Armrest width adjustable to 24" Depth: 17.5" or 19.75" (19.75" not available on 3wh models)
Batteries	U1, 40Amp, Grp 22



1700 SERIES

3-WHEEL DT/TA

4-WHEEL DT/TA

Weight Capacity	330 lbs. (S425), 350 lbs. (S700)
Wheel Options	12" pneumatic (S425), 13" pneumatic (S700)
Total Width	26"
Total Length	53"
Turning Radius	78" (S425), 98" (S700)
Max Speed	12 km/h (S425), 15 km/h (S700)
Seat-to-Floor Height	17.6" to 19.6" (S425), 17.7" to 19.7 (S700) - 3 seat positions only
Seating	Width: 19" Depth: 18" to 20"
Batteries	55Ah x 2 (S425), 75Ah x 2 (S700)



S-SERIES
PREMIUM SCOOTERS

S425

S700