



## ORDER FORM

February 2024

# Mygo Seat

Versatile therapy chair that grows with its user without sacrificing proximal pelvic positioning and support.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

## Mygo Seat

Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs

✓ OPTION	DESCRIPTION		RETAIL
<b>Seat Shell</b>			<b>\$</b>
Selection Required.			
<input type="checkbox"/> <b>117-610S</b>	<b>Size 1 Mygo Seat Shell</b>		<b>\$2,173</b>
117-610-1	Seat Base with Depth Adjustment		
117-610-2	Seat Back - Fixed		
117-610-3	Seat Back Reclining Mechanism		
117-610-4	Seat Base Hip Guides / Support		
117-610-5	Seat Base Medial Thigh Guides		
117-610-6	Seat Base Leg Adductors		
<input type="checkbox"/> <b>137-600S</b>	<b>Size 2 Mygo Seat Shell</b>		<b>\$2,408</b>
137-600-1	Seat Base with Depth Adjustment		
137-600-2	Seat Back - Fixed		
137-600-3	Seat Back Reclining Mechanism		
137-600-4	Seat Base Hip Guides / Support		
137-600-5	Seat Base Medial Thigh Guides		
137-600-6	Seat Base Leg Adductors		
<b>Seat Cover</b>			<b>\$</b>
Selection Required. Must match Seat Shell size.			
<input type="checkbox"/> 117-1851-08	Size 1 Mygo Seat Cover - Blue		\$597
<input type="checkbox"/> 117-1851-06	Size 1 Mygo Seat Cover - Gray		\$597
<input type="checkbox"/> 117-1851-07	Size 1 Mygo Seat Cover - Orange		\$597
<input type="checkbox"/> 117-1851-09	Size 1 Mygo Seat Cover - Pink		\$597
<input type="checkbox"/> 117-2851-08	Size 2 Mygo Seat Cover - Blue		\$597
<input type="checkbox"/> 117-2851-06	Size 2 Mygo Seat Cover - Gray		\$597
<input type="checkbox"/> 117-2851-07	Size 2 Mygo Seat Cover - Orange		\$597
<input type="checkbox"/> 117-2851-09	Size 2 Mygo Seat Cover - Pink		\$597
<b>Base</b>			<b>\$</b>
<input type="checkbox"/> 102-2753	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)		\$2,265
<input type="checkbox"/> 194-700	Mobility Base Interface Plate (Allows mounting Mygo Seat onto Sunrise Wheelchairs)		\$657
<b>Pelvic Fixation</b>			<b>\$</b>
Selection of either Pelvic Harness or Pelvic Cradle required.			
Hip Laterals require Pelvic Cradle			
<input type="checkbox"/> 117-801	Pelvic Harness & Brackets Small - Black		\$174
<input type="checkbox"/> 117-802	Pelvic Harness & Brackets Medium - Black		\$174
<input type="checkbox"/> 120-789	Pelvic Harness Spacer Pads - Black (Pair)		\$95
<input type="checkbox"/> 117-477	Pelvic Cradle Size 3 (old size 0)		\$425
<input type="checkbox"/> 117-478	Pelvic Cradle Size 4 (old size 1)		\$425
<input type="checkbox"/> 137-632	Hip Laterals (Pair)		\$357
<b>Dynamic Backrest Kit</b>			<b>\$</b>
<input type="checkbox"/> 137-661	Dynamic Backrest 100N Cylinder		\$492
<b>Armrest</b>			<b>\$</b>
<input type="checkbox"/> 117-1777	Size 1 Mygo Adjustable Armrests (Pair)		\$360
<input type="checkbox"/> 117-2777	Size 2 Mygo Adjustable Armrests (Pair)		\$360

## Mygo Seat

<input checked="" type="checkbox"/> OPTION	DESCRIPTION		RETAIL
--	-------------	--	--------

### Lower Leg Positioning \$

Footplate selection required with Sandals. Use wheelchair footplates if selection is not made here.

Sandal Selection Required with Sandal Raisers.

<input type="checkbox"/>	117-1776	Size 1 Mygo Angled Footplates (Pair)		\$360
<input type="checkbox"/>	117-2776	Size 2 Mygo Angled Footplates (Pair)		\$360
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)		\$269
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)		\$287
<input type="checkbox"/>	152-3600	Large Sandals w/ Straps (Pair)		\$310
<input type="checkbox"/>	117-787	Sandal Raiser 1" (Each)	QTY	\$105
<input type="checkbox"/>	117-765-05VIN			

### Harness \$

<input type="checkbox"/>	155-N001-L600	Mygo/SQSE Trunk Harness - Small		\$409
<input type="checkbox"/>	155-N001-L601	Mygo/SQSE Trunk Harness - Medium		\$409
<input type="checkbox"/>	155-N001-L602	Mygo/SQSE Trunk Harness - Large		\$409

### Lateral Support \$

Selection of Chest Support Requires Laterals. Not available with Complex Laterals.

<input type="checkbox"/>	117-837	Flip-Away Laterals with Cover (Pair)		\$549
<input type="checkbox"/>	137-634	Flip-Away Laterals Extended Height with Cover (Pair) <i>Required w/ Size 2 Seat</i>		\$674
<input type="checkbox"/>	117-803	Rigid Laterals with Cover (Pair)		\$324
<input type="checkbox"/>	137-633	Rigid Laterals Extended Height with Cover (Pair)		\$438
<input type="checkbox"/>	137-637	Flip-Away Complex Lateral Horizontal with Cover Left		\$364
<input type="checkbox"/>	137-638	Flip-Away Complex Lateral Horizontal with Cover Right		\$364
<input type="checkbox"/>	137-808-05	Chest Support - Black		\$201

### Headrest \$

Selection of Headrest Support Hardware (Contoured/Flat) Requires Cushion

Headrest Laterals only available with Flat Headrest

<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware		\$390
<input type="checkbox"/>	137-800-05	Contoured Headrest Cushion - Black		\$185
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware		\$243
<input type="checkbox"/>	117-850-05	Flat Headrest Cushion - Black		\$122
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with Covers - Black (Pair)		\$364

### Tray \$

Armrest selection required with Tray

<input type="checkbox"/>	117-646	Mygo Size 1 and 2 Black Tray		\$427
<input type="checkbox"/>	117-652-BLACK	Black Padded Tray (cover with tray)		\$622
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray		\$321

### Accessories \$

Shoulder Support is included on Mygo Size 2.

<input type="checkbox"/>	137-621	Shoulder Support Hardware for Mygo Size 1 Cushion selection required.		\$303
<input type="checkbox"/>	137-801-08	Shoulder Support Cushion - Blue		\$174
<input type="checkbox"/>	137-801-06	Shoulder Support Cushion - Gray		\$174
<input type="checkbox"/>	137-801-07	Shoulder Support Cushion - Orange		\$174
<input type="checkbox"/>	137-801-09	Shoulder Support Cushion - Pink		\$174

### Technical Specifications - Universal

Specification	Size 1	Size 2
Age (approximate)	3 - 10 years	8 - 14 years
Backrest Angle: Prone	10°	10°
Backrest Angle: Recline	25°	25°
Footplate: Abduction	8°	8°
Footplate: Adduction	12°	20°
Footplate: Plantar/Dorsiflexion	10° / 10°	10° / 10°

#### Size 2 Hi Lo Chassis (Excluding Push Bar)

Tilt in Space: Prone	10°	10°
Tilt in Space: Recline	25°	25°

### Technical Specifications - Imperial

Specification	Size 1	Size 2
User Height	41.3" - 59"	50" - 66"
User Weight	39.6 - 110 lbs	39.6 - 132 lbs
Seat Depth	10.6" - 16.5"	13.8" - 18.5"
Seat Width	8" - 12.8"	8.7" - 13.6"
Backrest Height	14.2" - 18.5"	18.1" - 22.4"
Seat to Sandal	8.5" - 13.8"	12.4" - 18.5"
Chest Width (Distance between laterals)	6.7" - 10.6"	6.7" - 10.6"
Armrest Height	6.3" - 8.3"	8.3" - 10.2"
Knee Width	3.5" - 4.3"	4.7" - 5.5"
Tray Size	21.6" x 18.9"	21.6" x 18.9"
Seat Unit Weight	22 lbs	31.9 lbs

#### Size 2 Hi Lo Chassis (Excluding Push Bar)

Top of Seat to Floor	13.4" - 25.8"	13.4" - 25.8"
Base Weight	26.5 lbs	26.5 lbs

### Technical Specifications - Metric

Specification	Size 1	Size 2
User Height	1050 mm - 1500 mm	1270 mm - 1680 mm
User Weight	18 - 50 kg	18 - 60 kg
Seat Depth	270 mm - 420 mm	350 mm - 470 mm
Seat Width	200 mm - 325 mm	220 mm - 345 mm
Backrest Height	360 mm - 470 mm	460 mm - 570 mm
Seat to Sandal	215 mm - 350 mm	315 mm - 470 mm
Chest Width (Distance between laterals)	170 mm - 270 mm	170 mm - 270 mm
Armrest Height	160 mm - 210 mm	210 mm - 260 mm
Knee Width	90 mm - 110 mm	120 mm - 140 mm
Tray Size	550 mm x 480 mm	550 mm x 480 mm
Seat Unit Weight	10 kgs	14.5 kgs

#### Size 2 Hi Lo Chassis (Excluding Push Bar)

Top of Seat to Floor	340 mm - 655 mm	340 mm - 655 mm
Base Weight	12 kgs	12 kgs


  
**SUNRISE MEDICAL.**
  
 Sunrise Medical Canada Inc.
   
 237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada

Customer Service: (800) 263-3390 Email: cscanada@sunmed.com www.sunrisemedical.ca

Specifications are subject to change without notice. © 2024 Sunrise Medical Canada Inc. MK-129849-CAN Rev (5)

