



ORDER FORM

Mygo Seat

EILK03

Versatile therapy chair that grows with its user without sacrificing proximal pelvic positioning and support.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Mygo Seat

Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs

✓ OPTION	DESCRIPTION		RETAIL
Seat Shell			\$
Selection Required.			
<input type="checkbox"/> 151BS24	Size 1 Mygo Seat Shell		
151BS26	Seat Base with Depth Adjustment		\$929
151ST88	Seat Back - Fixed		\$929
151ST90	Seat Back Reclining Mechanism		\$489
151GU3	Seat Base Hip Guides / Support		\$136
151KH19	Seat Base Medial Thigh Guides		\$169
151KB13	Seat Base Leg Adductors		\$136
<input type="checkbox"/> 151BS25	Size 2 Mygo Seat Shell		
151BS27	Seat Base with Depth Adjustment		\$1,031
151ST89	Seat Back - Fixed		\$1,031
151ST91	Seat Back Reclining Mechanism		\$489
151GU4	Seat Base Hip Guides / Support		\$136
151KH20	Seat Base Medial Thigh Guides		\$169
151KB14	Seat Base Leg Adductors		\$136
Seat Cover			\$
Selection Required. Must match Seat Shell size.			
<input type="checkbox"/> 151YJ06	Size 1 Mygo Seat Cover - Blue		\$606
<input type="checkbox"/> 151YJ07	Size 1 Mygo Seat Cover - Gray		\$606
<input type="checkbox"/> 151YJ08	Size 1 Mygo Seat Cover - Orange		\$606
<input type="checkbox"/> 151YJ09	Size 1 Mygo Seat Cover - Pink		\$606
<input type="checkbox"/> 151YJ10	Size 2 Mygo Seat Cover - Blue		\$606
<input type="checkbox"/> 151YJ11	Size 2 Mygo Seat Cover - Gray		\$606
<input type="checkbox"/> 151YJ12	Size 2 Mygo Seat Cover - Orange		\$606
<input type="checkbox"/> 151YJ13	Size 2 Mygo Seat Cover - Pink		\$606
Base			\$
<input type="checkbox"/> 151NB114	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)		\$2,296
<input type="checkbox"/> 151NB115	Mobility Base Interface Plate (Allows mounting Mygo Seat onto Sunrise Wheelchairs)		\$536
Pelvic Fixation			\$
Selection of either Pelvic Harness or Pelvic Cradle required.			
Hip Laterals require Pelvic Cradle			
<input type="checkbox"/> 1515W03	Pelvic Harness & Brackets Small - Black		\$167
<input type="checkbox"/> 1515W04	Pelvic Harness & Brackets Medium - Black		\$167
<input type="checkbox"/> 151NH108	Pelvic Harness Spacer Pads - Black (Pair)		\$94
<input type="checkbox"/> 1515W05	Pelvic Cradle Size 3 (old size 0)		\$463
<input type="checkbox"/> 1515W06	Pelvic Cradle Size 4 (old size 1)		\$514
<input type="checkbox"/> 151KC11	Hip Laterals (Pair)		\$345
Dynamic Backrest Kit			\$
<input type="checkbox"/> 151BT149	Dynamic Backrest 100N Cylinder		\$475
Armrest			\$
<input type="checkbox"/> 151AA170	Size 1 Mygo Adjustable Armrests (Pair)		\$339
<input type="checkbox"/> 151AA171	Size 2 Mygo Adjustable Armrests (Pair)		\$372

Mygo Seat				
✓ OPTION		DESCRIPTION		RETAIL
Lower Leg Positioning				\$
Footplate selection required with Sandals.				
Sandal Selection Required with Sandal Raisers.				
<input type="checkbox"/>	151F376	Size 1 Mygo Angled Footplates (Pair)		\$339
<input type="checkbox"/>	151F377	Size 2 Mygo Angled Footplates (Pair)		\$372
<input type="checkbox"/>	1518A107	Small Sandals w/ Straps (Pair)		\$278
<input type="checkbox"/>	1518A108	Medium Sandals w/ Straps (Pair)		\$304
<input type="checkbox"/>	1518A109	Large Sandals w/ Straps (Pair)		\$327
<input type="checkbox"/>	1518A110	QTY Sandal Raiser 1" (Each)		\$97
<input type="checkbox"/>	151XM16		Mygo Knee Pads (vinyl)	
Harness				\$
<input type="checkbox"/>	151NH105	Butterfly Harness (Small)		\$414
<input type="checkbox"/>	151NH106	Butterfly Harness (Medium)		\$414
<input type="checkbox"/>	151NH107	Butterfly Harness (Large)		\$414
Lateral Support				\$
Selection of Chest Support Requires Laterals. Not available with Complex Laterals.				
<input type="checkbox"/>	151KC12	Flip-Away Laterals with Cover (Pair)		\$591
<input type="checkbox"/>	151KC13	Flip-Away Laterals Extended Height with Cover (Pair) <i>Required w/ Size 2 Seat</i>		\$725
<input type="checkbox"/>	151KC14	Rigid Laterals with Cover (Pair)		\$338
<input type="checkbox"/>	151KC15	Rigid Laterals Extended Height with Cover (Pair)		\$474
<input type="checkbox"/>	151KC16	Flip-Away Complex Lateral Horizontal with Cover Left		\$378
<input type="checkbox"/>	151KC17	Flip-Away Complex Lateral Horizontal with Cover Right		\$378
<input type="checkbox"/>	151ZY25	Chest Support - Black		\$191
Headrest Support Hardware				\$
Selection of Headrest Support Hardware Requires Headrest Cushion				
<input type="checkbox"/>	151AW22	Contoured Headrest Support Hardware		\$393
<input type="checkbox"/>	151AW24	Flat Headrest Support Hardware		\$241
Headrest Cushion				\$
Headrest Laterals only available with Flat Headrest				
<input type="checkbox"/>	151Y027	Contoured Headrest Cushion - Black		\$180
<input type="checkbox"/>	151Y028	Flat Headrest Cushion - Black		\$123
<input type="checkbox"/>	151Y031	Flat Headrest Laterals with Covers - Black (Pair)		\$345
Tray				\$
Armrest selection required with Tray				
<input type="checkbox"/>	151A220	Mygo Size 1 and 2 Black Tray		\$406
<input type="checkbox"/>	151A221	Black Padded Tray (cover with tray)		\$633
<input type="checkbox"/>	1518A106	Grab Rail / Toy Bar for Tray		\$311
Accessories				\$
Shoulder Support is standard on Mygo Size 2.				
<input type="checkbox"/>	151AW25	Shoulder Support Hardware for Mygo Size 1 Cushion selection required.		\$287
<input type="checkbox"/>	151SS10	Shoulder Support Cushion - Blue		\$163
<input type="checkbox"/>	151SS11	Shoulder Support Cushion - Gray		\$163
<input type="checkbox"/>	151SS12	Shoulder Support Cushion - Orange		\$163
<input type="checkbox"/>	151SS9	Shoulder Support Cushion - Pink		\$163

Technical Specifications - Universal

Specification	Size 1	Size 2
Age (approximate)	3 - 10 years	8 - 14 years
Backrest Angle: Prone	10°	10°
Backrest Angle: Recline	25°	25°
Footplate: Abduction	8°	8°
Footplate: Adduction	12°	20°
Footplate: Plantar/Dorsiflexion	10° / 10°	10° / 10°

Size 2 Hi Lo Chassis (Excluding Push Bar)

Tilt in Space: Prone	10°	10°
Tilt in Space: Recline	25°	25°

Technical Specifications - Imperial

Specification	Size 1	Size 2
User Height	41.3" - 59"	50" - 66"
User Weight	39.6 - 110 lbs	39.6 - 132 lbs
Seat Depth	10.6" - 16.5"	13.8" - 18.5"
Seat Width	8" - 12.8"	8.7" - 13.6"
Backrest Height	14.2" - 18.5"	18.1" - 22.4"
Seat to Sandal	8.5" - 13.8"	12.4" - 18.5"
Chest Width (Distance between laterals)	6.7" - 10.6"	6.7" - 10.6"
Armrest Height	6.3" - 8.3"	8.3" - 10.2"
Knee Width	3.5" - 4.3"	4.7" - 5.5"
Tray Size	21.6" x 18.9"	21.6" x 18.9"
Seat Unit Weight	22 lbs	31.9 lbs

Size 2 Hi Lo Chassis (Excluding Push Bar)

Top of Seat to Floor	13.4" - 25.8"	13.4" - 25.8"
Base Weight	26.5 lbs	26.5 lbs

Technical Specifications - Metric

Specification	Size 1	Size 2
User Height	1050 mm - 1500 mm	1270 mm - 1680 mm
User Weight	18 - 50 kg	18 - 60 kg
Seat Depth	270 mm - 420 mm	350 mm - 470 mm
Seat Width	200 mm - 325 mm	220 mm - 345 mm
Backrest Height	360 mm - 470 mm	460 mm - 570 mm
Seat to Sandal	215 mm - 350 mm	315 mm - 470 mm
Chest Width (Distance between laterals)	170 mm - 270 mm	170 mm - 270 mm
Armrest Height	160 mm - 210 mm	210 mm - 260 mm
Knee Width	90 mm - 110 mm	120 mm - 140 mm
Tray Size	550 mm x 480 mm	550 mm x 480 mm
Seat Unit Weight	10 kgs	14.5 kgs

Size 2 Hi Lo Chassis (Excluding Push Bar)

Top of Seat to Floor	340 mm - 655 mm	340 mm - 655 mm
Base Weight	12 kgs	12 kgs



Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8, Canada



MK-129849-CAN Rev. 8 © 07/2025

