



Customer				Ship To			
Account				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City		Prov	
Buyer				PC	Phone		
Mark For				Email / Fax			

Please submit orders to: cscanada@sunmed.com

Seating Specifications

Typical WC Weight w/ Batteries Starting at 320 lbs. Weights vary depending on options


Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	17", 18", 19"

Drive Base Specifications

Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	300 lbs.	Drive Motors	4-pole 6mph
V6 Base Size: All Terrain	Width 28" Length 41.4" Turning Radius 23.5"	Drive Wheels	14" All Terrain Drive Wheels 14" Hybrid Drive Wheels
V6 Base Size: Hybrid	Width 27.3" Length 41.5" Turning Radius 23.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24 or 34
		Color Accents	9 colors
		Ground Clearance	3."

Client Measurements

a)	Top of head to seat	in
b)	Top of shoulder to seat	in
c)	Underarm to seat	in
d)	Elbow to seat	in
e)	Top of leg to seat	in
f)	Shoulder width	in
g)	Chest width	in
i)	Hip width	in
j)	Actual sitting depth	in
m)	Foot length	in
n)	Knee to foot	in
o)	Seat width	in



Notes

✓ OPTION	DESCRIPTION	RETAIL
SEATING OPTIONS		\$
<input type="checkbox"/> 600SF4	REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC
<input type="checkbox"/> 600SF5	MPS Seat and Back - Captain style seating	NC
<input type="checkbox"/> 600SF6	MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC
BASE		
Magic 360 Base - Includes Transit and lights End Item Number: EIMM		\$
<input type="checkbox"/> 600MW61	Magic Mobility - 360 No Power Seating Options	\$ 26,350.00
<input type="checkbox"/> 600MW62	Magic Mobility - 360 Single Power Seating Options	\$ 26,350.00
<input type="checkbox"/> 600MW63	Magic Mobility - 360 Multi Power Seating Options	\$ 26,350.00
<input type="checkbox"/> 600MW64	Magic Mobility - 360 No Power Seating Options - MPS (Captain)	\$ 26,350.00
<input type="checkbox"/> 600MW65	Magic Mobility - 360 Single Power Seating Options - MPS (Captain)	\$ 26,350.00
<input type="checkbox"/> 600MW66	Magic Mobility - 360 Multi Power Seating Options - MPS (Captain)	\$ 26,350.00
Drive Wheels		\$
<input type="checkbox"/> 600DW80	All Terrain Drive Wheels and Tires	NC
<input type="checkbox"/> 600DW82	Hybrid Drive Wheels and Tires - Solid	\$ 62.00
Optional Drive Tire Options		\$
<input type="checkbox"/> 6008A90	Tire Sealant for 2 Pneumatic Drive Tires	\$ 50.00
Base Accessories		\$
<input type="checkbox"/> 6008A91	Jack for Drive Tire	\$ 268.00
Frame Color		\$
<input type="checkbox"/> 600C185	Midnight Black	NC
<input type="checkbox"/> 600C186	Lucent Lime	NC
<input type="checkbox"/> 600C187	Sunset Red	NC
<input type="checkbox"/> 600C188	Horizon Blue	NC
<input type="checkbox"/> 600C189	Dusk Purple	NC
<input type="checkbox"/> 600C190	Vivid Pink	NC
<input type="checkbox"/> 600C191	Desert Orange	NC
<input type="checkbox"/> 600C192	Camouflage - Woods	NC
<input type="checkbox"/> 600C193	Camouflage - Desert	NC
<input type="checkbox"/> 600C171	Carbon Fiber	NC
<input type="checkbox"/> 600C194	American Flag	NC

Seat To Floor Height	\$
-----------------------------	-----------

360 Seat to Floor Heights start at 17" +/- 0.25"

<input type="checkbox"/> 600HT44	Lowest Seat To Floor Height	NC
<input type="checkbox"/> 600HT45	Seat To Floor Height +1"	NC
<input type="checkbox"/> 600HT46	Seat To Floor Height +2"	NC

Batteries	\$
------------------	-----------

<input type="checkbox"/> 600AG25	Group 24 Gel Battery (Pair) - Installed	\$ 1,184.00
----------------------------------	-----------------------------------------	-------------

Battery Charger	\$
------------------------	-----------

<input type="checkbox"/> 600A14	10A Battery Charger	NC
---------------------------------	---------------------	----

SEATING Seat Functions	\$
-------------------------------	-----------

<input checked="" type="checkbox"/> OPTION	DESCRIPTION	RETAIL
--------------------------------------------	-------------	--------

No Power Seat Options	\$
------------------------------	-----------

<input type="checkbox"/> 600ST87	No Power Seating System with Power Functionality	\$ 9,190.00
----------------------------------	--------------------------------------------------	-------------

Power Tilt	\$
-------------------	-----------

Utilized without Power Recline.

<input type="checkbox"/> 600PA64	50° Power Tilt in Space - 340 lb. weight capacity	\$ 10,085.00
----------------------------------	---------------------------------------------------	--------------

Power Seat Elevate	\$
---------------------------	-----------

Only available when Power Tilt or Power Tilt and Recline is selected.

<input type="checkbox"/> 600P16	12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 8,785.00
---------------------------------	--------------------------------------------------------	-------------

Power Recline	\$
----------------------	-----------

Utilized without Power Tilt. Requires No Power Seating with Power Functionality.

<input type="checkbox"/> 600PA65	150° Power Recline Only - 300 lb. weight capacity	\$ 7,178.00
----------------------------------	---------------------------------------------------	-------------

Power Tilt and Recline	\$
-------------------------------	-----------

<input type="checkbox"/> 600PA66	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 15,926.00
----------------------------------	---------------------------------------------------------	--------------

SEATING: SEAT SIZES	
----------------------------	--

<input checked="" type="checkbox"/> OPTION	DESCRIPTION	RETAIL
--------------------------------------------	-------------	--------

MPS (Captain) Seat Cushion (width x depth)	\$
---------------------------------------------------	-----------

MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.

<input type="checkbox"/> 600SI24	16" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI25	18" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI26	18" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI27	20" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI28	20" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI29	22" x 20" Vinyl Seat	\$ 608.00

REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes \$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/> 600PN13	Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/> 600PN14	Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 612.00
<input type="checkbox"/> 600PN15	Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 612.00

REHAB Seat Width

- 13"
 14"
 15"
 16"
 17"
 18"
 19"
 20"
 21"
 22"
 23"

REHAB Seat Depth

- 13"
 14"
 15"
 16"
 17"
 18"
 19"
 20"
 21"
 22"
 23"
 24"

REHAB Seating Cushion Options - \$

Note cushion size will be determined by Rehab Seat size selected above

<input type="checkbox"/> 600CU59	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 647.00
<input type="checkbox"/> 600CU60	JAY® Ion® (Size 22" Wide)	\$ 688.00
<input type="checkbox"/> 600CU61	JAY® Union® (Sizes 16" - 21" Wide)	\$ 753.00
<input type="checkbox"/> 600CU62	JAY® Union® (Size 22" Wide)	\$ 794.00
<input type="checkbox"/> 600CU64	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 1,030.00
<input type="checkbox"/> 600CU65	JAY® Fusion® (Size 22" Wide)	\$ 1,170.00

Lateral Thigh Support Pads and Hardware (must pick pad an \$

Lateral Thigh Support Pads

<input type="checkbox"/> 600KF22	Lateral Thigh Support Pad Short - Right (5"x5")	\$ 74.00
<input type="checkbox"/> 600KF23	Lateral Thigh Support Pad Short - Left (5"x5")	\$ 74.00

Lateral Thigh Support Mounts

<input type="checkbox"/> 600KF26	Lateral Thigh Removable Mount - Right	\$ 279.00
<input type="checkbox"/> 600KF27	Lateral Thigh Removable Mount - Left	\$ 320.00

SEATING: SEATBACK OPTIONS

MPS Backrest Options - Comes standard with Manual Recline \$

Backrest Width

<input type="checkbox"/> 600BN46	Standard Backrest Width - 16-20" Wide	NC
<input type="checkbox"/> 600BN47	Wide Backrest Width - 22-23" Wide	\$ 445.00

Backrest Height

<input type="checkbox"/> 600B80	Standard Backrest Height - 22" Tall	NC
<input type="checkbox"/> 600B81	Tall Backrest Height - 26" Tall	\$ 733.00

Backrest Material

<input type="checkbox"/> 600BU106	Cloth Backrest	NC
<input type="checkbox"/> 600BU107	Vinyl Backrest	\$ 404.00

MPS Back Accessories \$

<input type="checkbox"/> 600PH33	Push Handle - Needed for O2 holder, backpack, and cane holder	\$ 208.00
<input type="checkbox"/> 6009A13	Accessory Mount - Needed for aftermarket headrest and/or attendant control	\$ 157.00

REHAB Seating Backrest Options \$

Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.

JAY J3 Back Contour

<input type="checkbox"/> 600JP01	J3@ PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 899.00
<input type="checkbox"/> 600JP02	J3@ PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 1,912.00
<input type="checkbox"/> 600JP11	J3@ PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 971.00
<input type="checkbox"/> 600JP13	J3@ PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 1,052.00

JAY@ J3@ Back Shell Options and Height Ranges

<input type="checkbox"/> 600JR14	16.5" tall shell (Upper Thoracic Short)	NC
<input type="checkbox"/> 600JR15	18" tall shell (Upper Thoracic Med)	NC
<input type="checkbox"/> 600JR16	19.5" tall shell (Upper Thoracic Tall)	NC
<input type="checkbox"/> 600JR18	18" tall shell (Shoulder Height) for 22" wide	NC
<input type="checkbox"/> 600JR17	21" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR27	22.5" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR28	24" tall shell (Shoulder Height)	NC

REHAB Flip Forward Backrest Option \$

<input type="checkbox"/> 600D23	REHAB Flip Forward Backrest Option	#N/A
---------------------------------	-------------------------------------------	------

JAY@ Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only \$

JAY@ Lateral Thoracic Support Pads

<input type="checkbox"/> 600JW35	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$ 74.00
<input type="checkbox"/> 600JW38	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$ 74.00
<input type="checkbox"/> 600JW37	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$ 74.00
<input type="checkbox"/> 600JW40	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$ 74.00

JAY@ Lateral Thoracic Swing-Away Mount

<input type="checkbox"/> 600JW41	Lateral Thoracic Swing Away Bracket - Right	\$ 329.00
<input type="checkbox"/> 600JW42	Lateral Thoracic Swing Away Bracket - Left	\$ 329.00

JAY@ Anterior Thoracic Support \$

<input type="checkbox"/> 600ZR10	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 224.00
<input type="checkbox"/> 600ZR11	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 224.00
<input type="checkbox"/> 600K10	JAY@ Anterior Thoracic Support - Medium ATS8302	\$ 317.00
<input type="checkbox"/> 600K11	JAY@ Anterior Thoracic Support - Large ATS8303	\$ 317.00

Pelvic Positioning Options - Select One \$

Pelvic Positioning

<input type="checkbox"/> 600A848	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
<input type="checkbox"/> 600A849	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$ 101.00
<input type="checkbox"/> 600A850	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 157.00
<input type="checkbox"/> 600A851	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 268.00

Headrest Options - MPS Style \$

<input type="checkbox"/>	600NU12	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 544.00
<input type="checkbox"/>	600NU13	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 544.00

Headrest Options - REHAB Style \$

Headrest Pad

<input type="checkbox"/>	600HE39	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$ 371.00
<input type="checkbox"/>	600HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 438.00

Headrest Pad

<input type="checkbox"/>	600MH33	Whitmyer® Linx Headrest Mount with Quick Release	\$ 243.00
<input type="checkbox"/>	600MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 294.00

Armrest Style Options \$

<input type="checkbox"/>	600EG53	Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 1,762.00
<input type="checkbox"/>	600EG54	Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 2,501.00
<input type="checkbox"/>	600EG55	Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 3,120.00

Armrest Pad Options \$

<input type="checkbox"/>	600AB80	4" x 16" Tapered Pads - Right	\$ 146.00
<input type="checkbox"/>	600AE80	4" x 16" Tapered Pads - Left	\$ 146.00
<input type="checkbox"/>	600AB81	10" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE81	10" Waterfall Armrest Pad - Left	NC
<input type="checkbox"/>	600AB82	14" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE82	14" Waterfall Armrest Pad - Left	NC

LEGREST OPTIONS

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

Center Mount Fixed Legrest \$

Center Mount Fixed Legrest

<input type="checkbox"/>	600H210	Center Mount Flip-Up Footplate	\$ 601.00
--------------------------	---------	--------------------------------	-----------

Center Mount Fixed Legrest - Calfpads

<input type="checkbox"/>	600DI59	Fixed Center Mount Calf Pads	\$ 687.00
--------------------------	---------	------------------------------	-----------

Center Mount Power Elevating Legrest \$

Center Mount Power Elevating Legrest

<input type="checkbox"/>	600H213	Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 5,500.00
--------------------------	---------	-----------------------------------------------------------------------	-------------

Center Mount Power Elevating Legrest Extension

<input type="checkbox"/>	600E42	Center Mount Power Elevating Extension - Standard - Footplate range 5"-15"	NC
<input type="checkbox"/>	600E43	Center Mount Power Elevating Extension - Short - Footplate range 5"-12.5"	NC

Center Mount Power Elevating Legrest Footplates

<input type="checkbox"/>	600F117	Center Mount Power Elevating Footplate - Standard	NC
<input type="checkbox"/>	600F118	Center Mount Power Elevating Footplate - Short	NC

Swing-Away Legrest \$

Swing-Away Legrest

<input type="checkbox"/> 600H211	Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$ 864.00
----------------------------------	----------------------------------------------------------------------------------------------	-----------

S/A Non-Power Positioning Options

<input type="checkbox"/> 600LS15	Heel and Side Support - Attaches to the footplate stem and footplate	\$ 870.00
<input type="checkbox"/> 600LS16	Foot Strap - 1.5" wide strap	\$ 118.00

Swing-Away Power Elevating Legrest \$

Swing-Away Power Elevating Legrest

<input type="checkbox"/> 600H214	Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate	\$ 6,004.00
----------------------------------	----------------------------------------------------------------------------------------------	-------------

S/A Power Legrest Positioning Options

<input type="checkbox"/> 600DI58	2" Metal Heel Guards for 2 Piece Footplates	\$ 1,134.00
<input type="checkbox"/> 600HL15	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 168.00

ACCESSORIES \$

<input type="checkbox"/> 6008A93	Waterproof Wheelchair Cover	\$ 219.00
<input type="checkbox"/> 6008A94	Waterproof Joystick Cover	\$ 68.00
<input type="checkbox"/> 6008A95	Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 50.00
<input type="checkbox"/> 600A221	Deluxe Backpack	\$ 142.00
<input type="checkbox"/> 600A615	Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 283.00
<input type="checkbox"/> 600A719	Articulating Oxygen Holder for MPS Seating System	\$ 1,470.00
<input type="checkbox"/> 600A612	Cane Holder for MPS Seating System	\$ 943.00
<input type="checkbox"/> 6008A96	Luggage Rack for 360 - Attaches to base behind the chair (10.5" x 13")	\$ 492.00
<input type="checkbox"/> 6008A97	Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 791.00
<input type="checkbox"/> 600A613	Cell Phone Holder - Adjustable to fit most cell phones	\$ 438.00
<input type="checkbox"/> 600A614	Fishing Rod Holder	\$ 438.00

Be Adaptive Products

<input type="checkbox"/> 6008A98	Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 2,209.00
<input type="checkbox"/> 6008A99	Be Adaptive Trigger Mechanism	\$ 1,224.00
<input type="checkbox"/> 6008A100	Be Adaptive Limited Mobility Fishing Mount	\$ 2,209.00

ELECTRONICS

✓ OPTION	DESCRIPTION	RETAIL
Joystick		\$
<input type="checkbox"/> 600J121	Standard LED Joystick (w/o Screen Display)	NC
<input type="checkbox"/> 600J116	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
<input type="checkbox"/> 600J118	High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display	\$ 746.00
Joystick Mount		\$
<input type="checkbox"/> 600JM34	Fixed on Armrest - Depth adjustable 8.25"	NC
<input type="checkbox"/> 600JM45	Swing-Away Joystick Mount - Includes latch to lock in the forward position	\$ 836.00
<input type="checkbox"/> 600JM46	Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)	\$ 2,020.00
Joystick Mount - Location		\$
<input type="checkbox"/> 600HM16	Right Hand Armrest	NC
<input type="checkbox"/> 600HM17	Left Hand Armrest	NC
Display Language Selection R-net		\$
<input type="checkbox"/> 600L1	English	NC
<input type="checkbox"/> 600L3	French	NC
<input type="checkbox"/> 600L7	Spanish	NC
Joystick Handle Type		\$
Comes standard with PG Cone-Shaped Joystick Handle		
<input type="checkbox"/> 600JH13	Ball Handle	\$ 108.00
<input type="checkbox"/> 600JH26	Bodypoint Grooved Mushroom Joystick Handle	\$ 236.00
<input type="checkbox"/> 600JH24	Bodypoint – Rubber Dome Joystick Handle	\$ 252.00
<input type="checkbox"/> 600JH25	Bodypoint – U-Shaped Joystick Handle (4 inch)	\$ 252.00
Expandable Power Module		\$
Must select both items		
<input type="checkbox"/> 600EX05	Expandable 120A Controller - R-net	\$ 1,083.00
<input type="checkbox"/> 600EX06	Harness for 120A Expandable Controller	\$ 757.00
Power Seat Function		\$
<input type="checkbox"/> 600EO37	THRU DRIVE CNTRL 1 ACTUATOR- Must be selected if ordering 1 actuator function	\$ 1,422.00
<input type="checkbox"/> 600EO38	THRU DRIVE CNTRL 2+ ACTUATORS- Must be selected if ordering 2 or more actuator functions	\$ 2,860.00
Power Seat Control Option		\$
<input type="checkbox"/> 600QG11	Ctrl+5 Switch Box with Free Assignable Functions	\$ 1,002.00
<input type="checkbox"/> 600QG12	Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)	NC
Ctrl+5 Switch Type		\$
<input type="checkbox"/> 600EY22	Push Button	NC
<input type="checkbox"/> 600EY23	Toggle Switches	\$ 351.00
Power Seat Control Option - Mounting location		\$

<input type="checkbox"/>	600EY19	CTRL+5 MNT ARMREST RIGHT	NC
<input type="checkbox"/>	600EY20	CTRL+5 MNT ARMREST LEFT	NC
Attendant Controls			\$
<input type="checkbox"/>	600A05	Dual Control Module for Attendant	\$ 2,424.00
Attendant Control Mounting Location			\$
<input type="checkbox"/>	600DN5	Right Back Cane	NC
<input type="checkbox"/>	600DN6	Left Back Cane	NC
R-net Omni-2 - 3.5" LCD Display w/ Built-In Infrared & Blueto			\$
<input type="checkbox"/>	600ON14	R-net Omni-2 Only Input Control Module	\$ 3,558.00
<input type="checkbox"/>	600ON15	R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices	NC
R-net Omni-2 -Mounting Location			\$
<input type="checkbox"/>	600DS23	Right Hand Armrest	NC
<input type="checkbox"/>	600DS24	Left Hand Armrest	NC
Head Control Packages			\$
<input type="checkbox"/>	600EU25	Switch-It DualPro® Proportional Head Array - Facial Pads	\$ 6,953.00
<input type="checkbox"/>	600EU26	Switch-It DualPro® Proportional Head Array - Spot Pads	\$ 6,953.00