

MWD Crossover Compact Electric Wheelchair

EIPW60



Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		Province	
Buyer				PC	Phone		
Mark For				Email			

Please submit orders to: cscanada@sunmed.com

Seating Specifications

Typical WC Weight w/ Batteries Starting at 320 lbs. Weights vary depending on options

Power Seating 50° of CG Tilt / Power Recline / 12" Seat Elevator

Seat Type MPS and Rehab Style

Seat Size Range Width: 13-23" | Depth: 13-24"

STF Height Range 17", 18", 19"

Drive Base Specifications

Drive Wheel Position Mid Wheel Drive

Drive Control R-Net PM 120A controller

User Weight Capacity 350 lbs.

Drive Motors 4-pole 6mph

Base Size: Off-Road Width 26" | Length 39" | Turning Radius 22.5"

Drive Wheels 14" Off-Road Drive Wheels | 14" Urban Drive Wheels

Base Size: Urban Width 24" | Length 39" | Turning Radius 22.5"

Tires / Casters Pneumatic or Solid inserts

Charger Dual Mode


Battery Size Group 24

Ground Clearance 3"

Color Accents 10 colors

Client Measurements

a)	Top of head to seat	_____	in
b)	Top of shoulder to seat	_____	in
c)	Underarm to seat	_____	in
d)	Elbow to seat	_____	in
e)	Top of leg to seat	_____	in
f)	Shoulder width	_____	in
g)	Chest width	_____	in
i)	Hip width	_____	in
j)	Actual sitting depth	_____	in
m)	Foot length	_____	in
n)	Knee to foot	_____	in
o)	Seat width	_____	in



Notes

✓ OPTION		DESCRIPTION	RETAIL
Seating Options			\$
☐ 600SF4		REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC
☐ 600SF5		MPS Seat and Back - Captain style seating	NC
☐ 600SF6		MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC
Base			
Magic 360 Base - Includes Transit and lights		End Item Number: EIPW60	\$
☐ 600MW61		Magic Mobility - 360 No Power Seating Options	\$ 27,547.00
☐ 600MW62		Magic Mobility - 360 Single Power Seating Options	\$ 27,547.00
☐ 600MW63		Magic Mobility - 360 Multi Power Seating Options	\$ 27,547.00
☐ 600MW64		Magic Mobility - 360 No Power Seating Options - MPS (Captain)	\$ 27,547.00
☐ 600MW65		Magic Mobility - 360 Single Power Seating Options - MPS (Captain)	\$ 27,547.00
☐ 600MW66		Magic Mobility - 360 Multi Power Seating Options - MPS (Captain)	\$ 27,547.00
Drive Wheels			\$
☐ 600DW80		Off-Road Drive Wheels and Tires	NC
☐ 600DW82		Urban Drive Wheels and Tires - Solid	\$ 65.00
Optional Drive Tire Options			\$
☐ 6008A90		Tire Sealant for 2 Pneumatic Drive Tires	\$ 53.00
Base Accessories			\$
☐ 6008A91		Jack for Drive Tire	\$ 280.00
Frame Color			\$
☐ 600C185		Midnight Black	NC
☐ 600C186		Lucent Lime	NC
☐ 600C187		Sunset Red	NC
☐ 600C188		Horizon Blue	NC
☐ 600C189		Dusk Purple	NC
☐ 600C190		Vivid Pink	NC
☐ 600C191		Desert Orange	NC
☐ 600C192		Camouflage - Woods	NC
☐ 600C193		Camouflage - Desert	NC
☐ 600C171		Carbon Fiber	NC
☐ 600C194		American Flag	NC
Seat To Floor Height			\$
Seat to Floor Heights start at 17" +/- 0.25"			
☐ 600HT44		Lowest Seat To Floor Height	NC
☐ 600HT45		Seat To Floor Height +1"	NC
☐ 600HT46		Seat To Floor Height +2"	NC
Batteries			\$
☐ 600AG25		Group 24 Gel Battery (Pair) - Installed	\$ 1,238.00
Battery Charger			\$
☐ 600BJ28		10A Battery Charger	NC
☐ 600BJ38		12A Battery Charger	\$ 716.00

Seating: Seat Functions

✓ OPTION	DESCRIPTION	RETAIL
No Power Seat Options		\$
☐ 600ST87	No Power Seating System with Power Functionality	\$ 9,608.00
Power Tilt		\$
Utilized without Power Recline.		
☐ 600PA64	50° Power Tilt in Space	\$ 10,543.00
Power Seat Elevate		\$
Only available when Power Tilt or Power Tilt and Recline is selected.		
☐ 600PI16	12" Adjustable Seat Elevator	\$ 9,185.00
Power Recline		\$
Utilized without Power Tilt. Requires No Power Seating with Power Functionality.		
☐ 600PA65	150° Power Recline Only	\$ 7,505.00
Power Tilt and Recline		\$
☐ 600PA66	50° Tilt & 150° Power Recline	\$ 16,650.00

Seating: Seat Sizes

MPS (Captain) Seat Cushion (width x depth)		\$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.		
☐ 600SI24	16" x 18" Vinyl Seat	NC
☐ 600SI25	18" x 18" Vinyl Seat	NC
☐ 600SI26	18" x 20" Vinyl Seat	NC
☐ 600SI27	20" x 18" Vinyl Seat	NC
☐ 600SI28	20" x 20" Vinyl Seat	NC
☐ 600SI29	22" x 20" Vinyl Seat	\$ 636.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes		\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

☐ 600PN13	Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
☐ 600PN14	Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 640.00
☐ 600PN15	Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 640.00

REHAB Seat Width

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options		\$
-------------------------------	--	----

Note cushion size will be determined by Rehab Seat size selected above

☐ 600CU59	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 677.00
☐ 600CU60	JAY® Ion® (Size 22" Wide)	\$ 719.00
☐ 600CU61	JAY® Union® (Sizes 16" - 21" Wide)	\$ 787.00
☐ 600CU62	JAY® Union® (Size 22" Wide)	\$ 830.00
☐ 600CU64	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 1,076.00
☐ 600CU65	JAY® Fusion® (Size 22" Wide)	\$ 1,224.00

Seating: Seat Sizes			
✓ OPTION		DESCRIPTION	RETAIL
Support Pads and Hardware (must pick pad and mount)			\$
Lateral Thigh Support Pads			
☐ 600KF22		Lateral Thigh Support Pad Short - Right (5x5)	\$ 77.00
☐ 600KF23		Lateral Thigh Support Pad Short - Left (5x5)	\$ 77.00
Lateral Thigh Support Mounts			
☐ 600KF26		Lateral Thigh Removable Mount - Right	\$ 335.00
☐ 600KF27		Lateral Thigh Removable Mount - Left	\$ 335.00
Seating: Seatback Options			
MPS Backrest Options - Comes standard with Manual Recline			\$
Backrest Width			
☐ 600BN46		Standard Backrest Width - 16-20" Wide	NC
☐ 600BN47		Wide Backrest Width - 22-23" Wide	\$ 466.00
Backrest Height			
☐ 600B80		Standard Backrest Height - 22" Tall	NC
☐ 600B81		Tall Backrest Height - 26" Tall	\$ 766.00
Backrest Material			
☐ 600BU106		Cloth Backrest	NC
☐ 600BU107		Vinyl Backrest	\$ 422.00
MPS Back Accessories			\$
☐ 600PH33		Push Handle - Needed for O2 holder	\$ 217.00
☐ 6009A13		Accessory Mount - Needed for aftermarket headrest	\$ 164.00
REHAB Seating Backrest Options			\$
Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
JAY J3 Back Contour			
☐ 600JP01		J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 939.00
☐ 600JP02		J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 1,999.00
☐ 600JP11		J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 1,016.00
☐ 600JP13		J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 1,100.00
REHAB Seating Backrest Options			\$
JAY® J3® Back Shell Options and Height Ranges			
☐ 600JR41		15" tall shell (Upper Thoracic Med) for 22" wide	NC
☐ 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
☐ 600JR15		18" tall shell (Upper Thoracic Med)	NC
☐ 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
☐ 600JR18		18" tall shell (Shoulder Height) for 22" wide	NC
☐ 600JR17		21" tall shell (Shoulder Height)	NC
☐ 600JR27		22.5" tall shell (Shoulder Height)	NC
☐ 600JR28		24" tall shell (Shoulder Height)	NC
REHAB Flip Forward Backrest Option			\$
☐ 600OD23		REHAB Flip Forward Backrest Option	\$ 220.00

Seating: Seatback Options			
✓ OPTION		DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
JAY® Lateral Thoracic Support Pads			
☐ 600JW35		Lateral Thoracic (L x H) 4" x 5" Pad - Right	\$ 77.00
☐ 600JW38		Lateral Thoracic (L x H) 4" x 5" Pad - Left	\$ 77.00
☐ 600JW37		Lateral Thoracic (L x H) 6" x 4" Pad - Right	\$ 77.00
☐ 600JW40		Lateral Thoracic (L x H) 6" x 4" Pad - Left	\$ 77.00
JAY® Lateral Thoracic Swing-Away Mount			
☐ 600JW41		Lateral Thoracic Swing Away Bracket - Right	\$ 344.00
☐ 600JW42		Lateral Thoracic Swing Away Bracket - Left	\$ 344.00
JAY® Anterior Thoracic Support			\$
☐ 600ZR10		Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 234.00
☐ 600ZR11		Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 234.00
☐ 600K10		JAY® Anterior Thoracic Support - Medium ATS8302	\$ 332.00
☐ 600K11		JAY® Anterior Thoracic Support - Large ATS8303	\$ 332.00
Pelvic Positioning Options - Select One			\$
☐ 600A848		Positioning Belt, Auto Style Buckle, 2" x 60"	NC
☐ 600A849		Positioning Belt, Aircraft Buckle, Padded 2" x 50"	\$ 106.00
☐ 600A850		Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 164.00
☐ 600A851		Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 280.00
Headrest Options - MPS Style			\$
☐ 600NU12		Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 569.00
☐ 600NU13		Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 569.00
Headrest Options - REHAB Style			\$
Headrest Pad			
☐ 600HE39		Whitmyer® Plush Headrest Pad, 10"	\$ 388.00
☐ 600HE46		Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 458.00
Headrest Mount			
☐ 600MH33		Whitmyer® Linx Headrest Mount with Quick Release	\$ 254.00
☐ 600MH34		Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 307.00
Armrest Style Options			\$
☐ 600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 1,842.00
☐ 600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 2,615.00
☐ 600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 3,262.00
Armrest Pad Options			\$
☐ 600AB81		10" Waterfall Armrest Pad - Right	NC
☐ 600AE81		10" Waterfall Armrest Pad - Left	NC
☐ 600AB82		14" Waterfall Armrest Pad - Right	NC
☐ 600AE82		14" Waterfall Armrest Pad - Left	NC

Legrest Options			
✓ OPTION		DESCRIPTION	RETAIL
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
☐ 600H210		Center Mount (14.5" - 17" Leg Length) with Flip-Up Footplate (11" x 10")	\$ 628.00
Center Mount Fixed Legrest - Accessories			
☐ 600DI59		Fixed Center Mount Calf Pads	\$ 718.00
☐ 600LS16		Foot Strap - 1.5" wide strap	\$ 124.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest			
☐ 600H213		Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 5,750.00
Center Mount Power Elevating Legrest Extension			
☐ 600E42		Center Mount Power Elevating Extension - Standard (15" - 17" Leg Length)	NC
☐ 600E43		Center Mount Power Elevating Extension - Short (10" - 15" Leg Length)	NC
Center Mount Power Elevating Legrest Footplates			
☐ 600F117		Center Mount Power Elevating Footplate - Standard (5" x 12" each)	NC
☐ 600F118		Center Mount Power Elevating Footplate - Short (4.5" x 7" each)	NC
Swing-Away Legrest			\$
Swing-Away Legrest			
☐ 600H211		Swing-Away Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	\$ 903.00
Swing-Away Power Elevating Legrest			\$
Swing-Away Power Elevating Legrest			
☐ 600H214		Swing-Away Power Elevating Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	\$ 6,277.00
S/A Power Legrest Positioning Options			
☐ 600HL15		Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 176.00
Accessories			\$
☐ 6008A93		Waterproof Wheelchair Cover	
☐ 6008A94		Waterproof Joystick Cover	\$ 71.00
☐ 6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 53.00
☐ 600A221		Deluxe Backpack	\$ 148.00
☐ 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 296.00
☐ 600A719		Articulating Oxygen Holder for MPS Seating System	\$ 1,537.00
☐ 600A612		Cane Holder for MPS Seating System - Includes push handles	\$ 986.00
☐ 6008A96		Luggage Rack for 360 - Attaches to base behind the chair (10.5" x 13")	\$ 514.00
☐ 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 827.00
☐ 600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$ 458.00
☐ 600A614		Fishing Rod Holder	\$ 458.00
Be Adaptive Products			
☐ 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 2,309.00
☐ 6008A99		Be Adaptive Trigger Mechanism	\$ 1,279.00
☐ 6008A100		Be Adaptive Limited Mobility Fishing Mount	\$ 2,309.00

Electronics									
✓ OPTION			DESCRIPTION					RETAIL	
Joystick									\$
☐ 600J121			Standard LED Joystick (w/o Screen Display)					NC	
☐ 600J116			Color Joystick (2.3" Color LCD Screen Display, CJSM1)					NC	
☐ 600J118			High Resolution Color Joystick (CJSM2+BT) 3.5" Color LCD Display					\$ 780.00	
Joystick Mount									\$
☐ 600JM34			Fixed on Armrest - Depth adjustable 8.25"					NC	
☐ 600JM45			Swing-Away Joystick Mount - Includes latch to lock in the forward position					\$ 874.00	
☐ 600JM46			Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)					\$ 2,112.00	
Joystick Mount - Location									\$
☐ 600HM16			Right Hand Armrest					NC	
☐ 600HM17			Left Hand Armrest					NC	
Display Language Selection R-net									\$
☐ 600L1	English	NC	☐ 600L3	French	NC	☐ 600L7	Spanish	NC	
Joystick Handle Type									\$
Comes standard with PG Cone-Shaped Joystick Handle									
☐ 600JH13			Ball Handle					\$ 113.00	
☐ 600JH26			Bodypoint Grooved Mushroom Joystick Handle					\$ 247.00	
☐ 600JH24			Bodypoint – Rubber Dome Joystick Handle					\$ 264.00	
☐ 600JH25			Bodypoint – U-Shaped Joystick Handle (4 inch)					\$ 264.00	
Expandable Power Module									\$
Must select both items									
☐ 600EX05			Expandable 120A Controller - R-net					\$ 1,132.00	
☐ 600EX06			Harness for 120A Expandable Controller					\$ 791.00	
Power Seat Function									\$
☐ 600EO37			Thru Drive Control 1 Actuator - Must be selected if ordering 1 actuator function					\$ 1,486.00	
☐ 600EO38			Thru Drive Control 2+ Actuators - Must be selected if ordering 2 or more actuator functions					\$ 2,990.00	
Power Seat Control Option									\$
☐ 600QG11			Ctrl+5 Switch Box with Free Assignable Functions					\$ 1,048.00	
☐ 600QG12			Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)					NC	
Ctrl+5 Switch Type									\$
☐ 600EY22			Push Button					NC	
☐ 600EY23			Toggle Switches					NC	
Power Seat Control Option - Ctrl+5 Mounting location									\$
☐ 600EY19	Right Armrest	NC	☐ 600EY20	Left Armrest	NC				
Attendant Controls									\$
☐ 6000A05			Dual Control Module for Attendant - Includes accessory mount and push handle					\$ 2,534.00	
Attendant Control Mounting Location									\$
☐ 600DN5	Right Back Cane	NC	☐ 600DN6	Left Back Cane	NC				
R-net Omni2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth									\$
☐ 600ON14			R-net Omni2 Only Input Control Module					\$ 3,719.00	
☐ 600ON15			R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices					NC	
R-net Omni2 - Mounting Location									\$
☐ 600DS23	Right Hand Armrest	NC	☐ 600DS24	Left Hand Armrest	NC				
Head Control Packages									\$
☐ 600EU25			Switch-It DualPro® Proportional Head Array - Facial Pads					\$ 7,269.00	
☐ 600EU26			Switch-It DualPro® Proportional Head Array - Spot Pads					\$ 7,269.00	