

# FORTRESS

# S410



Fortress is proud to introduce the new S-410 Scooter!

NEWEST MEMBER OF THE  
**S-SERIES**  
PREMIUM SCOOTER RANGE

An internationally recognized name for reliability, performance and durability - the Fortress 4-wheel scooter is ideal for everyday use indoors and out. With the smallest turning radius in the Fortress lineup, riders will enjoy a comfortably-sized scooter without sacrificing easy access to their surroundings.

The new S-410 is available in Standard and High Performance models. The S-410 Std is perfect for daily use with a 10 km/hour top speed and 300 lb capacity. The High Performance option offers an increased weight capacity (up to 350 lbs) and additional traction when necessary thanks to a 120-amp controller that provides more torque and speeds up to 10.5 km/hour.

The S-410 features the large Fortress basket that provides riders with a convenient space to carry essential items. Support from the ergonomic padded seat makes for a smooth and comfortable ride.



## TECHNICAL DATA

	STD Model	HP Model
Maximum speed	10 km/h	10.5 km/h
Maximum user weight	300 lbs.	350 lbs.
Charger	8 amp	8 amp
Range	34 km	34 km
All around LED lighting	Yes	Yes
Suspension	Yes	Yes
Tires	10"	10"
Overall width /length	48"x 25.4"	48"x 25.4"
Seat width x depth	19" x 19"	19" x 19"
Seat height range	17.5" - 19.5"	17.5" - 19.5"
Turning radius	58.2"	58.2"
Battery group	Group 22	Group 22
Controller	90 amp	120 amp
Color options	Black, Blue, Red	

 **SUNRISE**  
MEDICAL

Sunrise Medical Canada Inc.

237 Romina Drive, Unit 3 · Concord, ON L4K 4V3

Telephone: 1-800-263-3390

Fax: 1-800-561-5834

E-Mail: [cscanada@sunmed.com](mailto:cscanada@sunmed.com)

[www.sunrisemedical.ca](http://www.sunrisemedical.ca)

Rev.A