



ZIPPIE® Sphynx™

March 2022



EIZ22

Mark For: _____ Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ PV: _____ PC: _____ Dealer Phone: _____ Fax: _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ <p style="text-align: center;">ADDITIONAL SHIPPING INFORMATION</p> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ PV: _____ PC: _____ Ship To Phone: _____ Fax: _____
--	--

ZIPPIE Sphynx

Each ZIPPIE Sphynx unit includes the following components:

- Silver & Black Folding Frame
- 10°, 20°, 30° Seat Tilt Settings
- Quick Adjust Recline Backrest
- Foot Operated Wheel Lock
- 12" Rear Wheels & 6" Casters (Solid)
- Quick Release Rear Wheels
- Adjustable Stroller Handle
- Height & Angle Adjustable Footplates
- Seat Sling with Depth Adjustability

Sphynx Model

Selection Required. Default Tilt setting is 10°.

	Item Code	Description	Retail
<input type="checkbox"/>	149WM103	ZIPPIE Sphynx (Includes upholstery, Transit optional)	\$3,435.00
<input type="checkbox"/>	149WM104	ZIPPIE Sphynx without Upholstery (Transit optional)	\$3,150.00
<input type="checkbox"/>	149WM111	ZIPPIE Sphynx Package (Includes Transit, Upholstery, Canopy, Underseat Basket, Rain Cover)	\$4,150.00

Transit System

Transit System includes Transit Brackets, Headrest Extension required

	Item Code	Description	Retail
<input type="checkbox"/>	149TO44	Transit System	NC
<input type="checkbox"/>	149TO40	OMIT transit (only available on transit optional Sphynx models)	NC

Pre-Set Seat Tilt

	Item Code	Description	Retail
<input type="checkbox"/>	149TA15	10° Tilt	NC
<input type="checkbox"/>	149TA16	20° Tilt	NC
<input type="checkbox"/>	149TA17	30° Tilt	NC

Seat Width

Selection Required. Seat width will have an adjustable Seat Depth range

	Item Code	Description	Retail
<input type="checkbox"/>	149W1	10" Seat Width (Adjustable Seat Depth of 10" to 15")	NC
<input type="checkbox"/>	149W3	12" Seat Width (Adjustable Seat Depth of 10" to 15")	NC
<input type="checkbox"/>	149W5	14" Seat Width (Adjustable Seat Depth of 11" to 16")	NC
<input type="checkbox"/>	149W7	16" Seat Width (Adjustable Seat Depth of 13" to 18")	NC
<input type="checkbox"/>	149W9	18" Seat Width (Adjustable Seat Depth of 14" to 19")	NC

Seat Depth

Selection Required. Seat depth is adjustable within the range stated above for the Seat Width

Item Code	Description	Retail
<input type="checkbox"/> 149SD3	10" Seat Depth (Only available on 10" and 12" Seat Widths)	NC
<input type="checkbox"/> 149SD4	11" Seat Depth (Only available on 10", 12" and 14" Seat Widths)	NC
<input type="checkbox"/> 149SD5	12" Seat Depth (Only available on 10", 12" and 14" Seat Widths)	NC
<input type="checkbox"/> 149SD6	13" Seat Depth (N/A on 18" Seat Width)	NC
<input type="checkbox"/> 149SD7	14" Seat Depth	NC
<input type="checkbox"/> 149SD8	15" Seat Depth	NC
<input type="checkbox"/> 149SD9	16" Seat Depth (N/A on 10" or 12" Seat Widths)	NC
<input type="checkbox"/> 149SD10	17" Seat Depth (Only available on 16" and 18" Seat Widths)	NC
<input type="checkbox"/> 149SD11	18" Seat Depth (Only available on 16" and 18" Seat Widths)	NC
<input type="checkbox"/> 149SD12	19" Seat Depth (Only available on 18" Seat Width)	NC

Seat-To-Footplate Range (Lower Leg Length)

Selection Required. Reference table below for achievable Seat-To-Footplate distances relative to the Seat Width & Depth.

Item Code	Description	Retail
<input type="checkbox"/> 149LO64	Extra Short Seat-To-Footplate Range	NC
<input type="checkbox"/> 149LO61	Short Seat-To-Footplate Range	NC
<input type="checkbox"/> 149LO62	Medium Seat-To-Footplate Range	NC
<input type="checkbox"/> 149LO63	Long Seat-To-Footplate Range	NC

Seat Width	Seat-to-Footplate Range Option	Achievable Seat-to-Footplate Distances*									
		Seat Depth									
		10"	11"	12"	13"	14"	15"	16"	17"	18"	19"
10" 12"	Extra Short	6"-12"	6"-11"	6"-10.5"	6"-10"	5"-10"	5"-10"	N/A	N/A	N/A	N/A
	Short	12"-15.5"	11"-15"	10.5"-15"	10"-14.5"	10"-14.5"	10"-14.5"	N/A	N/A	N/A	N/A
	Medium	16"-17.5"	15"-17.5"	14.5"-17"	14"-17.5"	14"-17.5"	13.5"-16.5"	N/A	N/A	N/A	N/A
	Long	18"-19"	16.5"-19"	16"-19"	16"-19"	16"-19"	15.5"-19"	N/A	N/A	N/A	N/A
14"	Extra Short	N/A	6"-11"	6"-10.5"	6"-10"	5"-10"	5"-10"	5"-10"	N/A	N/A	N/A
	Short	N/A	11"-15"	10.5"-15"	10"-14.5"	10"-14.5"	10"-14.5"	10"-14.5"	N/A	N/A	N/A
	Medium	N/A	15"-17.5"	14.5"-17"	14"-17.5"	14"-17.5"	13.5"-16.5"	13"-16"	N/A	N/A	N/A
	Long	N/A	16.5"-19"	16"-19"	16"-19"	16"-19"	15.5"-19"	15"-19"	N/A	N/A	N/A
16"	Extra Short	N/A	N/A	N/A	6"-10"	5"-10"	5"-10"	5"-10"	6"-10"	7"-10"	N/A
	Short	N/A	N/A	N/A	10"-14.5"	10"-14.5"	10"-14.5"	10"-14.5"	10"-14"	10"-14"	N/A
	Medium	N/A	N/A	N/A	14"-17.5"	14"-17.5"	13.5"-16.5"	13"-16"	13"-15.5"	13"-15.5"	N/A
	Long	N/A	N/A	N/A	16"-19"	16"-19"	15.5"-19"	15"-19"	15"-19"	14.5"-19"	N/A
18"	Extra Short	N/A	N/A	N/A	N/A	10"-14.5"	10"-14"	10"-14"	10"-14"	10"-13.5"	10"-11.5"
	Short	N/A	N/A	N/A	N/A	12"-15"	11.5"-15"	11.5"-14.5"	11.5"-14.5"	11.5"-14"	11.5"-13.5"
	Medium	N/A	N/A	N/A	N/A	14"-17"	13.5"-16.5"	13"-16"	13"-15.5"	13"-15.5"	13"-15.5"
	Long	N/A	N/A	N/A	N/A	15.5"-19"	15.5"-19"	15"-19"	15"-19"	14.5"-19"	14.5"-19"

*NOTE: All seat-to-footplate distances based on 10 deg tilt setting (seat to floor angle).

Sphynx Upholstery for 149WM103 Sphynx Model

Selection Required for Upholstery-Included Sphynx Model (149WM103). Includes seat and back cushions with tension adjustable straps and seat sling.

Item Code	Description	Retail
<input type="checkbox"/> 149UC37	Sphynx Upholstery - Black	STD
<input type="checkbox"/> 149UC49	Sphynx Upholstery - Red	NC
<input type="checkbox"/> 149UC50	Sphynx Upholstery - Blue	NC
<input type="checkbox"/> 149UC51	Sphynx Upholstery - Green	NC
<input type="checkbox"/> 149UC52	Sphynx Upholstery - Pink	NC
<input type="checkbox"/> 149UC53	OMIT upholstery	NC

JAY Options

Required to attach JAY order forms for bundled products

Item Code	Description	JAY Part #:	Retail
<input type="checkbox"/> 149JC20	JAY Cushion		See OF
<input type="checkbox"/> 149JB02	JAY Backrest		See OF

Towel Bars

Item Code	Description	Retail
<input type="checkbox"/> 149QT02	2" Towel Bar	NC
<input type="checkbox"/> 149QT01	4" Towel Bar	NC

Head Supports

Item Code	Description	Retail
<input type="checkbox"/> 149HE75	Headrest Extension <i>(Required with transit option)</i>	\$220.00
<input type="checkbox"/> 149ZH25	Padded Headwings <i>(May interfere with full torso vest)</i>	\$250.00

Trunk & Arm Positioning

Item Code	Description	Retail
<input type="checkbox"/> 149KI7	H-Harness <i>(Requires 3-Point Pelvic Belt and Headrest Ext., N/A with Full Torso Support Vest)</i>	\$152.00
<input type="checkbox"/> 149JY06	JAY Center Opening Dynamic <i>(Requires Headrest Ext. & Select harness size below)</i>	\$238.00
<input type="checkbox"/> 149JZ01	Small	NC
<input type="checkbox"/> 149JZ02	Medium	NC
<input type="checkbox"/> 149JZ03	Large	NC
<input type="checkbox"/> 149JZ04	X-Large	NC
<input type="checkbox"/> 149KI9	Full Torso Support Vest <i>(Requires Headrest Ext., N/A with H-Harness or Lateral Supports)</i>	\$326.00
<input type="checkbox"/> 149KI8	Adjustable Lateral Support & Scoli Strap <i>(Requires Armrests, N/A with Full Torso Vest)</i>	\$211.00
<input type="checkbox"/> 149AA1	Flip-Up Armrests (Angle Adjustable)	\$326.00
<input type="checkbox"/> 149A206	Upper Extremity Support Tray - Grey <i>(Requires armrests)</i>	\$286.00
<input type="checkbox"/> 149A207	Upper Extremity Support Tray - Clear <i>(Requires armrests)</i>	\$343.00

Pelvic Positioning

	Item Code	Description		Retail
<input type="checkbox"/>	149ZZ29	3-Point Pelvic Belt <i>(Requires H-Harness, N/A with Medial Thigh Support)</i>		\$125.00
<input type="checkbox"/>	149ZZ11	2-Point positioning belt w/ 4.5" pad		\$121.00
<input type="checkbox"/>	149ZZ24	V-Contour 4-Point Belt- Small Pad		\$227.00
<input type="checkbox"/>	149ZZ23	V-Contour 4-Point Belt- Med Pad		\$227.00
<input type="checkbox"/>	149KP05	Hip Guides - Inset Mount (Pair) <i>(Reduces seat width by 1" per side)</i>		\$250.00
<input type="checkbox"/>	149KP06	Hip Guides - In-Line Mount (Pair)		\$250.00
<input type="checkbox"/>	149CU30	Cushion 1" Foam Riser		\$125.00
<input type="checkbox"/>	149CU12	Wedge Cushion Insert		\$150.00
<input type="checkbox"/>	149KH11	Medial Thigh Support <i>(includes hardware, N/A with H-Harness & 3-Point Pelvic Belt)</i>		\$181.00


Foot & Leg Positioning

	Item Code	Description		Retail
<input checked="" type="checkbox"/>	149F26	Angle Adjustable Aluminum Footplates (Pair) Required		\$169.00
<input type="checkbox"/>	149HL5	Leg Strap		NC
<input type="checkbox"/>	149A54	Extra Leg Strap		\$60.00
<input type="checkbox"/>	149F304	Foot & Ankle Positioners - Small		\$250.00
<input type="checkbox"/>	149F305	Foot & Ankle Positioners - Medium		\$250.00
<input type="checkbox"/>	149F306	Foot & Ankle Positioners - Large		\$250.00
<input type="checkbox"/>	149F307	Foot & Ankle Positioners - XL		\$250.00

Accessories

	Item Code	Description		Retail
<input type="checkbox"/>	149LF16	Battery Utility Tray <i>(limit 50lbs, includes 1 utility hook)</i>		\$575.00
<input type="checkbox"/>	149LF17	Ventilator Tray <i>(Vertical mount)</i>		\$520.00
<input type="checkbox"/>	149LF6	Utility Hooks & ZIPPPIE Bag		\$94.00
<input type="checkbox"/>	149HE76	Headrest Cover - Canopy <i>(Includes hardware)</i>		\$250.00
<input type="checkbox"/>	149A468	Underseat Storage Basket <i>(N/A with Battery Tray)</i>		\$95.00
<input type="checkbox"/>	149LF5	Cup Holder		\$79.00
<input type="checkbox"/>	149A466	Curb Tipper		\$100.00
<input type="checkbox"/>	149A467	Rain Cover <i>(Requires Headrest Cover - Canopy)</i>		\$249.00
<input type="checkbox"/>	149A1	ZIPPPIE Backpack		\$108.00
<input type="checkbox"/>	149A469	Tool Kit		\$42.00

Child's Measurements (Not Required)



Child's Name: _____

Age: _____

Diagnosis: _____

Height: _____

Weight: _____

1). Top of Head to Seat: _____

7). Width of Head: _____

2). Top of Shoulder to Seat*: _____

8). Shoulder Width: _____

3). Under Arm to Seat: _____

9). Chest Width: _____

4). Trunk Depth: _____

10). Hip Width: _____

5). Actual Sitting Depth*: _____

11). Left Lower Leg Length*: _____

6). Foot Length: _____

12). Right Lower Leg Length*: _____

*Critical Dimension

