



A truly transportable tilt wheelchair, specifically designed for families on the move.

EIZ22



Customer				Ship To			
Company				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City		Prov	
Buyer				PC		Phone	
Mark For				Email / Fax			

INCLUDED COMPONENTS

- ✓ Silver & Black Folding Frame
- ✓ 10°, 20°, 30° Seat Tilt Settings
- ✓ Quick Adjust Recline Backrest
- ✓ Foot Operated Wheel Lock
- ✓ 12" Rear Wheels & 6" Casters (Solid)
- ✓ Quick Release Rear Wheels
- ✓ Adjustable Stroller Handle
- ✓ Seat Sling with Depth Adjustability

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

BASE

Zippie Sphynx		\$
<input type="checkbox"/> 149WM103	Zippie Sphynx - Includes upholstery, Transit optional	\$ 3,229.00
<input type="checkbox"/> 149WM104	Zippie Sphynx without Upholstery - Transit optional	\$ 3,521.00
<input type="checkbox"/> 149WM111	Zippie Sphynx Package - Includes Transit, Upholstery, Canopy, Underseat Basket, Rain Cover	\$ 4,254.00

FRAME

Transit System \$

Includes: Transit Brackets, and Headrest Extension when ordering Sphynx upholstery.

<input type="checkbox"/> 149TO44	Transit System	NC
<input type="checkbox"/> 149TO40	OMIT Transit - Only available on Transit Optional Sphynx Base	NC

Seat Tilt \$

The adjustable static seat tilt setting can be changed between 10°, 20°, or 30°.

<input type="checkbox"/> 149TA15	10° Seat Tilt	NC
<input type="checkbox"/> 149TA16	20° Seat Tilt	NC
<input type="checkbox"/> 149TA17	30° Seat Tilt	NC

✓ OPTION	DESCRIPTION	RETAIL	✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------	----------	-------------	--------

FRAME

Seat Width \$

Selection Required. Seat Width will have an adjustable Seat Depth range.

<input type="checkbox"/> 149W1	10" Seat Width - Adj Seat Depth of 10" to 15"	NC	<input type="checkbox"/> 149W7	16" Seat Width - Adj Seat Depth of 13" to 18"	NC
<input type="checkbox"/> 149W3	12" Seat Width - Adj Seat Depth of 10" to 15"	NC	<input type="checkbox"/> 149W9	18" Seat Width - Adj Seat Depth of 14" to 19"	NC
<input type="checkbox"/> 149W5	14" Seat Width - Adj Seat Depth of 11" to 16"	NC			

Seat Depth \$

Selection Required. Seat Depth is adjustable within the range stated above for Seat Width.

<input type="checkbox"/> 149SD3	10" Seat Depth - Only ava on 10" 12" SW	NC	<input type="checkbox"/> 149SD8	15" Seat Depth	NC
<input type="checkbox"/> 149SD4	11" Seat Depth - Only ava on 10" 12" 14" SW	NC	<input type="checkbox"/> 149SD9	16" Seat Depth - NA on 10" 12" SW	NC
<input type="checkbox"/> 149SD5	12" Seat Depth - Only ava on 10" 12" 14" SW	NC	<input type="checkbox"/> 149SD10	17" Seat Depth - Only ava on 16" 18" SW	NC
<input type="checkbox"/> 149SD6	13" Seat Depth - NA on 18" SW	NC	<input type="checkbox"/> 149SD11	18" Seat Depth - Only ava on 16" 18" SW	NC
<input type="checkbox"/> 149SD7	14" Seat Depth	NC	<input type="checkbox"/> 149SD12	19" Seat Depth - Only ava on 18" SW	NC

Front Seat Height Guide											
Tilt Angle	Seat Depth										
	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	
10°	18"	18"	18"	18.5"	18.5"	18.5"	19"	19"	19"	19.5"	
20°	17.5"	17.5"	18"	18.5"	18.5"	19"	19.5"	19.5"	20"	20.5"	
30°	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"	

FOOTREST

Seat-To-Footplate Range (Lower Leg Length) \$

Selection Required. Reference table below for achievable Seat-To-Footplate distances relative to the Seat Width and Depth.

<input type="checkbox"/> 149L064	Extra Short Seat-To-Footplate Range	NC
<input type="checkbox"/> 149L061	Short Seat-To-Footplate Range	NC
<input type="checkbox"/> 149L062	Medium Seat-To-Footplate Range	NC
<input type="checkbox"/> 149L063	Long Seat-To-Footplate Range	NC

Achievable Seat-To-Footplate Distances											
Seat Width	Seat-to-Footplate Range Option	Seat Depth									
		10"	11"	12"	13"	14"	15"	16"	17"	18"	19"
10"	Extra Short	6" - 12"	6" - 11"	6" - 10.5"	6" - 10"	5" - 10"	5" - 10"	NA	NA	NA	NA
	Short	12" - 15.5"	11" - 15"	10.5" - 15"	10" - 14.5"	10" - 14.5"	10" - 14.5"	NA	NA	NA	NA
12"	Medium	16" - 17.5"	15" - 17.5"	14.5" - 17"	14" - 17.5"	14" - 17.5"	13.5" - 16.5"	NA	NA	NA	NA
	Long	18" - 19"	16.5" - 19"	16" - 19"	16" - 19"	16" - 19"	15.5" - 19"	NA	NA	NA	NA
14"	Extra Short	NA	6" - 11"	6" - 10.5"	6" - 10"	5" - 10"	5" - 10"	5" - 10"	NA	NA	NA
	Short	NA	11" - 15"	10.5" - 15"	10" - 14.5"	10" - 14.5"	10" - 14.5"	10" - 14.5"	NA	NA	NA
	Medium	NA	15" - 17.5"	14.5" - 17"	14" - 17.5"	14" - 17.5"	13.5" - 16.5"	13" - 16"	NA	NA	NA
	Long	NA	16.5" - 19"	16" - 19"	16" - 19"	16" - 19"	15.5" - 19"	15" - 19"	NA	NA	NA
16"	Extra Short	NA	NA	NA	6" - 10"	5" - 10"	5" - 10"	5" - 10"	6" - 10"	7" - 10"	NA
	Short	NA	NA	NA	10" - 14.5"	10" - 14.5"	10" - 14.5"	10" - 14.5"	10" - 14"	10" - 14"	NA
	Medium	NA	NA	NA	14" - 17.5"	14" - 17.5"	13.5" - 16.5"	13" - 16"	13" - 15.5"	13" - 15.5"	NA
	Long	NA	NA	NA	16" - 19"	16" - 19"	15.5" - 19"	15" - 19"	15" - 19"	14.5" - 19"	NA
18"	Extra Short	NA	NA	NA	NA	10" - 14.5"	10" - 14"	10" - 14"	10" - 14"	10" - 13.5"	10" - 11.5"
	Short	NA	NA	NA	NA	12" - 15"	11.5" - 15"	11.5" - 14.5"	14.5"	11.5" - 14"	11.5" - 13.5"
	Medium	NA	NA	NA	NA	14" - 17"	13.5" - 16.5"	13" - 16"	13" - 15.5"	13" - 15.5"	13" - 15.5"
	Long	NA	NA	NA	NA	15.5" - 19"	15.5" - 19"	15" - 19"	15" - 19"	14.5" - 19"	14.5" - 19"

NOTE: All Seat-To-Footplate distances based on 10° Tilt Setting (Seat-To-Floor Angle).

✓ OPTION	DESCRIPTION	RETAIL	✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------	----------	-------------	--------

UPHOLSTERY

Upholstery for 149WM103 Sphynx Base \$

Selection Required for Upholstery-Included Sphynx Base (149WM103). Includes: Seat and back cushions with tension adjustable straps and seat sling.

<input type="checkbox"/> 149UC37	Sphynx Upholstery - Black	STD	<input type="checkbox"/> 149UC51	Sphynx Upholstery - Green	NC
<input type="checkbox"/> 149UC49	Sphynx Upholstery - Red	NC	<input type="checkbox"/> 149UC52	Sphynx Upholstery - Pink	NC
<input type="checkbox"/> 149UC50	Sphynx Upholstery - Blue	NC	<input type="checkbox"/> 149UC53	OMIT Upholstery	NC

SEATING

JAY Options \$

Required to attach JAY Order Forms for bundled products.

<input type="checkbox"/> 149JC20	JAY Cushion	JAY Part #:		See OF
<input type="checkbox"/> 149JB02	JAY Backrest	JAY Part #:		See OF

POSITIONING

Towel Bar \$

Default Towel Bar is 2" for more compact folded package. Select 4" Towel Bar to accommodate posterior contoured back supports.

<input type="checkbox"/> 149QT02		2" Towel Bar	STD
<input type="checkbox"/> 149QT01		4" Towel Bar	NC

Head Supports \$

<input type="checkbox"/> 149HE75		Headrest Extension - Required for Transit Option & Sphynx Upholstery	\$ 226.00
<input type="checkbox"/> 149ZH25		Padded Headwings - May interfere with Full Torso Support Vest	\$ 256.00

Trunk & Arm Positioning \$

<input type="checkbox"/> 149KI7		H-Harness - Requires 3-Point Pelvic Belt & Headrest Extension. NA with Full Torso Support Vest.	\$ 156.00	
<input type="checkbox"/> 149JY06		JAY Center Opening Dynamic - Requires Headrest Extension. Select harness size below.	\$ 244.00	
Select harness size:	<input type="checkbox"/> 149JZ01 Small	<input type="checkbox"/> 149JZ02 Medium	<input type="checkbox"/> 149JZ03 Large	<input type="checkbox"/> 149JZ04 X-Large
<input type="checkbox"/> 149KI9		Full Torso Support Vest - Requires Headrest Extension. NA with H-Harness or Lateral Supports.	\$ 334.00	
<input type="checkbox"/> 149KI8		Adjustable Lateral Support & Scoli Strap - Requires Armrests. NA with Full Torso Support Vest.	\$ 216.00	
<input type="checkbox"/> 149AA1		Flip-Up Armrests (Angle Adjustable)	\$ 334.00	
<input type="checkbox"/> 149A206		Upper Extremity Support Tray - Grey (Requires Armrests)	\$ 293.00	
<input type="checkbox"/> 149A207		Upper Extremity Support Tray - Clear (Requires Armrests)	\$ 352.00	

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

POSITIONING

Pelvic Positioning \$

<input type="checkbox"/> 149ZZ29	3-Point Pelvic Belt - Requires H-Harness. NA with Medial Thigh Support.	\$ 128.00
<input type="checkbox"/> 149ZZ11	2-Point Positioning Belt with 4.5" Pad	\$ 124.00
<input type="checkbox"/> 149ZZ24	V-Contour 4-Point Belt - Small Pad	\$ 233.00
<input type="checkbox"/> 149ZZ23	V-Contour 4-Point Belt - Medium Pad	\$ 233.00
<input type="checkbox"/> 149KP05	Hip Guides - Inset Mount (Pair) - Reduces seat width by 1" per side	\$ 256.00
<input type="checkbox"/> 149KP06	Hip Guides - In-Line Mount (Pair)	\$ 256.00
<input type="checkbox"/> 149CU30	Cushion 1" Foam Riser	\$ 128.00
<input type="checkbox"/> 149CU12	Wedge Cushion Insert	\$ 154.00
<input type="checkbox"/> 149KH11	Medial Thigh Support - Includes hardware. NA with H-Harness & 3-Point Pelvic Belt.	\$ 186.00

Foot & Leg Positioning \$

<input type="checkbox"/> 149F26	Angle Adjustable Aluminum Footplates (Pair) - Required	\$ 173.00
<input type="checkbox"/> 149HL5	Leg Strap	NC
<input type="checkbox"/> 149A54	Extra Leg Strap	\$ 62.00
<input type="checkbox"/> 149F304	Foot & Ankle Positioners - Small	\$ 256.00
<input type="checkbox"/> 149F305	Foot & Ankle Positioners - Medium	\$ 256.00
<input type="checkbox"/> 149F306	Foot & Ankle Positioners - Large	\$ 256.00
<input type="checkbox"/> 149F307	Foot & Ankle Positioners - X-Large	\$ 256.00

ACCESSORIES \$

<input type="checkbox"/> 149LF16	Battery Utility Tray - Limit 50lbs, includes 1 utility hook	\$ 589.00
<input type="checkbox"/> 149LF17	Ventilator Tray - Vertical mount	\$ 533.00
<input type="checkbox"/> 149LF6	Utility Hooks & Zippie Bag	\$ 96.00
<input type="checkbox"/> 149HE76	Headrest Cover - Canopy - Includes Hardware	\$ 256.00
<input type="checkbox"/> 149A468	Underseat Storage Basket - NA with Battery Tray	\$ 97.00
<input type="checkbox"/> 149LF5	Cup Holder	\$ 81.00
<input type="checkbox"/> 149A466	Curb Tipper	\$ 103.00
<input type="checkbox"/> 149A467	Rain Cover - Requires Headrest Cover - Canopy	\$ 255.00
<input type="checkbox"/> 149A1	Zippie Backpack	\$ 111.00
<input type="checkbox"/> 149A469	Tool Kit	\$ 43.00

MEASUREMENTS - NOT REQUIRED

Name			1. Top of Head to Seat		7. Width of Head	
Age			2. Top of Shoulder to Seat		8. Shoulder Width	
Diagnosis			3. Undarm to Seat		9. Chest Width	
Height			4. Trunk Depth		10. Hip Width	
Weight			5. Actual Sitting Depth		11. Left Lower Leg Length	
			6. Foot Length		12. Right Lower Leg Length	

Please use inches (") for unit of measure