



ORDER FORM

PAL Classroom Seat

The perfect balance of postural support and classroom practicality.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

PAL Classroom Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 77 lbs | Size 3 Weight Capacity: 88 lbs | Size 4 Weight Capacity: 110 lbs

✓ OPTION	DESCRIPTION	RETAIL
Seat		\$
Selection required. Includes Lap Strap.		
<input type="checkbox"/>	135-1600-03 Size 1 PAL Classroom Seat - Blue	\$1,059
<input type="checkbox"/>	135-1600-01 Size 1 PAL Classroom Seat - Green	\$1,059
<input type="checkbox"/>	135-2600-03 Size 2 PAL Classroom Seat - Blue	\$1,193
<input type="checkbox"/>	135-2600-01 Size 2 PAL Classroom Seat - Green	\$1,193
<input type="checkbox"/>	135-3600-03 Size 3 PAL Classroom Seat - Blue	\$1,328
<input type="checkbox"/>	135-3600-01 Size 3 PAL Classroom Seat - Green	\$1,328
<input type="checkbox"/>	135-4600-03 Size 4 PAL Classroom Seat - Blue	\$1,590
<input type="checkbox"/>	135-4600-01 Size 4 PAL Classroom Seat - Green	\$1,590
Side Pad		\$
Must match seat size.		
<input type="checkbox"/>	135-1635-03 Size 1 PAL Side Pads 1" Thick - Blue (Pair)	\$137
<input type="checkbox"/>	135-1635-01 Size 1 PAL Side Pads 1" Thick - Green (Pair)	\$137
<input type="checkbox"/>	135-2635-03 Size 2 PAL Side Pads 1" Thick - Blue (Pair)	\$137
<input type="checkbox"/>	135-2635-01 Size 2 PAL Side Pads 1" Thick - Green (Pair)	\$137
<input type="checkbox"/>	135-3635-03 Size 3 and 4 PAL Side Pads 1" Thick - Blue (Pair)	\$137
<input type="checkbox"/>	135-3635-01 Size 3 and 4 PAL Side Pads 1" Thick - Green (Pair)	\$137
<input type="checkbox"/>	135-1654-03 Size 1 PAL Side Pads 2" Thick - Blue (Pair)	\$137
<input type="checkbox"/>	135-1654-01 Size 1 PAL Side Pads 2" Thick - Green (Pair)	\$137
<input type="checkbox"/>	135-2654-03 Size 2 PAL Side Pads 2" Thick - Blue (Pair)	\$137
<input type="checkbox"/>	135-2654-01 Size 2 PAL Side Pads 2" Thick - Green (Pair)	\$137
<input type="checkbox"/>	135-3654-03 Size 3 and 4 PAL Side Pads 2" Thick - Blue (Pair)	\$137
<input type="checkbox"/>	135-3654-01 Size 3 and 4 PAL Side Pads 2" Thick - Green (Pair)	\$137
Ramped Seat Cushion		\$
Must match seat size. Additional option.		
<input type="checkbox"/>	135-1651-03 Size 1 PAL Ramped Seat Cushion - Blue	\$92
<input type="checkbox"/>	135-1651-01 Size 1 PAL Ramped Seat Cushion - Green	\$92
<input type="checkbox"/>	135-2651-03 Size 2 PAL Ramped Seat Cushion - Blue	\$92
<input type="checkbox"/>	135-2651-01 Size 2 PAL Ramped Seat Cushion - Green	\$92
<input type="checkbox"/>	135-3651-03 Size 3 PAL Ramped Seat Cushion - Blue	\$92
<input type="checkbox"/>	135-3651-01 Size 3 PAL Ramped Seat Cushion - Green	\$92
<input type="checkbox"/>	135-4651-03 Size 4 PAL Ramped Seat Cushion - Blue	\$92
<input type="checkbox"/>	135-4651-01 Size 4 PAL Ramped Seat Cushion - Green	\$92
Tray and Tray Pad		\$
Must match seat size.		
<input type="checkbox"/>	135-1620 Size 1 PAL Tray	\$249
<input type="checkbox"/>	135-2620 Size 2 PAL Tray	\$249
<input type="checkbox"/>	135-3620 Size 3 PAL Tray	\$249
<input type="checkbox"/>	135-4620 Size 4 PAL Tray	\$249
<input type="checkbox"/>	135-1656 Size 1 PAL Tray Pad	\$114
<input type="checkbox"/>	135-2656 Size 2 PAL Tray Pad	\$114
<input type="checkbox"/>	135-3656 Size 3 PAL Tray Pad	\$114
<input type="checkbox"/>	135-4656 Size 4 PAL Tray Pad	\$114

PAL Classroom Seat

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

Pommel		\$
--------	--	----

Size 1 Pommel suitable for Size 1 seat, Size 2 Pommel for Size 2, 3 and 4 seats.

<input type="checkbox"/>	135-1602-03	Size 1 PAL Pommel - Blue	\$182
<input type="checkbox"/>	135-1602-01	Size 1 PAL Pommel - Green	\$182
<input type="checkbox"/>	135-3602-03	Size 2 PAL Pommel - Blue	\$182
<input type="checkbox"/>	135-3602-01	Size 2 PAL Pommel - Green	\$182

Footrest		\$
----------	--	----

Must match seat size

<input type="checkbox"/>	135-1603	Size 1 PAL Footrest	\$139
<input type="checkbox"/>	135-2603	Size 2 PAL Footrest	\$139
<input type="checkbox"/>	135-3603	Size 3 PAL Footrest	\$139
<input type="checkbox"/>	135-4603	Size 4 PAL Footrest	\$139

Rocker / Stabilizer		\$
---------------------	--	----

Not available with PAL Ski with Casters

<input type="checkbox"/>	135-1604	Size 1 PAL Rocker/Stabilizer (Pair)	\$157
<input type="checkbox"/>	135-2604	Size 2 PAL Rocker/Stabilizer (Pair)	\$157
<input type="checkbox"/>	135-3604	Size 3 PAL Rocker/Stabilizer (Pair)	\$157
<input type="checkbox"/>	135-4604	Size 4 PAL Rocker/Stabilizer (Pair)	\$157

Caster		\$
--------	--	----

Not available with Rocker/Stabilizer

<input type="checkbox"/>	135-1610	Size 1 PAL Ski with Casters (Pair)	\$273
<input type="checkbox"/>	135-2610	Size 2 PAL Ski with Casters (Pair)	\$273
<input type="checkbox"/>	135-3610	Size 3 PAL Ski with Casters (Pair)	\$273
<input type="checkbox"/>	135-4610	Size 4 PAL Ski with Casters (Pair)	\$273



Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	27.6" - 37.4"	35.4" - 41.3"	39.4" - 51.2"	47" - 63"
Max User Weight	66 lbs	77 lbs	88 lbs	110 lbs
Seat to Footrest	5.5" - 9.4"	6.5" - 11.2"	9.1" - 13.8"	11.2" - 16"
Seat Depth	7.1" - 9.4"	8.7" - 11.4"	9.8" - 13"	12" - 16"
Backrest Height	10.4"	12.4"	15"	17.7"
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	5.9" - 9.8"	6.7" - 10.6"	9.1" - 13"	11.1" - 15.2"
Seat Weight	9.6 lbs	13.2 lbs	17.9 lbs	22 lbs
Tray Size	20"W x 18"D	21.3"W x 19.7"D	23.6"W x 21.7"D	28.7"W x 23.6"D

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	700 - 950 mm	900 - 1050 mm	1000 - 1300 mm	1200 - 1600 mm
Max User Weight	30 kg	35 kgs	40 kgs	50 kgs
Seat to Footrest	40 - 240 mm	165 - 285 mm	230 - 350 mm	280 - 400 mm
Seat Depth	80 - 240 mm	220 - 290 mm	250 - 330 mm	300 - 400 mm
Backrest Height	265 mm	315mm	380 mm	450 mm
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	50 - 250 mm	170 - 270 mm	230 - 330 mm	280 - 380 mm
Seat Weight	4.35kg	6 kgs	8.15 kgs	10 kgs
Tray Size	508mm W x 457mm D	541mm W x 500mm D	599mm W x 551mm D	729mmW x 599mmD



Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada

Customer Service: (800) 263-3390 E-mail: cscanada@sunmed.com www.sunrisemedical.ca

Specifications subject to change without notice. © 2024 Sunrise Medical Canada Inc. MK-129852-CAN Rev 7

