

Spex SuperHigh Contour Cushion

February 2024



MARK FOR: _____

Date: _____

Dealer Account: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ PV: _____ PC: _____

Dealer Phone: _____ Fax: _____

Confirm via: Fax E-mail

SUBMITTING FOR: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____


Address: _____

Ship to City: _____ PV: _____ PC: _____

Ship to Phone: _____ Fax: _____

PACKAGE INCLUDES:

1. Outer Cover
2. Incontinence Cover
3. Cushion Foam Unit
4. Positioning Pad Kit



1. SELECT SIZE

✓	WIDTH	LENGTH	ITEM CODE	\$
		10"	1105-1010-OSS	\$1,628.00
	10"	12"	1105-1012-OSS	\$1,628.00
		XX	1105-10XX-OSS	\$1,628.00
		10"	1105-1110-OSS	\$1,628.00
	11"	12"	1105-1112-OSS	\$1,628.00
		XX	1105-11XX-OSS	\$1,628.00
		12"	1105-1212-OSS	\$1,628.00
	12"	14"	1105-1214-OSS	\$1,628.00
		XX	1105-12XX-OSS	\$1,628.00
		12"	1105-1312-OSS	\$1,628.00
	13"	14"	1105-1314-OSS	\$1,628.00
		XX	1105-13XX-OSS	\$1,628.00
		14"	1105-1414-OSS	\$1,628.00
		16"	1105-1416-OSS	\$1,628.00
	14"	18"	1105-1418-OSS	\$1,628.00
		XX	1105-14XX-OSS	\$1,628.00
		14"	1105-1514-OSS	\$1,628.00
	15"	16"	1105-1516-OSS	\$1,680.00
		XX	1105-15XX-OSS	\$1,680.00
		16"	1105-1616-OSS	\$1,680.00
		18"	1105-1618-OSS	\$1,680.00
	16"	20"	1105-1620-OSS	\$1,680.00
		XX	1105-16XX-OSS	\$1,680.00
		17"	1105-1717-OSS	\$1,680.00
		18"	1105-1718-OSS	\$1,680.00
	17"	XX	1105-17XX-OSS	\$1,680.00
		18"	1105-1818-OSS	\$1,910.00
	18"	20"	1105-1820-OSS	\$1,910.00
		XX	1105-18XX-OSS	\$1,910.00
		18"	1105-1918-OSS	\$1,910.00
	19"	20"	1105-1920-OSS	\$1,910.00
		XX	1105-19XX-OSS	\$1,910.00
		20"	1105-2020-OSS	\$1,910.00
	20"	22"	1105-2022-OSS	\$1,910.00
		XX	1105-20XX-OSS	\$1,910.00
	Up to 20"	XX	1105-2929-OSS	By Quote Only
	Up to 25"	XX	1105-2939-OSS	By Quote Only

2. SELECT COLOUR NC

Your colour selection applies to any extra covers ordered


<input type="radio"/> Black (default)	<input type="radio"/> Granite 1111-0491-090	<input type="radio"/> Calypso 1111-0491-028
<input type="radio"/> Chilli 1111-0491-022	<input type="radio"/> Mango 1111-0491-030	<input type="radio"/> Ocean 1111-0491-060
<input type="radio"/> Leaf 1209-2491-054	<input type="radio"/> Bumblebee 1209-2491-044	<input type="radio"/> Royal 1209-2491-077

3. CONFIGURE (tick any configurations that apply)


✓ ITEM	IMAGE	NOTES	\$
<input type="radio"/> Leg Length Discrepancy 1111-0292-000		T = ____mm Cutout on Left /Right (circle one)	\$331.00
<input type="radio"/> Front Edge Cut-Back for Leg Contractures 1111-0317-000		U = ____mm	\$331.00
<input type="radio"/> Remove Medial Thigh Support 1111-0335-000		NC	
<input type="radio"/> Specify Medial Thigh Support Dimensions 1111-0294-000		K = ____mm L = ____mm M = ____mm N = ____mm	\$568.00
<input type="radio"/> Specify Lateral Thigh Support Dimensions 1111-0295-000		O = ____mm P = ____mm Q = ____mm	\$568.00
<input type="radio"/> Increase Rear Half Height to Full Lateral Thigh Support 1111-0304-000			\$568.00
<input type="radio"/> Additional Pressure Relief, Extra SunMate Overlay 1111-0391-000			\$568.00

4. SELECT ACCESSORIES Input your cushion dimensions

✓ ACCESSORY	IMAGE	\$
<input type="radio"/> Extra Outer Cover (A) 1125-____-6SS		\$201.00
<input type="radio"/> Extra Incontinence Cover (B) 1135-____-7SS		\$275.00
<input type="radio"/> Rigidizer Base Panel (can be inserted inside cover for a rigid surface) 1110-XXXX-200		\$338.00
<input type="radio"/> Extra Positioning Kit 1111-0273-000		\$147.00
<input type="radio"/> 5° Wedge 1209-2724-300		\$147.00



LIFE-SIZE LATERAL TEMPLATE



BROAD 6
230 x 175mm

NARROW 4
115 x 240mm

BROAD 5
195 x 165mm

NARROW 3
115 x 210mm

BROAD 4
160 x 155mm

NARROW 2
95 x 180mm

BROAD 3
140 x 150mm

NARROW 1
95 x 145mm

BROAD 2
125 x 140mm

NARROW 0
85 x 115mm

BROAD 1
105 x 130mm

Spex® provides the tools to make decisions for faster results, including this lateral pad size guide to scale. When deciding which lateral support pad style and size will suit your clients, place this page against their torso to make an informed decision on the spot!