

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ PV: _____ PC: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation E-mail: _____
 Confirm Via: Fax E-mail

Submitting for:

Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ PV: _____ PC: _____
 Ship To Phone: () _____ Fax: () _____

✓ Option # Definition Price ✓ Option # Definition Price

1. FRAME

Base Model			
<input type="checkbox"/>	196WM94	Quickie QX (Fixed Axle) 300 lbs.	\$1,299.00

Base Model			
<input type="checkbox"/>	196TO38	Transit System WC19	NC

Frame Style			
<i>See Seat Height Matrix on page 2</i>			
<input type="checkbox"/>	196FS9	Standard Frame	NC
<input type="checkbox"/>	196FS30	Hemi Frame	NC

2. SIZE / SEATING

Frame Width			
<input type="checkbox"/>	196W5	14"	NC
<input type="checkbox"/>	196W6	15"	NC
<input type="checkbox"/>	196W7	16"	NC
<input type="checkbox"/>	196W8	17"	NC
<input type="checkbox"/>	196W9	18"	NC
<input type="checkbox"/>	196W10	19"	NC
<input type="checkbox"/>	196W11	20"	NC

Seat Depth			
<i>Frame Depths equals Seat Depths.</i>			
<input type="checkbox"/>	196SD7	14"	NC
<input type="checkbox"/>	196SD8	15"	NC
<input type="checkbox"/>	196SD9	16"	NC
<input type="checkbox"/>	196SD10	17"	NC
<input type="checkbox"/>	196SD11	18"	NC
<input type="checkbox"/>	196SD71	19"	NC
<input type="checkbox"/>	196SD13	20"	NC

Frame Width - Built 4 Me Options **NEW!**

Not available with EXO standard or Tension adjustable upholstery, Not available with Ballistic Nylon or 3DX Tension adjustable upholstery.

<input type="checkbox"/>	196W59	B4ME - 21" Width	\$284.00
<input type="checkbox"/>	196W60	B4ME - 22" Width	\$284.00

Seat Sling			
<input type="checkbox"/>	196SL1	Standard Seat Sling	NC
<input type="checkbox"/>	196SL6	Omit Seat Sling	NC



Quickie QX, Ultra Lightweight Folding Wheelchair



Quickie QX featuring the Jay X2™ cushion

A1 is client responsibility. A2 is clinical justification.

3. CASTERS / SEAT HEIGHT

Front Seat Height

See Front Seat Height Matrix. 12.5 & 13.0 only available w/ 60° hanger
12.5" & 13" only available w/ Alum Angle Adj Footplate

<input type="checkbox"/>	196SH34	12.5"	NC	<input type="checkbox"/>	196SH14	16.0"	NC
<input type="checkbox"/>	196SH35	13.0"	NC	<input type="checkbox"/>	196SH15	16.5"	NC
<input type="checkbox"/>	196SH40	13.5"	NC	<input type="checkbox"/>	196SH16	17.0"	NC
<input type="checkbox"/>	196SH41	14.0"	NC	<input type="checkbox"/>	196SH17	17.5"	NC
<input type="checkbox"/>	196SH42	14.5"	NC	<input type="checkbox"/>	196SH18	18.0"	NC
<input type="checkbox"/>	196SH12	15.0"	NC	<input type="checkbox"/>	196SH20	18.5"	NC
<input type="checkbox"/>	196SH13	15.5"	NC	<input type="checkbox"/>	196SH21	19.0"	NC
				<input type="checkbox"/>	196SH23	19.5"	NC

Caster Wheels

<input type="checkbox"/>	196CW3	4" Low Profile Polyurethane	NC
<input type="checkbox"/>	196CW44	4" x 1.25" Semi Pneumatic	NC
<input type="checkbox"/>	196CW5	5" Low Profile Polyurethane	NC
<input type="checkbox"/>	196CW7	5" x 1.5" Semi Pneumatic	\$95.00
<input type="checkbox"/>	196CW34	6" Polyurethane	NC
<input type="checkbox"/>	196CW35	6" x 1.5" Semi Pneumatic	\$95.00
<input type="checkbox"/>	196CW32	8" Polyurethane	NC
<input type="checkbox"/>	196CW39	8" x 1.5" Semi Pneumatic	\$134.00

Caster Fork

<input type="checkbox"/>	196CF71	3.5" (standard stem bolt only)	NC
<input type="checkbox"/>	196CF12	5.25"	NC
<input type="checkbox"/>	196CF17	6"	NC
<input type="checkbox"/>	196CF23	7"	NC

Caster Stem Bolt

<input type="checkbox"/>	196CS1	Standard	NC
<input type="checkbox"/>	196CS5	3/4"	NC
<input type="checkbox"/>	196CS9	1 1/2"	NC

Front Seat Height Matrix

Front Seat Height Matrix - Standard Stem Bolt								
Caster	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	14.5	15.0	12.5	13.0	16.0	17.0	14.0	15.0
4 x 1.25 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 Low Poly	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 x 1.5 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
6 Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
8 Poly	N/A	N/A	N/A	N/A	19.0	19.0	17.0	17.0
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.0	16.5	17.0
Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
4 x 1.25 Semi Pneu	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
5 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
5 x 1.5 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
6 Poly	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
6 x 1.5 Semi Pneu	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
8 Poly	19.5	19.5	17.5	18.0	19.5	19.5	17.5	19.0
8 x 1.5 Semi Pneu	18.5	19.5	16.5	17.5	19.0	19.5	17.0	18.5
Seat squeeze and seat depth can affect FSH availability For 3/4" and 1 1/2" Stem bolt matrix see page 5								

4. REAR WHEELS / REAR SEAT HEIGHTS

Rear Seat Height

See Rear Seat Height Matrix. 12.5" & 18.5" RSH not available with fixed axle.

<input type="checkbox"/>	196RH9	12"	NC
<input type="checkbox"/>	196RH11	13"	NC
<input type="checkbox"/>	196RH12	13.5"	NC
<input type="checkbox"/>	196RH13	14"	NC
<input type="checkbox"/>	196RH14	14.5"	NC
<input type="checkbox"/>	196RH15	15"	NC
<input type="checkbox"/>	196RH16	15.5"	NC
<input type="checkbox"/>	196RH17	16"	NC
<input type="checkbox"/>	196RH18	16.5"	NC
<input type="checkbox"/>	196RH19	17"	NC
<input type="checkbox"/>	196RH20	17.5"	NC
<input type="checkbox"/>	196RH21	18"	NC
<input type="checkbox"/>	196RH21	19"	NC

Rear Wheel Size

See Rear Tire & Wheel Matrix for availability

<input type="checkbox"/>	196RW3	20" Mag Wheel	NC
<input type="checkbox"/>	196RW8	22" Mag Wheel	NC
<input type="checkbox"/>	196RW1	24" Mag Wheel	NC
<input type="checkbox"/>	196RW103	20" Mag Wheel 5 Spoke (Xcore)	NC
<input type="checkbox"/>	196RW104	22" Mag Wheel 5 Spoke (Xcore)	NC
<input type="checkbox"/>	196RW102	24" Mag Wheel 5 Spoke (Xcore)	NC
<input type="checkbox"/>	196RW25	Omit	(\$44.00)

Rear Wheel Tire

See Rear Tire & Wheel Matrix for availability

<input type="checkbox"/>	196RT1		Pneumatic	NC
<input type="checkbox"/>	196RT2	A1	Pneumatic w/ Airless Insert	\$228.00
<input type="checkbox"/>	196RT3		Full Poly	NC
<input type="checkbox"/>	196RT39	A1	High-Pressure Clincher Black	\$96.00

Handrim

Natural Fit only available w/ 22, 24, 25 Lite Spoke. Plastic N/A w/ 25, 26"

<input type="checkbox"/>	196HR1	Aluminum Anodized	NC
<input type="checkbox"/>	196HR4	Plastic	\$134.00

Seat Squeeze Matrix

Frame Depth	Max Seat Squeeze
14"	2.0"
15" - 16"	2.5"
17" - 18"	2.5"
19" - 20"	3.0"

Rear Seat Height Matrix

Rear Wheel Size	Available Rear Seat Heights					
	Pneumatic Tire		Full Poly Tire		HP Clincher Tire	
	Min Squeeze	Max Squeeze	Min Squeeze	Max Squeeze	Min Squeeze	Max Squeeze
20"	12" - 16"	12" - 16"	12" - 16"	12" - 16"	N/A	N/A
22"	13.5" - 17.5"	13" - 17"	13" - 17"	13" - 17"	13" - 17"	13" - 17"
24"	14" - 18"	13.5" - 16.5"	13.5" - 17.5"	13.5" - 16.5"	13.5" - 17.5"	13.5" - 16.5"
26"	15" - 19"	15" - 17"	N/A	N/A	15" - 19"	14.5" - 16.5"
12.5" and 18.5" RSH not available with QX						

A1 is client responsibility. A2 is clinical justification.

4. REAR WHEELS / REAR SEAT HEIGHTS (CONT)

Rear Tire & Wheel Matrix

Tire	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Pneumatic	20,22,24	20,22,24,26	20,22,24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	20,22,24
Full Poly	20,22,24	20,22,24	20,22,24
High-Pressure Clincher Blk	N/A	22,24,26	N/A
20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm			

One Arm Drive

Available only with 24" MBL Lite spoke wheels. Std Aluminum black handrim. 1AD option is only available on 14" - 19" frame widths.

<input type="checkbox"/>	196AL16	A2	One Arm Drive Adj. Axle Plate	\$1,276.00
	196HM1	A2	<input type="checkbox"/> Right Hand Mount	NC
	196HM2	A2	<input type="checkbox"/> Left Hand Mount	NC

Axle

<input type="checkbox"/>	196AX1		Quick Release Steel	NC
<input type="checkbox"/>	196AX11	A2	Quad Quick Release Steel N/A on Spoke	\$168.00
<input type="checkbox"/>	196AX45	A2	Threaded Axle	\$138.00

5. WHEEL LOCKS

Wheel Lock

<input type="checkbox"/>	196WL1		High Mount Push	NC
<input type="checkbox"/>	196WL2		High Mount Pull	NC
<input type="checkbox"/>	196WL30	A1	Compact	\$150.00
<input type="checkbox"/>	196WL24	A1	Unilateral Wheel Lock Right	\$295.00
<input type="checkbox"/>	196WL23	A1	Unilateral Wheel Lock Left	\$295.00
<input type="checkbox"/>	196WL10		Omit Wheel Lock	NC
<input type="checkbox"/>	196WL12		Do Not Mount	NC

Extension Handle

Ext. handles only available with high mount push or pull wheel locks

<input type="checkbox"/>	196EH7	E0961	Extension Handle (Pair)	\$63.00
--------------------------	--------	-------	-------------------------	---------

6. HANGERS AND FOOTRESTS

Front Hangers

Footplate mounts with Extension Tube using clamping hardware

<input type="checkbox"/>	196H145		60° Swing In/Out Hanger	NC
<input type="checkbox"/>	196H139		70° Swing In/Out Hanger	NC
<input type="checkbox"/>	196H141		80° Swing In/Out Hanger	NC
<input type="checkbox"/>	196H159	A2	Articulating Hanger N/A w/ Alum Angle Adj Front Mount footplates.	\$485.00
<input type="checkbox"/>	196H160	A2	Elevating Hanger - Ext Mount N/A w/ Alum Angle Adj Front Mount footplates.	\$313.00
<input type="checkbox"/>	196H161	A2	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$313.00
<input type="checkbox"/>	196H144		Omit Hanger & Footplate	NC

Extension Tubes

See Lower Leg Length Matrix for range

<input type="checkbox"/>	196E1		Short	NC
<input type="checkbox"/>	196E2		Medium	NC
<input type="checkbox"/>	196E5		Long	NC

Footplate

Heel loop not included, see Footplate Option section below

<input type="checkbox"/>	196F1		Composite Footplate	NC
<input type="checkbox"/>	196F271		Composite Angle Adj-(16"-22"w)	\$162.00
<input type="checkbox"/>	196F3	A2	Aluminum Angle Adj. Footplate	\$152.00
<input type="checkbox"/>	196F5		Aluminum Only available w/Elevating Legrest	NC

Hanger & Footplate Options

<input type="checkbox"/>	196HL12		Adjustable Heel Loop	NC
<input type="checkbox"/>	1968A02	A1	Gel Cover For Hangers (Pair)	\$150.00

7. BACKREST / SEATING

Back Type

<input type="checkbox"/>	196BT10		Standard	NC
<input type="checkbox"/>	196BT100		Angle Adjustable Back (-15° to +3°)	\$319.00

Backrest Height

<input type="checkbox"/>	196B5		Med (13"-17")	NC
<input type="checkbox"/>	196B9		Tall (16"-20")	NC

Push Handles

Integral Push Handles, Folding Height Adjustable & Omit 8° for Std Backrest, Straight for Angle Adj. Backrest

<input type="checkbox"/>	196PH6		Integral Push Handles	NC
<input type="checkbox"/>	196PH5	A1	Fold-Down Push Handles Only available with Adjustable Upholstery. Used in combination w/Angle Adjust only avail in 16"-20" back height	\$156.00
<input type="checkbox"/>	196PH17	A1	Folding Push Handle with Ht. Adj. Back	\$286.00
<input type="checkbox"/>	196PH2		Omit Push Handles	NC

A1 is client responsibility. A2 is clinical justification.

Customer Service: (800) 263-3390 Fax: (800) 561-5834 E-mail: cscanada@sunmed.com www.sunrisemedical.ca
Specifications are subject to change without notice. © 2019 Sunrise Medical Canada Inc. MK-100379-CAN REV (2).

7. BACKREST / SEATING (CONT)

Back Upholstery			
<input type="checkbox"/>	196BU1	Standard Upholstery	NC
<input type="checkbox"/>	196BU70	Backrest Tension Adj Ballistic Nylon	\$175.00
<input type="checkbox"/>	196BU14	Omit Back Upholstery	(\$21.00)

8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away			
<i>Some limitations when used with the Angle Adjustable Back in the rear most setting</i>			
<i>Please circle a setting, will default to 11" if no selection is made</i>			
<input type="checkbox"/>	196AA18	Padded Swing-Away STD (9" to 13")	\$256.00
	Setting	<input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13"	

Armrest - Height Adjustable			
<i>Flip Back HA Classic N/A w / 14"-15" depth</i>			
<input type="checkbox"/>	196AA104	Single Post Ht. Adj - Std (9.75"-14.75")	\$313.00
<input type="checkbox"/>	196AA105	Single Post Ht Adj - Low (8.25"-12.25")	\$313.00
<input type="checkbox"/>	196AA65	Flip Back Ht Adj LITE (10"-15")	\$331.00
<input type="checkbox"/>	196AA103	Dual Post Flip Back Ht Adj	\$313.00
<input type="checkbox"/>	196AA159	Flip Back Ht Adj Classic (11"-15")	\$351.00
<input type="checkbox"/>	196AA14	Omit Armrest	NC

Armpad - For Height Adjustable			
<input type="checkbox"/>	196AM23	Desk Length (10")	NC
<input type="checkbox"/>	196AM24	Full Length - Full (14")	NC

Armrest Options			
<input type="checkbox"/>	1968A01	Gel Cover For Swing Away Arms (pair)	\$150.00

Side Guards			
<input type="checkbox"/>	196SG33	Plastic - Standard (8.5" High)	\$173.00
<input type="checkbox"/>	196SG29	Plastic - Low (6.5" High)	\$173.00

9. KEY OPTIONS

Anti Tips			
<input type="checkbox"/>	196AT1	Anti Tips (Pair)	\$105.00
<input type="checkbox"/>	196A86	Anti Tips (Pair) Low Profile	\$111.00

Pack Accessories			
<input type="checkbox"/>	196A1	A1 Backpack - Black/Charcoal Adult	\$103.00
<input type="checkbox"/>	1961A01	A1 Seat Pouch Black	\$77.00

Positioning Belt			
<input type="checkbox"/>	196A141	1.5" Auto Buckle w/ Clamps	\$90.00
<input type="checkbox"/>	196A804	Positioning Belt Slide Clamps	\$71.00

Spoke Guards			
<input type="checkbox"/>	1964A01	Quickie Black Spoke Guard	\$145.00
<input type="checkbox"/>	1964A02	Quickie Clear Spoke Guard	\$145.00

Miscellaneous			
<input type="checkbox"/>	1968A03	A1 Caddy	\$83.00
<input type="checkbox"/>	196A2100	A1 Crutch Holder	\$106.00

10. COLOR

Matte Frame Colors			
<input type="checkbox"/>	196C116	Matte Black	NC
<input type="checkbox"/>	196C150	Matte Electric Blue	NC

Additional Options Through Built for Me			
<i>To inquire about any unique customization for Quickie products, contact: Phone: (800) 263-3390, Fax: (800) 561 5834</i>			



A1 is client responsibility. A2 is clinical justification.

Reference page 12 of the AADL Wheelchair APL for options and catalogue numbers for A wheelchairs.

Customer Service: (800) 263-3390 Fax: (800) 561-5834 E-mail: cscanada@sunmed.com www.sunrisemedical.ca
Specifications are subject to change without notice. © 2019 Sunrise Medical Canada Inc. MK-100379-CAN REV (2).

Front Seat Height Matrix Cont.

Front Seat Height Matrix - 3/4" Stem								
Caster	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	14.0	15.0
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	14.5	15.5
5 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5
6 Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	17.0	17.0
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.5	19.5	16.5	17.0
Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
4 x 1.25 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
5 Low Poly	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0
5 x 1.5 Semi Pneu	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0
6 Poly	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5
6 x 1.5 Semi Pneu	18.0	19.5	16.5	17.5	19.0	19.5	17.0	18.5
8 Poly	N/A	N/A	18.0	18.5	N/A	N/A	18.0	19.5
8 x 1.5 Semi Pneu	19.5	19.5	17.5	18.5	N/A	N/A	18.0	19.5

Seat squeeze and seat depth can affect FSH availability

Front Seat Height Matrix - 1 1/2" Stem								
Caster	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5
5 Low Poly	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0
6 Poly	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5
4 x 1.25 Semi Pneu	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5
5 Low Poly	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0
5 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0
6 Poly	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
6 x 1.5 Semi Pneu	19.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
8 Poly	N/A	N/A	19.0	19.5	N/A	N/A	19.0	19.5
8 x 1.5 Semi Pneu	N/A	N/A	18.0	19.0	N/A	N/A	18.5	19.5

Seat squeeze and seat depth can affect FSH availability

Lower Leg Length Range Matrix

Extension Mount Hanger Leg Lengths					
Hanger	Ext Tube	Seat-to-Footrest Range			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	10.25" - 14.75"	10.25" - 14.75"	9.25" - 13.25"	N/A
	Medium	11.00" - 16.00"	11.00" - 16.00"	9.25" - 14.75"	N/A
	Long	12.50" - 16.00"	12.50" - 16.00"	10.75" - 16.25"	N/A
70° Swing In/Out	Short	11.00" - 15.75"	11.00" - 15.75"	10.00" - 14.50"	N/A
	Medium	12.00" - 16.50"	12.00" - 16.50"	10.75" - 16.50"	N/A
	Long	13.75" - 16.50"	13.75" - 16.50"	12.25" - 16.50"	N/A
80° Swing In/Out	Short	10.50" - 15.25"	10.50" - 15.25"	10.50" - 15.25"	N/A
	Medium	11.75" - 16.50"	11.75" - 16.50"	11.75" - 16.50"	N/A
	Long	13.50" - 16.50"	13.50" - 16.50"	11.75" - 16.50"	N/A

All values are ±0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights (see "Min Allowable Front Seat Heights" matrix for minimum FSH information).



Min Allowable Front Seat Heights (for Extension Mount Hangers)

Hanger	Ext Tube	Min Allowable FSH			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	14.0	13.5	12.5	N/A
	Medium	14.5	14.5	12.5	N/A
	Long	16.0	16.0	13.5	N/A
70° Swing In/Out	Short	14.0	14.0	13.0	N/A
	Medium	15.0	15.0	13.5	N/A
	Long	16.5	16.5	15.0	N/A
80° Swing In/Out	Short	14.0	14.0	13.5	N/A
	Medium	14.5	14.5	14.5	N/A
	Long	16.5	16.5	15.5	N/A

ALR & ELR Hanger Leg Lengths

Hanger	Extension	Seat Sling to Footrest		Min Allowable FSH
		Min	Max	
ALR	Short	9 1/2	15 1/2	11 1/2
	Med	11	17 1/2	13
	Long	13	19 1/2	15
ELR Ext. Mount	Short	9	15	11
	Med	9	17	11
	Long	10	19	12
ELR Front Mount	No Ext (top hole)	4 1/2	8 1/2	6 1/2
	Short	8	11	10
	Med	10	13	12
	Long	12	15	14

Rear Wheel & Handrim Matrix

Handrim	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Aluminum Anodized	20,22,24	20,22,24,26	20,22,24
Plastic Coat	20,22,24	22,24	20,22,24
Natural Fit - Standard Grip	24	22,24,26	22,24
Natural Fit - Super Grip	24	22,24,26	22,24



A1 is client responsibility. A2 is clinical justification.

Customer Service: (800) 263-3390 Fax: (800) 561-5834 E-mail: cscanada@sunmed.com www.sunrisemedical.ca
 Specifications are subject to change without notice. © 2019 Sunrise Medical Canada Inc. MK-100379-CAN REV (2).