



# TIGAFX

» The rigid chair that folds



# TIGA fx - For people on the move!

A folding chair that is lightweight and rigid.

Airline travel? Car travel? No problem!

The TIGA fx folds into a briefcase style shape, ideal for storage in tight spaces. It fits\* into aircraft overhead storage lockers and the boot of your car!

The Q-Lock system ensures all the benefits of a folding chair without compromising strength. The system is simple, clean and streamlined enough to meet all your travel needs. It's created individually for you, fully customised to your measurements, needs and personal preferences.

## Key Features


- » Aluminium 7020 customised frame
- » Q Lock folding frame and backrest
- » Integrated fixed COG or adjustable COG available
- » Ergonomic seat and back optional
- » Fold flat sideguards
- » Integrated footrest with carbon fibre cover option
- » Wide accessories range
- » Choice of anodised accent colours
- » Airtech upholstery system
- » Max user weight 125 kg




\*Fits most commercial aircraft and standard sized cars



Q-lock system for easy folding of front frame and backrest



"Fold flat" folding sideguards



Ergonomic seat and back optimises user posture to increase stability, facilitate propulsion, and help skin protection.



## FrontWheel

Enables you to travel easily across a wide range of terrains by keeping your front castor wheels off the floor at all times. Attaches to TIGA, TIGA fx and Hi Lite.