



MANUAL WHEELCHAIR

REFERENCE GUIDE

A COMPREHENSIVE
GUIDE TO SELECTING
THE RIGHT QUICKIE
MANUAL WHEELCHAIR



LIVE
WITH
OUT
LIMITS



		RIGID					TILT-IN-SPACE			FOLDING						
		K0005					E1161			K0005				K0007		
		Nitrum®	Nitrum® Hybrid	5"	QRi	GP / GPV	GPVA	IRIS / IRIS SE	SR45	ACCESS	QS5 X SA	QS5 X FF	QUICKIE 2	QXi	M6	
Specifications subject to change without notice.																
Design	Frame Style	Rigid - Open	Rigid - Dual Tube	Rigid - Open	Rigid - Open	Rigid - Traditional	Rigid - Swing-Away	Rotation in Space	Rotation in Space	Posterior Tilt	Folding	Folding	Folding	Folding	Folding	
	Standard Weight(2)	18.3 lbs (Std. Config.) 10.8 lbs (Transport)	20.3 lbs (Std. Config.) 12.8 lbs (Transport)	21.5 lbs (Std. Config.) 12.9 lbs (Transport)	22 lbs (Std. Config.) 13.4 lbs (Transport)	GP: 23.75 lbs GPV: 23.6 lbs	26 lbs	38 lbs (w/footrests, 18x18)	39 lbs	43 lbs (Std. Config.) 32 lbs (Transport)	22.3 lbs (Std. Config.) 16.5 lbs (Transport)	21.8 lbs (Std. Config.) 16.1 lbs (Transport)	24.5 lbs (Std. Config.) 13.6 lbs (Transport)	28.8 lbs	53 lbs (w/o footrests)	
	User Weight Capacity(4)	275 lbs	300 lbs	265 lbs	265 lbs	250 lbs(1)	250 lbs(1)	300 lbs (Std) 400 lbs (HD)	265 lbs	300 lbs	300 lbs	300 lbs	300 lbs (Std) 350 lbs (HD)	300 lbs	650 lbs	
Adjustability	Seat Width	12" to 20"	12" to 20"	12" to 20" (1)	12" to 20"	12" to 22" (1)	12" to 22" (1)	14" to 22"	14" to 20"	14" to 22"	13" to 20"	13" to 20"	12" to 22" (1)	14" to 20"	22" to 30"	
	Seat Depth	12" to 20"	12" to 20"	12" to 20" (1)	12" to 20"	12" to 18" (1) (GP) 12" to 22" (1) (GPV)	12" to 22"	14" to 22"	14" to 22"	14" to 22"	13" to 20"	13" to 20"	12" to 20" (1)	14" to 20"	18" to 22"	
	Back Angle Adjustment	Fore (-) 18° to Aft (+) 10°	Fore (-) 18° to Aft (+) 10°	Fore (-) 17° to Aft (+) 8°	Fore (-) 17° to Aft (+) 8°	Fore (-) 8° to Aft (+) 4°	Fore (-) 8° to Aft (+) 4°	Fore (-) 5° to Aft (+) 30°	Fore (-) 5° to Aft (+) 30°	-30° to 5° Recl. back: 0° to 50°	Fore (-) 10° to Aft (+) 15°	Fore (-) 10° to Aft (+) 15°	Fore (-) 3° to Aft (+) 15°	Fore (-) 5° to Aft (+) 14°	N/A	
	Camber	0°, 3°, 6°	0°, 3°, 6°	0°, 3°, 6°	0° to 3°	0° to 11°	0° to 11°	0°	0°	0°	0°	0°, 2°, 4°	0°, 2°, 4°	0°, 2°, 4°	0°, 3°	0°
	Tension Adj. Upholstery	Yes	Yes	Yes	Yes	Yes	Yes	No	No	Yes (Nylon, EXO)	Yes (Nylon, 3DX, EXO)	Yes (Nylon, 3DX, EXO)	Yes	Yes	No	
Clinical Applications	Front Seat-to-Floor Ht.	16" to 21"	16" to 21"	16" to 21"	15" to 20"	15.5" to 23"	15.5" to 23"	12.5" to 19"	14" to 20.5"	12" to 20"	14" to 22"	14" to 22"	13.5" to 22"	12.5" to 19.5"	17" to 20"	
	Rear Seat-to-Floor Ht.	13" to 20"	13" to 20"	13" to 20"	14" to 20"	14.5" to 20.5" (GP) 14.5" to 20.6" (GPV)	14.5" to 20.5"	12.5" to 19"	14" to 20.5"	12" to 20"	14" to 20"	14" to 20"	13" to 22"	12" to 19"	17" to 20"	
	Footrest / Frame Angle	75°, 80°, 85°	75°, 80°, 85°	75°, 80°, 85°	75°, 80°, 85°	70°, 80°, 85°	Swing-Away	60°, 70°, 80°, 90°, Elevating	70°, 80°, Elevating	60°, 70°, 80°, ELR, ALR	60°, 70°, 80°, ELR, ALR	75°, 85° 1" Inset	60°, 70°, 80°, 90° Swing-Away	60°, 70°, 80°, Elevating	Swing-Away	
	Fold-Down Back	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Optional	Yes	Optional, Half-Fold	Optional, Half-Fold	No	No	No	
	Adjustable Back Height	Adjust. 10" to 20" Fixed: 8" to 19"	Adjust. 10" to 20" Fixed: 8" to 19"	Adj. 10" to 20"	Adj. 14" to 20"	Adj. 8" to 20"	Adj. 8" to 20"	Adj. 15" to 21" or 18" to 24"	Adj. 15" to 21" or 18" to 24"	15" to 21" 18" to 24"	Adj. 10" to 20"	Adj. 10" to 20"	9" to 21"	13" to 17" or 16" to 20"	16.75" to 20.13"	
	Center of Gravity (COG)	0" to 4"	0" to 4"	-1.5" to 4"	1" to 3"	0.5" to 3.5"	0.5" to 3.5"	N/A	N/A	N/A	1" to 3"	1" to 3"	1" to 3.5"	1" to 3"	N/A	
	Tilt Range	-	-	-	-	-	-	40° or 55°	45°	-5 to 20° or 0 to 25°	-	-	-	-	-	
Lifestyle Options	Power Assist(4,5)	R90 / Xtender	R90 / Xtender	R90 / Xtender	R90	Xtender	Xtender	R20	R20	No	R90	R90	R90 / R20 / Xtender	R20	R20	
	Front Suspension	Yes	Yes	Yes	Yes	Built-4-Me	Built-4-Me	Built-4-Me	Built-4-Me	Built-4-Me	Yes	Yes	Yes	Built-4-Me	No	
	Rear Suspension	Yes	No	Yes	No	No	No	No	No	No	No	No	No	No	No	
	Transit	Yes(3)	Yes(3)	Yes(3)	Yes(3)	No	No	Yes(3)	Yes	Yes	Yes	Yes	Yes(3)	Yes	No	
	Warranty	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	Frame: Ltd.Lifetime Parts: 1 Year	

(1) Additional options available through QUICKIE Built-4-Me.

(2) Transport Weight = all items that can be removed without tools to depict transporting into or out of a vehicle. Standard Config. = 16" x 16" chair with wheel locks.

(3) Conforms with ANSI/RESNA WC-19 standards. Please refer to QUICKIE & ZIPPY WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at www.SunriseMedical.com.

(4) WARNING: Addition of options can limit the maximum allowable user weight. Refer to order form for options available on each chair.

(5) Check chair compatibility on the R90 product page.



IMPORTANT NOTE: The color options shown are for reference only. They are not intended as accurate representations. Refer to the QUICKIE Color Brochure or a QUICKIE Color Ring for representations of QUICKIE colors.


Kolorfusion™ available on select QUICKIE chairs.




Our Built-4-Me Team is dedicated to creating product customizations beyond what is listed on our product order forms. Experts in the engineering and manufacturing of custom rehab products will listen to your specific requests and build a chair as individual as your fingerprint!

JAY[®] Seating Selection

JAY CUSHIONS




Basic
E2601, E2602



Soft Combi™ P
E2605, E2606



Zip®
E2607



GO™
E2603, E2604



Easy®
E2607, E2608




Ion®
E2607, E2608



Union®
E2607, E2608



X2®
E2622



J2®
E2622, E2623, E2624, E2625



J2® Deep Contour
E2622, E2623, E2624, E2625



Fusion®
E2622, E2623



Fusion® CryoFluid™
E2622, E2623

JAY BACKS



GO™
E2611, E2612



J3® Carbon Fiber
E2613



J3® PA
E2613, E2614



J3® PL
E2615



J3® PD, PDL, PDC
E2620