



ORDER FORM			
Customer Dealer Information		Shipping Address	Identical To Billing Address
Dealer		Client	
Client		Address	
Address			
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Information	
Purchase Order No.:		Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):		User Disability:	
		User Weight:	kg User Height
Build4me			

Serial No:

RGK Octane Sub4

SBAXZ010

TIGA Sub4

\$9,328.00

Max User Weight 110kg (242lbs)

Please Note: All RGK chairs are hand built and have a manufacturing measurement tolerance of +/- 5mm

Product Style

XA070

Open Frame

No Cost

XA080

Closed Frame (Non certified Sub4 option)

No Cost

CAD Approval

XW009

Quote only (No CAD For Approval)

XW010

Order CAD Drawing For Approval
(2 Customer Amendments ONLY)

\$300.00

XW011

Order (No CAD Required)

XW012

Quote & CAD For Approval
(2 Customer Amendments ONLY)

\$550.00

Key Measurements

XY010

Seat Width (Inside of sideguards)

mm (Min - 250mm / Max - 440mm)

XY020

Seat Length

mm (Min - 250mm / Max - 440mm)

XY030

Extra Seat Length (Provides additional leg support and hand grip)

mm (Max - 440mm Total Seat)

No Cost

XY050

Front Frame Height

mm (Min - 380mm / Max - 520mm)

XY060

Rear Frame Height

mm (Min - 360mm / Max - 500mm)

XY070

Frame Length

mm (Min = XY020 + XY030 +140mm / Max - 625mm)

XY072

Castor Position



mm Optimise Measurements

XA020

Ergonomic Seat

mm Specify Length - (Standard = 120mm)

\$125.00




Front Frame

XY080	Seat Taper Design	(When combined with XY091, dimension is in between the centre of the main bends)	No Cost
	mm		
XY091	Parallel Taper Design	(When combined with XY080, dimension must be smaller by 20mm)	No Cost
	mm		


Footrest Measurements

XY100	Round Design	(Castor wheel size 4", max height achievable 110mm)	Std
XY105	Square Design	(Castor wheel size 4", max height achievable 130mm)	No Cost
XY115	Front Footrest Height		
	mm		
XY120	Rear Footrest Height		
	mm	STD = 20mm lower than front height	
XY125	Footrest Length	Maximum Sub4 length = 140mm (certified weight) (160mm and over requires heavy duty upgrade XB150)	
	mm		



Footrest Style

XB015	Tubular Aluminium Footrest Sub4 - Height Adjustable 15.88 x 1.6mm	Std
XB010	Tubular Aluminium Footrest Non Sub4 - Height Adjustable 16 x 3mm	No Cost
XB200	High Level Footrest	\$125.00
	mm Indicate front footrest height	



Footrest Cover

XB045	Lightweight Carbon Cover	(With weight saving cut-out)	Std
	No Cut-out	(Non certified Sub4 Option)	
	Footplate Width Sub4 certified	(STD if no option selected)	
	Footplate Width Standard	(Optimised plate width) (Non certified Sub4 Option)	



Castor Wheels

XC135	Aluminium Centre 4" x 1.25" Silver Hub / Black Tyre (Std) Black Hub / Black Tyre Purple Hub / Black Tyre Green Hub / Black Tyre Red Hub / Black Tyre Blue Hub / Black Tyre	Std
SBAXC155	Forged Carbon Centre 4" x 1.4" 	No Cost
XC140	Aluminium Centre 5" x 1.4" (Silver Hub/Black Tyre) 	\$239.00
XC145	Aluminium Centre 6" x 1.4" (Silver Hub/Black Tyre) 	\$319.00

Castor Forks

XC190	Performance Castor Fork (Std) Note: Colour follows accent	Std
XC192	High Trail Castor Fork Note: Greater distance between contact patch and steering axis reduces flutter at higher speeds. Colour follows accent	No Cost

Castor Pin Width

XY330	Castor Pin Width Optimised Measurements (Minimum width for clear castor rotation and maximum wheelbase length)
-------	--




Wheel Lock Options

SBAXD036	Compact Scissor with Lightweight Fixing	No Cost
XD020	Underseat Scissor 	No Cost
XD050	Push to Lock 	No Cost
SBAXD085	Eagle Wheel Locks 	\$227.00
SBAXD086	Eagle Wheel Locks with L/W Fixing  (Only with fixed centre of gravity)	\$273.00
	Black (Std)	
	Purple	
	Red	
	Blue	
XD010	Omit Wheel Locks	

Sideguards

XE055	Screw on Lightweight Carbon Level with Tyre Individual Above Tyre (Follow Tyre - Max 25mm) mm		No Cost
XE065	Screw on Lightweight Carbon with Fender (Fender centre of tyre) 10mm Above Tyre (Std) Individual Above Tyre (Min 10mm) mm		No Cost

Backrest

XG010	Fixed Backrest (Non Folding)		Std
SBAXG020	Folding Backrest with 2 Position Snap Lock & Angle Adjust		\$737.00
XY230	Backrest Height (Max 400mm) mm		
XY360	Backrest Angle -5 -3 0 +3 +5 Individual (Max -10/+20 Deg) Deg		
XG100	Ergonomic Backrest Vertical Height mm Rearward Position mm		\$150.00
XG110	Inward Taper (Inside Of Backrest Posts - Specify Width Or Per Side) mm		\$150.00
XG120	Outward Flare (Inside Of Backrest Posts - Specify Width Or Per Side) mm		\$150.00

Camber Bar

XY390	Camber Angle 1 2 3		
SBAXH100	Single Post Fixed Centre of Gravity (25-125mm) mm		Std
SBAXH130	Square Back Adjustable Centre of Gravity (25-125mm) mm		\$892.00
SBAXH140	Round Back Adjustable Centre of Gravity (60-125mm) mm		\$892.00

Push Handles

XG070	Screw in Push Handles		\$232.00
XG080	Flip Down Push Handles		\$454.00

Anti-Tip

XJ100	Single Sport Anti-Tip		\$519.00
-------	-----------------------	---	----------

Frame Finish

SBAXN030	Brushed Frame	No Cost
SBAXN050	Polished Frame	No Cost
XN010	Powder Coat Frame 	No Cost
Paint Code		

Accent Colour

Accent colour includes, forks & castor caps

Black
Red
Purple
Silver
Blue

Seating Options

XP045	Sub4 Seat Sling Strap Adjustable	Std
XR050	Carbon Fibre Seat Plate 	\$494.00
mm (Std XY020 - 20mm)		

Sub4 Backrest Upholstery Options

XP035	Airtech Sub4 backrest strap adjustable	Std
-------	--	-----

Non Sub4 Backrest Upholstery Options

XP050	Airtech backrest strap adjustable 	No Cost
XP070	Backrest upgrade to leather corners (Black)  (Only with airtech upholstery XP050)	\$125.00
XP110	Backrest with zip pocket  (Only with airtech upholstery XP050)	\$131.00
XP170	Individual backrest embroidery (Please supply JPEG)  (Only with airtech upholstery XP050)	\$247.00

Backrest Upholstery Stitching Colour

Silver	Red
Black	Orange
Blue	Purple
Individual	

RGK Logo Colour

White (Std)	Red
Black	Orange
Blue	Purple
Individual	

Rear Wheel Size									
XY420 24" Wheels					XY430 25" Wheels				
XY460 Rear Wheel Gap - 15 20 25									
Rear Wheel Style									
XK075	Spinergy BLXL 18 Blade Spoke (N/A in Silver Hub)								Std
SBAXK150	Spinergy CLX 18 Carbon Blade Spoke (N/A in Silver Hub)								\$1,599.00
SBAXK040	Supply Without Wheels	Wheel Type:			Tyre Type			NC	
XK226	Titanium 100mm quick release axle pins								
Spinergy Spoke Colour									
		Black		Red		Purple			
		White		Blue		Orange			
		Yellow		Pink		Green			
XK180	Individual Spoke Design			Inside		Outside		\$219.00	
Spinergy Hub / Rim Colour									
XK190	Black Hub (Std)								Std
XK200	Silver Hub								\$131.00
XK201	Red Hubs & Rims								\$421.00
XK202	Purple Hubs & Rims								\$421.00
XK203	Blue Hubs & Rims								\$421.00
XK204	Gold Hubs & Rims								\$555.00
Handrims Style									
XL010	Aluminium hard anodised 6 tab								Std
XL160	Ellipse 6 Tab								\$650.00
XL020	Tetra Grip 6 Tab								\$412.00
XL100	Para Grip 6 Tab								\$412.00
XL040	Titanium 6 Tab								\$459.00
XL035	Natural Fit LT								\$765.00
XL170	Carbolife Curve L (BLK)								\$1,435.00
XL150	Carbolife Gekko -	Black	Silver		\$1,448.00				
Tyre Style									
XM010	Kenda -	Grey	Blue	Red	(Grey standard if no option selected)				Std
XM060	Schwalbe One (Lightest Option)								\$137.00
XM020	Schwalbe Marathon Plus								\$101.00
XM040	SHOX Solid (26" N/A with Fusion Wheels)								\$165.00
XM130	Kenda Black								\$166.00
XM050	Kenda K885 mountain bike tyre (24",25" only)								\$41.00
Inner Tube									
XY500	Schraeder (Car Valve) - Standard if no option selected								No Cost
XY490	Presta (Sports Valve)								No Cost
Wheel Build Up Options/Accessories									
XY480	Narrow Set Handrim								
XY485	Wide Set Handrim								
XL090	Remove Excess Mounting Tabs								\$150.00
Accessories:									
XK145	Black Gloss Spokeguards			Fitted		Supply Non - Fitted		\$655.00	

Cushion

XP080	Cushion Foam	\$125.00
	Soft (CMS040)	
	Medium (CMS050)	
	50mm	
	75mm	
	100mm	

Cushion Cover

XP085	Nylon	Std
XP090	Airlite	\$101.00
XP100	Leather	\$262.00

Extra Frame Options

XR060	Welded calf strap mounting loops	(Specify Height from Floor)	\$159.00
	Optimise Measurement		



Accessories

XR080	Frontwheel Assembly (Only with Tubular and angle adjustable Footrest) - Select Colour:	\$1,085.00
	Black (Standard)	
	Blue	
	Silver	
	Purple	
	Red	
XP120	Frame protectors Neoprene zip fasten (Black)	\$87.00
XP420	Frame protectors Nylon velcro fasten (Black)	\$83.00
XR010	Seat pouch velcro attach (Black)	\$61.00
XR200	Rucksack with RGK logo (Black)	\$156.00
XR090	Tool Kit	No Cost
XS065	Sub4 Calf Strap	No Cost

Spec Notes



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8