

Hi·Lite

» Flexible, rigid and
lightweight

RGK



Hi Lite - Maximum efficiency and driving comfort

Its 25mm frame tubing combined with its clean frame design provides the Hi Lite with a classic and elegant look.

Available in open and closed frames, its adjustable centre of gravity provides confidence to manoeuvre easily throughout your daily life. **Feel the comfort and efficiency of the rigid frame fully tailored to your body shape and personal needs.**

Key Features

- » TITANIUM rigid frame
- » Elegant look with 25mm frame tubing
- » Efficient and comfortable drive
- » Available in open and closed frame
- » Fix COG (integrated in the frame) or adjustable
- » Integrated fixed backrest or folding back available

- » Ergonomic seating and back optional
- » Wide accessories range
- » Airtech tension adjustable seat and back upholstery
- » Max. user weight 125 kg

Open Frame




Closed Frame



Customised castor fork position (height, width and angle). Low ride fork



Integrated adjustable lateral supports. Additional stability for those users that require more lateral support than provided by the upholstery



Ergonomic seat and back optimises user seated posture improving stability, facilitating propulsion and reducing skin breakdown risk



Airtech tension adjustable upholstery (Std)