

Mark For:

Submitting for: **Quote** **Order**

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ PV: _____ PC: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ PV: _____ PC: _____
 Ship To Phone: () _____ Fax: () _____

Standard Features

- High back captain seat with: sliding, seat, back recline, fully adjustable flip-back armrests (width, angle, depth adjustable)
- 8 Amp off-board Charger
- Deluxe LED light package with: front & rear lights, turn signals, 4-way flashers, automatic brake lights. Energy-efficient, long lasting
- Deluxe extra-large front basket
- Front and rear suspension (Front suspension is independent)
- Maximum Speed of 15km/h (9.3 mph)
- 350 lbs. Weight Capacity
- Large Battery compartment fits Grp22 batteries (batteries not included)
- Easy-adjust delta tiller with auto leveling headlight and dual finger/thumb throttle controls
- 13" wheels w/ pneumatic tires- Left and right rearview mirrors



Base Model

Model	Part Number	MSRP	Qty
Fortress S700 4-Wheel, Luxury Performance Scooter, Carbon Grey Metallic	22300010	\$5,816.00	

Other Options/Accessories

Option/Accessory	Part number	MSRP	Qty
Cane/crutch holder - rear bumper mount	22219095	\$192.00	
Flat-free (foam filled) tires, set of 4	22219102	\$613.00	
Group 22NF BrightWay battery (each) (sealed lead acid)	BW-12550	\$460.00	
Group 24 BrightWay battery 70Ah AGM	BW-12750	\$501.00	
Rear storage bag	22219092	\$207.00	
Safety Flag	SP-000 228 229	\$61.00	