

LECKEY



MyWay

All the support kids need to do it
their way

Supporting Children

I'm ready to go, MyWay!

MyWay is the leader of the pack when it comes to maximizing opportunities for independent movement and giving children the freedom to explore their world. The genius behind this is **MyWay's** revolutionary supportive harness and minimalist open-front, hands-free design.

UNIQUE CORSET-LIKE HARNESS FOR OPTIMAL STEPPING PATTERN

- Provides core trunk stability required to develop head and limb control
- Hug-like effect which increases sensory feedback, dampens involuntary movement and increases confidence and chance of succeeding

OPEN-FRONT DESIGN FOR EXPLORATION AND PARTICIPATION

- Being hands-free, children can easily touch, feel, learn, discover, take part, play—the choice is theirs
- Open design allows therapists to get close to the child and work with their lower limbs while they are in **MyWay**, facilitating guided movement and maximizing the potential for therapeutic benefits



EASIER TRANSFERS

- Direct transfers eliminate the need for additional transfer slings and make the most out of therapy time
- Delivers superior comfort as the harness does not slip into the armpits



Freedom to Explore

Meet Gemma

After Gemma's first experience of MyWay she achieved independent stepping which she hadn't been able to do in any other walking device.

“This is the best piece of equipment I have seen in a long time! Gemma is having fun while exercising”

[Gemma's Mom]



There's always more to MyWay

Multi-functional casters, lightweight frame, prone, gross and fine height adjustment combine to help initiate stepping and control movement.

Frame Size	1	2	3
Age (approx.)	1-5	4-11	9-16
User Shoulder Height (to floor)	25"-36"	34"-47"	42"-59"
User Saddle Height (inside leg)	8.5"-17.5"	16"-24"	21.5"-32"
Max. User Weight	55 lbs.	110 lbs.	176 lbs.
Frame Weight	26.4 lbs. Weight for size 1 includes frame + upper support unit	14.3 lbs.	20.9 lbs.
Upper Support Unit Weight	Size 1 frame & upper support unit do not separate	23 lbs.	24 lbs.
Frame Width	22.5"	26.5"	28.5"
Frame Length	28.5"	33.5"	43.5"

Harness Size	1	2	3	4	5	6
Age (approx.)	1-4	2-6	5-9	8-12	2-15	Up to 16
Waist at Navel Level	18"-22"	19.5"-23.5"	21.5"-25.5"	23.5"-27.5"	27.5"-31.5"	31.5"-35.5"
Max. User Weight	110 lbs.	110 lbs.	110 lbs.	176 lbs.	176 lbs.	176 lbs.
Color	Purple	Green	Blue	Red	Orange	Grey