



Serial No:



ORDER FORM

Customer Dealer Information	Shipping Address	Identical To Billing Address
Dealer	Client	
Client	Address	
Address		
E-Mail:	E-Mail:	
Phone:	Phone:	
Order Information	User Information	
Purchase Order No.:	Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):	User Disability:	
	User Weight: kg	User Height
Build4me		

Serial No:

Product Style

HLTXZ010	HILITE		\$6,855.00
XA080	Closed design		
XA090	Open design	Closed Design	Open Design



Max User Weight 125kg (275lbs)

Please Note: All RGK chairs are hand built and have a manufacturing measurement tolerance of +/- 5mm

CAD Approval

XW009	Quote only (No CAD For Approval)	
XW010	Order CAD Drawing For Approval <i>(Includes 2 CAD Amendments)</i>	\$300.00
XW011	Order (No CAD Required)	
XW012	Quote & CAD For Approval <i>(Includes 2 CAD Amendments)</i>	\$550.00

XY010	Seat Width (Inside of sideguards)	
	mm (Max - 500mm)	
XY020	Seat Length	
	mm (Max - 500mm)	
XY030	Extra Seat Length (Provides additional leg support and hand grip)	
	mm	
XY050	Front Frame Height	
	mm (Min - 380mm / Max - 550mm)	
XY060	Rear Frame Height	
	mm (Min - 360mm / Max - 550mm)	
XY070	Frame Length	
	mm (Min = XY020 + XY030 +140mm)	
XY072	Castor Position	
	mm Optimise Measurements TX forward castor position	
XA020	Ergonomic Seat	\$125.00
	mm Specify Length - (Standard = 120mm)	



Front Frame

XY080	Seat Taper Design	XY091	Parallel Taper Design
	Specify width		Specify width
	mm		mm
	(When combined with XY091, dimension is inbetween the centre of the main bends)		(When combined with XY080, dimensions must be smaller by 20mm)

Footrest Design

XY100	Round Design (Castor wheel size 4", max height achievable 110mm)	No Cost
XY105	Square Design (Castor wheel size 4", max height achievable 110mm)	No Cost
XB200	High Level Footrest (Only with XB010, XB080 & XB150)	\$207.00

Footrest Individual Measurements

XY115	Front Footrest Height (Minimum 50mm) mm
XY120	Rear Footrest Height (Minimum 50mm) mm STD = 20mm lower than front height
XY125	Footrest Length (Std = 100mm) mm (160mm and over requires heavy duty upgrade XB150)

Footrest Style and Cover

XB010	Tubular aluminium footrest (Height adjustable)	No Cost
XB020	Tubular Titanium footrest (Height adjustable)	No Cost
XB030	Flip up plastic footrest (Angle/height adjustable, with closed frame only)	\$363.00
XB070	Angle adjustable aluminium footrest (Angle/height adjustable)	\$213.00
XB150	Heavy Duty Tubular Footrest (Height adjustable, for Freewheel use)	\$185.00
HLTXB080	Tubular footrest welded into frame (No height adjustment)	\$268.00

Footrest Cover

XB050	Carbon fibre cover (Not with XB070 & XB030)	\$255.00
XB170	Aluminium cover (Not with XB070 & XB030)	\$101.00
XB040	Plastic cover (Not with XB070 & XB030)	\$101.00
XB165	Rubber grip tape (3M safety walk) - (not with XB030)	\$89.00

Castor Wheels

XC060	Soft roll 4" x 1.3" (Black only)	No Cost
XC110	Aluminium centre 4" x 1" (Silver hub and black urethane tyre only)	\$239.00
XC130	Aluminium centre 4" x 1.4" (Indicate Colour) Silver hub & black tyre Black hub & black tyre Red hub & black tyre Gold Hub & black tyre Green hub & black tyre Blue hub & black tyre Purple hub & black tyre	\$239.00
XC070	Soft roll 5" x 1.3" (Black only)	No Cost
XC100	Pneumatic Aluminium Centre 5" x 1" (Grey only)	\$255.00
XC120	Aluminium Centre 5" x 1" (Silver hub and black urethane tyre only)	\$239.00
XC140	Aluminium Centre 5" x 1.4" (Silver hub and black urethane tyre only)	\$239.00
XC080	Soft roll 6" x 1.3" (Black hub / grey tyre only)	No Cost
XC145	Aluminium centre 6" x 1.4" (Silver hub and black urethane tyre only)	\$319.00

Castor Forks

XC190	Performance Castor Fork (Std) Note: Colour follows accent	Std
XC192	High Trail Castor Fork Note: Greater distance between contact patch and steering axis reduces flutter at higher speeds. Colour follows accent	No Cost
XC210	Froglegs Suspension Castor Fork - Black - Silver	\$652.00

Castor Pin Width

XY330	Castor Pin Width
	mm Optimise measurements

Wheel Lock Options

XD020	Underseat Scissor	No Cost
XD050	Push To Lock	No Cost
XD030	Compact Scissor	\$274.00
XD046	Out front push to lock	\$421.00
XD035	Compact scissor with lightweight fixing (Only with fixed centre of gravity)	\$363.00
XD085	Eagle wheel locks (Indicate colour)	\$482.00
XD086	Eagle WL with L/W fixing (Indicate colour) (Only with fixed centre of gravity)	\$523.00
	Black Purple Red Blue	
XD010	Omit wheel locks	

Sideguards

Aluminium options

XE030	Screw on aluminium	\$337.00
XE080	Lift out aluminium	\$548.00
HLTXE130	Fold flat aluminium	\$882.00
XE180	Flip back aluminium	\$405.00
XE040	Screw on aluminium with fender	\$287.00
XE090	Lift out aluminium with fender	\$596.00
HLTXE140	Fold flat aluminium with fender	\$760.00
XE190	Flip back aluminium with fender	\$344.00

Carbon Fibre options

XE050	Screw on carbon fibre	\$748.00
XE100	Lift out carbon fibre	\$789.00
HLTXE150	Fold flat carbon fibre	\$1,121.00
XE200	Flip back carbon fibre	\$618.00
XE060	Screw on carbon fibre with fender	\$981.00
XE110	Lift out carbon fibre with fender	\$1,166.00
HLTXE160	Fold flat carbon fibre with fender	\$1,362.00
XE210	Flip back carbon fibre with fender	\$768.00

ABS Plastic options

XE070	Screw on ABS plastic	\$172.00
XE120	Lift out ABS plastic	\$502.00
HLTXE170	Fold flat ABS plastic	\$690.00

Titanium options

HLTXE220	Welded with 6mm bead	\$771.00
----------	----------------------	----------

Sideguard measurements

Optimised Measurements

25mm above tyre (Follow tyre)	Individual above tyre (Follow tyre)
Level with tyre (Follow tyre)	mm
Sideguard length mm	Sideguard height mm

Sideguard measurements

Fender Measurements

10mm above tyre (Standard)	Individual above tyre (Min = 10mm)
Centre of tyre (Standard if no option selected)	mm
Outside of tyre	

Armrest options

XF010	Lift out armrest with plastic sideguard (Indicate fixed height)	\$610.00
XF020	Lift out armrest with aluminium sideguard (Indicate fixed height)	\$481.00

Height (Top of seat rail to top of armrest) mm

Backrest Height

XY230	Fixed height backrest (Maximum 400mm)	XY245	Adjustable height backrest (Maximum 400mm) (+/- 25MM)	\$255.00
	Specify height mm		Specify height mm	

Backrest Design

XG010	Fixed backrest (Non-folding)	HLTXG020	Folding backrest with 2 position snap lock (Black)	\$761.00
-------	---------------------------------	----------	---	----------

Ergonomic Options

XY360	Backrest angle -5 -3 0 +3 +5 Individual Deg (Max -10/+10 Deg)	XG100	Ergonomic backrest Vertical height mm Rearward position mm	\$150.00
XG110	Inward taper (Inside of backrest posts) Specify width or per side mm			\$150.00
XG120	Outward flare (Inside of backrest posts) Specify width or per side mm			\$150.00

Push Handle Options

XG070	Screw in push handles	\$232.00
XG080	Flip down push handles (N/A with XY245)	\$454.00
XG085	Flip down push handles height adjustable (N/A with XY245)	\$742.00
XG097	Integrated push handle (Quickie design) (N/A with XY245)	\$696.00
XG090	Stroller push handles height adjustable Height from floor to top of handles mm Optimise (1000mm from floor)	\$454.00

Backrest Options

XG040	Heavy duty backrest upgrade (Recommended for hand cycling) (Not with XY245)	\$213.00
XG130	Extra deep square brace bar (Recommended for Jay backs) (STD-90mm)	\$118.00

ILSA Backrest Options

XG140	Integrated lateral support	\$493.00
	A-Left mm A-Right mm	
	B-Left mm B-Right mm	
	C-Left mm C-Right mm	
	D-Left mm D-Right mm	

Camber Bar

Adjustable Centre Of Gravity

XH110	Square back adjustable centre of gravity	No Cost
	(25-150mm)	

(Max seat width 500mm)

Specify measurement

mm

XH120	Round back adjustable centre of gravity	No Cost
	(60-150mm)	

Specify measurement

mm

Fixed Centre Of Gravity

HLTXH090	V Style Fixed centre of gravity	\$675.00
	(25-150mm)	
	Specify measurement	

mm

Camber Angle

XY390	Camber angle
	0 1 2 3 5

Camber Bar Options

XH040	Handbike camber bar fixed	\$412.00
-------	---------------------------	----------

Camber angle (Handbike camber bar)

0 5

Castor height from floor for handbike fitting
(Standard 30mm if not included)

30 50

Anti-Tip Options

XJ010	Swingaway anti tip LEFT side	\$172.00
	(Only with square back frame. Not with XH040)	

XJ020	Swingaway anti tip RIGHT side	\$172.00
	(Only with square back frame. Not with XH040)	

XJ100	Single sport anti tip	\$519.00
	(Centre mounted as standard)	

Rear Wheel Size								
XY420	24" Wheels		XY430	25" Wheels		XY440	26" Wheels	
Rear Wheel Gap								
XY460	Rear Wheel Gap							
	15	20	25	Individual	mm			
Rear Wheel Options								
XK010	Fusion 16 Spoke (Silver hub, black spokes)						Std	
XK320	Spinergy XLX 12 Spoke, cross/radial spoked						\$1,099.00	
XK050	Spinergy LX 12 Spoke						\$1,030.00	
XK060	Spinergy Spox 18 Spoke						\$685.00	
XK150	Spinergy CLX 18 Carbon blade spoke (N/A in silver hub)						\$3,233.00	
XK070	Spinergy BLXL 18 blade spoke (N/A in silver hub)						\$1,099.00	
XK040	Supply without Wheels		Wheel Type		Tyre Type	NC		
Spinergy Spoke Colour options								
Black spokes		White spokes		Yellow spokes				
Red spokes		Blue spokes		Pink spokes				
Purple spokes		Orange spokes		Green spokes				
XK180	Individual spoke design (Specify colour)			Inside	Outside	\$219.00		
Spinergy Hub Colour options								
XK190	Black Hub (Std option)						Std	
XK200	Silver Hub (n/a with BLXL/CLX)						\$131.00	
XK201	Red Hub/Rims (LX, XLX & BLXL only)						\$421.00	
XK202	Purple Hub/Rims (LX, XLX & BLXL only)						\$421.00	
XK203	Blue Hub/Rims (LX, XLX & BLXL only)						\$421.00	
XK204	Gold Hub/Rims (LX, XLX & BLXL only)						\$555.00	
Offroad Wheels								
XK300	Offroad wheel set 22" with handrim, MTB tyre & axle pins (Replaces 24" & 25" wheel with XM010)						\$961.00	
XK310	Offroad wheel set 24" with handrim, MTB tyre & axle pins (Replaces 26" wheel with XM010)						\$961.00	

Handrim Position / Options

XY480	Narrow set handrim	No Cost
XY485	Wide set handrim	No Cost
XL090	Remove excess mounting tabs (Recommended for small wheel gaps)	\$150.00

Handrim Styles

XL010	Aluminium hard anodised 6 tab (Standard if no option selected)	Std
XL040	Titanium 6 tab	\$459.00
XL160	Elipse 6 tab (replacement for surge)	\$650.00
XL035	Natural Fit LT	\$765.00
XL020	Tetra grip 6 tab	\$412.00
XL150	Carbolife Gekko BLACK SILVER	\$1,448.00
XL100	Para grip 6 tab	\$412.00
XL170	Carbolife Curve L (Black)	\$1,435.00

Tyre Options

XM010	Kenda Grey Blue Red (Grey standard if no option selected)	Std
XM060	Schwalbe one (24", 25" only)	\$137.00
XM130	Kenda black	\$29.00
XM050	Kenda K885 mountain bike tyre (24", 25" only)	\$41.00
XM020	Schwalbe marathon plus	\$101.00
XM040	Shox solid airless (Black) (26" N/A with Fusion Wheels)	\$165.00

Inner Tube Valve Options

XY500	Schraeder inner tube valve (Standard if no option selected)	No Cost
XY490	Presta inner tube valve	No Cost

Wheel Accessories

XK220	Titanium 130mm quick release axle pins (Pair) (N/A with XK010 & XK070)	\$165.00
XK240	Quad release axle nuts (Pair)	\$248.00
XK145	Black gloss spokeguards Fitted Supply non - fitted	\$655.00

Finish Options

XN010 Powder coated frame

No Cost

Frame

Gloss

Black gloss 103
White gloss 100
Grey gloss 132
Cream gloss 135
Red gloss 108
Yellow gloss 120
Forest green gloss 143
Green gloss 107
Dark blue gloss 151
Royal blue gloss 154
Blue gloss 174
Duckegg Blue gloss 141
Pink gloss 188
Orange gloss 180
Purple gloss 157

Frame

Matt

Black matt 105
White matt 101
Grey matt 131
Light grey matt 133
Navy blue matt 152

Frame

Special

Chrome mirror 300
Sparkle silver 170
Metallic black matt 104
Metallic grey matt 130
Black sandpaper 201
Grey sandpaper 202
White sandpaper 200

XN150

Powder coated components

No Cost

(Frame colour if not specified)

Specify component colour code from above list

HLTXN030

Bead Blast

(All parts black 105 unless specified)

\$582.00

XN165

Bead Blast components

\$258.00

HLTXN050

Polished Frame

(All parts black 103 unless specified)

\$1,350.00

XN170

Polished Components

\$219.00

Accent Colour

Note: Accent colour includes forks, clamps & castor caps.

Black accent

Silver accent

Red accent

Blue accent

Purple accent

Upholstery

Seat upholstery

XP040 Seat sling strap adjustable (Std)

XR050 Carbon fibre seat plate

Specify length of plate (mm) (Std XY020 - 20mm)

\$494.00

Backrest upholstery

XP050 Airtech backrest strap adjustable

XP070 Backrest upgrade to leather corners (Black) \$125.00

Backrest stitching colour

Silver

Black

Blue

Individual:

Red

Orange

Purple

RGK Logo colour

White (std)

Green

Blue

Individual:

Red

Orange

Purple

Model Logo colour

White (std)

Green

Blue

Individual:

Red

Orange

Purple

Backrest Accessories

XP110

Backrest with zip pocket (N/A with sub4 backrest)

\$131.00

XP170

Individual backrest embroidery (Please provide JPEG)

\$247.00

Cushion

XP080	Cushion foam (Select density and size below)				\$125.00
	Soft (CMS040)	Medium (CMS050)			
	50mm	75mm	100mm		

Cushion cover upgrades

XP085	Nylon cushion cover (Black)	Std
XP090	Air lite cushion cover (Black)	\$101.00
XP100	Leather cushion cover (Black)	\$262.00

Frame Options

XR060	Welded calf strap mounting loops (To locate calf strap on front frame) Specify height from floor	\$159.00
	Optimise measurements	

Accessories

XR080	Frontwheel assembly - Specify colour: (Only with tubular and angle adjustable style footrest)	Black	Blue		\$1,085.00
		Silver	Purple	Red	
XP120	Frame protectors neoprene zip fasten (Black)				\$87.00
XP420	Frame protectors nylon velcro fasten (Black)				\$83.00
XR010	Seat pouch velcro attach (Black)				\$61.00
XR200	Rucksack with RGK logo (Black)				\$156.00
XR090	Tool kit				Std
XS060	Calf strap				Std



1000 Creditstone Road, Unit 2, Concord, ON L4K4 4P8