

**GEMINO 20  
GEMINO 20 M  
GEMINO 20 S**

February 2024



Product sold with backrest in Canada.


**ORDER FORM**
**Invoice address or customer number:**

 Customer number:  
 Name:  
 Street:  
 Town: Postal code:  
 Country:  
 E-mail:  
 Tel:

**Delivery address:**  Deliver to Invoice address

 Customer number:  
 Name:  
 Street:  
 Town: Postal code:  
 Country:  
 E-mail:  
 Tel:

**Order information**
 Order  Quote  
 Order date:  
 Your reference:  
 Reimbursement code:

**Customer Service: [cscanada@sunmed.com](mailto:cscanada@sunmed.com)**
**Base price \$633.00**
**DIMENSIONS**


|                                     | Gemino 20   | Gemino 20 M | Gemino 20 S |
|-------------------------------------|-------------|-------------|-------------|
| <b>A</b> Total height               | 78-100 cm   | 69-88 cm    | 65-77 cm    |
| <b>B</b> Total length               | 65 cm       | 65 cm       | 68 cm       |
| <b>C</b> Total width                | 60 cm       | 60 cm       | 55 cm       |
| <b>D</b> Width between push handles | 47 cm       | 47 cm       | 42 cm       |
| <b>E</b> Seat height                | 62 cm       | 55 cm       | 48 cm       |
| Height folded                       | 80 cm       | 76 cm       | 69 cm       |
| Length folded                       | 65 cm       | 65 cm       | 68 cm       |
| Width folded                        | 23 cm       | 23 cm       | 23 cm       |
| Turning circle                      | 84 cm       | 84 cm       | 79 cm       |
| <b>F</b> Front / rear wheels        | 20 x 3.6 cm | 20 x 3.6 cm | 20 x 3.6 cm |
| Total weight without basket         | 7.1 kg      | 7 kg        | 6.6 kg      |
| <b>G</b> Weight basket, standard    | 0.3 kg      | 0.3 kg      | 0.25 kg     |
| Max. user weight                    | 150 kg      | 130 kg      | 125 kg      |
| Max. weight in basket               | 5 kg        | 5 kg        | 5 kg        |
| Recommended user height             | 150-200 cm  | 135-170 cm  | 125-165 cm  |

ORDER FORM GEMINO 20 / GEMINO 20 M / GEMINO 20 S

| SIZES AND COLOURS                          |                   | Price \$ | Quantity | Total |
|--|-------------------|----------|----------|-------|
| Gemino 20 + HC-9012222 Adj. Comfort Back   |                   |          |          |       |
| <input type="checkbox"/> HC-7120000        | Gemino 20, grey   | 633.00   |          |       |
| Gemino 20 M + HC-9012222 Adj. Comfort Back |                   |          |          |       |
| <input type="checkbox"/> HC-7120200        | Gemino 20 M, grey | 633.00   |          |       |
| Gemino 20 S + HC-9012222 Adj. Comfort Back |                   |          |          |       |
| <input type="checkbox"/> HC-7120100        | Gemino 20 S, grey | 633.00   |          |       |

| ACCESSORIES   |   | Price \$ | Quantity | Total |
|---|---|----------|----------|-------|
| Bag for Gemino 20 S. <i>Not available for the 20 and 20M.</i> |   |          |          |       |
| <input type="checkbox"/> HC-7104007                           | Bag, advanced with lid and handles, black/black     | 117.00   |          |       |
| Other accessories   |   |          |          |       |
| <input type="checkbox"/> HC-7109011                           | Crutch holder/ Cane holder                          | 97.00    |          |       |
| <input type="checkbox"/> HC-7109140                           | Table tray  | 97.00    |          |       |
| <input type="checkbox"/> HC-7109060                           | Umbrella, black                                     | 240.00   |          |       |
| <input type="checkbox"/> HC-7109050                           | Basket for oxygen cylinder - Gemino 20, Gemino 20 M | 240.00   |          |       |
| <input type="checkbox"/> HC-7109051                           | Basket for oxygen cylinder - Gemino 20 S            | 240.00   |          |       |
| <input type="checkbox"/> HC-7109080                           | Slowdown brake, set                                 | 230.00   |          |       |
| <input type="checkbox"/> HC-7109095                           | Light   | 117.00   |          |       |
| <input type="checkbox"/> HC-7109100                           | Infusion holder                                     | 240.00   |          |       |

| GEMINO 20 / GEMINO 20 M / GEMINO 20 S |  | Price \$ |
|---------------------------------------|--|----------|
| Base price                            |  |          |
| Option price                          |  |          |
| Total price                           |  |          |



Sunrise Medical  
 237 Romina Drive, Unit 3  
 Concord, ON  
 L4K 4V3, Canada  
 Tel. +1.800.263.3390  
 Fax: +1.800.561.5834  
 E-mail: [cscanada@sunmed.com](mailto:cscanada@sunmed.com)  
[www.sunrisemedical.ca](http://www.sunrisemedical.ca)

# GEMINO ACCESSORIES



## BAG -GEMINO 20 S



7104007 Bag, advanced, black/black **\$117.00**

## OTHER ACCESSORIES



// 7109011 **\$97.00**  
Crutch holder



// 7109140 **\$97.00**  
Table tray



// 7109095 **\$117.00**  
Light



// 7109060 **\$240.00**  
Umbrella, black



// 7109050 & 7109051 **\$240.00**  
Basket for oxygen cylinder\*



// 7109100 **\$240.00**  
Infusion holder