

GEMINO 20
GEMINO 20 M
GEMINO 20 S

April 2025



Product sold with backrest in Canada.


ORDER FORM
Invoice address or customer number:

 Customer number:
 Name:
 Street:
 Town: Postal code:
 Country:
 E-mail:
 Tel:

Delivery address: Deliver to Invoice address

 Customer number:
 Name:
 Street:
 Town: Postal code:
 Country:
 E-mail:
 Tel:

Order information
 Order Quote
 Order date:
 Your reference:
 Reimbursement code:

 Customer Service: cscanada@sunmed.com
Base price \$642.00
DIMENSIONS


	Gemino 20	Gemino 20 M	Gemino 20 S
A Total height	78-100 cm	69-88 cm	65-77 cm
B Total length	65 cm	65 cm	68 cm
C Total width	60 cm	60 cm	55 cm
D Width between push handles	47 cm	47 cm	42 cm
E Seat height	62 cm	55 cm	48 cm
Height folded	80 cm	76 cm	69 cm
Length folded	65 cm	65 cm	68 cm
Width folded	23 cm	23 cm	23 cm
Turning circle	84 cm	84 cm	79 cm
F Front / rear wheels	20 x 3.6 cm	20 x 3.6 cm	20 x 3.6 cm
Total weight without basket	7.1 kg	7 kg	6.6 kg
G Weight basket, standard	0.3 kg	0.3 kg	0.25 kg
Max. user weight	150 kg	130 kg	125 kg
Max. weight in basket	5 kg	5 kg	5 kg
Recommended user height	150-200 cm	135-170 cm	125-165 cm

ORDER FORM GEMINO 20 / GEMINO 20 M / GEMINO 20 S

SIZES AND COLOURS		Price \$	Quantity	Total
Gemino 20 + HC-9012222 Adj. Comfort Back				
<input type="checkbox"/> HC-7120000	Gemino 20, grey	642.00		
Gemino 20 M + HC-9012222 Adj. Comfort Back				
<input type="checkbox"/> HC-7120200	Gemino 20 M, grey	642.00		
Gemino 20 S + HC-9012222 Adj. Comfort Back				
<input type="checkbox"/> HC-7120100	Gemino 20 S, grey	642.00		

ACCESSORIES		Price \$	Quantity	Total
Bag for Gemino 20 S. <i>Not available for the 20 and 20M.</i>				
<input type="checkbox"/> HC-7104007	Bag, advanced with lid and handles, black/black	119.00		
Other accessories				
<input type="checkbox"/> HC-7109011	Crutch holder/ Cane holder	98.00		
<input type="checkbox"/> HC-7109140	Table tray	98.00		
<input type="checkbox"/> HC-7109060	Umbrella, black	244.00		
<input type="checkbox"/> HC-7109050	Basket for oxygen cylinder - Gemino 20, Gemino 20 M	244.00		
<input type="checkbox"/> HC-7109051	Basket for oxygen cylinder - Gemino 20 S	244.00		
<input type="checkbox"/> HC-7109080	Slowdown brake, set	233.00		
<input type="checkbox"/> HC-7109095	Light	119.00		
<input type="checkbox"/> HC-7109100	Infusion holder	244.00		

GEMINO 20 / GEMINO 20 M / GEMINO 20 S		Price \$
Base price		
Option price		
Total price		

GEMINO ACCESSORIES



BAG -GEMINO 20 S



7104007 Bag, advanced, black/black **\$119.00**

OTHER ACCESSORIES



// 7109011 **\$98.00**
Crutch holder



// 7109140 **\$98.00**
Table tray



// 7109095 **\$119.00**
Light



// 7109060 **\$244.00**
Umbrella, black



// 7109050 & 7109051 **\$244.00**
Basket for oxygen cylinder*



// 7109100 **\$244.00**
Infusion holder



Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2 Concord,
ON
L4K 4P8, Canada
Tel. +1.800.263.3390
Fax: +1.800.561.5834
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca